DE

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Abbreviated instructions for use

Dear customer,

Thank you for choosing a quality product from Dentaurum.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaurum.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurum GmbH & Co. KG | Turnstr. 31 | 75228 Ispringen | Germany

2. General device description, intended purpose and characteristics, technical description

The tioLogic® TWINFIT implant system contains specially coordinated instruments, abutments, and accessory components for the prosthetic restoration of the tioLogic® TWINFIT implants. Only original components of the tioLogic® TWINFIT implant system should be combined in accordance with the Instructions for use.

The abutment is selected by the dentist/dental technician to suit the individual patient's needs.

The diameter of the implant and the gingival thickness determine the choice of abutment.

Selection aids with conical and platform connector geometry for the abutment series S, M and L are available to help choose the right abutment on the model (gingival height and angulation).

3. Area of application

Ball and tioLOC abutments tioLogic® TWINFIT are prosthetic components for implant-retained, gingiva-supported restorations on tioLogic® TWINFIT implants.

4. Medical indications/contraindications

4.1 Ball abutments tioLogic® TWINFIT

Indications

- Implant-retained, gingiva-supported coverdenture
- Incorporation into an existing coverdenture as temporary restoration or into a partial denture framework is possible
- Prosthetic restoration on at least two implants

Contraindications.

- Uneven number of implants per jaw
- Unfavorable positioning of implants which prevent a tangential rotation axis
- Restorations with mixed retention (tooth/implant)
- In combination with other retention elements
- Divergence of more than 20° between the implants
- Heavy mechanical strain with ø 3.3 mm implants

4.2 tioLOC abutments tioLogic® TWINFIT

Indications

- Implant-retained, gingiva-supported coverdenture with low profile
- Incorporation into an existing coverdenture is possible
- Prosthetic restoration on two or four implant posts

Contraindications

- Uneven number of implants per jaw
- Unfavorable positioning of implants which prevent a tangential rotation axis
- Restorations with mixed retention (tooth/implant)
- In combination with other retention elements
- Implant divergences of more than 10° when using tioLogic® TWINFIT implants ø 3.3 mm are contraindicated
- Intermaxillary space should be at least 4.0 mm

5. Information on use, processing, fitting and insertion

5.1 Ball abutment tioLogic® TWINFIT

Ball abutments with platform connector geometry are available for the S, M and L series of abutments. The gingival heights are 1.5 mm, 3.0 mm and 4.5 mm; this refers to the distance between the upper edge of the implant and the lower edge of the socket for the PentaGrip insertion key tioLogic® TWINFIT ISO shank of the ball abutment. This lower edge should be approx. 1.0 mm above the gingival line. To achieve optimum retention, all ball abutments should be positioned at the same level. The ball is 2.25 mm in diameter.

The components of the ball abutment may not be modified, for example, by grinding. Only the matrix withdrawal force can be adjusted from 2 N up to 12 N by activating the inner matrix. The ball abutment is placed with the PentaGrip insertion key tioLogic® TWINFIT ISO shank.

Tightening torque

- Ball abutment on a model: manually
- Ball abutment intra-orally: 35 Ncm

5.2 tioLOC abutment tioLogic® TWINFIT

tioLOC abutments with platform connector geometry are available for the S, M and L series of abutments. The gingival heights are 1.5 mm, 3.0 mm, and 4.5 mm, referring to the cylindrical area of the tioLOC abutment up to the lower edge of the retention element. The upper edge of the tioLOC abutments should be approx. 1.0 mm above the gingival line. The head with retainer is identical for all tioLOC abutments. tioLOC abutment components may not be modified, for example, by grinding.

- Implant divergences of more than 10° when using tioLogic® TWINFIT implants ø 3.3 mm are contraindicated
- Intermaxillary space should be at least 4.0 mm

For more information on processing and on the matrix withdrawal forces for tioLOC abutments TWINFIT, see www.dentaurum.com (Prosthetics Manual tioLogic® TWINFIT REF 989-913-20).

The tioLOC abutment is inserted with the insertion key tioLOC.

Tightening torque

- tioLOC abutment on the model: manually
- tioLOC abutment intra-orally: 30 Ncm

6. Composition

- Ball abutment: titanium alloy (Ti6Al4V)
- Laboratory implant: aluminum
- tioLOC abutment: titanium alloy (Ti6Al4V)
- 4Base tioLOC cap: (Ti6Al4V)
- PentaGrip insertion key - ISO shank 4Base/ball abutment: 1.4197
- Insertion key tioLOC abutment: 1.4197

Information about the composition is included in the material list for our implantology products (REF 989-801-05), or see www.dentaurum.com. The SSCP is available at <https://ec.europa.eu/tools/eudamed>.

7. Safety information

- The product should not be used if there is a known allergic reaction to one or more of the material components.

- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. All articles used intra-orally should therefore be secured against swallowing and/or aspiration.

- All serious incidents arising from the use of the product should be reported to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.

- Different types of alloy in the oral cavity can lead to galvanic reactions.

- If the packaging shows serious damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.

- Only use tioLogic® TWINFIT prosthetic components in combination with tioLogic® TWINFIT implants.

8. Storage and shelf life

Store in a dry place.

9. Information for single use products

Ball and tioLOC abutments are designed for single use only. Reconditioning of a ball or tioLOC abutment (recycling) or reuse on patients are not permitted since it can then no longer be guaranteed that the article can be reprocessed safely or can function safely.

10. Information for multiple use products

The PentaGrip insertion key and insertion key tioLOC abutment are intended for multiple use. The products must be disinfected, cleaned and sterilized before being used again. Instruments and accessories must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels. Additional information can be found at www.dentaurum.com (Processing Instructions Prosthetic Components, Instruments and Accessories REF 989-801-09).

11. Information on disinfection/cleaning/sterilizing

tioLogic® TWINFIT prosthetic components are delivered non-sterile. They are for single use on one patient only. They must be cleaned and disinfected before being used on a patient. Dentaurum generally recommends that the product be sterilized in addition. Only sterilized prosthetic components may be used if bleeding occurs.

Additional information can be found at www.dentaurum.com (Processing Instructions Prosthetic Components, Instruments and Accessories REF 989-801-09).

12. Product overview/scope of delivery, incl. REF

Please refer to the tioLogic® TWINFIT product catalog for the tioLogic® TWINFIT product range.

13. Quality information

Dentaurum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaurum has no influence on the processing on site.

14. Explanation of symbols

 Please refer to the label. Additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00)

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