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## Abbreviated instructions for use

**Dear customer,**

Thank you for choosing a quality product from Dentaurum.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at [www.dentaurum.com](#), even though you may frequently use the same product.

#### 1. Manufacturer

Dentaurum GmbH & Co. KG | Turnstr. 31 | 75228 Ispringen | Germany

#### 2. General device description, intended purpose and characteristics, technical description

The tiologic® TWINFIT implant system contains specially coordinated instruments, abutments, and accessory components for the prosthetic restoration of the tiologic® TWINFIT implants. Only original components of the tiologic® TWINFIT implant

system should be combined in accordance with the Instructions for use.

The abutment is selected by the dentist/dental technician to suit the individual patient's needs.

The diameter of the implant and the gingival thickness determine the choice of abutment. Selection aids with conical and platform connector geometry for the abutment series S, M and L are available to help choose the right abutment on the model (gingival height and angulation).

#### 3. Area of application

The titanium abutments tiologic® TWINFIT are prosthetic components for single tooth and bridge restorations on tiologic® TWINFIT implants.

#### 4. Medical indications/contraindications

**Titanium abutments tiologic® TWINFIT**

Indications

- Individual crown and bridgework (cemented, individually screw-retained)
- Abutment on implant Ø 3.3 mm only in the area of the upper lateral incisors and lower incisors

Contraindications.

- Primary bracing of abutments
- Correction of angulation exceeding 20°, or exceeding 15° in the case of implants Ø 3.3 mm
- Cast-on technique

#### 5. Information on use, processing, fitting and insertion

**Titanium abutment tiologic® TWINFIT**

All titanium abutments tiologic® TWINFIT are available in conical and platform connector geometry for the abutment series S, M and L.

The straight and angled (S 15°, M 20°, L 20°) titanium abutments can be slightly altered to adapt to the gingival line (straight 1.5 mm, 2.5 mm, 4.0 mm/angled 1.5 mm, 3.0 mm). These abutments have an exactly defined crown margin and an integrated rotational stop.

The angled titanium abutments also have different gingival heights (labial/ palatal). The universal titanium abutments are anatomically shaped and can be prepared as required for the individual case.

For digital manufacture of individual single tooth and bridge restorations, CAD/CAM titanium blocks (PreForms) are available for one-piece abutments and CAD/CAM titanium bases for hybrid abutments.

The prosthetic screw (not included in delivery) serves only to help the dental technician fix the abutment on the working cast. A new AnoTite screw should be used for the final fitting in the mouth.

**Tightening torque**

- Prosthetic screw, on model: manually
- AnoTite screw intra-orally: 30 Ncm

**Attention: Do not confuse the prosthetic screw with the AnoTite screw. The prosthetic screw is colored silver and the AnoTite screw is anodized gold.**

#### 6. Composition

- Titanium abutments: titanium alloy (Ti6Al4V)

- Prosthetic screws:
  - AnoTite screw: titanium alloy (Ti6Al4V)
  - Prosthetic screw: titanium alloy (Ti6Al4V)

- Hex keys:
  - Hexagon key - ISO shank 1.3: 1.4197
  - Hexagon key - ratchet 1.3: 1.4197

- Hexagonal ball end key - ISO shank 1.3: 1.4305

Information about the composition is included in the material list for our implantology products (REF 989-801-05), or see [www.dentaurum.com](#). The SSCP is available at <https://ec.europa.eu/tools/eudamed>.

#### 7. Safety information

- The product should not be used if there is a known allergic reaction to one or more of the material components.

- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. All articles used intra-orally should therefore be secured against swallowing and/or aspiration.

- All serious incidents arising from the use of the product should be reported to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.

- Different types of alloy in the oral cavity can lead to galvanic reactions.

- If the packaging shows serious damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.

- Only use tiologic® TWINFIT prosthetic components in combination with tiologic® TWINFIT implants.

#### 8. Storage and shelf life

Store in a dry place.

#### 9. Information for single use products

Titanium abutments are designed for single use only. Reconditioning of a titanium abutment (recycling) or reuse on patients are not permitted since it can then no longer be guaranteed that the article can be reprocessed safely or can function safely.

#### 10. Information for multiple use products

The hex key 1.3 and the hexagonal ball end key 1.3 are intended for multiple use. The product must be disinfected, cleaned and sterilized before being used again. Instruments and accessories must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels. Additional information can be found at [www.dentaurum.com](#) (Processing Instructions Prosthetic Components, Instruments and Accessories REF 989-801-09).

#### 11. Information on disinfection/cleaning/sterilizing

tiologic® TWINFIT prosthetic components are delivered non-sterile. They are for single use on one patient only. They must be cleaned and disinfected before being used on a patient. Dentaurum generally recommends that the product be sterilized in addition. Only sterilized prosthetic components may be used if bleeding occurs.

Additional information can be found at [www.dentaurum.com](#) (Processing Instructions Prosthetic Components, Instruments and Accessories REF 989-801-09).

#### 12. Product overview/scope of delivery, incl. REF

Please refer to the tiologic® TWINFIT product catalog for the tiologic® TWINFIT product range.

#### 13. Quality information

Dentaurum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaurum has no influence on the processing on site.

#### 14. Explanation of symbols

 Please refer to the label. Additional information can be found at [www.dentaurum.com](#) (Explanation of symbols REF 989-313-00).

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