
Gebrauchsanweisung | Instructions for use |
Mode d'emploi | Modo de empleo | Modalità d'uso

CE 0483

tomas®-abutments
tomas®-transfer cap
tomas®-laboratory pins

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Dear customer

Thank you for choosing a quality product from Dentaureum.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit.

Instructions for use cannot describe every eventuality and possible application. Should you have any questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaureum.com, even though you may frequently use the same product.

1. Manufacturer

Dentaureum GmbH & Co. KG | Turnstr. 31 | 75228 Ispringen | Germany

2. General description


A temporary, reliable skeletal anchorage option is created for orthodontic treatment with the aid of an endosteally anchored mini-implant (tomas[®]-pin). The tomas[®]-abutment/tomas[®]-transfer cap/tomas[®]laboratory pin extend the range of tomas[®]-pin indications. With these abutments and laboratory elements as well as further Dentaureum components that are not included in the scope of delivery (see www.dentaureum.com) various orthodontic appliances can be fabricated, e.g. for skeletal anchored distalization, mesialization, intrusion, rapid palatal expansion.

All abutments, except for the tomas[®]-abutment EP, connect with the hexagon of the tomas[®]-pin SD and the tomas[®]-pin EP via a snap mechanism. The snap mechanism guarantees the secure fit of the tomas[®]-abutment on the tomas[®]-pin head and the easy removal of the tomas[®]-abutment from the hexagon in vertical direction.

The tomas[®]-abutment universal is a one-piece construction. The abutment wings have indents lengthways and crossways to insert wires (ø 1.1 mm). The other tomas[®]-abutments with the snap mechanism consist of three pieces (see Fig. 1). The base body of the abutment is the connection to the tomas[®]-pin. The upper part has a groove that fits a ring with a welded coupling element (tube or wire). Due to the welded lid, the coupling elements can be turned 360°. This way, the appliance can be positioned in any direction independent from the position of the tomas[®]-pin. Rotating forces from the appliance or tooth movements are not transferred onto the tomas[®]-pin.

The tomas[®]-abutment EP consists of a loop with a welded tube. This abutment can only be used with the tomas[®]-pin EP.

Fig. 1

The tomas®-transfer caps are made of plastic and their main purpose is to transfer the oral status onto the model. They fit on all tomas®-pin heads. The tomas®-transfer caps are vertically plugged onto the hexagon of the tomas®-pin. The hold is achieved by the fit only. Together with the tomas®-pin EP, the tomas®-transfer caps can be used as a support for temporary crowns.


The tomas®-laboratory pins are used to fabricate skeletal anchored appliances in the laboratory. Analogous to the tomas®-pin heads there are two versions of the tomas®-laboratory pin. The tomas®-laboratory pin SD has a head with cross slot, the tomas®-laboratory pin EP has a mushroom-shaped head. The tomas®-laboratory pins are provided with retentions that guarantee secure fixing in the model material.

3. Product overview and indications

The tomas®-abutments are designed to fabricate a wide variety of skeletal anchored appliances. The following product overview lists some of the possible indications. Further information can be found on the Dentaurem website (www.dentaurem.com).

- Direct coupling means that the orthodontic appliance is directly connected to the tomas®-pin via the tomas®-abutment. The tomas®-pin serves as counter bearing for the force exerted by the appliance.
- Indirect coupling means that the orthodontic appliance is fixated on the teeth, just as in the conventional techniques. The tooth or the group of teeth that serve as counter bearing are coupled with the tomas®-pin via the tomas®-abutment using a sufficiently dimensioned wire.

Description	REF	Quantity	Figure	Indication
tomas®-abutment universal	302-025-00	2 pieces		Indirect coupling. Direct coupling with various appliances for skeletal anchorage.
tomas®-abutment tube 1.1	302-025-11	1 piece		Direct coupling with various appliances for bi- or unilateral mesialization, distalization, intrusion. The tomas®-abutment wire 6 and the tomas®-abutment wire 12 can also be used.
tomas®-abutment tube 1.5	302-025-15	1 piece		Use as hybrid RPE appliance or skeletal RPE appliance. Direct coupling with hyrax® click or hyrax®.
tomas®-abutment double tube 1.5	302-025-30	1 piece		Direct coupling of various appliances for bilateral or unilateral mesialization, distalization. In combination with tomas®-abutment U-wire for tomas®-pin placed in median position in the palate.

Description	REF	Quantity	Figure	Indication
tomas®-abutment wire 6	302-025-06	1 piece		Indirect coupling. Direct coupling with various appliances for unilateral mesialization, distalization, intrusion. The tomas®-abutment tube 1.1 can also be used.
tomas®-abutment U-wire	302-025-21	1 piece		Direct coupling of various appliances for bilateral or unilateral mesialization, distalization. In combination with tomas®-abutment double tube 1.5 for tomas®-pins placed in median position in the palate.
tomas®-abutments median	302-025-51	1 pair		Combination of the tomas®-abutment U-wire (REF 302-025-21) and tomas®-abutment double tube 1.5 (REF 302-025-30). This combination is meant for two tomas®-pins which have been positioned sagittally in the region of the palatal suture (median insertion). The distance between the tomas®-pins must be at least 6 mm.
tomas®-abutment plain	302-026-00	1 piece		To weld individual elements. Anchoring temporary crowns.

Description	REF	Quantity	Figure	Indication
tomas®-abutment EP	302-027-00	2 piece		Direct coupling between tomas®-pin EP and amda® (advanced molar distalization appliance). Direct coupling with various appliances for bi- or unilateral mesialization, distalization, intrusion.
tomas®-transfer cap	302-028-01	2 pieces		Transferring oral status to a model. Anchoring temporary crowns.
tomas®-laboratory pin EP	302-029-01	1 piece		Laboratory analog for the tomas®-pin EP to manufacture skeletal anchored orthodontic appliances in the laboratory.
tomas®-laboratory pin SD	302-030-01	1 piece		Laboratory analog for the tomas®-pin SD to manufacture skeletal anchored orthodontic appliances in the laboratory.

4. Precautions

Before use, operators should ensure that they have carefully read and followed the instructions for use supplied.

We recommend operators complete a relevant tomas® training course prior to use, as the instructions for use can only cover some of the many areas of application.

With tomas®-abutments / tomas®-transfer caps / tomas®-laboratory pins only use original tomas® components according to the instructions for use. The tomas®-abutments and tomas®-transfer caps were specially developed for the tomas®-pin.

If the tomas®-abutments are used in the mouth as single components without an

appliance, they need to be fixed to a safety thread (dental floss or similar). All tomas®-abutments and tomas®-transfer caps are provided with loops or undercuts for this purpose.

If wires or similar elements are soldered onto tomas®-abutments, these may deform. A perfect fit can no longer be ensured. Consequently, avoid soldering.

5. Application

5.1 tomas®-abutment

The tomas®-abutments are coupled to the tomas®-pin via a snap mechanism. This form-fitting connection uses the hexagon of the tomas®-pin. When sliding the tomas®-abutment on the head of the tomas®-pin it is important to observe the position of the hexagon. The tomas®-abutment snaps into place by using the finger to press lightly. Due to its construction it can be somewhat difficult to slide the abutment on the first time. The tomas®-abutment is positioned correctly if the hexagon and the groove of the tomas®-pin are completely covered (see Fig. 2). To release the snap mechanism the tomas®-abutment is pulled vertically from the tomas®-pin.

Fig. 2

On the tomas®-abutment universal the abutment wings have indents lengthways and crossways to insert wires (ø 1.1 mm). For temporary fixation or first fixation, you can insert a wire in the length groove. The wire can be permanently fixed by welding, but also with adhesive or ligatures. If the lower groove lengthways is used to fixate the wire, the abutment cannot be used with the tomas®-pin EP.

Put the tomas®-abutments with welded attachments (tubes, wires) on the tomas®-pin or tomas®-laboratory pin according to the instructions and then fabricate orthodontic appliance.

5.2 tomas®-abutment EP



The tomas®-abutment EP can only be used with the tomas®-pin EP. The tube (ø 1.05 mm) enables the connection to an orthodontic appliance. The loop of the tomas®-abutment EP is slid over the mushroom-shaped head of the tomas®-pin EP. The reciprocal force of the appliance pushes the loop under the mushroom-shaped head and anchors the abutment to the pin. As additional safety measure, the head of the tomas®-pin EP and the tomas®-abutment EP can be surrounded with adhesive.

5.3 tomas®-transfer cap

To transfer the clinical position of the tomas®-pin to a model, slide the tomas®-transfer cap over the head of the tomas®-pin. In the case of the tomas®-pin SD, the tomas®-transfer cap should be positioned so that the retention wings follow the course of a slot. There are additional markings of the slots for positioning on the bottom of the tomas®-transfer cap (see Fig. 3).

The tomas®-transfer cap is anchored by alignment on the hexagon of the tomas®-pin and the upper cylindrical part of the gingival collar. The cap is positioned correctly if the lower border of the tomas®-transfer cap is aligned with the lower border of the cylinder (see Fig. 4). If parts of the cylinder or the entire cylinder lie beneath the gingiva, the tomas®-transfer cap has to be shortened. If an impression of two tomas®-pins next to each other is required, the retention wings of the tomas®-transfer caps should not touch one another. The wings can be shortened by milling.

Before taking the impression, make sure the tomas®-transfer cap is fitted correctly on the tomas®-pin. Transferring the position of the tomas®-pin correctly can only be guaranteed with a distortion-resistant impression tray and a silicone or polyether as impression material. Using alginate or plastic impression trays can lead to inaccuracies. To fix the tomas®-transfer cap securely in the impression, the cap has to be surrounded with impression material before taking the impression of the jaw. After the impression material has hardened completely, remove the impression tray from the mouth carefully. It is important to observe the insertion direction of the tomas®-pins.

5.4 tomas®-laboratory pin

The tomas®-laboratory pin is the laboratory analog for the respective tomas®-pin. Before making the model, make sure the tomas®-transfer cap is securely fixed within the impression material. Then stick the tomas®-laboratory pin into the tomas®-transfer cap. With the tomas®-laboratory pin EP, the direction is irrelevant.

If you want to reproduce the exact slot position of the tomas®-pin SD, you can use the arrows on the base of the tomas®-transfer caps as a guide to align the laboratory analog. These markings show the course of the slots and serve as guide when placing the tomas®-laboratory pins SD (see Fig. 4 + 5).

The tomas®-laboratory pin must be fixed securely within the impression cap. In order for the tomas®-laboratory pin not to loosen from the tomas®-transfer cap during the filling process of the impression material, it is important to fix the pins with wax. Then, make the model.

6. Composition

The composition is included in the materials list; please refer to the catalog or www.dentaurum.com.

7. Safety instructions

The product should not be used if there is a known allergic reaction to one or more of the material components.

8. Information for single use products

The tomas®-abutments, tomas®-transfer caps und tomas®-laboratory pins are designed for single use only. Reconditioning of products that have been inserted previously or reuse on patients is not permitted.

9. Quality information

Dentaurum ensures faultless quality of its products. These recommendations are based upon Dentaurum's own experiences. The user is responsible for the processing of the products. Responsibility for failures cannot be taken by Dentaurum as we have no influence on the processing on site.

10. Explanation of symbols used on the label

Please refer to the label. Additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00).

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