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
## EN

**Instructions for use**

**Dear customer,**  
Thank you for choosing a quality product from Dentaurum. It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current instructions for use supplied with the product and stored in the internet at [www.dentaurum.com](http://www.dentaurum.com), even though you may frequently use the same product.

1. **Manufacturer**  
Dentaurum GmbH & Co. KG
I Turnstr. 31 | 75228 Ispringen | Germany

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Do not loosen these 2 screws, as the ratchet will lose its torque function.

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2. **General device description, intended purpose and characteristics, technical description**

The torque ratchet with 10 – 50 N cm torque control and block function for higher torque is a dental instrument for the manual preparation of the implant site, for tightening and loosening screws for implants and prosthetic components and for inserting implants.

3. **Area of application**

Dental instrument for tightening and loosening screws for implants and prosthetic components and for inserting implants

4. **Medical indications/contraindications**

**Indications**

The torque ratchet is indicated for manual thread tapping and insertion of implants in the jawbone. It can also be used to fix the AnoTite screw and the prosthetic screw (according to the drilling protocol for implants from Dentaurum)

**Contraindications.**

All further indications not mentioned here.

5. **Installation / Assembly**

The torque ratchet is a precision instrument that can be disassembled. It is not sterile on delivery. After each use the torque ratchet must be disassembled, disinfected, cleaned, lubricated and sterilized according to the Instructions for use to ensure proper functioning. Read the Instructions for use prior to handling. Handling and use of the product are carried out without direct contact from the manufacturer and remain the responsibility of the user. The dental professional is solely responsible for any damages that may occur. The proper functioning of the torque ratchet should be checked before each use to guarantee torque precision. The torque ratchet should be used only by personnel trained in dentistry.

6. **Information on use, processing, initial implementation and handling**

To set the correct torque, the torque adjustment screw ⊗ is turned clockwise to the required torque line on the scale sleeve ⊗. The torque ratche is additionally provided with a blocking function. To set the blocking function, turn the adjustment screw to the "∞" symbol (Fig. 1). Do not turn too tightly. For storage, turn the torque adjustment screw back until the spring is as relaxed as possible (approx. 10 Nm).

The pressure point for exact torque release is at the head of the torque adjustment screw ⊗, see arrow. When the set torque has been reached, the scale sleeve ⊗ will bend around the axis in the ratchet head ⊗ (Fig. 2). The release is audible and perceptible. After the torque release, DO NOT apply more pressure as this could damage the ratchet. When you let go of the torque adjustment screw ⊗, the ratchet returns to its initial position.

The word 'IN' on the ratchet head shows that the ratchet is in the correct position for tightening. The word 'OUT' stands for loosening the torque.

7. **Composition**

Sleeve, ratchet latch, ratchet head	1.4197
Springs	1.4310
Washer	PEEK (polyether ether ketone)
Remaining components	1.4305

Information about the composition is included in the material list for our implantology products (REF 989-801-05), or see [www.dentaurum.com](http://www.dentaurum.com). The SSCP is available at <https://ec.europa.eu/tools/eudamed>.

8. **Safety information**

All components of the torque ratchet must be checked before use for wear and to ensure they function correctly. In the event of damage, malfunction or risk of malfunction, the ratchet must be returned to the manufacturer. To ensure maximum precision, the single components of the torque ratchet are coordinated and adjusted during manufacturing. During cleaning and assembling of the torque ratchet ensure that no components from a different torque ratchet are used, as the tolerance values can then no longer be ensured to be correct. The replacement parts cannot be purchased separately. Store the torque ratchet with the spring relaxed at approx. 10 Nm. The instrument should not be sterilized in the packaging provided by the manufacturer.

Exceeding the torque specified by Dentaurum can cause mechanical damage to components, to the implants, and destruction of bone structures. If the packaging shows serious damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary. The block function mode should be used with extreme caution. After use the valve must be reset to standard torque to prevent subsequent faults. The precision of the instrument is ± 10,0 % on leaving the factory.

9. **Storage and shelf life**

**Service life**

3 years from the manufacturer supply date.

**Storage**

After sterilization, the ratchet must be stored with the spring relaxed in a dry, dust-free, ventilated place that is free of corrosive vapors.

10. **Information for multiple use products**

The torque ratchet is intended for multiple use. It must be disinfected, cleaned, lubricated and sterilized before being used again. For more details, see the section "Information on disinfection, cleaning and sterilization".

11. **Notes on maintenance**

11.1. **Disassembly**

Before cleaning (regardless of the selected cleaning method), the torque ratchet must be dismantled into the individual parts. This can be done without tools. Completely unscrew the torque adjustment screw ⊗, and remove the spring ⊗ and the ratchet head ⊗ with threaded rod.

Take care not to lose the plastic washer ⊗ as this would have a negative impact on the instrument's precision. (The plastic washer needs only to be removed if there is visible contamination. It can be pulled off if necessary and replaced after cleaning.)

**Remove ratchet wheel**

Pull back the pin ⊗ in the direction of the arrow using your thumb and index finger and remove the ratchet wheel ⊗ (Fig. 4 – 6).

11.2. **Assembly**

To assemble the torque ratchet correctly, connect the components in the following order: first pull back the pin ⊗ as described above and insert the ratchet wheel ⊗.

**Caution:** To avoid confusion, the ratchet wheel ⊗ can only be inserted from one side. If several torque ratchets are in use, do not interchange the individual parts. Each individual part belongs to one instrument.

11.3. **Lubricating point (▲)**

Lubricate the areas marked with the "drop" symbol lightly with maintenance oil for instruments.

Ensure that only instrument oils (paraffinic white oil without corrosion inhibitor or other additives) are used, which – depending on the maximum sterilization temperature used – are approved for steam sterilization and are certified as biocompatible. The oil should be used sparingly.

Reassemble the ratchet and perform a function test.