DE

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Abbreviated instructions for use

Dear customer,

Thank you for choosing a quality product from Dentaurum.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. In case of questions or ideas, please contact your local representative. As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaurum.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurum GmbH & Co. KG
I Turnstr. 31
I 75228 Ispringen
I Germany

2. General device description, intended purpose and characteristics, technical description

The tiologic® TWINFIT implant system contains specially coordinated instruments, abutments, and accessory components for the prosthetic restoration of the tiologic® TWINFIT implants. Only original components of the tiologic® TWINFIT implant system should be combined in accordance with the Instructions for use.

The abutment is selected by the dentist/dental technician to suit the individual patient's needs.

The diameter of the implant and the gingival thickness determine the choice of abutment. Selection aids with conical and platform connector geometry for the abutment series S, M, and L are available to help choose the right abutment on the model (gingival

height and angulation).

3. Area of application

The 4Base abutments tiologic® TWINFIT are prosthetic components for bridge and bar restorations on tiologic® TWINFIT implants.

4. Medical indications / contraindications

4Base abutments tiologic® TWINFIT

Indications

- Anchorage of implant-supported, splinted bridges and bars in the maxilla in combination with at least 6 tiologic® TWINFIT implants and, in the mandible, in combination with at least 4 or more tiologic® TWINFIT implants.
- Contraindications.**
 - For biomechanical reasons we recommend that the following angulations are not exceeded:

4Base abutment with 0°:	-20° - 20°
4Base abutment with 20°:	0° - 40°
4Base abutment with 30°:	10° - 50°
 - Restorations with mixed retention (tooth/4Base)
 - Combination with other anchorage elements
 - Grinding and shortening 4Base abutments

5. Information on use, processing, fitting and insertion

5.1. 4Base abutments tiologic® TWINFIT

4Base abutments with conical and platform connector geometry are available for the S, M and L series of abutments. 4Base abutments are provided in 3 angulations: 0°, 20° and 30°. The seating surface of the 4Base caps (ø 4.6 mm) and the cone (20°) of 4Base abutments are always identical for the S, M and L series of abutments, so that only one size of impression posts, closure caps etc. is required.

The gingival heights of the straight 4Base abutments are 1.5 mm, 2.5 mm and 4.0 mm, that of the angulated 4Base abutments is 3.2 mm.

In case of steep inclination of the implants, splinting the denture is absolutely essential. When fitting the angulated 4Base abutments, keep in mind that the abutments could have a larger diameter than the implants. As part of the abutments may be below the bone line, bone may have to be removed in this region, if required, so that the abutment sits flush on the implant. Prerequisite for the successful use of the 4Base system is the best possible accurate angular position of the implant. The more accurate this angle can be maintained, the easier the prosthetic treatment, as the crowns of the abutments are then positioned parallel to one another. The mounted insertion guide facilitates the insertion of the angulated 4Base abutments. At the same time, they secure the pre-assembled AnoTite screw in the 4Base abutment. The insertion guide can be bent sideways to tighten the AnoTite screw.

The 4Base abutments should **NOT** be ground or shortened. When handling the 4Base abutments, ensure that the occlusal screw is not damaged by tweezers etc.

Every construction must be checked for passive fit prior to placing it intra-orally (Sheffield test). This involves placing the construction on the 4Base abutment and only one prosthetic screw (not included in delivery) is fixed in the 4Base abutment. If the construction is raised and a gap is created between the construction and 4Base abutment, stresses are present and must be eliminated.

5.2 4Base impression post, open and closed tiologic® TWINFIT

To achieve stable apposition of the gingiva with 4Base abutments, it is possible to take an open or closed impression on a 4Base abutment that is fitted in the mouth. Since the 4Base impression posts have no rotational security, they may only be used for bridge and bar constructions.

5.3 4Base closure screw tiologic® TWINFIT

The 4Base closure screw is screwed directly onto the 4Base abutment after impression-taking and protects the surface for the caps on the 4Base abutments.

5.4 4Base plastic cap tiologic® TWINFIT

The plastic cap, that burns residue-free, can be individually adapted to the space available. The processing instructions given by the investment material manufacturer must be adhered to for investing and preheating. After casting, the casting ring must be cooled to room temperature to allow vesting and to achieve stress-free cast objects. The cast object should never be blasted with abrasive agents since these can damage the inner configuration. Vesting is carried out either by steam cleaning or an ultrasonic bath.

5.5.4Base titanium cap tiologic® TWINFIT – lasering

The 4Base titanium cap is lasered with prefabricated bar segments made of titanium (e.g. rematitan®, Dentaurum). The prefabricated titanium bars are shortened and fixed in the correct horizontal axis for laser-welding with a parallelometer / fixator for the laser technique (e.g. Paralas, Dentaurum). The bar components and titanium caps are then lasered together. The only filler material which can be used is pure titanium wire (e.g. Dentaurum). Should further questions arise on lasering and processing titanium, please call the Dentaurum Hotline at: + 49 72 31 / 803 -410.

5.6 4Base titanium cap tiologic® TWINFIT – adhesive technique

The 4Base titanium cap is available for the passive fit technique (adhesive technique). Once a cast bar has been produced, it is bonded intra-orally to the titanium adhesive base. This helps compensate possible inaccuracies of fit that could create tensions.

5.7 4Base scan cap titanium tiologic® TWINFIT

The surface of the 4Base scan cap titanium is optimized for digital capture, both intra-orally and on the model without scanning spray. For this purpose, the 4Base scan cap is screwed onto the 4Base abutment and scanned. The matching process of the scan data with the STL data takes place according to the specifications of the respective software manufacturer. The straight 4Base abutments (0°) are screwed into the tiologic® TWINFIT implant with the PentaGrip insertion key tiologic® TWINFIT – ISO shank. The angulated 4Base abutments (20° and 30°) have a rotationally secure fit in the tiologic® TWINFIT implants and are fixed with the AnoTite screw 9.0 mm. A new AnoTite screw should be used for the final fitting in the mouth.

To fix bridge and bar constructions on 4Base abutments, special prosthetic screws and AnoTite screws are available.

The special AnoTite screws with a reduced screw head height can be used for limited occlusal spaces. Such screws are not included in delivery.

Screws with hexagon screw head geometry are tightened manually with the hex key ratchet 1.3 or mechanically with the hex key ISO shank 1.3.

Tightening torque

- 4Base abutments, straight on model: by hand
- 4Base abutments, straight intra-orally: 35 Nm
- Prosthetic screw on model: manually
- AnoTite screw intra-orally: 30 Nm
- Prosthetic screw for 4Base caps on model: by hand
- AnoTite screw for 4Base caps intra-orally: 25 Nm
- Sure-grip screw impression post on 4Base abutment intra-orally: 15 Nm
- Closure screw on 4Base abutment intra-orally: 15 Nm

Attention: Do not confuse the prosthetic screw with the AnoTite screw. The prosthetic screw serves only to help the dental technician fix the abutment on the working cast. The prosthetic screw is colored silver and the AnoTite screw is anodized gold.

6. Composition

- 4Base abutments: titanium alloy (Ti6Al4V)

- 4Base caps:

- 4Base plastic cap: polycarbonate
- 4Base titanium cap, laser-weldable: titanium alloy (Ti6Al4V)
- 4Base titanium cap, adhesive technique: titanium alloy (Ti6Al4V)
- Closure screw: titanium alloy (Ti6Al4V)
- 4Base impression post, open: titanium alloy (Ti6Al4V) / screw: 1.4305
- 4Base impression post, closed:
 - impression post and screw: titanium alloy (Ti6Al4V) / impression cap: POM
- 4Base scan cap titanium: titanium alloy (Ti6Al4V)
- 4Base tiolOC cap: titanium alloy (Ti6Al4V)
- 4Base insertion guide: PEEK (polyether ether ketone)

- 4Base laboratory implant:** titanium alloy (Ti6Al4V)

Prosthetic screws:

- AnoTite screw: titanium alloy (Ti6Al4V)
- Prosthetic screw: titanium alloy (Ti6Al4V)
- Prosthetic screw for 4Base cap: titanium alloy (Ti6Al4V)
- AnoTite screw for 4Base cap: titanium alloy (Ti6Al4V)
- PentaGrip insertion key - ISO shank 4Base/ball abutment:** 1.4197

- Hex key:**

- Hexagon key - ISO shank 1.3: 1.4197
- Hexagon key - ratchet 1.3: 1.4197
- Hexagonal ball end key - ISO shaft 1.3: 1.4305

Information about the composition is included in the material list for our implantology products (REF 989-801-05), or see www.dentaurum.com. The SSCP is available at https://ec.europa.eu/tools/eudamed.

7. Safety information

- The product should not be used if there is a known allergic reaction to one or more of the material components.
- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. All articles used intra-orally should therefore be secured against swallowing and/or aspiration.
- All serious incidents arising from the use of the product should be reported to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.
- Different types of alloy in the oral cavity can lead to galvanic reactions.
- If the packaging shows serious damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.
- Only use tiologic® TWINFIT prosthetic components in combination with tiologic® TWINFIT implants.

8. Storage and shelf life

Store in a dry place.

9. Information for single use products

4Base abutments and caps are designed for single use only. Reconditioning of a 4Base abutment or cap (recycling) or reuse on patients are not permitted since it can then no longer be guaranteed that the article can be reprocessed safely or can function safely.

10. Information for multiple use products

The PentaGrip insertion key, the hex key 1.3, the hexagonal ball end key and the 4Base scan cap titanium are intended for multiple use. The product must be disinfected, cleaned and sterilized before being used again. Instruments and accessories must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels. Additional information can be found at www.dentaurum.com (Processing Instructions Prosthetic Components, Instruments and Accessories REF 989-801-09).

11. Information on disinfection/cleaning/sterilizing

tiologic® TWINFIT prosthetic components are delivered non-sterile. They are for single use on one patient only. They must be cleaned and disinfected before being used on a patient. Dentaurum generally recommends that the product be sterilized in addition. Only sterilized prosthetic components may be used if bleeding occurs.

Additional information can be found at www.dentaurum.com (Processing Instructions Prosthetic Components, Instruments and Accessories REF 989-801-09).

12. Product overview/scope of delivery, incl. REF

Please refer to the tiologic® TWINFIT product catalog for the tiologic® TWINFIT product range.

13. Quality information

Dentaurum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaurum has no influence on the processing on site.

14. Explanation of symbols

📘 Please refer to the label. Additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-3

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