DE

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Instructions for use

Dear customer,

Thank you for choosing a quality product from Dentaurem.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaurem.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurem GmbH & Co. KG | Turnstr. 31 | 75228 Ispringen | Germany

2. General device description, intended purpose and characteristics, technical description

The osteotomes are available in angled form with a concave tip and can be used in the anterior and posterior region of the jaw. They can be used to work the bone by pushing, turning or knocking. Osteotomes are used to manually form the available bone by condensing it laterally or apically.

Depending on the bone situation, osteotomes are used in ascending order after the marking drill and the ø 2,0 mm depth drill. The osteotomes are available in various diameters. The diameter is indicated on the osteotome. A sequence of lines on the tip of the instrument indicates different heights.

The number of osteotomes used to prepare the implant site depends on the following factors:

- Width of the alveolar ridge
- Implant diameter
- Thickness and structure of the cortical bone

3. Area of application

Osteotomes are used to condense or displace local bone in the soft maxilla and mandible (bone quality 3 and 4) in order to achieve sufficient primary stability of tiologic® ST / tioLogic® TWINFIT implants.

Osteotomes can also be used to elevate the sinus floor.

4. Medical indications/contraindications Indications

Osteotomes are used before inserting tiologic® ST / tioLogic® TWINFIT implants.

Contraindications.

General contraindications:

General contraindications for dental surgery procedures apply. These include: reduced immunodeficiency / steroid treatment / blood coagulation disorders / uncontrolled endocrine diseases / rheumatic disorders / bone system diseases / cirrhosis of the liver / drug, alcohol or tobacco abuse / depression, psychopathic disorders / poor patient compliance / chronic underlying inflammatory diseases / incomplete physical growth of patient.

Local / personal contraindications:

osteomyelitis / radiotherapy in the head region / recurring mucosal diseases / temporomandibular joint dysfunctions / parafunctions / lack of vertical or horizontal bone availability, jaw defects, inadequate bone quality / poor oral hygiene.

5. Information on use, processing, fitting and insertion

The osteotomy product range is supplied exclusively to doctors, dentists and dental technicians. It should only be used by doctors or dentists who are familiar with dental implantology, including diagnosis, preoperative planning, surgical technique and prosthetic treatment.

Before use, operators should ensure that they have carefully read and understood all tiologic® ST / tioLogic® TWINFIT Instructions for use / manuals. These abbreviated instructions do not replace the Instructions for use / manuals for the tiologic® ST / tioLogic® TWINFIT implant system. As the instructions and manuals cannot provide all information for immediate use, we strongly recommend that, before using the system, implantologists attend a tiologic® ST / tioLogic® TWINFIT implant system training course offered by Dentaurem to learn the correct techniques.

- Please refer to the **Product Catalog** and the **Surgery Manuals** for information on the selection of components for the clinical procedure.

- Please refer to the **Product Catalog** and the **Prosthetic Manuals** for information on the selection of components for the prosthetic procedure.

Before using this product, the implantologist must examine the patient thoroughly and explain the procedure in detail. Dentaurem recommends full clinical, radiological, photographic and statistical documentation.

The implantologist should ensure the products cannot be aspirated during intra-oral use. Improper surgical and prosthetic procedures may damage the implant or lead to bone loss.

When using the instruments, proper safety clearance to the remaining teeth and tissues must be maintained. Drilling or working the bone with osteotomes beyond the planned depth due to an incorrect calculation of the drill depth in relation to the radiograph can lead to permanent nerve or tissue damage. A 1.5 mm safety clearance to the nervus mandibularis or nervus alveolaris inferior must be maintained. Smoothing the bone edge will reduce the available bone height.

Not all components are available in every country.

6. Composition

- Osteotomes, angled, concave: 1.4021
- Osteotome tray: Radel® (PPSU) & silicone

Information about the composition is included in the material list for our implantology products (REF 989-801-05), or see www.dentaurem.com. The SSCP is available at https://ec.europa.eu/tools/eudamed.

7. Safety information

- Using third party instruments and components may impair safety and function of the tiologic® ST / tioLogic® TWINFIT implant systems.
- The product should not be used if there is a known allergic reaction to one or more of the material components.
- Different types of alloy in the oral cavity can lead to galvanic reactions.
- All serious incidents arising from the use of the product should be reported to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.
- If the packaging shows serious damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.
- Only use tiologic® ST surgical and prosthetic components with tiologic® ST implants!
- Only use tiologic® TWINFIT surgical and prosthetic components with tiologic® TWINFIT implants!
- All osteotomes have length markings on the instrument tip and laser markings indicating color and diameter.

8. Warnings

Warning: Beware of perforating the maxillary sinus!

9. Storage and shelf life

Store in a dry place.

10. Information on disinfection/cleaning/sterilizing

All osteotomes are intended for multiple use and are supplied non-sterile. All osteotomes must be disinfected, cleaned and sterilized before every use. Instruments and accessories must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels. Additional information can be found at www.dentaurem.com (Processing Instructions Prosthetic Components, Instruments and Accessories REF 989-801-09).

Important: To guarantee adequate cleaning conform with Robert Koch Institute regulations, the osteotomes must be cleaned separately and not in the osteotome tray. The osteotome tray serves only to hold the osteotome during sterilization and then to store the osteotome. For sterilization, the loaded osteotome tray is packed in disposable sterilization packaging dedicated to the sterilization process in accordance with ISO 11607-1/2.

11. Technical data

All osteotomes have length markings on the instrument tip and laser markings indicating color and diameter.

12. Product overview /scope of delivery, incl. REF

350-103-00	Osteotome tray, not equipped
350-170-00	Osteotome ø 2,0 – 2,5 mm
350-172-00	Osteotome ø 2,6 – 3,1 mm
350-174-00	Osteotome ø 3,0 – 3,5 mm
350-176-00	Osteotome ø 3,5 – 4,0 mm
350-178-00	Osteotome ø 4,0 – 4,6 mm
350-180-00	Osteotome ø 4,8 – 5,3 mm

13. Quality information

Dentaurem ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaurem has no influence on the processing on site.

14. Explanation of symbols

 Please refer to the label. Additional information can be found at www.dentaurem.com (Explanation of symbols REF 989-313-00).

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