DE
----

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## Instructions for use

**Dear customer,**

Thank you for choosing a quality product from Dentaum.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit.

Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaum.com, even though you may frequently use the same product.

**1. Manufacturer**

Dentaum GmbH & Co. KG | Turnstr. 31 | 75228 Ispringen | Germany

**2. General device description, intended purpose and characteristics, technical description**

CITO mini<sup>®</sup> implants are designed for insertion in the endosteal region of the maxilla or mandible. The one-piece implants are fitted with a ball ø 1.8 mm, on top of which matrices of the system are fixated once they have been polymerized into the patient’s dental prosthesis. The CITO mini<sup>®</sup> implant system contains specially coordinated instruments, abutments and accessories for insertion of the implants and for the dental prosthetic restoration. Only original components of the CITO mini<sup>®</sup> implant system should be combined in accordance with the Instructions for use.

**3. Area of application**

CITO mini<sup>®</sup> implant types can be used in both the mandible and maxilla for surgical immediate implantation, delayed immediate implantation and delayed implantation using either the one-stage or two-stage technique.

**4. Medical indications / contraindications**

**Indications**

CITO mini<sup>®</sup> implant types are used to stabilize prostheses in the edentulous maxilla and mandible. The indication for an implant should be made with consideration of possible advantages and disadvantages, as well as risks of implant treatment and alternative treatments.

With all implantological cases, the implant diameter and length of the CITO mini<sup>®</sup> implant types should be in proportion to the prosthetic restoration.

When choosing the number of implants to be inserted, it is important to consider the bone quality and whether the implants will be loaded immediately.

To fixate full dentures, Dentaum recommends at least four implants in the mandible and six implants in the maxilla\* .

The insertion of CITO mini<sup>®</sup> implants in soft bone is not indicated.

For restorations in the maxilla, use CITO mini<sup>®</sup> implants ø 2.2 mm or ø 2.5 mm. When positioning the implants, pay attention to keeping an adequate distance between the implants and to ensuring that the static force on the bone is optimally applied\*\*.

Immediate implant loading is only possible if a torque of 35 Ncm was reached during insertion. Torque should not exceed 40 Ncm during insertion.

\* see Planning Scheme CITO mini<sup>®</sup> (989-502-45) or www.dentaum.com

\*\* The size of the matrix should also be considered

**Contraindications**

**General contraindications:**

General contraindications for dental surgery procedures apply. These include: reduced immunodeficiency / steroid treatment / blood coagulation disorders / uncontrolled endocrine diseases / rheumatic disorders / bone system diseases / cirrhosis of the liver / drug, alcohol or tobacco abuse / depression, psychopathic disorders / poor patient compliance / chronic underlying inflammatory diseases / incomplete physical growth of patient / insertion of CITO mini<sup>®</sup> implants into soft bone is not indicated.

**Local / personal contraindications:**

Osteomyelitis / radiotherapy in the head region / recurring mucosal diseases / tem-poromandibular joint dysfunctions / parafuncions / lack of vertical or horizontal bone availability, jaw defects, inadequate bone quality / poor oral hygiene.

**5. Information on use, processing, fitting and insertion**

The CITO mini<sup>®</sup> product range is supplied exclusively to doctors, dentists and dental technicians. It should only be used by doctors, dentists or dental technicians who are familiar with dental implantology, including diagnosis, preoperative planning, surgical techniques and prosthetic restorations.

Before use, dental professionals should ensure that they have carefully read and understood the full CITO mini<sup>®</sup> Instructions for use / manuals. As the instructions and manuals cannot provide all information for immediate use, we strongly recommend that, before using the system, implantologists attend a CITO mini<sup>®</sup> implant system training course offered by Dentaum to learn the correct techniques.

■ Refer to the order form and the Surgery + Prosthetics Manual (989-360-20) for information on the selection of components for the clinical and prosthetic procedure.

Before using this product, the implantologist must examine the patient thoroughly and explain the procedure in detail. Dentaum recommends full clinical, radiological, photographic and statistical documentation.

The CITO mini<sup>®</sup> implant system components can be documented, e.g. in the patient file or PatientPass, using the additional labels.

The implantologist should ensure the products cannot be aspirated during intra-oral use.

**Not all components are available in every country.**

**6. Composition**

■ **Titanium grade 5 (titanium alloy)** DIN EN ISO 5832-3

<b>Chemical composition (% by mass)</b>	<b>Al</b>	<b>5.5 % – 6.75 %</b>
	<b>V</b>	<b>3.5 % – 4.5 %</b>
	<b>Fe</b>	<b>0.3 % max.</b>
	<b>C</b>	<b>0.08 % max.</b>
	<b>N</b>	<b>0.05 % max.</b>
	<b>H</b>	<b>0.015 % max.</b>
	<b>O</b>	<b>0.2 % max.</b>
	<b>Ti</b>	<b>Rest</b>
<b>Physical and mechanical properties</b>	<b>0.2 % yield strength</b>	<b>780 MPa</b>
	<b>Tensile strength</b>	<b>860 MPa</b>
	<b>Elongation at rupture</b>	<b>10 % min.</b>

Information about the composition is included in the material list for our implantology products (REF 989-801-05), or see www.dentaum.com. The SSCP is available at https://ec.europa.eu/tools/eudamed.

**7. Adverse reactions**

CITO mini<sup>®</sup> implant types are not known to have any side effects or produce any irritation.

**8. Safety information**

■ The product should not be used if there is a known allergic reaction to one or more of the material components.

■ If the packaging shows serious damage, the product must be checked for integrity and cleanliness before use, and if necessary must be disposed of.

■ Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. All articles used intra-orally should therefore be secured against swallowing and/or aspiration.

■ All serious incidents arising from the use of the product should be reported to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.

■ Different types of alloy in the oral cavity can lead to galvanic reactions.

■ Only use CITO mini<sup>®</sup> prosthetic components in combination with CITO mini<sup>®</sup> implants!

**9. Warnings**

Though placement of dental implants has a high rate of success and implants have a long durability, successful treatment cannot be guaranteed. The implantologist should note and document any problematic cases and inform the manufacturer.

**Indications**

CITO mini<sup>®</sup> implant types are used to stabilize prostheses in the edentulous maxilla and mandible. The indication for an implant should be made with consideration of possible advantages and disadvantages, as well as risks of implant treatment and alternative treatments.

With all implantological cases, the implant diameter and length of the CITO mini<sup>®</sup> implant types should be in proportion to the prosthetic restoration.

When choosing the number of implants to be inserted, it is important to consider the bone quality and whether the implants will be loaded immediately.

To fixate full dentures, Dentaum recommends at least four implants in the mandible and six implants in the maxilla\* .

The insertion of CITO mini<sup>®</sup> implants in soft bone is not indicated.
For restorations in the maxilla, use CITO mini<sup>®</sup> implants ø 2.2 mm or ø 2.5 mm. When positioning the implants, pay attention to keeping an adequate distance between the implants and to ensuring that the static force on the bone is optimally applied\*\*.

Immediate implant loading is only possible if a torque of 35 Ncm was reached during insertion. Torque should not exceed 40 Ncm during insertion.
\* see Planning Scheme CITO mini<sup>®</sup> (989-502-45) or www.dentaum.com

Dentaum guarantees the sterility of the CITO mini<sup>®</sup> implants until the expiry date on the packaging, provided that the original packaging is undamaged. No guarantee can be given for the sterility after the expiration date and the CITO mini<sup>®</sup> implants may no longer be used on patients after this date. If the sterile packaging is damaged, the product may not be used or re-sterilized. **Packaging that has been opened may not be re-sterilized.** A CITO mini<sup>®</sup> implant that has been removed from the packaging but has not been inserted, may also no longer be used or re-sterilized.

All other components of the CITO mini<sup>®</sup> implant system are supplied non-sterile.

**12. Information for single use products**

CITO mini<sup>®</sup> implants are designed for single use only. Reconditioning of a CITO mini<sup>®</sup> implant (recycling) and its reuse on a patient are not permitted.

**13. Information for multiple use products**

Rotating instruments and many accessories are intended for multiple use. The product must be disinfected, cleaned and sterilized before being used again. This applies in particular to components that can be disassembled, e. g. the torque wrench. Rotary instruments – used with proper care and provided that they are not damaged or contaminated – can be reused in dense bone 15 to 20 times; any further reuse or the use of damaged and / or contaminated instruments must be avoided. The dental professional bears responsibility. No liability is accepted if these instructions are disregarded.

Instruments and accessories must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels. Additional information can be found at www.dentaum.com (Processing Instructions Prosthetic Components, Instruments and Accessories REF 989-801-09).

**14. Information on disinfection, cleaning and sterilizing**

Additional information can be found at www.dentaum.com (Processing Instructions Prosthetic Components, Instruments and Accessories REF 989-801-09).

**15. Notes on compliance / Conformity declaration**

Development, clinical testing, production and quality control of the CITO mini<sup>®</sup> prod-uct range are completed in accordance with the Medical Device Directive 93/42/EEC. Sections 9 and 10 of our General Terms of Delivery and Payment apply with regard to warranty or liability – unless otherwise provided in this quick reference guide. Warranty and liability are rendered void in particular if the products are not used by the dental professional or a third party in accordance with the Instructions for use; this also applies if the CITO mini<sup>®</sup> product range is used in combination with products from other manufacturers which have not been specifically recommended for use by Dentaum. Dentaum has no control over processing and application of the product. These are the sole responsibility of the dental professional.

Technical advice (oral and written) is based on the scientific and technical knowl- edge available at the time the product is put on the market. It does not release the dental professional from the responsibility of personally checking the suitability of the products for the intended indications and applications. Advice is only given as a non-binding recommendation, which cannot be assumed to provide any form of assurance or guarantee.

All products are subject to continuous development based on current scientific knowl- edge and we reserve the right to make changes in the construction, design or material of the products.

**16. Product overview / scope of delivery, incl. REF**

Please refer to the CITO mini<sup>®</sup> order form for the CITO mini<sup>®</sup> product range.

**17. Quality information**

Dentaum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaum has no influence on the processing on site.

**18. Explanation of symbols**

 Please refer to the label. Additional information can be found at www.dentaum.com (Explanation of symbols REF 989-313-00).

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