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Stand der Information: 05/10

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Dear Customer

Thank you for choosing this high-quality product from Dentaaurum. In order to achieve the best results with this product, it is important to carefully study and follow these instructions for use.

The written instructions cannot cover all possible situations that can occur during use. For this reason, our hotline is available to answer any other questions or concerns that may arise.

Due to constant developments, we recommend, even when using the same products regularly, that you study the enclosed current instructions for use or refer to our website at www.dentaaurum.de.

1. Manufacturer

Dentaaurum GmbH & Co. KG
Turnstraße 31 · 75228 Ispringen · Germany

2. With reference to our quality standards

Dentaaurum ensures a faultless quality of the products that we manufacture. These recommendations are based upon our own experiences. The user himself is responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence on the processing on site.

3. Range of indications

The Herbst IV is a further development of the standard Herbst I hinge. It has a ball hinge that allows the patient greater freedom of movement, increasing comfort and tolerance during treatment. Retention using snap rings (C-clips) facilitates fitting by the orthodontist and replaces the relatively awkward insertion of screws.

4. Fabricating the Herbst IV appliance

Pour the impressions, which include the fitted, waxed out bands on the first premolars and molars, using Class IV dental stone, e.g. Rapidur® (REF 164-450-00). Then mount the models with plaster in the orthodontic Split Fixator (REF 072-004-00) to the correct relationship using the activating bite registered by the dentist. Bend remaloy® wire, half-round, e.g. \varnothing 1.50 x 0.75 mm (REF 528-155-00), or round, e.g. \varnothing 1.00 mm (REF 528-100-00), to form a lingual arch and then weld and/or solder the wire to the bands. Fit a palatal arch in the upper jaw to provide stabilization, e.g. remanium® GOSHGARIAN palatal bar or round 1.0 mm remaloy® wire (REF 528-100-00). Weld and/ or solder the base of the retention hinges to the bands, e.g. using the Master 3000 (REF 079-000-00).

Another option is to fabricate CrCo splints instead of using bands. The bases of the retention hinges should not be cast on in this technique but laser welded for reasons of accuracy and biocompatibility. Ensure that the screws are at the correct angle to the axis between the upper and lower jaw. It is important that the retention hinges operate friction free to guarantee perfect functioning of the telescope mechanism. Friction can damage the appliance – in particular the bands.

Fitting the telescope tubes and rods:

First, shorten the telescope tubes according to the activating bite and then fit the telescope rods, which should be approx. 2 mm longer. If the telescope rods are too short, they could come out of the telescope tube if the patient's mouth is opened too wide. If the telescope rods are too long, they could cause mucosal irritation (risk of injury!). The correct length should be checked again intraorally and adjusted if necessary.

Note:

We recommend using “double bands” to minimize the risk of shear forces damaging the bands. A second band is welded over the fitted band, e.g. using the Master 3000 (REF 079-000-00). Generally a band that is 2 to 3 sizes larger will fit over the first band. Pre-welded “double bands” can be ordered from Dentaaurum.

Fitting intraorally:

The bands are fitted intraorally using glass ionomer cement, e.g. Glas Tec (REF 163-540-00). After the cement has set, the prepared retention hinges are fitted on the bases and secured with a snap ring, C-clip. When the hinges are fitted correctly on the base, it should be possible to insert the C-clip into the two slots of the ball unit so that the ends are clipped in at both sides. We recommend using Weingart universal pliers (REF 003-120-00) or How pliers (REF 003-111-00) to facilitate gripping the C-clips.

The first recall appointment should be one week following fitting of the appliance!

If necessary, the Herbst IV appliance can be activated using various spacers, refer to Product range.

5. Safety instructions

Ensure that the C-clip is fitted parallel to the slot of the ball hinge! If the C-clip is fitted at an angle it can be bent and become loose in the mouth. Always check that the snap rings, C-clips, are properly seated after fitting the Herbst IV appliance and at each recall appointment. If the C-clips are repeatedly fitted and removed, they lose their spring tension and should be replaced!

6. Product range

Herbst IV set REF 607-115-00 1 Set
with ball and socket hinge

For welding or soldering to bands or cast splints. Ball and socket hinge for great lateral movement, with a „C“ clip to hold the rods in place. No screws required. Greater patient comfort.

Contents:

- 2 x Telescope tube with ball joint
- 2 x Telescope rod with ball and socket hinge
- 2 x Base for ball and socket hinge
- 4 x Bases
- 10 x C-clip
- 1 x Instructions for use

Snap rings (C-Clips) REF 607-116-00 10 pieces
for Herbst IV

Bases REF 607-118-00 10 pieces
for Herbst IV

Telescope tube REF 607-125-00 1 piece
with ball and socket hinge
for Herbst IV

Telescope rod with ball REF 607-128-00 1 piece
and socket hinge
for Herbst IV

Spacers
for Herbst I/II/IV
Length: 1 mm REF 607-103-00 10 pieces
Length: 2 mm REF 607-104-00 10 pieces
Length: 3 mm REF 607-105-00 10 pieces
Length: 4 mm REF 607-107-00 10 pieces
Length: 5 mm REF 607-112-00 10 pieces

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Dentaurum GmbH & Co. KG