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Dentaurum GmbH & Co. KG
Turnstraße 31 · 75228 Ispringen · Germany

2. *tomas*[®]-pin product range

Description ¹	REF	Thread design	Sterile	Non-sterile
<i>tomas</i> [®] -pin 06	302-006-11	Self-tapping	●	
<i>tomas</i> [®] -pin 08	302-008-00	Self-tapping	●	
<i>tomas</i> [®] -pin 10	302-010-00	Self-tapping	●	
<i>tomas</i> [®] -pin N 06	302-006-10 <p>302-006-30</p>	Self-tapping		●
<i>tomas</i> [®] -pin N 08	302-008-10 <p>302-008-30</p>	Self-tapping		●
<i>tomas</i> [®] -pin N 10	302-010-10 <p>302-010-30</p>	Self-tapping		●
<i>tomas</i> [®] -pin SD 06	302-106-00	Self-drilling	●	
<i>tomas</i> [®] -pin SD 08	302-108-00	Self-drilling	●	
<i>tomas</i> [®] -pin SD 10	302-110-00	Self-drilling	●	
<i>tomas</i> [®] -pin SD-N 06	302-106-10 <p>302-106-30</p>	Self-drilling		●
<i>tomas</i> [®] -pin SD-N 08	302-108-10 <p>302-108-30</p>	Self-drilling		●
<i>tomas</i> [®] -pin SD-N 10	302-110-10 <p>302-110-30</p>	Self-drilling		●

¹ N = Non-sterile; SD = Self-drilling

3. Short description

The *tomas*[®]-pin is an endosteally anchored mini screw that provides temporary, reliable anchorage for orthodontic treatment. The head of the *tomas*[®]-pin can be used with various orthodontic appliances, according to each indication, in order to achieve or support the required tooth movement.

tomas[®] is a system which consists of various matching components to insert the mini screw. The *tomas*[®]-pin is made from titanium Grade 5 according to the ASTM (American Society for Testing and Materials).

4. Indications and contra-indications

The *tomas*[®]-pin is a temporary orthodontic anchor that can be used for the following treatments:

- Space closure in class-I-occlusion
- Distalization and mesialization of teeth
- Alignment of molars
- Intrusion of teeth
- Sliding mechanics in Class II
- Preventing the protrusion of incisors
- Oligodontia
- Changing the tooth position in preprosthetic treatment

General contra-indications

The *tomas*[®]-pin shall not be used if the patient has a history of immune deficiency, steroid therapy, problems with blood clotting, uncontrolled endocrinological disease, rheumatic disease, bone disease, cirrhosis of the liver or any other acute disease.

Local contra-indications:

The *tomas*[®]-pin shall not be used if the patient suffers from osteomyelitis, receives radio therapy of the head, has receding gingival disease or unsatisfactory oral hygiene. Furthermore the *tomas*[®]-pin must not be used if there is insufficient bone structure or possible bone defects in the area in which the screw is to be inserted.

5. Availability, safety issues and guarantees

The *tomas*[®]-pin shall only be inserted by orthodontists, dentists periodontists, oral surgeons and facial surgeons. The user shall ensure to have read and understood the instructions for use and the *tomas*[®] insertion manual (REF 989-631-20) supplied with the product.

It is recommended that the user takes part in a *tomas*[®] training course as only a few of the many possible applications are described in the written instructions for use.

Only original *tomas*[®] components shall be used according to the instructions for use.

Each patient must be examined and informed about the *tomas*[®]-pins before use.

The *tomas*[®]-pins are designed for single use only. Recycling or reuse of *tomas*[®]-pins on the patient is not permitted.

6. Form of presentation

tomas[®]-pins are supplied sterile and non-sterile; refer to the product range in Section 2.

All the other components apart from the *tomas*[®]-drill 1.1 SP, long pilot drill (REF 302-103-11), *tomas*[®]-drill SD 1.1 SP, short pilot drill (REF 302-103-40) and the *tomas*[®]-punch (see the *tomas*[®] insertion manual, REF 989-631-20) are delivered non-sterile and shall be sterilized prior to use.

6.1. Sterile pins (*tomas*[®]-pin and *tomas*[®]-pin SD)

The *tomas*[®]-pins are delivered sterile, wrapped in blister foil packaging in a glass ampoule. The packaging shall be removed just prior to insertion.

Dentaurum guarantees the sterility of the *tomas*[®]-pins in their undamaged, original packaging, up to the expiration date stated on the outside. Once the time limit has been exceeded, the sterility can no longer be guaranteed and the *tomas*[®]-pin shall no longer be used on a patient.

If the sterile packaging is damaged, do not re-sterilize. **Once opened, components shall not be re-sterilized!** *tomas*[®]-pins that have been removed from the packaging but not inserted shall also not be reused or re-sterilized.

6.2. Non-sterile pins (*tomas*[®]-pin N and *tomas*[®]-pin SD-N)

These *tomas*[®]-pins shall be cleaned, disinfected and sterilized before use.

Non-sterile *tomas*[®]-pins are cleaned, disinfected and packed germ-free and particle-free by the manufacturer. Prior to clinical use, non-sterile *tomas*[®]-pins shall again be cleaned and disinfected and then sterilized by the operator. Throughout the entire procedure, from removing the pin from the packaging to insertion, the *tomas*[®]-pin thread shall not be touched or contaminated in any way. Wear sterile gloves and always grip the head of the *tomas*[®]-pin with sterile tweezers. *tomas*[®]-pins shall be cleaned in an ultrasonic unit. A 5 % n-propanol solution (made up of 70 % n-propanol and deionised water) shall be used as cleaning agent. This shall be degassed for 5 min (timed on the ultrasonic unit) before the actual cleaning process.

Remove the *tomas*[®]-pin from the packaging and place it in the basket of the ultrasonic cleaner. Clean it in the ultrasonic unit for 8 min at 25 °C/77 °F. Then remove the cleaning basket. Allow the cleaning liquid to drain from the *tomas*[®]-pin and the basket.

The pin shall be disinfected with 70 % n-propanol. Insert the cleaning basket for 20 min in the disinfectant (room temperature). Then allow the *tomas*[®]-pin to dry for 15 min at room temperature.

The cleaning agent and disinfectant shall be regularly renewed according to the manufacturer's instructions.

The dry *tomas*[®]-pin can be sterilized separately or together with other *tomas*[®]-pin components in the *tomas*[®]-tray/-tool tray. Place the *tomas*[®]-pin in one of the empty inserts in the tray for sterilization.

Shrink-wrap each *tomas*[®]-pin in a similar sized foil bag for separate sterilization. Pins shall be sterilized in an autoclave. Adhere strictly to the autoclave manufacturer's instructions for use. Steam sterilization shall be completed at a temperature of 134 °C/273 °F with a minimum hold time of 25 min.

7. Insertion of the *tomas*[®]-pin

7.1. General instructions

When inserting the *tomas*[®]-pin, ensure that all the hygienic measures required for invasive surgery are completed, e.g. sterile working area, sterile gloves, face mask etc.

Reliable functioning of the *tomas*[®]-pin depends on rigid anchorage in the bone (primary stability) and placing the head in the region of the attached gingiva (alveolar gingiva). When using the *tomas*[®]-pin as an anchor, ensure that the head and surrounding tissue are not subjected to any detrimental mechanical effects (e.g. movement of the mucosa, effect of bands and/or tongue, manipulation).

The *tomas*[®]-pin can be inserted in the maxilla buccally or on the alveolar ridge. In the mandible the pin shall always be inserted buccally and on the centre of the alveolar ridge. According to the latest findings, insertion of the *tomas*[®]-pin on the lingual of the mandible or near extraction wounds, tooth follicles or deciduous teeth is not advisable.

7.2. Procedure (general)

- Select the insertion area. **Immediate load force required at 90 degrees or**

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