Gebrauchsanweisung | Instructions for use
Mode d'emploi | Modo de empleo
Modalità d'uso

D
DENTAURUM

- Informationen zu Produkten finden Sie unter www.dentaurum.com
- For more information on our products, please visit www.dentaurum.com
- Vous trouverez toutes les informations sur nos produits sur www.dentaurum.com
- Para más información sobre nuestros productos, visite www.dentaurum.com
- Informazioni sui prodotti sono disponibili nel sito www.dentaurum.com

Stand der Information | Date of information | Mise à jour |
Fecha de la información | Data dell'informazione: 2023-12

Änderungen vorbehalten | Subject to modifications |
Sous réserve de modifications | Reservado el derecho de modificación |
Con riserva di apportare modifiche

CE 0483

D
DENTAURUM
Turnstr. 31 | 75228 Ispringen | Germany | Tel. +49 72 31 / 803-0 | Fax +49 72 31 / 803-295
www.dentaurum.com | info@dentaurum.com

Printed by Dentaurum Germany 2023-12/BR1-2

989-609-50

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EN

Dear customer,

Thank you for choosing a quality product from Dentaurum. Please read these instructions for use to ensure you use this product in a safe and simple way and that you and your patients gain full benefit. Should you have any questions or ideas, please contact your local representative.

1. Manufacturer

Dentaurum GmbH & Co. KG | Turnstr. 31 | 75228 Ispringen | Germany

2. General product description

discovery® pearl is a ceramic bracket manufactured using the CIM process (Ceramic Injection Molding). The bracket is made from ultrapure polycrystalline aluminum oxide. discovery® pearl ceramic brackets are also offered as a product variant with a positioning guide for the standard systems.

3. Intended purpose

For the orthodontic treatment of tooth misalignments.

4. Intended users

All products are intended for dental professionals only.

4.1 Informing the patient

Inform the patient about the possible consequences and risks of bracket treatment. This applies especially in the case of teeth with enamel that already shows signs of prior damage (cracks, crazing, splitting, etc.).

4.2 Conditioning teeth

Check each tooth thoroughly for any existing damage to the enamel to avoid damaging the tooth during debonding.

4.3 Selection of bracket adhesive

discovery® pearl brackets can be bonded using all commercially available bracket adhesives.

4.4 Removing and handling

Grip the discovery® pearl brackets mesially-distally or on the bar of the positioning guide using crimping tweezers to remove them from the packaging. The base of the discovery® pearl brackets should not be touched with your fingers or contaminated in any way, as this could otherwise have a detrimental effect on the bond.

4.4.1 Procedure for discovery® pearl brackets without positioning guide

Grip the discovery® pearl brackets mesially-distally using crimping tweezers (e.g. REF 025-277-00) to remove them from the packaging and handle according to the following steps. The base of the brackets should not be touched with your fingers or contaminated in any way, as this could otherwise have a detrimental effect on the bond.

4.4.2 Procedure for discovery® pearl brackets with positioning guide

Grip the discovery® pearl brackets mesially-distally or on the bar of the positioning guide using crimping tweezers to remove them from the packaging. The base of the discovery® pearl brackets should not be touched with your fingers or contaminated in any way, as this could otherwise have a detrimental effect on the bond.

4.5 Using the bracket adhesive

Prepare the teeth and the bracket adhesive according to the manufacturer's instructions. Apply a small amount of adhesive to the center of the bracket base, enough to wet the entire bracket base. This guarantees an optimum bond.

4.6 Indirect bonding technique

discovery® pearl brackets are bonded on the plaster or the plastic model with the bracket adhesive. The adhesive layer should not be removed afterwards. Thermoformed splints or silicone transfers can be used as a transfer matrix. Indirect bonding techniques using a sugar solution or other water-soluble laboratory adhesives that later have to be completely removed are not suitable. This would destroy the mechanical retention of discovery® pearl and it could no longer be guaranteed that the bracket will hold on the tooth.

4.7 Placing and aligning brackets

After applying the adhesive, place, align and press the bracket in the usual way onto the tooth. Remove completely any excess adhesive before it has cured.

4.8 Ligating the archwire

We recommend using Tensic®, rematitan® LITE or the translucent pearl archwires at the beginning of treatment. Once the adhesive has cured and the position guides removed, ligate the archwires as usual with wire or plastic ligatures (e.g. Dentalastics® pearl REF 774-563-01).

4.9 Bracket debonding

First, remove the archwire. For debonding the discovery® pearl brackets, we recommend the special discovery® pearl debonding instrument, REF 019-001-00, which grips the body of the bracket directly at the base. The bracket is gently removed from the bracket adhesive by slightly tipping it sideways in the direction of the slot and archwire. Alternatively, the straight bracket removal pliers (REF 004-346-00) can be used for debonding the discovery® pearl brackets. Place the blades of the bracket removal pliers mesio-distally on the bracket base, directly on the surface of the tooth. The bracket is gently removed from the tooth by slightly twisting the pliers.

Note:

Despite their small, delicate shape there is no risk of fracture with discovery® pearl

brackets provided they are used correctly. Use of other instruments can result in bracket fracture. If the bracket fractures during incorrect debonding, we recommend carefully removing the fractured segments and the adhesive from the enamel using bracket removing pliers (e.g. REF 003-349-00) or a scaler (e.g. REF 027-349-00). **Avoid damaging the enamel, danger of enamel chipping. Use an aspirator and wear protective glasses.**

5. Composition

discovery® pearl brackets are made of the following material: Al₂O₃. The composition is included in the orthodontic materials list; please refer to the catalog or www.dentaurum.com. The SSCP is available at <https://ec.europa.eu/tools/eudamed> and www.dentaurum.com.

6. Safety information

The product should not be used if there is a known allergic reaction to one or more of the material components. Should the dental professional and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the dental professional and/or the patient is resident are informed accordingly.

7. Storage and shelf life

Store discovery® pearl brackets in the closed packaging. If the packaging shows serious damage, the product must be checked for integrity and cleanliness before use, and if necessary must be disposed of.

8. Information for single use products

The product is intended for single use only. Reconditioning of a product that has already been used (recycling) and its reuse on a patient is not permitted. Micro cracks may result from loading during treatment and in particular bracket debonding that could cause fractures in the ceramic if reused!

9. Notes on identification system (positioning guides)

9.1 discovery® pearl brackets without positioning guide

The brackets can be identified and matched to the correct tooth by means of the FDI marking, e.g. 21, on their base.

9.2 discovery® pearl brackets with positioning guide

- The color of the positioning guide, the tooth symbol and the FDI-marking identify the tooth/bracket combination.
- Maxilla: right = red
left = green
- Mandible: right = yellow
left = blue

- The occlusal-gingival distance can be determined by the occlusal indicator. The distance between the occlusal end and the center of the slot is 6 mm.
- The tooth symbol always points from occlusal to distal.
- Removal of position guides
Once the adhesive has cured, grip the mid-section of the positioning guides using tweezers or Weingart pliers and remove.

10. Further information

- Technical information such as torque, angulation, size etc. can be found on the relevant product label or in the product catalog.
- There are no particular precautions to be taken when disposing of discovery® pearl brackets after treatment.
- Should the dental professional and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the dental professional and/or the patient is resident are informed accordingly..

11. Notes on disposal

When disposing of the product, please adhere to the national regulations that are valid for your region and - if available - the details outlined in the safety data sheets.

12. Product range

For the complete product range of discovery® pearl brackets, please see our Orthodontics Catalog.

13. Quality information

Dentaurum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaurum has no influence on the processing on the site.

14. Explanation of symbols used on the label

 Please refer to the label. These instructions for use and additional information can be found at www.dentaurum.com (explanation of symbols REF 989-313-00).

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