Platorit®

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Instructions for use

Mode d'emploi

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Instructions for use Platorit®

Special investment material for precious metals.

Platorit® is a gypsum- and graphite free investment material. It can be used for high melting gold-platinum ceramic bonding alloys as well as all other precious metals.

Platorit® must only be used with the Dentaaurum mixing fluid. The thermal expansion can be controlled by diluting the mixing fluid with distilled water.

- Caution:** Investment contains silica.
Avoid inhalation of dust.
- Muffel preparation:** Platorit® can be used with or without a metal muffle ring.
The metal ring must always be lined with one layer Kera-Vlies®.
- Mixing fluid:** The mixing fluid is used undiluted for high melting refined alloys.
With standard precious metal alloys or for metal-on-metal procedures, the mixing fluid can be used at a 90% concentration (9 parts mixing fluid : 1 part distilled water).
- Mixing proportions:** 150 g (1 packet à 150 g) + 25 ml
300 g (2 packets à 150 g) + 50 ml
450 g (3 packets à 150 g) + 75 ml
- Separation:** Apply Dentaaurum separating fluid, Lubrofilm® REF 112-050-00/
Lubrofilm® plus REF 112-053-00 and remove the surplus.
- Mixing:** Mechanical agitation under vacuum: 60 seconds
Manual stirring in a bowl: at least 60 seconds
- Vacuum investing:** A vacuum investor should be used to achieve greater investment density.
- Working time:** 5-6 minutes. Warm mixing fluid and room temperatures as well as soiled bowls shorten the working time.

DIN EN ISO 15912**Type 1, class 1**

Thermal expansion: 1.35%*

* Liquid 100 %

Setting time:	40 minutes		
Preheating:	<p>Attention: Always place muffle in the burnout furnace with the sprue former down and preheat at a rate of 5 °C/41 °F per minute to 250 °C/482 °F. Holding time 60 minutes at 250 °C/482 °F. Thereafter, slowly heat to 850 °C/1562 °F for high melting gold-platinum ceramic bonding alloys and to 750 °C/1382 °F for standard recious metal alloys. The final temperature depends upon the metal being cast. Hold at final temperature for 30 min (or longer, depending on the muffle size). At the end of the week, always dewax the muffles after investing.</p>		
Deflasking:	Cool the muffle to room temperature. Push the investment out of the ring and separate it from the casting.		
Mixing fluid:	Protect from frost; preferably order supplies before the onset of winter. Store in a cool place. Always close the bottle immediately after use.		
Special note:	Platorit® should never be used with gypsum-containing bowls or scoops. Always use separate measuring flasks an spatulas.		
Ordering:	40 x 150 g packets (6 kg)	REF 106-620-00	
	1000 ml mixing fluid	REF 106-625-00	
	1200 ml mixing liquid are sufficient for 6 kg powder.		
Shelf life:	24 months provided it is stored under dry conditions.		
Lubrofilm®/ Lubrofilm® plus	Wax surface tension reducing medium Lubrofilm®		
	100 ml spray bottle	REF 112-050-00	
	100 ml spray bottle	REF 112-053-00	
	1000 ml refill bottle	REF 112-051-00	
Example:	Batch number 030762	03	07
		month	year
		of production	

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