
# Castorit<sup>®</sup>-super

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Instructions for use

Mode d'emploi

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## Dear Customer

Thank you for having chosen a quality product from Dentaureum.

In order to use this product at its best for you and your patients, it is important to study and follow these directions for use carefully.

The written instructions cannot cover all eventualities during operation. For this reason our Hotline is available to answer any other questions and ideas that may arise.

Due to a constant development we recommend, even when you use the same products frequently, to study the enclosed latest issue of the instructions for use or refer to our website at [www.dentaureum.de](http://www.dentaureum.de).

# Instruction for use Castorit®-super

## Investment Type I, class 1

**Castorit®-super** is an investment free of gypsum and graphite.

## Safety Tips

**Investments contain quartz. Avoid inhalation of dust as this may cause lung diseases (silicosis / lung cancer). Wearing of face mask type FFP 2 - EN 149:2001 is recommended. Open pouch with scissors and avoid dust development when pouring into mixing bowl. Rinse empty pouch with water before disposal.**

**Remove dust from working area only with a wet cloth.**

**To avoid dust when devesting place ring into water until it is completely soaked after ring has reached room temperature.**

**When sandblasting use suction with fine dust filter.**

## Availability

### Castorit®-super powder

6 kg box with 100 x 60 g bags

**REF 106-700-00**

6 kg box with 40 x 150 g bags

**REF 106-701-00**

**Castorit®-super mixing liquid A** 1 l  
(ready-to-use) for standard gold alloys

**REF 106-702-00**

**Castorit®-super mixing liquid B** 1 l  
(ready-to-use) for platinum-gold alloys, palladium and low gold alloys

**REF 106-703-00**

**Castorit®-super mixing liquid C** 1 l  
(ready to use) for palladium/gold and palladium/copper alloys ①

**REF 106-704-00**

1.4 l mixing liquid is sufficient for 6 kg investment.

Kera-Vlies® ring liner 1.0 x 50 mm, 25 m

**REF 127-250-00**

## Storage and Transport Instructions

### Investment Mixing Liquid:

To prevent errors in mixing, we recommend the use of the mixing liquid specified for each group of alloys. To be used at room temperature (18 – 22 °C/64 – 72 °F).

Temperature influences the rate of expansion. Do not subject the liquid to frost! Take care with deliveries in winter. The liquid must not be used once thawed! Mixing liquid is sensitive to low temperature and cannot be delivered in frosty conditions. Please order in good time to avoid delays due to climate conditions. DO NOT STORE BELOW 5 °C/41 °F. Avoid direct exposure to sunlight. Close container immediately after use.

Shelf life exceeds 24 months when in sealed container.

**Investment Powder:** Store in a cool dry place. Reseal bags tightly once opened.

Shelf life exceeds 36 months when in sealed container.

## Using

Alloy	Range of Application	Investment material	Mixing Liquid	Special measures
Standard gold alloys (gold content 75 % +). Semi-precious alloys with a gold content of 50 % +.	Crown and bridge-work, inlays	Castorit®-super	<b>A</b>	Line the casting ring with a 1.0 mm thick, <b>moistened</b> , wool-based, liner (Kera-Vlies®).
Au/Pt – Au/Pd/Cu bonding alloys. Pd/Cu/Au – Pd/Ag alloys and bonding alloys. Semi-precious alloys with a high Ag or Pd content.	Crown and bridge-work	Castorit®-super	<b>B</b>	Place the 1.0 mm thick, <b>moistened</b> "Kera-Vlies®" in the casting ring.
① Pd/Ag Alloys Pd/Cu	Telescopic crowns	Castorit®-super	② <b>C+</b>	Place the 1.0 m thick, <b>moistened</b> "Kera-Vlies®" in the casting ring. More viscous mix, 150 g : 34 ml.
			③ <b>B++</b>	Hygroscopic investing. Line the ring with 1.0 mm thick, <b>moistened</b> "Kera-Vlies®".

+

++ Increased expansion

① High palladium content semi-precious alloys, used for telescopic crowns etc., require an investment material with a very high rate of expansion.

Avoid producing non stress-relieved plastic caps, as this would impede the setting expansion of the investment material. The expansion can be increased noticeably by:

- (a) using wax or stress-relieved plastic caps (cut a groove),
- (b) employing the hygroscopic investing method ③ or
- (c) using the mixing fluid C with altered mixing ratio ②.

### Further increase of expansion

Change from mixing liquid **A** to **B** or use Castorit®-super C with mixing liquid **C**.

### Reducing the expansion – increased friction –

By diluting the mixing liquids **A** and **B**.

Castorit®-super + mixing liquid **A** diluted with 10 – 20 % distilled water.

Castorit®-super + mixing liquid **B** diluted with 10 – 20 % distilled water.

**Caution:** Never mix **Castorit®-super** in plaster contaminated mixing bowls or scoops. Use separate mixing bowls and spatulas.

### Preparing the casting ring

The metal casting ring should be lined with Kera-Vlies®, 1 mm thick. The preproportioned packages of 60 and 150 g are suitable for the standard casting ring sizes.

### Mixing ratio

The prepackaged sizes of 60 and 150 g correspond to the standard casting ring sizes so that there is no waste of material during use.

60 g (1 bag with 60 g)	+	14 ml
150 g (1 bag with 150 g)	+	35 ml
300 g (2 bags with 150 g)	+	70 ml
450 g (3 bags with 150 g)	+	105 ml

### Reducing surface tension

Use Dentaurem's Lubrofilm® (100 ml), REF 112-050-00 to reduce surface tension. Blow dry prior to investing.

### Mixing Time

Vacuum mixing: 60 seconds

Hand mixing in a bowl: 60 seconds minimum

Vacuum investing is recommended for all castings. For heavy flowing alloys, please apply air vents.

## Working time

4 – 5 minutes. Temperature of the mixing fluid 18 – 22 °C/64 – 72 °F. If a warmer fluid or powder is used, the working time is shorter. Colder fluid or duplicating material lengthens the setting time. It is advisable to store the mixing fluid in a refrigerator on very hot days but do not freeze.

Mixing speed and shape of the mixing bowl influence the mixture and therefore the rate of expansion. Ideal mixing speed approx. 360 – 400 min<sup>-1</sup> (Airvac).

First add the liquid to the mixing bowl, then the powder, stir vigorously by hand. Ensure the mixing bowl is clean! No plaster residues!

**Warning:** The investment will set rapidly after the working time has elapsed. Setting temperature is about 60 °C/140 °F. Do not alter mixing ratio when room temperature is higher. Keep mixing liquid in refrigerator.

## Setting time

40 minutes

**Attention:** Always place muffle in the burnout furnace with the sprue former down.

## Preheating

### Procedures for Castorit®-super

1. Place muffle in cold burn out furnace and preheat at a rate of 5 °C/41 °F per minute to 250 °C/482 °F. Holding time 60 minutes at 250 °C/482 °F. Thereafter, heat at a constant rate to the final temperature.

Final temperatures:

750 °C/1382 °F Standard precious metal alloys.

850 – 900°C/1562 – 1652 °F Pt/Au ceramic bonding alloys, low golds and palladium.

Hold at final temperature for 30 minutes (depending upon the size of the casting ring).

When working over the weekend, dewax immediately after the investment material has set (see above).

2. **Transfer technique:** This can only be used for the following sizes of casting ring:

1: ø 30 mm

3: ø 48 mm

6: ø 65 mm

After 40 min, when the investment has set, place the casting ring in a furnace preheated to 250 °C/482 °F. Hold at temperature for 1 hour and then place into a furnace preheated to 500 °C/932 °F. Heat soak for 60 minutes, then place into a furnace preheated to the final temperature (700, 850, 900 °C/1292, 1562, 1652 °F, see 1). Heat soak for 30 minutes and cast.

## Manipulation of fit

An additional 0.2–0.3 % linear increase in the maximum expansion is possible with the following procedure. Place the casting ring in the cold furnace, preheat to 250 °C/482 °F, hold for 1 hour, cool to room temperature, then heat at a constant rate to the final temperature.

**Warning:** too quick preheating (less than 2 hours) will cause a reduction in expansion – i. e. too tight castings.

Build ups in wax allow a higher setting expansion of the investment than those made of pattern resin. Build ups made of resin result in tighter castings.

## Casting process

Vacuum-pressure casting or centrifugal casting

## Devesting

Cool the casting ring to room temperature. Push the investment out of the casting ring. Remove all investment from the casting using instruments and brushes.

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