Paratherm

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hereby declares that the design and construction of the laboratory equipment described below, including the version marketed by us, comply with the basic regulations governing safety and health as stated in the EEC guidelines. This declaration will become invalid if the laboratory equipment is modified or altered in any way without our prior consent.

Description of the unit	Surveying, blocking-out and waxing unit	
REF	094-270-00	
Type	Paratherm	
Start with unit No.	32-080	
EEC guidelines	89/336/EEC	Electromagnetic compatibility
	92/031/EEC	Modified guideline
	73/023/EEC	Electrical equipment used within certain voltage limits
	93/068/EEC	Modifications of the above guidelines
Applied unified standards	EN 55011	
	EN 50082	
	EN 61010	
Date/manufacture's signature	14.07.2005
Signatory		- p. p. Dipl.-Ing.(FH) K.Merkle - Production Manager Mechanic

DENTAUURUM
GmbH & Co. KG
Turnstraße 31 - 75228 Ispringen - Germany
Telefon +49 7231 / 803 - 0 - Fax +49 7231 / 803 - 295

1.1 General information

Thank you for having chosen a Dentaureum unit. Please study these instructions for use carefully so that you may enjoy the benefits of using this Paratherm unit.

2 Safety instructions

⚠ Caution Read through the instructions for installing and operating the unit carefully before switching it on.

2.1 Disposal

When a Dentaureum product is finally no longer in use the disposal regulations of that particular country apply. Dentaureum or the dental trade is available to answer questions regarding the correct disposal of any specific product.

Note:

This symbol indicates, products marked with it should not be disposed off together with garbage of households. The legislator disallows commercial customers to return electronic waste via municipal collection points. Further information is available through Dentaureum or the dental trade.

2.2 Delivery

- 1 Paratherm unit with surveying / blocking-out arm
- 1 Analyzing rod (A)
- 1 Graphite holder (B)
- 10 Replacement graphite markers (C)
- 4 Undercut gauges 0.25 / 0.35 / 0.50 / 0.75 mm (D)
- 2 Thermal blocking-out instruments 0°, 2° (E)
- 1 Plastic cover
- 1 Fuse T 0.25 A
- 1 Instructions for use with EEC-Declaration of conformity

3 Areas of use

Cast partial denture, combination techniques, conical crowns, positioning of snap attachments and attachments.

4 Installation

Set the necessary voltage on the underside of the unit (110 / 220 / 240 V). Then connect the unit to the mains supply. Stand the Paratherm on a level working area.

5 Operating the Paratherm

Mount the master model on the model table already available or supplied with the unit. Loosen the knurled screw ② and raise the sliding rods to the necessary level. Fasten with retainer ring ①. Grip the instrument eg. analyzing rod in the surveying handpiece ④ and push the selected blocking-out instrument into the thermal blocking-out handpiece ⑤.

Place the model table ⑥ in the most convenient position on the work surface ⑦. Loosen the locking screws ③, so that the arms may be moved freely. Switch on the thermal blocking-out handpiece ⑧ and set the desired temperature (maximum 160°C / 320°F) on the controller ⑨.

Caution! Only change the blocking-out instruments ⑩ whilst they are still warm.

5.1 Electric wax knife

The knife is connected by means of the socket on the rear side of the unit. Only operate the wax-knife with an instrument in the place.

Set the desired temperature on the potentiometer ⑪. Maximum temperature 280°C / 536°F.


Caution! Only change the knife whilst it is still warm.

5.2 Locking the surveying arm (Positioning attachments)

Rotate the surveying handpiece ④ forwards of the centre of the work surface and lock into the position with screws ② ③.

Grip the wire aid in the surveying arm and fix the level using ③. With the model already clamped into position, guide the model table onto the wire aid. Correct the level.

When ordering spare parts please quote the unit number (type plate or rear side of the unit).


6 Spare parts

1 Thermal blocking-out instrument	0°	REF 094-275-00	(E)
1 Thermal blocking-out instrument	2°	REF 094-277-00	(E)
1 Undercut gauge	0.25 mm	REF 094-283-00	(D)
1 Undercut gauge	0.35 mm	REF 094-284-00	(D)
1 Undercut gauge	0.50 mm	REF 094-285-00	(D)
1 Undercut gauge	0.75 mm	REF 094-286-00	(D)
1 Analyzing rod		REF 094-287-00	(A)
1 Graphite holder		REF 094-288-00	(B)
10 Replacement graphite markers		REF 094-289-00	(C)

7 Recommended accessories

(Not included in original kit REF 094-270-00)

1 Thermal blocking-out instrument	1°	REF 094-276-00	(E)
1 Thermal blocking-out instrument	3°	REF 094-278-00	(E)
1 Thermal blocking-out instrument	4°	REF 094-279-00	(E)
1 Thermal blocking-out instrument	6°	REF 094-280-00	(E)
1 Thermal blocking-out instrument	8°	REF 094-281-00	(E)

Model table **REF 086-340-00**

Electric wax knife **REF 094-272-00**

Handpiece, cable and 1 wax knife
In order to regulate the temperature a transformer and potentiometer have been built into the Paratherm (REF 094-270-00)

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