



Thermomat

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Electric wax knife with plug-in transformer

This economical modelling unit, specially developed for dental applications, will occupy little space on any work bench having a wall socket.

Wax knife and points

Manufactured from a special alloy which prevents heat distortion of the knife up to 400 °C/752 °F.

Transformer

Heavy duty plug-in transformer, adjustable to 280 °C/536 °F. Extremely simple wax liquefaction. Mains indicator. 230 volt or 110 volt operation. Dimensions: 72 x 95 x 42 mm.

Handpiece

Insulated, anatomically-shaped plastic handle with practical, non rotating socket for the wax knife or modelling points. Flexible lead (Length 95 cm). Length without knife 130 mm, ø 13 mm.

Ordering

Handpiece with lead and wax knife, adjustable transformer 230 V **REF 089-620-00**

Handpiece with lead and wax knife, adjustable transformer 110 V, with plug converter **REF 089-611-00**

Recommended accessories

Wax knife, small (supplied with unit) **REF 089-602-00**

Modelling point, fine, A **REF 089-603-00**

Modelling point, coarse, B **REF 089-604-00**

Modelling spatula, C **REF 089-606-00**

Wax knife, D **REF 089-605-00**

Instructions

Connect to correct voltage (230 V or 110 V). Insert the transformer into a plug next to your work area. For 110 V, use the plug converter. Connect the pin of the handpiece with the plug-in transformer. Turn on the unit by pushing down the switch. The green LED is illuminated. Set the desired temperature with the temperature regulator. Exchange the modelling tips only at lowest possible handpiece temperatures. Remove the tip with an appropriate tool. When the new tip is inserted, make sure it locks into position.

Attention: Do not touch wax knife and tip of instrument. Hot tip may cause burns.

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