
Gebrauchsanweisung | Instructions for use |
Mode d'emploi | Modo de empleo | Modalità d'uso

CE 0483

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Dear customer,

Thank you for choosing a quality product from Dentaureum. Please read these Instructions for use to ensure you use this product in a safe and simple way and that you and your patients gain full benefit.

Should you have any questions or ideas, please contact your local representative.

1. Manufacturer

Dentaureum GmbH & Co. KG | Turnstr. 31 | 75228 Ispringen | Germany

2. General device description

amda® is a skeletal anchorage device for the distalization and/or mesialization of upper molars. The device is based on a sliding mechanism. The active element is a spring encapsulated in a telescope element (amda® telescope ; Fig. 1, ❶). The maximum spring travel is 4 mm and the maximum spring force is 4.7 N (approx. 470 g). The anterior and posterior ends of the telescopic tube are fitted with stop screws (identical in construction to the tomas®-stop screw, REF 302-013-01). The amda® telescope and the palatal sheath on the molars are connected with a removable connector (amda® connector; Fig. 1, ❷ and ❸). The amda® connector is available in 3 lengths, each for the right and the left side (see amda® product overview).

The amda® telescopes can be pushed onto the amda® palatal arch (Fig. 1, ❹). The stop screws are used to fix the telescopes in the desired position. The anterior arch of the amda® palatal arch is equipped with 2 tomas®-abutment EP (REF 302-027-00; Fig. 1, ❺). The abutments are used to connect amda® to the tomas®-pin EP (not included in the delivery; Fig. 1, ❻), achieving the desired skeletal anchorage. Alternatively, the anterior part of the amda® palatal arch can be embedded in a Nance acrylic pad.

3. Intended purpose

Skeletal anchorage device for the distalization of upper molars.

4. Intended users

All products are intended for dental professionals only.

5. Indications

amda® can be used for bi- or unilateral distalization or mesialization of upper molars depending on how it is fitted. Find more details in these Instructions for use.


Fig. 1: amda® components

6. Contraindications

The amda® appliance must not be used for patients with a nickel allergy. Warnings are found on the respective product labels. Also see 9. Composition.

7. Patient target group

All patients for whom a fixed appliance can be placed as part of orthodontic treatment.

8. Processing steps

8.1 Adaptation and incorporation

To adapt and incorporate amda®, follow these steps:

- Inserting 2 tomas®-pin EP
- Fabricating a model
- Adapt amda® to model
- Inserting and activating amda®
- Monitoring clinical progress
- Removing amda®

8.1.1 Inserting of the tomas®-pin EP

amda® is skeletally anchored in the palate by 2 tomas®-pin EP. The choice of tomas®-pin EP length depends on the bone availability and the gingival thickness. Most cases can be treated with 8 mm tomas®-pin EP (REF 302-208-00). The tomas®-pin EP is available in two more lengths: 10 mm (REF 302-210-00) and 6 mm (REF 302-206-00).

The bone availability can be measured using X-rays. The midline and the third palatine ruga can help to determine the correct insertion point. Trace an imaginary transverse line between the distal contact points of the upper cuspids. If the cuspids are nonexistent or mispositioned, use the third palatine ruga as orientation. The tomas®-pin EP should be inserted 3 to 6 millimeters posterior to the incisive foramen. Position the two tomas®-pin EP 3 to 4 millimeters left and right of the midline, as parallel as possible and in the same anterior-posterior plane. The pin heads should sit at the same height. Insert the two tomas®-pin EP as perpendicular as possible to the plane of the anterior palate. Hold the insertion instrument at an inclined angle to the occlusal plane during insertion. Align the tips of the tomas®-pin EP with the posterior region of the root tips of the upper incisors (Fig. 2). Base the selection of tomas®-pin EP length and insertion angle on the bone availability and the space needed for a later retraction of the upper incisors. We recommend the longer version of the tomas®-driver (REF 302-004-60) for the insertion. Find more details on inserting tomas®-pin EP in the tomas®-pin Instructions for use and in the tomas® manual.


Fig. 2: Insertion angle of the tomas®-pin EP (8 mm) at the anterior palate

8.1.2 Fabricating a model

Adapt bands (e. g. dentaform® Snap or dentaform®) with conventional palatal sheaths to the upper first molars. Put the selected bands on the molars. Put 2 tomas®-transfer cap (REF 302-028-00) on the tomas®-pin EP. The caps may require trimming (see tomas® manual).

Transferring the position of the tomas®-pin EP and the bands correctly can only be guaranteed with a distortion-resistant impression tray and a silicone or polyether impression material. Using alginate or plastic impression trays can lead to inaccuracies. To fix the tomas®-transfer caps and bands securely in the impression, the transfer cap and the palatal sheaths must be surrounded with a light bodied impression material. Then, make the impression of the jaw. After the impression material has hardened completely, remove the impression tray from the mouth carefully. It is important to observe the insertion path of the tomas®-pin EP. The tomas®-transfer caps and the bands should be in the impression. Check if the position is correct and adjust, if necessary. If the parts can be repositioned accurately, the impression does not have to be repeated.

After disinfecting the impression, you can start making the working cast. Before pouring the model material into the impression, put the tomas®-laboratory pin EP (REF 302-029-01) into the tomas®-transfer caps. The tomas®-laboratory pins must sit tightly in the transfer caps. In order for the tomas®-laboratory pins not to loosen from the tomas®-transfer caps during the impression process, it is important to fix the pins with wax. If necessary, block out the bands with wax. Then, make the model.

After the model material has hardened, carefully remove the impression from the model. The bands should sit on the model. If this is not the case, put the bands back in the correct position. Remove the tomas®-transfer caps.

8.1.3 Adapting amda® to model

8.1.3.1 Adapting the amda® palatal arch

First, remove both amda® telescopes from the amda® palatal arch. Loosen the stop screws and slide the amda® telescopes off distally. Hang the eyelets of the tomas®-abutment EP over the heads of the tomas®-laboratory pin EP.

The anterior arch of the amda® palatal arch should be parallel to the arch of the anterior palate. Position the long ends of the amda® palatal arch parallel to the occlusal plane at the height of the trifurcation of the first molars (Fig. 3). This way, the width of the amda® palatal arch is 6

adjusted to the width of the palate. Bend the anterior, arched part of the amda® palatal arch using Aderer 3-Prong pliers Maxi (REF 003-202-00). Always use two pliers to make adjustments and bends. Do not bend or adjust with your hands. Do not bend the posterior, straight part of the amda® palatal arch. Depending on the length of the palate and the required distalization distance of the molars (plus a 5 mm buffer), the amda® palatal arch must be shortened. Once amda® has been placed in the patient's mouth, you can accurately measure and shorten the posterior part of the amda® palatal arch.

Fig. 3: Location of the amda® palatal arch

8.1.3.2 Aligning the amda® telescopes

Slide the amda® telescopes with the amda® connector onto the amda® palatal arch. The larger telescope tube should face posterior. Move the amda® telescopes on the amda® palatal arch until the fixing elements of the amda® connectors are at the height of the palatal sheaths of the molar bands (Fig. 4). Align the fixing screws of the amda® telescopes towards occlusal. Tighten the posterior screw so that the amda® telescope can no longer move on the amda® palatal arch. This screw remains tightened for now.

8.1.3.3 Adjusting the amda® connectors

First, check the vertical distance between the amda® telescope and the molar palatal sheath. The amda® connector element that fits into the palatal sheath should vertically be at the height of the palatal sheath and horizontally parallel to the tube of the sheath. The amda® starter set is delivered with mounted amda® connectors (10 mm). If the distance between the parallel part of the amda® palatal arch and the palatal sheath is larger or smaller, you can use the 13 mm or 7 mm amda® connector. They are included in the starter set. For fine corrections the amda® connector can be removed from the retaining tube on the amda® telescope. To achieve the desired height and angulation of the insertion element, bend the S-shaped part of the amda® connector with two pliers. To achieve the desired rotation/derotation and uprighting of

Fig. 4: The anterior-posterior orientation of the amda® telescope

the molars, bend the insertion element in the desired direction with two pliers. Depending on the type and scope of the correction, insert the amda® connector back into the amda® telescope and check the position. At this point (or later), adjust or correct the amda® connector to upright or rotate the molars.

After the amda® connectors have been aligned, insert the connectors back into the retaining tube on the amda® telescope (Fig. 5). The distalization appliance can now be fitted into the patient's mouth.

Fig. 5: amda® ready to be fitted into in the patient's mouth

8.1.4 Inserting and activating amda®

8.1.4.1 Inserting amda®

First, cement the molar bands. First, cement the molar bands. After the cement has hardened and the cement excess has been removed, amda® can be fitted.

Start by inserting the amda® connectors in the palatal sheath on the molars. Then, loosen the posterior stop screws and move the amda® palatal arch until the eyelets of the tomas®-abutment EP are on top of the mushroom-shaped heads of the tomas®-pin EP. Align the abutments so that they are parallel to each other at the anterior part of the amda® palatal arch. Then, pull the amda® palatal arch back towards anterior until the eyelets of the tomas®-abutment EP sit tightly below the

mushroom-shaped heads of the tomas®-pin EP. Tighten the distal screws of the amda® telescopes.

Check the position and orientation of the amda® palatal arch and the amda® telescopes. If fine corrections are necessary, remove the appliance from the mouth. Perform the extraoral corrections as described before. During fine corrections, check the posterior length of the amda® palatal arch. It should be 5 mm longer than the actual distalization distance. Do not limit the mobility of the soft palate or traumatize the latter. If the amda® palatal arch needs to be shortened, mark the corresponding spots. Shorten the arch extra-orally. It can already be shortened in the laboratory. To prevent injuries, surround the distal ends of the amda® palatal arch with a drop of light-curing adhesive.

Then, insert amda® again. Start with the amda® connectors and insert them in the palatal sheaths on the molar bands. Loosen the posterior screws of the amda® telescopes. Then, place the tomas®-abutment EP on the tomas®-pin EP. Move the palatal arch towards anterior until the abutments sit tightly below the mushroom-shaped heads of the tomas®-pin EP. Tighten the posterior screws of the amda® telescopes.

We recommend covering the tomas®-pin EP and the tomas®-abutment EP with a light-curing adhesive, so that these elements are firmly attached to the amda® palatal arch.

8.1.4.2 Activating amda®

Activate amda® one side at a time. Position mosquito forceps in front of the anterior stop screws on the amda® palatal arch. This makes activation easier. Depending on the desired initial force (see Fig. 6), move the anterior stop screw with the inner telescopic tube using mosquito forceps towards posterior and tighten the screw. Then, loosen the posterior stop screw and remove the mosquito forceps. This step creates the distalization force and activates the first side.

Repeat this procedure for the other side. If needed, you can adjust a different amount in force on each side. The maximum spring travel is 4 mm and the maximum spring force is 4.7 N (approx. 470 g). The force that acts upon the tooth is determined by different factors (among others friction of the amda® telescopes). The force is always lower than the maximal spring force and can be measured with appropriate instruments (e.g. measuring gauge REF 040-712-00).

Double-check one last time.

8.2 Monitoring clinical progress

Do the first checkup a few days after fitting the appliance into the patient's mouth. Do regular checkups every 4 to 6 weeks. Generally, no movement can be seen in the first weeks. Check if the amda® palatal arch remains firmly attached to thetomas®-pin EP and the anterior stop screws.

As distalization progresses, the force of the spring in the amda® telescope decreases (Fig. 6). The reduced force is clearly shown by a larger distance between the anterior stop screw and the telescope and/or by a larger distance between the anterior and the posterior stop screw (Fig. 7). If the spring in the amda® telescope needs to be activated again, loosen the anterior stop screw and fix mosquito forceps in front of it. Then, push the screw towards posterior using the mosquito forceps until the desired force is achieved (Fig. 6). Afterwards, tighten the anterior stop screw again.


Fig. 6: Force-travel diagram of an amda® telescope

8.2.1 Retracting maxillary anterior teeth

If the desired distalization is achieved, amda® can be used for indirect skeletal anchorage of the molars, e.g. to retract the upper incisors. Loosen the anterior stop screw, push it towards posterior using mosquito forceps and tighten it again. Then, tighten the distal stop screws. Next, bond brackets to the upper jaw teeth. After the leveling phase, the upper incisors can be retracted.

8.2.2 Retention phase

If the desired distalization of the molars and retraction of the upper incisors is achieved, amda® can be used for retention. To do this, tighten both the anterior and the posterior stop screws. If, after the molars have been distalized and the upper incisors have been retracted, the molars show a mesial rotation, they can be derotated during this phase. To do this, bend the amda® connector extra-orally with two pliers in the desired direction and remount them.

Alternatively, the amda® can be replaced with a conventional transpalatal arch (TPA). The TPA can be used for the rotation of the molars and then for retention.


Fig. 7: After distalization

8.3 Removing amda®

Once the treatment is finalized, amda® can be removed from the mouth. To do this, loosen all stop screws from the amda® telescopes and pull the amda® connectors out of the palatal sheaths. Then, loosen the tomas®- abutment EP from the tomas®-pin EP. If the abutments were secured with a light-curing adhesive, remove the adhesive with Weingart pliers first. The appliance can now be taken out of the mouth. Depending on further treatment plans, the tomas®-pin EP can remain in situ or be removed (see tomas®-pin Instructions for use).

8.4 Additional amda® indications

8.4.1 Unilateral distalization with amda®

amda® can be used for unilateral distalization. To do this, prepare and incorporate the appliance as described before. Activate the side to be distalized as described. On the other side, mount the amda® telescope on the amda® palatal arch in a way that it does not exert any force. To fix it, tighten up both stop screws.

8.4.2 Mesialization with amda®

amda® was designed for distalization. However, the appliance can also be used for bi- or unilateral mesialization (Fig. 8). In addition, amda® can be used for the distalization of one jaw side and mesialization of the other side.


Fig. 8: Mesialization with amda®

For mesialization, turn the amda® telescope 180° and place it on the amda® palatal arch again. The stop screw with the thinner telescope tube is then on the posterior end. The right amda® connector must be replaced by the left connector and vice versa. Follow the same steps to adapt, incorporate, monitor, and remove amda® as above. To activate mesialization, tighten the posterior stop screw after activating it and loosen the anterior stop screw (now on the larger telescopic tube).

9. Composition

amda® is made of the following materials: DIN 1.4301, DIN 1.4305 and DIN 1.4310. The material composition is included in the orthodontics materials list; please refer to the catalog or www.dentaurum.com. The SSCP is available at <https://ec.europa.eu/tools/eudamed> and www.dentaurum.com.

10. Safety information

The force of the amda® telescope must be adjusted so that no molar root resorption can occur. The product should not be used if there is a known allergic reaction to one or more of the material components. If the packaging shows serious damage, the product must be checked for integrity and cleanliness before use, and if necessary must be disposed of.

11. Information for single use products

The product is intended for single use only with the exception of the hexagon socket keys. Reconditioning (recycling) of the product and its reuse is not permitted, as perfect functionality of the product cannot be guaranteed. The Processing Instructions Instruments and Accessories (REF 989-801-09) are available at www.dentaurum.com.

12. Further information

Should the dental professional and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the dental professional and/or the patient is resident are informed accordingly.

13. Notes on disposal

When disposing of the product, please adhere to the national regulations that are valid for your region and - if available - the details outlined in the safety data sheets.

14. Product range

For the complete product range of amda®, please see our orthodontics catalog or www.dentaurum.com.

15. Quality information

Dentaurum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as we, Dentaurum, have no influence on the processing on site.

16. Explanation of symbols used on the label

Please refer to the label. These Instructions for use and additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00).

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