## remanium® 2000+

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### EN

## remanium® 2000+

**CoCrMoW - Dental alloy based on cobalt for metal-ceramic bonding, type 5 acc. to DIN EN ISO 22674 / DIN EN ISO 9693 for melting and casting**

### Composition (% by mass)

Co	Cr	Mo	W	Si
63.0	23.0	7.0	5.0	1.5

Additional elements < 1%: Mn, N

This product contains cobalt.

Free from nickel, beryllium, lead, cadmium and iron.

### Technical data

Yield strength*	R <sub>p,0,2</sub>	700 MPa
Tensile strength*	R <sub>m</sub>	900 MPa
Hardness*	H	340 HV10
Elongation at rupture*	A <sub>s</sub>	4.5%
Modulus of elasticity*	E	200 GPa
Density	ρ	8.6 g/cm <sup>3</sup>
Solidus temperature	T <sub>s</sub>	1,290 °C / 2354 °F
Liquidus temperature	T <sub>l</sub>	1,415 °C / 2579 °F
Pouring temperature (recommended)	T <sub>g</sub>	1,515 °C / 2759 °F
CTE (25 °C - 500 °C / 77 °F - 932 °F)	α	14.0 x 10 <sup>-6</sup> K <sup>-1</sup>
Weight per casting	m	approx. 6 g

\* Variations of ± 10 % are possible depending on the batch.

### Dear customer,

Thank you for choosing a quality product from Dentaurum. It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative. As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at [www.dentaurum.com](http://www.dentaurum.com), even though you may frequently use the same product.

### 1. Manufacturer

Dentaurum GmbH & Co. KG | Turnstraße 31 | 75228 Ispringen | Germany

### 2. General product description

Non-precious alloy for the crown and bridge technique.

### 3. Area of application

remanium® 2000+ is a cobalt-based cast alloy. It is particularly suitable for the production of crowns and bridges, double crowns, implant-supported superstructures, adhesive bridges, secondary parts for cast partial dentures (combination work) and for metal ceramics.

### 4. Wax-up

Construction must be in line with dental regulations in order to fulfil mechanical requirements. Use only materials (e.g. wax, plastic) that burn without residue. Detailed information can be found in the processing instructions.

### 5. Investing

Investment materials for non-precious alloys with high overall expansion must be used, such as e.g. rema® CC, Castorit® all speed, Castorit®-super C, Trivest and rema® TT or rema® CAD/Vest (Dentaurum).

### 6. Melting and casting

Only suitable ceramic-based melting crucibles should be used for remanium® alloys. Care must be taken to ensure that for each alloy separate melting crucibles are used which are not contaminated and/or damaged. The use of new material is a prerequisite for ensuring clear batch traceability. Suitable melting and casting processes are high-frequency melting with vacuum pressure casting or centrifugal casting machines, arc melting with vacuum pressure casting or open melting. After the breakup of the melt, cast off immediately to avoid overheating of the alloy. Following the casting process, the casting ring is carefully removed and should cool slowly in the air.

### 7. Finishing

The finishing of the framework and the surfaces should be carried out by means of finely-toothed hard metal burs, ceramic-bonded grinding points or diamond grinding points.

### 8. Welding

Suitable welding processes such as laser welding and TIG impulse welding produce a solder-free, mechanically strong and corrosion-resistant joint. It is important to observe the framework design, the surface structures, the welding sequence and the welding parameters of each individual laser machine. Suitable welding wires are available as filler material for all remanium® alloys, e.g. CoCr welding wire 0.25 mm (REF 528-215-10), 0.35 mm (REF 528-210-10) and 0.5 mm (REF 528-200-10).

### 9. Soldering

Soldering should be avoided if possible to reduce the number of materials. If, however, it is necessary to solder, please use a suitable solder and flux which are appropriate for the composition and melting temperature of

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