



Processing Instructions ceraMotion® Me

Veneering ceramic for dental alloys



Contents

BASIC LINE / INDIVIDUAL LINE / TOUCH UP

| | |
|-----------------------|---|
| Classification | 4 |
| Framework design | 5 |
| Framework preparation | 6 |
| Firing check | 7 |

BASIC LINE

| | |
|---------------------------------|----|
| Paste opaque | 8 |
| Powder opaque | 10 |
| Layering scheme: Basic layering | 11 |
| Layering | 12 |
| Firing | 14 |
| Correction layering | 15 |
| Processing | 17 |
| Glaze firing | 18 |
| Completion | 19 |

INDIVIDUAL LINE

| | |
|---|----|
| Framework preparation with ceramic shoulder | 20 |
| Opaque | 21 |
| Layering scheme: Individual layering | 23 |

| | |
|-------------------------------------|----|
| Mixing table for shoulder materials | 24 |
| Shoulder | 26 |
| Layering | 28 |
| Correction layering and completion | 33 |
| Completion | 34 |

TOUCH UP

| | |
|-------------------|----|
| Creative firing | 36 |
| Additional firing | 37 |
| Repair | 38 |

BASIC LINE / INDIVIDUAL LINE / TOUCH UP

| | |
|--|-------|
| Firing chart | 40 |
| Physicochemical data | 41 |
| Product overview | 42 |
| Furnace cleaning | 44 |
| Processing tips for ceraMotion® Me ceramic | 45–61 |


Classification

ceraMotion® Me is a veneering ceramic Type I Class 1 (according to DIN EN ISO 6872) for veneering frameworks made of precious metal alloys or non-precious metal alloys.

Indication

Allocation of ceramic to framework material

ceraMotion® veneering ceramics



ceraMotion® Me is suitable for veneer dental alloys with a thermal expansion of $13.9 \cdot 10^{-6} \text{ K}^{-1}$ to $15.1 \cdot 10^{-6} \text{ K}^{-1}$ (25 °C to 500 °C).

ceraMotion® Me must not be used for veneering frameworks made of high-performance ceramics (Al_2O_3 , ZrO_2), titanium/titanium alloys, alloys containing more than 30% silver, or dental alloys outside the specified CTE range.

If there are known intolerances to any of the ingredients, ceraMotion® Me must not be used.

*CTE – coefficient of thermal expansion (10⁻⁶ K⁻¹, 25 °C to 500 °C).

Framework design

The framework replicates the tooth in an anatomical form reduced in size; it is necessary to avoid corners and edges in the framework. The veneering ceramic can only be fired at a maximum layer thickness of 2 mm.




Fig. 1: Framework design, anterior crown




Fig. 2: Framework design, posterior tooth crown

Framework preparation

Please observe the manufacturer's instructions for finishing, blasting and oxidation firing.

Finish remanium® alloys with a cross-cut tungsten carbide bur, blast with Al_2O_3 (125 μm) and clean. Oxidation firing is not necessary for remanium®(Fig. 3), but is recommended as cleaning and control firing.

Further information on framework preparation for non-precious alloys can be found online at www.remanium-kompodium.de



Fig. 3: Framework

Note:

Dental alloys containing zinc (Zn) must be treated in an acid bath for 5 – 10 min after oxidation according to the manufacturer's instructions.



Fig. 4: Correctly fired firing sample



Fig. 5: Underfired firing sample

Firing check

To adjust the firing temperature of your furnace, we recommend carrying out a test firing. This is the only way to find the right firing procedure for your furnace.

To prepare the test sample, mix Transpa material T with Me Standard Modelling Liquid (REF 254-001-10).

Carry out the first dentin firing. Place the firing sample on platinum foil and not on firing cotton, as otherwise there is a risk of clouding.

The furnace temperature is correct if the firing sample comes out of the furnace clear, translucent, and with sharp edges (see Fig. 4).

If the temperature is too high, the sample is very shiny and no longer has any sharp edges. If the final temperature is too low, the sample is milky white (see Fig. 5).

Please lower or raise the final temperature accordingly in increments of 10 °C. Then fire another sample.



Fig. 6: Application with a brush



Fig. 7: Paste Opaque after the first firing

Paste Opaque

Paste Opaque:

Paste Opaque is suitable for all precious and non-precious bonding alloys as well as galvano gold.

Apply Paste Opaque evenly to the framework; no wash firing is required (observe the instructions of the alloy manufacturer).

Note:

Gently mix Paste Opaque in the jar with a glass or agate spatula before use. The paste should have a creamy consistency. To restore the correct consistency after mixing, Paste Liquid (REF 254-006-02) may be added in very small amounts.

Avoid contact between Paste Opaque and water; clean the brush with Paste Liquid.



Fig. 8: Second application of Paste Opaque



Fig. 9: Paste Opaque fully fired

| | Start temperature °C | Drying time (min) | Heating rate °C /min | Vacuum start °C | Vacuum end °C | Firing temperature °C | Holding time |
|--------------------|----------------------|-------------------|----------------------|-----------------|---------------|-----------------------|------------------------|
| Me Primer CAD/CAM | 500 | 8 | 75 | 500 | 980* | 980* | 1 min (with vacuum) |
| Paste Opaque 1 + 2 | 500 | 8 | 75 | 500 | 950/980* | 950/980* | 1 min (with vacuum) |

This replaces the first firing Paste Opaque if ceraMotion® Me Primer CAD/CAM is used.

* For non-precious alloys, fire Paste Opaque 1 or ceraMotion® Me Primer CAD/CAM at 980 °C.

Powder Opaque

Opaque can be used on all indicated dental alloys. Apply Opaque evenly to the framework in 1 or 2 layers; no wash firing is required.



Fig. 10: First application of Opaque




Fig. 11: Opaque fully fired

Note: Mix Opaque with Powder BOL Liquid (REF 254-008-10) to a creamy consistency. Opaque can also be sprayed on using a spray-on system.

| | Start temperature °C | Drying time (min) | Heating rate °C /min | Vacuum start °C | Vacuum end °C | Firing temperature °C | Holding time |
|---------------------|----------------------|-------------------|----------------------|-----------------|---------------|-----------------------|---------------------|
| Powder Opaque 1 + 2 | 500 | 6 | 75 | 500 | 930 | 930 | 1 min (with vacuum) |

Layering scheme: Basic layering



- Opaque
- Dentin
- Incisal

Fig. 12: Basic layering scheme

Layering

Build up the complete anatomical tooth shape with Dentin, then cut back the Dentin in the incisal third. We recommend the mixing liquid Me Standard Modelling Liquid (REF 254-001-10).

Note:

Stains / Body Stains can be mixed into the ceramic materials at up to 10 %.



Fig. 13: Complete anatomical shape



Fig. 14: Cutback of the dentin in the incisal third



Fig. 15: Addition of incisal material



Fig. 16: Layering prior to the first firing

Layering

Addition of incisal material

Incisal allocation table:

| Dentin shade | Incisal Standard | Incisal Opal | Incisal Transpa |
|--|------------------|--------------|-----------------|
| A1, A2, B1 | I 1 | IO 1 | IT 1 |
| A3, A3.5, B2, B3, B4, C1, C2, C3, D2, D3, D4 | I 2 | IO 2 | IT 2 |
| A4, C4 | I 3 | IO 3 | IT 3 |

Note:

The layers should be somewhat over-contoured to compensate for shrinkage during sintering (Fig. 15 + 16). To control shrinkage on bridges, the layering should be separated interdentally down to the framework before the first dentin firing.

| | Start temperature °C | Drying time (min) | Heating rate °C /min | Vacuum start °C | Vacuum end °C | Firing temperature °C | Holding time | Long-term cooling |
|----------|----------------------|-------------------|----------------------|-----------------|---------------|-----------------------|--------------|-------------------|
| Dentin 1 | 500 | 6 | 55 | 500 | 870 | 870 | 2 min | – |

Firing

ceraMotion® Me was developed for rapid cooling.

This also applies to use with non-precious alloys.

The surface should appear shiny after firing.

For long-span bridges and restorations with large pontics, support pins are recommended in all crowns.

The values specified here are guidelines and must be adjusted individually according to the manufacturer's instructions and age-related deviations of the dental furnaces.

The firing charts assume furnaces that are regularly calibrated with fine silver.

All information has been compiled carefully, but is provided without guarantee.

Correction layering

Results after the first dentin firing and correction layering.



Fig. 17: Result after the first dentin firing



Fig. 18: Shape corrections with Dentin and Incisal after the first dentin firing

Correction layering



Note:

For bridge designs, first fill the inter-dental spaces and the basal surface of the pontic with Dentin.

Fig. 19: Shape corrections with Dentin and Incisal after the first dentin firing

| | Start temperature °C | Drying time (min) | Heating rate °C /min | Vacuum start °C | Vacuum end °C | Firing temperature °C | Holding time | Long-term cooling |
|----------|----------------------|-------------------|----------------------|-----------------|---------------|-----------------------|--------------|-------------------|
| Dentin 2 | 500 | 4 | 55 | 500 | 870 | 870 | 1 min | - |

Processing

Shape corrections and finishing

Make the desired shape corrections with suitable instruments.

Grind the entire surface evenly and clean thoroughly before the glaze firing.



Fig. 20: Completion

Glaze firing

Individual shade effects can be achieved by applying Stains/Body Stains to the surface (Fig 21). If required, apply Glaze (glaze material), mixed with Stains Liquid REF 254-010-02 to the entire restoration.



Fig. 21: Application of Stains/Glaze

| | Start temperature °C | Drying time (min) | Heating rate °C /min | Vacuum start °C | Vacuum end °C | Firing temperature °C | Holding time | Long-term cooling |
|---------------------|----------------------|-------------------|----------------------|-----------------|---------------|-----------------------|--------------|-------------------|
| Fixing Stains | 500 | 4 | 75 | 500 | 860 | 860 | 20 s | – |
| Shine without glaze | 500 | 4 | 75 | 500 | 870 | 870 | 1 min | – |
| Shine with glaze | 500 | 6 | 75 | 500 | 860 | 860 | 1 min | – |

Completion

The finished restoration after glaze firing.



Fig. 22: Finished restoration labial



Fig. 23: Finished restoration labial

Framework preparation with ceramic shoulder

Please observe the manufacturer's instructions for finishing, blasting and oxidation firing.

Finish remanium® alloys with a cross-cut tungsten carbide bur, blast with Al_2O_3 (125 μm) and clean. Oxidation firing is not necessary for remanium® (Fig. 1).

Further information on framework preparation for non-precious alloys can be found online at www.remanium-kompodium.de



Fig. 1: Reduced framework, blasted

Opaque



Fig. 2: Individual Paste Opaque



Fig. 3: Application with a brush

Opaque




Fig. 4: Paste Opaque fired with white band inserted



Fig. 5: Paste Opaque fired with orange effect inserted

Layering scheme: Individual layering



- Opaque
- Shoulder
- Base Dentin
- Dentin
- Incisal

Fig. 6: Individual layering system

Mixing table for shoulder materials

ceraMotion® Me offers four shoulder materials, divided into the shade groups A-B-C-D. The shoulder material “white” can be mixed accordingly to produce all shade graduations from A1 to D4, as specified in the mixing table. Adding the shoulder material “transparent” increases the translucency of the shoulder. Use Shoulder Liquid mixing liquid (REF 254-004-02).

| Tooth shade | A | B | C | D | white |
|-------------|------|------|------|------|-------|
| A1 | 50% | | | | 50% |
| A2 | 65% | | | | 35% |
| A3 | 70% | | | | 30% |
| A3.5 | 100% | | | | |
| A4 | 100% | | | | |
| B1 | | 35% | | | 65% |
| B2 | | 80% | | | 20% |
| B3 | | 90% | | | 10% |
| B4 | | 100% | | | |
| C1 | | | 50% | | 50% |
| C2 | | | 75% | | 25% |
| C3 | | | 85% | | 15% |
| C4 | | | 100% | | |
| D2 | | | | 60% | 40% |
| D3 | 60% | | | 30% | 10% |
| D4 | | | | 100% | |

Notes

Shoulder



Fig. 7: First application of Shoulder material



Fig. 7: First application of Shoulder material

| | Start temperature °C | Drying time (min) | Heating rate °C /min | Vacuum start °C | Vacuum end °C | Firing temperature °C | Holding time | Long-term cooling |
|----------------|----------------------|-------------------|----------------------|-----------------|---------------|-----------------------|--------------|-------------------|
| Shoulder 1 + 2 | 500 | 6 | 55 | 500 | 900 | 900 | 1 min | - |



Fig. 9: Second application of Shoulder material



Fig. 10: Fired shoulder

Layering

Build up the complete anatomical tooth shape with Dentin.



Fig. 11: Complete anatomical shape



Fig. 12: Cutback of the dentin in the incisal third



Fig. 13: Application of the Transpa margin



Fig. 14: Application of Dentin Modifier Fluo

Note: The individual layering shown is a suggestion and must be adapted to achieve the desired effects.

Layering



Fig. 15: White band added,
orange effect in the cervical area



Fig. 16: Alternate layering with I 2 and IO 2



Fig. 17: Cutback, application of Dentin Modifier Fluo orange, blended out to the incisal edge



Fig. 18: Addition of Dentin and Incisal Opal

Build-up



Fig. 19: Result after the first dentin firing



Fig. 20: Result after the first dentin firing

Correction layering and completion



Fig. 21: Layering with Dentin, Incisal and Transpa 1/1



Fig. 22: Processing and finishing of the surface

Completion



Fig. 23: Application of Stains/Glaze



Fig. 24: Finished restoration



Fig. 25: Finished restoration



Fig. 26: Finished restoration

Creative firing / additional firing:

a) Creative firing

- Natural-looking marginal and incisal ridges
- Individual glaze levels on a single restoration
- Correction and glaze firing in one step



Fig. 27: Application of Dentin/Incisal TU materials mixed with Modelling Liquid



Fig. 28: Finished restoration

Note: It can be carried out with or without glaze material. When using glaze material, first wetten the entire surface, apply the staining, then apply the Touch Up materials over it.

| | Start temperature °C | Drying time (min) | Heating rate °C /min | Vacuum start °C | Vacuum end °C | Final temperature °C | Holding time | Long-term cooling |
|---|----------------------|-------------------|----------------------|-----------------|---------------|----------------------|--------------|-------------------|
| Glaze and correction with or without glaze material | 500 | 6 | 75 | 500 | 860 | 860 | 20 s | – |

b) Additional firing

- Subsequently add occlusal contact points
- Basal relining



Fig. 29: Application of Touch Up for missing occlusion



Fig. 30: Basal relining

Note: It can be carried out with or without glaze material. When using glaze material, first wetten the entire surface, apply the staining, then apply the Touch Up materials over it.

Repair

- Repair of a worn restoration



Fig. 31: Preparing the framework



Fig. 32: Finished restoration after glaze firing

Note: Restorations that have been worn in the mouth must be dried in the preheating furnace. Clean the restoration; the surface must be roughened or blasted. Heat the restoration in the preheating furnace from room temperature to 600 °C at 5 °C/min; holding time 2–4 hours. After the holding time, remove the object from the furnace immediately.

Apply Opaque medium and fire using the specified firing procedure.

For repair work, mix 50% Touch Up with 50% Base Dentin, Dentin or Incisal.

| | Start temperature °C | Drying time (min) | Heating rate °C /min | Vacuum start °C | Vacuum end °C | Final temperature °C | Holding time | Long-term cooling |
|--|----------------------|-------------------|----------------------|-----------------|---------------|----------------------|--------------|-------------------|
| Touch Up Opaque | 500 | 8 | 75 | 500 | 860 | 860 | 20 s | – |
| Dentin firing for repair | 500 | 6 | 75 | 500 | 860 | 860 | 20 s | – |
| Glaze firing <u>without</u> glazing materials for repair | 500 | 4 | 75 | 500 | 860 | 860 | 20 s | – |
| Glaze firing <u>with</u> glazing materials for repair | 500 | 6 | 75 | 500 | 860 | 860 | 20 s | – |

Firing chart (universal)

| | Start temperature °C | Drying time (min) | Heating rate °C /min | Vacuum start °C | Vacuum end °C | Firing temperature °C | Holding time | Long-term cooling |
|--|----------------------|-------------------|----------------------|-----------------|---------------|-----------------------|---------------------|-------------------|
| Me Primer CAD/CAM | 500 | 8 | 75 | 500 | 980** | 980** | 1 min (with vacuum) | – |
| Paste Opaque 1 + 2 | 500 | 8 | 75 | 500 | 950/980** | 950/980** | 1 min (with vacuum) | – |
| Powder Opaque 1 + 2 | 500 | 6 | 75 | 500 | 930 | 930 | 1 min (with vacuum) | – |
| Shoulder 1 + 2 | 500 | 6 | 55 | 500 | 900 | 900 | 1 min | – |
| Dentin 1 | 500 | 6 | 55 | 500 | 870 | 870 | 2 min | – |
| Dentin 2 | 500 | 4 | 55 | 500 | 870 | 870 | 1 min | – |
| Correction material* | 500 | 4 | 75 | 500 | 860 | 860 | 20 s | – |
| Fixing Stains | 500 | 4 | 75 | 500 | 860 | 860 | 20 s | – |
| Shine without glaze | 500 | 4 | 75 | 500 | 870 | 870 | 1 min | – |
| Shine with glaze | 500 | 6 | 75 | 500 | 860 | 860 | 1 min | – |
| Touch Up Glaze and correction with or without glaze material | 500 | 6 | 75 | 500 | 860 | 860 | 20 s | – |

This replaces the first firing Paste Opaque if ceraMotion® Me Primer CAD/CAM is used.

* The correction material must be mixed 1:1 with Base Dentin, Dentin or Incisal.

** For non-precious alloys, fire Paste Opaque 1 or ceraMotion® Me Primer CAD/CAM at 980 °C.

Firing

ceraMotion® Me was developed for rapid cooling.

This also applies to use with non-precious alloys.

The surface should appear shiny after firing.

For long-span bridges and restorations with large pontics, support pins are recommended in all crowns.

Physicochemical data (according to DIN EN ISO 6872) ceraMotion® Me

| ceraMotion® Me | Type | Class | Coefficient of thermal expansion α (25 °C to 500 °C or T_g) in $1 \cdot 10^{-6} K^{-1}$ | Glass transition temperature T_g in °C | Chemical solubility L in $\mu g/cm^2$ | Flexural strength σ in MPa |
|----------------|------|-------|---|--|--|--------------------------------------|
| Primer | I | 1 | 12.1 | 610 | 40 | 120 |
| Opaque | I | 1 | 12.3 | 590 | 35 | 125 |
| Dentin | I | 1 | 12.5 | 550 | 35 | 95 |
| Incisal | I | 1 | 12.5 | 550 | 35 | 95 |
| Modifier | I | 1 | 12.5 | 550 | 35 | 95 |
| Touch Up | I | 1 | 12.3 | 520 | 30 | 100 |
| Glaze | I | 1 | 10.9 | 480 | 50 | 130 |
| Stains | I | 1 | 8.4 | 530 | 45 | 130 |

BASIC LINE / INDIVIDUAL LINE / TOUCH UP

Product overview

| | | |
|-----------------------------|-------|---|
| Opaque, Paste Opaque | O, PO | OA1- OD4 |
| Opaque Modifier | PO | gingival, orange |
| Paste Opaque Modifier | POM | gingival, orange |
| Primer CAD/CAM | PR | |
| Shoulder | SM | A, B, C, D, white, transpa |
| Gingival | G | 1, 2, 3, 4 |
| Incisal Value | IV | 1, 2, 3 |
| Value Modifier Bright | VM B | incisal opal, dentin opal |
| Incisal Modifier Value | IM V | opal pure, opal pink, opal violet |
| Transpa Modifier Value | TM V | orange, lemon, salmon |
| Transpa Value | TV | Transpa |
| Base Dentin | BD | A-D |
| Base Dentin Modifier | BDM | salmon, caramel, ochre, ivory, lemon, vanilla, brown |
| Dentin | D | A-D |
| Dentin Modifier Chroma | DM C | A, B, C, orange |
| Dentin Modifier Fluo | DM F | cream, yellow, orange |
| Incisal | I | 1, 2, 3 |
| Incisal Opal | IO | 1, 2, 3 |
| Incisal Transpa | IT | 1, 2, 3 |
| Transpa | T | Transpa |
| Incisal Modifier | IM | opal honey, opal white, opal blue, grey, opal grey |
| Chroma Concept Paste Opaque | CC PO | 1, (bleach), 2, 3, 4 |
| Chroma Concept Opaque | CC O | |
| Chroma Concept Dentin | CC D | 1 (bleach), 2 (bleach), 3, 4, 5, 6, 7, 8, 9, 10, 11, 12 |
| Chroma Concept Incisal | CC I | 1 (bleach) |
| Touch Up Paste Opaque | TU BD | light, medium, dark |
| Touch Up Base Dentin | TU BD | light, medium, dark |
| Touch Up Dentin | TU D | light, medium, dark |
| Touch Up Incisal | TU I | medium, opal, transpa |
| Touch Up Gingival | TU G | 2, 4 |
| Correction | C | Transpa |
| Glaze | GL | Transpa |
| Body Stains | B ST | A, B, C |
| Stains | ST | 1 white, 2 vanilla, 3 yellow, 4 orange, 5 pink, 6 purple, 7 blue, 8 grey, 9 olive green, 10 olive yellow, 11 medium brown, 12 red brown, 13 black |
| Liquids | | Standard Modelling Liquid, Modelling Liquid, Modelling Liquid +, Paste Liquid, Powder BOL Liquid, Shoulder Liquid, Stains Liquid, Contrast Marker |

Notes

Note

- The basis for fabricating good ceramics is the exact firing temperature of your ceramic furnace. We therefore recommend you regularly check the temperature controls on your furnace.
- The chambers should be cleaned from time to time to avoid contamination of the ceramic.
- Keep furnaces closed. To avoid moisture entering the chamber, always close the furnace after use and change to night mode where appropriate.
- For restorations temporarily worn in the mouth with ceraMotion® Me, the furnace cleaning instructions are crucial for successful further processing.

Furnace cleaning

Ceramic furnaces should be cleaned regularly to remove deposits on the inner walls of the firing chamber.

We generally recommend:

- Cleaning the furnace with carbon fire chips REF 260-317-00
- Also clean the firing tray
- Base temperature: 600 °C
- Drying time: 1 min
- Heating rate: 100 °C – 120 °C/min
- Final temperature: 1050 °C
- Holding time: 10 min

Carry out the firing program without vacuum. Observe the furnace manufacturer's instructions!

Processing tips for ceraMotion® Me ceramic

| NO. | PROCESSING TIPS FOR THE FOLLOWING PROBLEM |
|-----|---|
|-----|---|

- | | |
|----|---|
| 1 | Irregular shade of the framework surface after oxidation firing. |
| 2 | Metal frameworks are distorted. |
| 3 | Shades are too light and not sufficiently transparent. Ceramic is porous. |
| 4 | Ceramic surfaces too rough. |
| 5 | Surfaces are too smooth. Edges and contours become rounded. |
| 6 | Poor bond between the ceramic layers. |
| 7 | Glaze level insufficient. |
| 8 | Longitudinal cracks after firing. |
| 9 | Cracks after firing. |
| 10 | Cracks or microbubbles basally or on the shoulder after firing. |
| 11 | Cracks and chipping in the marginal region. |
| 12 | Chipping during dentin firing. |
| 13 | Delayed cracks in the ceramic. |
| 14 | Bubbles in the opaque. |
| 15 | Bubbles in high-gold alloy with high zinc content. |
| 16 | Bubbles originating from the framework. |
| 17 | Bubbles in the ceramic. |
| 18 | Bubbles during firing of temporarily worn bridges. |

Tips for ceraMotion® Me ceramic

| No. | Problem | Cause |
|-----|---|---|
| 1 | Irregular shade of the framework surface after oxidation firing. | <ul style="list-style-type: none"> • Unsuitable or contaminated abrasive tools. • Wrong or old casting crucible. <hr/> <ul style="list-style-type: none"> • Contamination due to incorrect blasting and cleaning. <hr/> <ul style="list-style-type: none"> • Porosities in the metal framework with inclusions that generate gas. Cause: <ul style="list-style-type: none"> – Too little casting metal – Old casting metal (reused too often) – 50:50 old-to-new rule not met – Incorrect positioning of sprues – Grinding in different directions during finishing – which creates overlaps, especially with precious alloys – unsuitable diamond grinding point for finishing metal frameworks |

| | Remedy |
|--|--|
| | <ul style="list-style-type: none"> • Use abrasive tools suitable for the type of alloy. Use separate abrasive tools for each type of alloy. |
| | <ul style="list-style-type: none"> • Use the casting crucible for one metal only. Use a new crucible if too much slag is formed. Only use ceramic crucibles. |
| | <ul style="list-style-type: none"> • Blast the framework at a 45° angle using a disposable blasting unit at 2–3 bar pressure. Ultrasonic cleaning with distilled water or steam blasting. |
| | <ul style="list-style-type: none"> • Observe the manufacturer's instructions. Different types of alloy (precious and non-precious alloys) require specific procedures (finishing, oxidation, pickling, etc.). |

Tips for ceraMotion® Me ceramic

| No. | Problem | Cause |
|-----|--|--|
| 1 | Irregular shade of the framework surface after oxidation firing. | <ul style="list-style-type: none"> <li data-bbox="468 303 908 328">• Insufficient grinding of the framework surface. <li data-bbox="468 405 626 430">• Soldered joints. <li data-bbox="468 478 896 530">• Carbon, hydrogen and/or oxygen enrichment in the alloy. |
| 2 | Metal frameworks are distorted. | <ul style="list-style-type: none"> <li data-bbox="468 552 735 577">• Final temperature too high. <li data-bbox="468 625 661 650">• Heat rate too high. <li data-bbox="468 670 666 695">• Connector too thin. |
| 3 | Shades are too light and not sufficiently transparent. Ceramic is porous. | <ul style="list-style-type: none"> <li data-bbox="468 718 789 743">• Preheating temperature too high. <li data-bbox="468 763 727 788">• Final temperature too low. <li data-bbox="468 808 758 833">• Vacuum switched on too late. <li data-bbox="468 852 984 877">• Insufficient or no vacuum during the ongoing program. <li data-bbox="468 897 789 922">• Moisture in the furnace chamber. <li data-bbox="468 942 1003 967">• Unsuitable and/or too thickly applied separating medium. |

| | Remedy |
|--|---|
| | <ul style="list-style-type: none"> • Rework the entire surface to be veneered. Reduce oxides, surface porosities and investment material residues. This also applies to milled frameworks or those fabricated by laser melting. |
| | <ul style="list-style-type: none"> • Thoroughly rework, pickle and blast the solder joints. |
| | <ul style="list-style-type: none"> • Observe the alloy manufacturer's processing instructions. Take the flame setting, casting temperature and crucible recommendations into account. |
| | <p>We recommend you run a test firing to find the right firing temperature for your furnace, as this is the only way to assess the correct firing procedure.</p> <p>Use Transpa T material, mixed with Modelling Liquid (REF 254-000-10) and carry out the first dentin firing.</p> <p>Place the firing sample on platinum foil and not on firing cotton, as otherwise there is a risk of clouding. The furnace temperature is correct if the firing sample comes out of the furnace clear and translucent, with sharp edges.</p> <p>Do not use baby oil or similar products as separating media.</p> |

Tips for ceraMotion® Me ceramic

| No. | Problem | Cause |
|-----|--|--|
| 3 | Shades are too light and not sufficiently transparent. Ceramic is porous. | <ul style="list-style-type: none"> • Use of metal spatulas during mixing. • Ceramic remixed with modeling liquid. • Contaminated brush water. |
| 4 | Ceramic surfaces too rough. | <ul style="list-style-type: none"> • Final temperature too low. |
| 5 | Surfaces are too smooth. Edges and contours become rounded. | <ul style="list-style-type: none"> • Final temperature too high. |
| 6 | Poor bond between the ceramic layers. | <ul style="list-style-type: none"> • Final temperature too low. • See point 12. • The ceramic was not evenly moist during layering and may have been dry. • The surface of the ceramic was contaminated prior to the subsequent firings. |

| Remedy | |
|--------|--|
| | <ul style="list-style-type: none"> • Use glass and/or agate spatulas for mixing to prevent metal abrasion. |
| | <ul style="list-style-type: none"> • Only use distilled water to remix the ceramic. |
| | <ul style="list-style-type: none"> • Replace brush water. |
| | <p>The correct firing temperature results in a shiny ceramic with sharp edges. If the ceramic is rough, the temperature is too low. Increase the temperature in 10 °C increments and fire a new test sample.</p> <ul style="list-style-type: none"> • Check the furnace temperature; see remedies 2, 3, 4 and 5. • Check the vacuum. |
| | <ul style="list-style-type: none"> • See remedy point 12. |
| | <ul style="list-style-type: none"> • Maintain a consistent moisture level during layering. |
| | <ul style="list-style-type: none"> • After grinding and prior to the subsequent firing, remove any grinding dust or other contaminants, such as grease or separating medium, and steam clean. |

Tips for ceraMotion® Me ceramic

| No. | Problem | Cause |
|-----|--|--|
| 7 | Glaze level insufficient. | <ul style="list-style-type: none"> Contamination of the ceramic surface by grinding dust or residues from silicone and rubber polishers. |
| 8 | Longitudinal cracks after firing. | <ul style="list-style-type: none"> Ceramic material not separated down to the opaque prior to the first dentin firing. |
| 9 | Cracks after firing. | <ul style="list-style-type: none"> Incorrect framework design. Framework finished with edges too sharp. Framework completely covered with ceramic. No oxidation firing carried out. Ceramic material not separated down to the opaque prior to the first dentin firing. Long-term cooling. CTE of the framework outside the specified range or very close to the limit. |

| Remedy | |
|--------|--|
| | <ul style="list-style-type: none"> • Clean well. |
| | <ul style="list-style-type: none"> • Increase the final temperature by 10 °C and fire again. |
| | <ul style="list-style-type: none"> • Mix the materials with Me Standard Modelling Liquid (REF 254-001-10). |
| | <ul style="list-style-type: none"> • Separate the build-up down to the opaque prior to the first firing in order to control shrinkage. |
| | <ul style="list-style-type: none"> • Model reduced tooth shapes, design the framework to be stable. |
| | <ul style="list-style-type: none"> • Round off edges using suitable burs. |
| | <ul style="list-style-type: none"> • Model collars or cooling surfaces on the framework such that heat can dissipate from the framework. |
| | <ul style="list-style-type: none"> • Oxidation firing according to the alloy manufacturer's instructions. |
| | <ul style="list-style-type: none"> • See point 8. |
| | <ul style="list-style-type: none"> • Fast cooling; opening the furnace immediately after the firing process, i.e., the ceramic furnace should be fully open within 15 seconds at the latest. Wetten cracks with Me glaze material (REF 252-270-02) and carry out a second dentin firing; no slow cooling. |
| | <ul style="list-style-type: none"> • CTE should be in the range 13.9 °C – 15.0 °C/25 °C – 500 °C. If the specification refers to the temperature range 25 °C – 600 °C, it may be slightly lower at 25 °C – 500 °C. |

Tips for ceraMotion® Me ceramic

| No. | Problem | Cause |
|-----|--|--|
| 10 | Cracks or microbubbles basally or on the shoulder after firing. | <ul style="list-style-type: none"> • Oil-based separating media. |
| 11 | Cracks and chipping in the marginal region. | <ul style="list-style-type: none"> • Corrections to the inside of the framework; coarse reworking in the event of fit inaccuracies. • Unfavorable stump preparations. • Trial wearing of the restoration without bonding (cement). • Coarse removal after framework try-in. • Framework too thin in the marginal region. |
| 12 | Chipping during dentin firing. | <ul style="list-style-type: none"> • Furnace start temperature too high. • Furnace opening too narrow. • Firing tray and pins too hot. • Predrying time too short. • The displayed temperature does not always reflect the actual conditions in the furnace chamber (depending on the positioning of the thermocouple and heat radiation). • Surface contaminated, possibly by separating medium or grinding dust, (acts like a separating layer). |

| Remedy |
|--|
| <ul style="list-style-type: none"> • Use separating media for low-fusing ceramics. |
| <ul style="list-style-type: none"> • Carry out framework try-ins in case of inaccurate impressions or unfavorable preparations. |
| <ul style="list-style-type: none"> • Defined preparation margins; chamfer preparations where necessary. |
| <ul style="list-style-type: none"> • Avoid trial wearing of non-cemented restorations. |
| <ul style="list-style-type: none"> • Use the shepherd's hook only in the interdental region. |
| <ul style="list-style-type: none"> • Framework must not be thinner than 0.3 mm. |
| <ul style="list-style-type: none"> • Lower the base temperature to 450 °C. |
| <ul style="list-style-type: none"> • Use cold firing trays and pins. |
| <ul style="list-style-type: none"> • Increase the predrying times for larger restorations. |
| <ul style="list-style-type: none"> • Do not place the object on the firing tray too early. |
| <ul style="list-style-type: none"> • Clean the surface thoroughly prior to application to ensure bonding. |

Tips for ceraMotion® Me ceramic

| No. | Problem | Cause |
|-----------|---------------------------------------|--|
| 13 | Delayed cracks in the ceramic. | <ul style="list-style-type: none">• Excessive blasting of the inner crown, possibly with too much pressure and the incorrect grain size.• Excessive localized steam cleaning. |
| 14 | Bubbles in the opaque. | <ul style="list-style-type: none">• Insufficient predrying of the opaque and/or too high preheating temperature. If the proportion of opaque liquid is too high and the predrying time is too short, or the preheating temperature is too high, the liquids escape explosively. This results in bubbles, cracking and impaired bonding (chipping).• Oil residues caused by blowing with compressed air.• Separating media, finger grease and cream residues on the surfaces to be veneered weaken the bond and cause bubbles and cracks in the opaque. |

| Remedy | |
|--------|---|
| | <ul style="list-style-type: none"> • Blast with 50 µm aluminum oxide or glass beads at a pressure below 2 bar. |
| | <ul style="list-style-type: none"> • Avoid thin areas on the framework – min. 0.3 mm. |
| | <ul style="list-style-type: none"> • Carefully steam clean the restoration. |
| | <ul style="list-style-type: none"> • Wetten cracks with Me glaze material (REF 252-270-02) and carry out a second dentin firing; no slow cooling. |
| | <ul style="list-style-type: none"> • Lower the preheating temperature to 450 °C. Ensure that both the firing tray and the pins are as close to room temperature as possible. Take account of the radiant heat from the furnace. Lift position fully lowered. Increase the predrying time to 1–2 min. Ensure that the object is placed in the furnace once the preheating temperature has been reached, i.e., the furnace must have cooled back down fully to the preheating temperature after a previous firing. |
| | <ul style="list-style-type: none"> • Clean the brush for paste opaque with Paste Liquid (REF 254-006-02). Do not allow contact with water. Adjust the opaque consistency with a small amount of Opaque Liquid. |
| | <ul style="list-style-type: none"> • Check compressed air system; avoid compressed air if necessary. |
| | <ul style="list-style-type: none"> • Handle the surfaces cleanly. |

Tips for ceraMotion® Me ceramic

| No. | Problem | Cause |
|-----|---|---|
| 15 | Bubbles in high-gold alloy with high zinc content. | <ul style="list-style-type: none">• No blasting or pickling. |
| 16 | Bubbles originating from the framework. | <ul style="list-style-type: none">• Incorrect framework conditioning: smeary finishing leads to unfavorable changes to the framework surface, especially with precious metal alloys.• Contamination: Severely contaminated furnaces (if also used for other ceramic systems, bonder firings), bonder firings for electroformed restorations, furnace soldering, stump investment materials, etc. <ul style="list-style-type: none">• Voids or porosities after casting. |

| | Remedy |
|--|--|
| | <ul style="list-style-type: none"> • Strictly follow the alloy manufacturer’s instructions regarding blasting or pickling. |
| | <ul style="list-style-type: none"> • Use only hard metal burs and grind in one direction. After finishing, blast the framework with aluminum oxide (125 µm–250 µm) in a disposable blasting unit at a pressure of 2–3 bar and at a 45° angle. Then steam clean. Oxidation firing according to the alloy manufacturer’s instructions. |
| | <ul style="list-style-type: none"> • Carry out furnace cleaning firing with carbon fibre pads more frequently. (REF 260-317-00). Furnace cleaning: <ul style="list-style-type: none"> – Also clean the firing tray – Base temperature: 600 °C – Drying time: 1 min – Heating rate: 100 °C – 120 °C/min – Final temperature: 1050 °C – Holding time: 10 min – Carry out the firing program without vacuum. – Observe the furnace manufacturer’s instructions! |
| | <ul style="list-style-type: none"> • Open the voids and porosities and laser-weld or solder them. |

Tips for ceraMotion® Me ceramic

| No. | Problem | Cause |
|-----|---|--|
| 17 | Bubbles in the ceramic. | <ul style="list-style-type: none">• Dirt particles incorporated during layering.• Separating media on the ceramic surface, poorly cleaned surfaces (grinding particles acting like a separating layer).• Ceramic repeatedly remixed with modeling liquid.• Bubbles originating from the framework (point 16).• Bubbles in the opaque (point 14). |
| 18 | Bubbles during firing of temporarily worn bridges. | <ul style="list-style-type: none">• Worn restorations were not dried carefully. |

| | Remedy |
|--|---|
| | <ul style="list-style-type: none"> • Cover the material (close the ceramic jars after use and protect the powder on the mixing plate). Clean the restorations under running water after each grinding operation. |
| | <ul style="list-style-type: none"> • Work cleanly. |
| | <ul style="list-style-type: none"> • Only use distilled water to remix the ceramic. |
| | <ul style="list-style-type: none"> • See point 16. |
| | <ul style="list-style-type: none"> • See point 14. |
| | <ul style="list-style-type: none"> • Clean the restoration. The surface must be roughened or blasted. Heat the restoration in the preheating furnace from room temperature to 600 °C at a rate of 5 °C per minute. Holding time at least 2–4 hours. Remove directly from the furnace and carry out correction firings. |

Notes



CONTACT
DENTAURUM

Dentaurum GmbH & Co. KG



Turnstr. 31
75228 Ispringen · Germany



info@dentaurum.com
www.dentaurum.com



+49 72 31/803-0



ONLINE SHOP
SHOP.DENTAURUM.COM