Estimado cliente

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EN

Dear customer,

Thank you for choosing a quality product from Dentaurum. Please read these Instructions for use to ensure you use this product in a safe and simple way and that you and your patients gain full benefit. In case of questions or ideas, please contact your local representative.

1. Manufacturer

Dentaurum GmbH & Co. KG
Turnstr. 31 | 75228 Ispringen | Germany

2. General device description

The hyrax® click screw without retention legs (Hygienic Rapid Expansion Screw, Figs. 1 – 3) is a special expansion screw for rapid palatal expansion. It has a screw locking feature to prevent unwanted reversing. The body of the RPE screw is welded onto individual retention legs, connectors or cap splints which have been either cast or produced by the SLM method in metal. The appliances are cemented or fixed to the tomas®-pin as usual. Follow the instructions from the cement manufacturer.

3. Intended purpose

Used for rapid palatal expansion with a fixed appliance.

4. Intended users

All products are intended for dental professionals only.

5. Indications

In case of transverse under-development of the upper jaw.

6. Contraindications

Stainless steel expansion screws must not be used in patients with a nickel and/or chromium allergy. Warnings are found on the respective product labels. Also see 9. Composition.

7. Patient target group

All patients for whom a fixed appliance with an expansion screw can be placed as part of orthodontic treatment.

8.1. Fabrication of a hyrax® click appliance with individualized retention legs and cap splints

The following steps should be followed to manufacture individualized retention legs and cap splints by casting or by using the SLM method (Figs. 4 – 5).

Individualized cap splints are fitted on the upper jaw first premolars and the first molars, each left and right (Figs. 4 – 7). The retention legs are connected to the hyrax® click screw in such a way that two weld seams form along the lateral groove (Fig. 7).

The hyrax® click screw should be “suspended” at least 1 mm – 2 mm above the palatal vault/gingiva (Fig. 6) and its dorsal side should be in line with the mesio-palatal cusps of the first molars. In order to avoid leverage, it must be situated parallel to the occlusal plane. The hyrax® click screw is positioned so that it may be turned in a dorsal direction (Figs. 6 – 8) in order to make the task of inserting the safety key (Fig. 10) and activating the screw easier.

8.2. Tips for working with hyrax® click screws without retention legs

For the connection to be durable, the individualized retention legs must be welded onto the groove of the hyrax® click screw body (Fig. 7), for example using the desktop Compact welding unit (REF 090-578-50). As an average value for NiCr, the initial laser setting on the unit is as follows: 240 volt, 4 ms, 0.8 mm. The use of a shielding gas is indispensable during the welding process.

When using or adding new material to reinforce the welding area, we recommend the NiCr welding wire, ø 0.50 mm (REF 528-220-00, Fig. 9).

9. Composition

The hyrax® click screws are made of the following materials: DIN 1.4301, DIN 1.4305 and DIN 1.4310. The material composition is included in the orthodontics materials list; please refer to the orthodontics catalog or www.dentaurum.com. The SSCP is available at https://ec.europa.eu/tools/eudamed and www.dentaurum.com.

10. Safety information

Always use only the blue safety key (Fig. 10) supplied to prevent swallowing.

The product should not be used if there is a known allergic reaction to one or more of the material components.

If the packaging shows serious damage, the product must be checked for integrity and cleanliness before use, and if necessary must be disposed of.

11. Information for single use products

The product is intended for single use only. Reconditioning (recycling) of the product and its reuse is not permitted, as perfect functionality of the product cannot be guaranteed.

12. Further information

Should the dental professional and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the dental professional and/or the patient is resident are informed accordingly.

13. Product range

For the complete range of hyrax® click screws, please see our orthodontics catalog or visit www.dentaurum.com.

14. Quality information

Dentaurum ensures faultless quality of its products. These recommendations are based upon Dentaurum's own experiences. The dental professional is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as Dentaurum has no influence on the processing on site.

15. Explanation of symbols used on the label

⚠ Please refer to the label. These Instructions for use and additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00).

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Information for the dentist/patient

Before the hyrax® click appliance is cemented into position, open the spindle (screw) by making 2 – 3 quarter turns in the direction of the arrow to test the mobility. Then turn back to the original position. The finished hyrax® click appliance is cemented into position in the usual manner. Follow the instructions from the cement manufacturer. Following this, the screw is turned four times 90° (quarter turn) using the safety key supplied for this purpose. One full 360° turn of the spindle = 0.8 mm expansion. Following instruction, the patient should complete two quarter turns (2 x 90°) in the morning and two quarter turns in the evening until the desired expansion is achieved. After each step of activation (quarter turn), a clear „click“ can be perceived. Patients: Please contact the dentist who is treating you if you have any further questions.

Note: Always use only the blue safety key supplied to prevent swallowing! The safety key is designed to take a maximum force of 35 N ± 3.5 kg. Please note that the force exerted on the safety key wire should always be exerted on the closed side of the wire clasp so that the clasp cannot loosen from the axis. Improper use may lead to overload, and consequently to failure of the safety key, thus potentially excluding any guarantee and liability. Regular visits are recommended to monitor progress. We recommend securing the hyrax® click appliance by fitting a ligature or blocking out with acrylic after the desired expansion is reached, in spite of the feature to prevent unwanted reversing. For stabilization, the appliance should remain in the mouth for approximately 2–3 months. A retention plate should be worn following removal of the appliance.

⚠ The product should not be used if there is a known allergic reaction to one or more of the material components.

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