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For long durability of the crucibles please observe the following recommendations:

- Use only dry, clean crucibles for casting
- Sandblast oxidized crucibles with AL_2O_3 (110–250 μm)
- Use original Dentaurem-casting cylinders only
- Place casting cylinders according to inserts onto the crucible (please refer to instructions for use of casting machines)
- Place crucible correctly into crucible support (please refer to instructions for use)
- Observe the melting times for the different ingot sizes

Note:

Condition of the electrode and adjustment of the distance to the ingot influence the melting of the metal, thus increasing the thermal exposure of the crucible.

We recommend the following:

- Check electrode distance and adjust correctly
- Check tip of the electrode
- Sandblast oxidized electrode

Extended melting times result in overheating and distortion of the crucible

- Remove crucible **immediately** after casting and quench in water
- Do not allow slag to cool on the crucible
- Do not knock on crucible with hammer
- Do not remove adhering residues of titanium by force, but rather by lightly knocking on the back of the **slag**
- Small residue amount of Titanium on the crucible may be removed with a rubber polisher
- Replace worn or distorted crucibles

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