



Preparation protocol.

Preparation protocol with CITO mini® preparation instruments.
Taking into account different bone qualities.

Adapt preparation protocol according to the indication and situation of each patient.
The insertion of CITO mini® implants in soft bone is not indicated.



		Soft bone quality	Medium bone quality	Hard bone quality
		 $\varnothing 1.8^3$ $\varnothing 2.2$ $\varnothing 2.5$	 $\varnothing 1.8^3$ $\varnothing 2.2$ $\varnothing 2.5$	 $\varnothing 1.8^3$ $\varnothing 2.2$ $\varnothing 2.5$
Marking drill ²			X X X	X X X
Stepped drill for $\varnothing 1.8$ mm ¹			min. 6.0 mm	X ⁴
Stepped drill for $\varnothing 2.2$ mm ¹		The insertion of CITO mini® implants in soft bone is not indicated.	min. 6.0 mm	X ⁴
Stepped drill for $\varnothing 2.5$ mm ¹			min. 6.0 mm	X ⁴

¹ The insertion depth/length of the drill depends on the desired primary stability and the bone quality. The depth markings/scaling help to avoid damaging neighboring structures (e.g. nerves) due to perforations that are too deep.

² Optional use.

³ Not intended for maxillary prostheses.

⁴ Should the bone be very hard, it may be necessary to use a stepped drill one size larger to ream the cortical bone.