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Dentaurum GmbH & Co. KG · Turnstr. 31 · 75228 Ispringen · Germany

2. General device description

Prosthetic abutments and components replace missing teeth, restore chewing function, and thereby improve quality of life. The implants of the tioLogic® implant family are used to fixate corresponding abutment components, depending on the indication. Prosthetic abutments and components serve to produce prosthetic meso-constructions, which are then restored with a prosthetic superstructure.

Abutment components					
Series of abutments*	S	M		L	
Implant diameter in mm	3.3	3.7	4.2	4.8	5.5

*all of the abutment components listed in these Instructions for use are available in the S, M and L series of abutments.

tioLogic® TWINFIT scan abutments, titanium, are intended for digital recording of the intra-oral implant position of the tioLogic® TWINFIT implants (axis inclination and position of the antirotational mechanism). This can be performed either on the implant in the mouth or on the laboratory implant on the working cast. The scan abutments, titanium, are supplied with a retaining screw. The tioLogic® TWINFIT CAD/CAM scan abutments, titanium, are reusable. The retaining screw is intended for single use and is marked on the blister packaging with .

Scan abutment			
Series of abutments	S	M	L
Length in mm	12.0		

tioLogic® TWINFIT laboratory implants for printed and cast models are used to create working casts. They precisely reproduce the position of the tioLogic® TWINFIT implants (axis inclination and position of the antirotational mechanism). The laboratory implants are provided together with a counter screw for model fixation. Processing and use take place exclusively in the (practice) laboratory or in the manufacturing center and not on the patient. This is not a medical device.

Laboratory implant			
Series of abutments	S	M	L
Length in mm	12.0		

tioLogic® TWINFIT CAD/CAM titanium blocks, PreForm are prefabricated blanks/components for the machined production of individual, single-piece abutments (and gingiva formers). The titanium blocks are designed for fixation in the Medentika PreFace abutment holder and subsequent computer-controlled machined production (CAM). The CAD/CAM titanium blocks are supplied with either a conical or platform connector geometry, each with an AniTite screw.

Titanium block			
Series of abutments	S	M	L
Length in mm	19.9		

tioLogic® TWINFIT CAD/CAM titanium bases serve as an adhesive base for customized, implant-supported dental restorations or for the production of two-part abutments. With almost all CAD/CAM titanium bases, the screw aperture can be angulated up to 20° in relation to the implant axis. The CAD/CAM titanium bases are supplied with either a conical or platform connector geometry, each with an AniTite screw. In addition, modeling sleeves are available specifically for CAD/CAM titanium bases.

Titanium base			
Series of abutments	S	M	L
Length in mm	3.5		

The modeling sleeves for CAD/CAM titanium bases are intended for the production of mesostructures and crown frameworks that are waxed up with a modeling wax or acrylic, or are ground with a milling cutter for acrylic to achieve the final shape. The wax-up is

then transferred to a framework using the casting technique and bonded to the CAD/CAM titanium base. Processing and use take place exclusively in the (practice) laboratory or in the manufacturing center and not on the patient. This is not a medical device.

tioLogic® TWINFIT CAD/CAM titanium bases VARIO serve as an adhesive base for customized, implant-supported dental restorations or for the production of two-part abutments. The shortenable bonding cylinder maximizes the bonding surface between the CAD/CAM titanium base VARIO and the ceramic mesostructure. With the CAD/CAM titanium bases VARIO, the screw aperture can be angulated up to 20° in relation to the implant axis. The CAD/CAM titanium bases VARIO are supplied with either a conical or platform connector geometry, each with an AniTite screw.

Titanium base VARIO	Conical			Platform		
Series of abutments	S	M	L	S	M	L
Gingival height in mm	1.6			0.6		
Length in mm	8.5					
Shortenable segments/length in mm	3 x 1.5					

CAD/CAM titanium adhesive bases are intended specifically for bonding multi-unit bridge restorations or full-arch restorations in the edentulous mandible or maxilla. The titanium adhesive bases feature a platform connector geometry and do not have an antirotational mechanism. For this reason they are therefore not approved for single restorations.

Titanium adhesive base			
Series of abutments	S	M	L
Length in mm	3.0		

The tioLogic® TWINFIT implant system contains specially coordinated abutments and components, including instruments, for the insertion and prosthetic restoration of the implants.

3. Intended purpose

Impression components are intended for three-dimensional recording and transfer of the intra-oral position of the implants.

tioLogic® TWINFIT laboratory implants for printed and cast models are designed to represent the position of the implant (axis inclination and position of the antirotational mechanism) in a working cast.

Abutments are intended to connect the implants with the prosthetic restoration.

Modeling sleeves for CAD/CAM titanium bases are designed to be processed as wax-ups for the production of cast mesostructures and crown frameworks on CAD/CAM titanium bases.

4. Intended users

The tioLogic® TWINFIT product range is supplied exclusively to qualified personnel such as surgeons, dentists and dental technicians. It should only be used by qualified personnel who are familiar with dental implantology, including diagnosis, preoperative planning, surgical technique and prosthetic restoration.

5. Medical indications / contraindications

5.1 Indications

5.1.1 tioLogic® TWINFIT scan abutments, titanium

- Localization of the exact implant position in the intra-oral scan
- Localization of the exact implant position in the model scan

5.1.2 tioLogic® TWINFIT laboratory implants for printed and cast models

- Precise reproduction of positions when printing models
- Precise reproduction of positions for cast models (dental stone)

5.1.3 tioLogic® TWINFIT CAD/CAM titanium blocks

- For the production of individual crown and bridge techniques as a connecting piece between the implant and the prosthetic abutment (cemented, individually screw-retained)
- Adhesive bond with ceramic mesostructures
- Abutment on implant \varnothing 3.3 mm only for replacement in the area of the upper lateral incisors and lower incisors

5.1.4 tioLogic® TWINFIT CAD/CAM titanium bases

- For the production of individual crown and bridge techniques as a connecting piece between the implant and the prosthetic abutment (cemented, individually screw-retained)
- Adhesive bond with ceramic mesostructures
- Angulation of the screw aperture to max. 20°
- Abutment on implant \varnothing 3.3 mm only for replacement in the area of the upper lateral incisors and lower incisors

5.1.5 tioLogic® TWINFIT CAD/CAM titanium bases VARIO

- For the production of individual crown and bridge techniques as a connecting piece between the implant and the prosthetic abutment (cemented, individually screw-retained)
- Adhesive bond with ceramic mesostructures
- Angulation of the screw aperture to max. 20°
- Abutment on implant \varnothing 3.3 mm only for replacement in the area of the upper lateral incisors and lower incisors

5.1.6 Modeling sleeves for CAD/CAM titanium bases

- Only use in conjunction with CAD/CAM titanium bases
- Wax-up for mesostructures and crown frameworks

5.1.7 tioLogic® TWINFIT CAD/CAM titanium adhesive bases

- Individual bridgework (cemented, individually screw-retained)

5.2 Contraindications

- In patients with a known intolerance to individual components.
- The contact surface between the implant and abutment must not be blasted or otherwise mechanically processed
- Insufficient oral hygiene
- Bruxism

5.2.1 tioLogic® TWINFIT scan abutments, titanium

- Use as restoration

5.2.2 tioLogic® TWINFIT laboratory implants for printed and cast models

- Use of laboratory analogs with defective or damaged connector geometry
- Use on patients

5.2.3 tioLogic® TWINFIT CAD/CAM titanium blocks

- Primary bracing of titanium blocks
- Correction of angulation exceeding 20°, or exceeding 15° in the case of implants \varnothing 3.3 mm
- Cast-on technique
- Direct laser welding
- Restorations with a length ratio to the length of the implant which exceeds 1:1.25
- Single crowns with a cantilever bridge unit
- Abutment on implant \varnothing 3.3 mm only for replacement in the area of the upper and lower posterior teeth

5.2.4 tioLogic® TWINFIT CAD/CAM titanium bases

- Primary bracing of titanium bases
- Correction of angulation exceeding 20°, or exceeding 15° in the case of implants \varnothing 3.3 mm
- Cast-on technique
- Direct laser welding
- Restorations with a length ratio to the length of the implant which exceeds 1:1.25
- Single crowns with a cantilever bridge unit
- Angulation of the screw aperture of more than 20°
- Abutment on implant \varnothing 3.3 mm only for replacement in the area of the upper and lower posterior teeth

5.2.5 tioLogic® TWINFIT CAD/CAM titanium bases VARIO

- Primary bracing of titanium bases
- Correction of angulation exceeding 20°, or exceeding 15° in the case of implants \varnothing 3.3 mm
- Cast-on technique
- Direct laser welding
- Restorations with a length ratio to the length of the implant which exceeds 1:1.25
- Single crowns with a cantilever bridge unit
- Angulation of the screw aperture of more than 20°
- Abutment on implant \varnothing 3.3 mm only for replacement in the area of the upper and lower posterior teeth

5.2.6 Modeling sleeves for CAD/CAM titanium bases

- Use on patients
- Use on other abutments than on CAD/CAM titanium bases
- Use as direct restoration (without cast)

5.2.7 tioLogic® TWINFIT CAD/CAM titanium adhesive bases

- Single crowns
- Primary bracing of titanium adhesive bases
- Correction of angulation exceeding 30°
- Cast-on technique
- Direct laser welding

6. Patient target group

All skeletally mature patients undergoing implant treatment. Please consider the indications / contraindications.

7. Information on use, processing, fitting and insertion

Prior to use, the dental professional ensures that they have carefully read and understood the tioLogic® TWINFIT Instructions for use.

We strongly recommend that, prior to using Dentaaurum products, dental professionals attend a Dentaaurum training course to familiarize themselves with the Dentaaurum products.

Carry out thorough radiological, psychological and physical examinations before implantation and inform patients of the process in detail. Dentaaurum recommends full clinical, radiological, photographic and statistical documentation.

The CAD/CAM abutments and components of the tioLogic® TWINFIT implant system can be documented, e.g., in the patient file or PatientPass, using the additional labels.

Precautions:

Before surgical implantation, it is important to carry out a careful patient anamnesis and, if required, a third-party anamnesis by the patient's doctor (GP/PCP) in order to clarify if 1) implantation may be difficult due to an anatomical situation, 2) there is the risk of a serious surgical problem or a general risk, 3) wound healing/osseointegration may be difficult, or 4) good hygiene and care of the implant, abutment and restoration may worsen. Do not fit dental implants if any of the above cases or a combination hereof is severe or uncontrollable.

Prosthetic restoration of the implant must not be performed until the soft tissue has healed without irritation. Furthermore, a radiographic check-up should be performed before the prosthetic restoration.

Important:

The dental professional should ensure the products cannot be aspirated during intra-oral use.

CAD/CAM data sets for 3shape™, dental wings and exocad can be downloaded free of charge at www.dentaaurum.de/cadcam. The data sets include all the information needed for all types of restorations as a complete package, i.e., for printed implant models, hybrid abutments, one-piece abutments, as well as bridges and bars.

7.1 tioLogic® TWINFIT CAD/CAM scan abutments, titanium

The tioLogic® TWINFIT scan abutments, titanium, are available in the S, M or L series of abutments. Place the scan abutment intra-orally on the tioLogic® TWINFIT implant or extra-orally on the laboratory implant in the master model (pay attention to the antirotational mechanism) and secure it with the hex key AF 1.3 and the enclosed red retaining screw. The long cylinder at the connection point of the scan abutment, titanium, serves for an easy fit into the (laboratory) implant and for exact axial detection of the position. The enclosed red retaining screw features a shortened thread, which only grips in the (laboratory) implant when the scan abutment has been inserted into the connection

point in the correct position. The surface of the scan abutment, titanium, is optimized for digital recording, both intra-orally as well as on the model without scanning spray. The scanning and matching process of the scan data within the respective CAD software is performed according to the specifications of the individual software manufacturer. The scan abutments, titanium, are intended for multiple use while maintaining the integrity of the scan body and the connection to the (laboratory) implant. The retaining screw is intended for single use and is marked on the blister packaging with .

Tightening torque	
Retaining screw with shortened thread on the model / in the mouth	15 Ncm


7.2 tioLogic® TWINFIT laboratory implants for printed and cast models

The laboratory implant can be used in both analog and digital workflows. The cavities of the laboratory implants are stored in the CAD/CAM libraries of the various software manufacturers to achieve an optimum fit and the correct positioning of the laboratory implants with the counter screw in 3D model printing. The parameters of the respective software must be checked and adjusted if necessary.

To ensure that the supplied counter screw can fixate the laboratory implant at the apical end in the printed model, the model base should be printed with a sufficient height **1**.

If a detachable mucosa is created, it should be lower than 1.0 mm below the interface to ensure the optimal locking and retaining function of the laboratory implant in the neck area of the printed model **2**. The 3 retention tabs should still be visible in the printed model for monitoring purposes.


When the model is printed with the planned cavities, the laboratory implant for printed and cast models can be inserted in accordance with the intended S, M or L series of abutments or 4Base. For this purpose, the laboratory implant is placed on the manual insertion key and fixed with the counter screw. The laboratory implant is then aligned in the cavity and pressed in. The counter screw **3** supplied with the laboratory implant can be used to fixate the laboratory implant in the end position from apical.

Laboratory implants for printed and cast models are intended for multiple uses in the laboratory, provided that the outer contour and the connection to the impression and abutment components remain intact.

Tightening torque	
Counter screw on the model	15 Ncm

7.3 tioLogic® TWINFIT CAD/CAM titanium blocks, PreForm


By using original Dentaurem CAD/CAM PreForm titanium blocks, it is possible to produce highly precise abutments for the prosthetic restoration that are specific to each patient, have a precise fit and can be directly screw-retained. The titanium blocks were designed for the Medentika® PreFace® abutment holders. These ensure a perfect connection to the tioLogic® TWINFIT implant

interface, both for in-house production and for production in a manufacturing center. CAD/CAM PreForm titanium blocks are available for all three series of abutments S, M and L with conical and platform connector geometry. Following selection of the respective indication in the data set in the corresponding CAD software, the design is made according to the specifications of the respective software manufacturer and according to the rules for dental prosthetics.

Tightening torque	
Prosthetic screw on the model	30 Ncm
AnoTite screw intra-orally	30 Ncm

To achieve optimum preloading of the AnoTite screw, the AnoTite screw must be retightened with the same torque approx. 5 minutes after insertion.

7.4 tioLogic® TWINFIT CAD/CAM titanium bases


The geometry of the tioLogic® TWINFIT CAD/CAM titanium base is designed to ensure a reliable, aesthetic bond with the ceramic mesostructure. tioLogic® TWINFIT CAD/CAM titanium bases are available in all three series of abutments S, M and L with conical and platform connector geometry. Following selection of the respective indication in the data set in the corresponding CAD software, the design is made according to the specifications of the

respective software manufacturer and according to the rules for dental prosthetics.

The angulated screw aperture offers the optimum solution for the complete segment that is occlusally screw-retained. The screw aperture can be inclined at an angle of up to 20° to the implant axis, allowing discrete emergence in the palatal area, particularly in the anterior region.

CAUTION:

Due to the higher gingival height in the conical design of the CAD/CAM titanium base, a recess in the area to be bonded was constructed for the angulated screw aperture. This creates free space for a maximum angulation of 20°. In order to achieve the correct angle when planning, the scanning surface on the scan abutment must be opposite the desired exit of the angulated screw aperture. The recess in the CAD/CAM titanium base (angulation alignment) **1** is flush with a cam from the inner connection C. In the case of the scan abutment, the upper scanning surface **2**

is located opposite a cam from the inner connection **3**. During scanning, care must be taken to determine where the exit of the screw aperture should be located. For angulated screw apertures, use the AnoTite screw for angulated screw apertures L 8.5 mm and the hex ball point key AF 1.3. Before bonding, the head of the prosthetic screw, which fixes the CAD/CAM titanium base, is covered with wax. The bonding surfaces of the ceramic sleeve and of the CAD/CAM titanium base are blasted with aluminum oxide (50 µm/2 bar). It is also advisable to use the silanization method to condition this bonding area on the CAD/CAM titanium base. Before bonding, ensure the surfaces are dry and free from grease. Adhesives must be used according to the manufacturer's instructions (e.g.: "PANAVIA™F2.0" by Kuraray Noritake or "Multilink® Hybrid Abutment" by Ivoclar Vivadent AG). Dentaurem GmbH & Co. KG gives no warranty for the shelf life of the meso/superstructure and the adhesive. After bonding, excess material is removed and the ceramic abutment is put into position (observing the antirotational mechanism). The all-ceramic crown is completed according to the manufacturer's instructions.


Important information:

- The wall thickness of the ceramic abutment must NOT be less than 0.5 mm.
- Preparation of a chamfer with angled inner edge and a minimum step of 0.5 mm.
- The contact surface between the implant and abutment must not be blasted or otherwise mechanically processed
- The height and diameter of the area to be bonded must not be reduced

Tightening torque	
Prosthetic screw on the model	30 Ncm
AnoTite screw intra-orally	30 Ncm
Prosthetic screw for angled screw aperture on the model	25 Ncm
AnoTite screw for angled screw aperture in the mouth	25 Ncm

To achieve optimum preloading of the AnoTite screw, the AnoTite screw must be retightened with the same torque approx. 5 minutes after insertion.

7.5 tioLogic® TWINFIT CAD/CAM titanium bases VARIO


The geometry of the tioLogic® TWINFIT CAD/CAM titanium base VARIO is designed to ensure a reliable, aesthetic bond with the ceramic mesostructure. tioLogic® TWINFIT CAD/CAM titanium bases VARIO are available in all three series of abutments S, M and L with conical and platform connector geometry. Following selection of the respective indication in the data set in the corresponding CAD software, the design is made according to the specifications of the respective software manufacturer and according to the rules for dental prosthetics.

The CAD/CAM titanium base VARIO features a height-adjustable bonding cylinder that, depending on the tooth geometry, maximizes the bonding surface between the CAD/CAM titanium base VARIO and the ceramic crown. For this purpose, three horizontal grooves are provided on the bonding cylinder that allow for precise adjustment of the overall height using a SiC separating disc (Supercut STM separating disc, REF 130-114-00).

In addition, the CAD/CAM titanium base VARIO features an adjustable fenestration that allows for angulation of the screw aperture by up to 20° relative to the implant axis. The individual segments of the fenestration can be removed using combined round-nose wire bending pliers (Jarabak Light Wire pliers, REF 003-125-00). The angulated screw aperture offers the optimum solution for the complete occlusally screw-retained segment and allows for a discrete emergence in the palatal area, particularly in the anterior region.

CAUTION:

Scan abutment titanium M, tiologic® TWINFIT, platform | CAD/CAM titanium base VARIO, tiologic® TWINFIT, conical


Due to the higher bonding cylinder of the CAD/CAM titanium base VARIO, a fenestration in the area to be bonded was constructed for the angulated screw aperture to create free space for a maximum angulation of 20°. After adapting the bonding cylinder to the required height, the fenestration is opened completely step-by-step.

It is possible to open the bonding cylinder on a pre-prepared surface ① for an angulated screw aperture. In order to achieve the correct angle when planning, the scanning surface ② on the scan abutment and the desired exit of the angulated screw aperture must be on the same side. The cam of the inner connection ③ is located on the opposite side in both cases.

CAUTION:

When digitally planning the ceramic crown in the CAD/CAM software, ensure that the screw aperture and the fenestration of the CAD/CAM titanium base VARIO are precisely aligned with each other.

To better facilitate the geometric adjustment of the CAD/CAM titanium base VARIO, it is recommended to screw it onto a laboratory implant using a prosthetic screw. Adjust the height of the bonding cylinder to match the specified cutting grooves taking into account the geometry of the ceramic crown using a SiC separating disc (Supercut STM separating disc, REF 130-114-00). Break out the fenestration in stages using combined flat, round-nose wire bending pliers with the smallest possible diameter in the anterior area. The flat, round-nose wire bending pliers - with the round nose on the inside - should be placed directly at the junction to the next segment. It is recommended to use Dentaureum GmbH & Co. KG Jarabak Light Wire pliers (REF 003-125-00). Then deburr the lateral edges of the fenestration down to the outer edge (see gold-colored line).

Check the deburred fenestration by inserting the head of the prosthetic screw for the angulated screw aperture laterally through the fenestration of the CAD/CAM titanium base VARIO.

Check that the screw aperture and the fenestration of the CAD/CAM titanium base VARIO are precisely aligned with each other. Fit the CAD/CAM titanium base VARIO to the crown and check whether the prosthetic screw for the angulated screw aperture can be inserted through the screw aperture.

Before bonding, cover the head of the prosthetic screw, which fixates the CAD/CAM titanium base VARIO, with wax. The bonding surfaces of the all-ceramic crown and of the CAD/CAM titanium base VARIO are blasted with aluminum oxide (50 µm/2 bar). It is also advisable to use the silanization method to condition this bonding area on the CAD/CAM titanium base VARIO. Before bonding, ensure the surfaces are dry and free from grease. Use adhesives according to the manufacturer's instructions (e.g.: "PANAVIA™F2.0" by Kuraray Noritake or "Multilink® Hybrid Abutment" by Ivoclar Vivadent AG) and position the all-ceramic crown. Dentaureum GmbH & Co. KG gives no warranty for the shelf life of the meso/superstructure and the adhesive. After bonding, remove the excess, loosen the prosthetic screw and remove the CAD/CAM titanium base VARIO from the laboratory implant. The all-ceramic crown is completed according to the manufacturer's instructions.

Important information:

- The wall thickness of the ceramic abutment must NOT be less than 0.5 mm.
- Preparation of a chamfer with angled inner edge and a minimum step of 0.5 mm.
- The contact surface between the implant and abutment must not be blasted or otherwise mechanically processed
- The bonding cylinder may only be shortened at the intended grooves using a SiC separating disc (Supercut STM separating disc, REF 130-114-00).
- The fenestration of the CAD/CAM titanium base VARIO must be opened in steps and completely with the aid of flat, round-nose pliers (Jarabak Light Wire pliers, REF 003-125-00). When breaking out, the pliers are placed directly at the transition to the separation line
- After complete opening of the fenestration, the outer edge of the fenestration is deburred using a hard metal bur (super hard metal bur, conical, round face; REF 123-593-00 super hard metal bur) up to the outer edge
- The diameter of the bonding cylinder must not be reduced

Tightening torque	
Prosthetic screw on the model	30 Ncm
AnoTite screw intra-orally	30 Ncm
Prosthetic screw for angled screw aperture on the model	25 Ncm
AnoTite screw for angled screw aperture in the mouth	25 Ncm

To achieve optimum preloading of the AnoTite screw, the AnoTite screw must be retightened with the same torque approx. 5 minutes after insertion.

7.6 Modeling sleeve for CAD/CAM titanium bases

The modeling sleeve can be ground with a milling cutter for acrylics and shaped into the desired form of the mesostructure or the crown framework with standard burn-out modeling waxes or acrylics.

The wax-up is then transferred to a framework using the casting technique and bonded to the CAD/CAM titanium base.

For additional information on bonding, see 7.4 CAD/CAM titanium base.

7.7 tioLogic® TWINFIT CAD/CAM titanium adhesive bases

CAD/CAM titanium adhesive bases are designed specifically for adhesion in the case of multi-unit bridge restorations or full-arch restorations in the edentulous mandible or maxilla. The titanium adhesive bases have a platform connector geometry for all series of abutments S, M and L that ensures the prosthesis sits optimally. There is no antirotational mechanism to enable maximum freedom when positioning the base on the implant. The cone prepared for bonding has retention grooves. It is 3.0 mm high and can compensate for divergences of up to 30°.

By using the bonding technique, the ceramic bridges, manufactured using CAD/CAM technology, can be fitted passively. Before bonding, the head of the prosthetic screw, which fixates the CAD/CAM titanium base, is covered with wax. The bonding surfaces of the ceramic sleeve and of the CAD/CAM titanium base are blasted with aluminum oxide (50 µm/2 bar). It is also advisable to use the silanization method to condition this bonding area on the CAD/CAM titanium base. Before bonding, ensure the surfaces are dry and free from grease. Adhesives must be used according to the manufacturer's instructions (e.g.: "PANAVIA™F2.0" by Kuraray Noritake or "Multilink® Hybrid Abutment" by Ivoclar Vivadent AG). Dentaureum GmbH & Co. KG gives no warranty for the shelf life of the meso/superstructure and the adhesive. After bonding, excess material is removed and the ceramic abutment is put into position (observing the antirotational mechanism). The all-ceramic crown is completed according to the manufacturer's instructions.

Important information:

- NOT approved for single restorations.
- The wall thickness of the ceramic abutment must NOT be less than 0.5 mm.
- Preparation of a chamfer with angled inner edge and a minimum step of 0.5 mm.
- The contact surface between the implant and abutment must not be blasted or otherwise mechanically processed
- The height and diameter of the cone prepared for bonding must not be reduced

Tightening torque	
Prosthetic screw on the model	30 Ncm
AnoTite screw intra-orally	30 Ncm

To achieve optimum preloading of the AnoTite screw, the AnoTite screw must be retightened with the same torque approx. 5 minutes after insertion.

7.8 Additional components

The prosthetic screw (not included in scope of delivery) serves to help the dental technician fixate the abutment on the working cast temporarily.

A new AnoTite screw should always be used for the final fitting and fixation of the abutment in the mouth at the prescribed torque.

Overview of additional components:

Torque ratchet
Hex key
Hex ball point key for angulated screw apertures
Insertion key – ratchet
PentaGrip insertion key – ISO shank
Adapter, ISO shank, hexagon/ratchet
Prosthetic screw on the model

Caution:

Do not confuse the prosthetic screw with the AnoTite screw. The prosthetic screw is colored silver and the AnoTite screw is anodized gold.

8. Technical specification

Composition

Product	Composition
tioLogic® TWINFIT CAD/CAM scan abutment, titanium	3.7165
tioLogic® TWINFIT laboratory implant, for printed and cast models	3.7165
tioLogic® TWINFIT CAD/CAM titanium blocks, PreForm	3.7165
tioLogic® TWINFIT CAD/CAM titanium bases	3.7165
tioLogic® TWINFIT CAD/CAM titanium bases VARIO	3.7165
Modeling sleeve for CAD/CAM titanium bases	PC
tioLogic® TWINFIT CAD/CAM titanium adhesive bases	3.7165
AnoTite screw	3.7165
AnoTite screw for angulated screw apertures	3.7165
Retaining screw	3.7165
Counter screw	3.7165

Information about the composition is included in the materials list for our implantology products (REF 989-801-05), or see www.dentaurum.com/ifu. The SSCP is available at <https://ec.europa.eu/tools/eudamed> and www.dentaurum.com.

9. Adverse reactions and residual risks

Possible complications and risks during prosthetic treatment (temporary adverse reactions) are:

Pain, swellings, difficulties speaking, inflammation of the gums.

Possible complications and risks during prosthetic treatment (permanent adverse reactions) are:

Fractured dentures, aesthetic problems.

10. Safety information and precautions

- Do not use the product if there is a known allergic reaction to one or more of the material components.
- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. Therefore, secure all articles used intra-orally against swallowing and/or aspiration.
- Different types of alloy in the oral cavity can lead to galvanic reactions.
- If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.
- All tioLogic® TWINFIT prosthetic components may only be used in combination with tioLogic® TWINFIT implants! Use of third-party components can negatively impact the function and safety. In such cases, Dentaureum GmbH & Co. KG accepts no liability for replacement products or guarantees.
- Observe the indications for the individual products. Consider all indications that are not listed explicitly to be contraindications.
- Disinfection, cleaning and sterilizing should not be done in the original packaging.
- Products must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.
- Technical advice (oral and written) is based on the scientific and technical knowledge available at the time the product is put on the market. It does not release the dental professional from the responsibility of personally checking the suitability of the products for the intended indications and applications. Advice is only given as a non-binding recommendation and cannot be assumed to provide any form of assurance, liability or guarantee.

Severe incidents:

Report all severe incidents arising from the use of the product to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.

Severe incidents are identified as malfunctions and a drop in performance or deterioration of the properties of a medical device which resulted or could have resulted in a direct or indirect temporary or permanent serious deterioration in a patient's state of health, in the death of a patient, a dental professional or other persons or as posing a serious danger to public health.

11. Storage and shelf life

Store the tioLogic® TWINFIT implants, abutments and components always in a dry, dark place at room temperature in the original packaging. If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.

12. Disposal

Adhere strictly to the rules that apply for your region and the details outlined in the safety data sheets, if available. Products to be disposed of should be treated as medical waste and be decontaminated in compliance with the relevant regulations for medical waste disposal.

13. Information for single use products

Prosthetic abutments and components are supplied non-sterile and are intended for single use on a single patient only. The product must be disinfected, cleaned and sterilized before first use.

Reconditioning of abutments and components (recycling) as well as reuse on patients are not permitted, as safe reconditioning and/or functional safety cannot be guaranteed and may lead to a loss in mechanical, chemical and/or biological properties. Reuse may cause infection. Do not use implants after the expiry date stated on the label.

14. Information for multiple use products

The scan abutments, titanium, are supplied non-sterile and are intended for multiple use. The product must be disinfected, cleaned and sterilized before being used initially and before being used again. Scan abutments, titanium, must be replaced if they cannot be clearly identified or if their function is impaired due, for example, to poor readability of the markings or labels.

15. Information on disinfection / cleaning / sterilizing

The tioLogic® TWINFIT prosthetic abutments and components are supplied non-sterile. They are for single use on one patient only, unless explicitly intended for multiple use. They must be cleaned and disinfected before and after being used on a patient (e.g., transfer to dental laboratory). Dentaaurum generally recommends that the product be sterilized in addition, as the use of non-sterile components may cause infections. Only sterilized prosthetic components may be used if bleeding occurs.

16. Product overview and scope of delivery

Please refer to the tioLogic® TWINFIT product catalog for the tioLogic® TWINFIT product range.

17. Quality information

Dentaaurum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaaurum has no influence on the processing on site.

18. Explanation of symbols

Please refer to the label. See chapter "Additional information".

19. Information on MRI safety

The products within these Instructions for use are MRI safe.

20. Additional information

Additional information can be found at www.dentaaurum.com/ifu

Instructions for use	REF
Screws	989-804-98
tioLogic® TWINFIT Additional components	989-805-88
Torque ratchet	989-719-00
Preparation instructions for the implant system components	989-801-07
Explanation of symbols used on the label	989-313-00
MRI safety of Dentaaurum products	989-804-42

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