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## Abbreviated instructions for use

**Dear customer,**

Thank you for choosing a quality product from Dentaum Implants.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit.

Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaum.com, even though you may frequently use the same product.

### 1. Manufacturer

Dentaum Implants GmbH | Turnstr. 31 | 75228 Ispringen | Germany

### 2. General description of product, intended use and characteristics, technical description

The tioLogic® TWINFIT implant system contains specially coordinated instruments, abutments, and accessory components for the prosthetic restoration of the tioLogic® TWINFIT implants. Only original components of the tioLogic® TWINFIT implant system should be combined in accordance with the Instructions for use / user manuals.

### 3. Area of application

For optimal gingiva management and for transferring the oral situation to the working cast.

### 4. Medical indications/contraindications

**Temporary abutment tioLogic® TWINFIT**

**Indications**

- Temporary restoration for up to 180 days in situ

**Contraindications**

- In situ longer than 180 days
- Primary bracing of abutments
- Single restoration with cantilever bridge unit
- Restorations with a length ratio to the length of the implant which exceeds 1:1.25

### 5. Information on use, processing, initial implementation and handling

**5.1 Closure screws tioLogic® TWINFIT**

After implant insertion and monitoring of the operation site, the wound can be closed by suturing. Ensure that the wound closure is saliva-proof and that there is good blood circulation.

Before wound closure, remove the closure screw from the packaging by pressing the hex key 1.3 into the head of the screw and pulling the screw off the silicone lid. The closure screw is then inserted into the implant and tightened. Closure screws should fit flush on the implant to ensure that bone tissue growth cannot penetrate into the implant. The closure screws are marked with S, M or L corresponding to the series of abutments.

Closure screws are intended for single use only.

**Tightening torque**

- Closure screw: manually or 15 Ncm

Should a closure screw lose its sterility during the operation, it is possible to purchase closure screws individually as a non-sterile product.

**Important: Closure screws purchased separately should be sterilized before they are inserted into the implant.**

**5.2 Gingiva former tioLogic® TWINFIT**

Gingiva formers, cylindrical or anatomically shaped, or 4Base abutments used directly – particularly gentle on the tissue – are available for the clinician to ensure optimal gingiva management. Gingiva former with platform connector geometry are available for the S, M and L series of abutments.

The anatomical gingiva formers are designed to form a wide gingival contour. Depending on the type of prosthetic restoration, this can make it easier for the clinician to fit the restoration.

The gingiva formers are selected according to the series of abutments, gingival height and insertion depth of the implant. They are available for the series of abutments S, M or L and in different gingival heights (laser-printed on gingiva former). Selection aids with platform connector geometry for the abutment series S, M and L are available to help choose the right gingiva former on the model (gingival height and angulation).

**Tightening torque**

Gingiva formers can also be used with open healing of the implant for specific indications and for preserving the soft tissue.

**IMPORTANT: The gingiva formers and 4Base abutments should be cleaned and sterilized before insertion in the implant.**

If a temporary restoration is fitted, the denture should be relieved during gingiva forming. The impression should not be taken until the tissue is completely free of inflammation.

**Tightening torque**

- Gingiva former: manually or 15 Ncm
- 4Base abutment: 35 Ncm
- 4Base closure screw: 15 Ncm

**5.3 Impression post tioLogic® TWINFIT**

The impression can be taken using either the open or closed technique. Relevant components with platform geometry for the abutment series S, M and L are available for both impression techniques.

**IMPORTANT: When positioning the impression post on the inserted implant, ensure it is seated without a gap. The impression component should NOT come into contact with the individual impression tray.**

In the case of removable restorations (restorations with 4Base, ball abutments, tioLOC) the impression can also be taken with other special impression components over the respective primary abutments.

Silicone or polyether impression materials are recommended for impression-taking due to their high precision and elastic recovery.

**Open impression technique:**

Impression posts are available in lengths of 10.0 mm and 14.0 mm with the corresponding screws to cater for different occlusal spaces. The impression posts are laser-printed with S, M or L on the retention surface and at the interface.

In order to make the abutment series more recognizable, the interface is additionally marked with dots:

- 1 dot is equivalent to abutment S
- 2 dots are equivalent to abutment M
- 3 dots are equivalent to abutment L

The screw is pushed down before fitting the impression post. This provides additional guidance when fitting the post. The inner connection is shorter with an open impression to ensure a compression-free impression even with divergent axes.

**Tightening torque**

- Screw impression post intra-orally: manually or 15 Ncm
- Screw impression post in the laboratory implant: manually or 15 Ncm

**Closed impression technique:**

Components for the closed impression technique include impression posts, retaining screws and impression caps. They are laser-printed or marked with the series of abutments S, M or L. The enclosed red retaining screw has a shortened thread which will only grip in the (laboratory) implant if the impression post has been inserted in the correct position into the connection. The corresponding impression aid S, M or L is fitted taking the vertical retention grooves into consideration until it perceptibly and audibly clicks into place.

**IMPORTANT: The impression caps are items for single use only and are NOT suitable for sterilization or multiple use.**

The design of the retention grooves ensures that they can be positioned without coming into contact with the adjacent teeth. The impression is taken according to standard criteria.

**Tightening torque**

- Retaining screw, intra-orally: manually, or 15 Ncm
- Retaining screw on laboratory implant: manually, or 15 Ncm

**5.4 Temporary abutment tioLogic® TWINFIT**

In aesthetically relevant areas, the peri-implant structures are retained with a temporary abutment.

It is possible to fit a long-term, non-functional immediate temporary restoration on implants if there is absolute primary stability and no recession of the implant site.

Temporary abutments with platform connector geometry are available for the S, M and L series of abutments.

It is supplied non-sterile and is made from high-strength plastic (PEEK), which can be quickly and easily customized. The temporary abutment can be faced directly with composite or fitted with a temporary crown or bridge. In both cases the abutment is secured intra-orally with the screw for the temporary abutment; the contours are marked and adjusted extra-orally. The implant should **NOT** be subjected to mechanical stress.

The restoration should only be shortened as far as the upper edge of the screw for the temporary abutment.

**Tightening torque**

- Screw for temporary abutment on the model; manually
- Screw for temporary abutment in the mouth: 15 Ncm

### 6. Composition

- Closure screws: titanium alloy (Ti6Al4V)
- Gingiva former: titanium alloy (Ti6Al4V)
- Impression posts:
  - open – impression post: titanium alloy (Ti6Al4V)/ screw: 1.4305 closed – impression post and screw: titanium alloy (Ti6Al4V)/ impression cap: POM
- Temporary abutment: PEEK (polyether ether ketone) /
- Retaining screw: titanium alloy (Ti6Al4V)
- Prosthetic screws:
  - AnoTite screw: titanium alloy (Ti6Al4V)
  - Prosthetic screw: titanium alloy (Ti6Al4V)
- Hex keys:
  - Hexagon key - ISO shank 1.3 (1.4197)
  - Hexagon key - ratchet 1.3 (1.4197)

### 7. Safety information

- The product should not be used if there is a known allergic reaction to one or more of the material components.
- Due to the small size, it is possible that an article be swallowed or aspirated. Aspiration can lead to difficulty in breathing or death due to asphyxiation. All articles used intra-orally should therefore be secured against swallowing and/or aspiration.
- All serious incidents arising from the use of the product should be reported to the manufacturer and the competent authority in the country in which the user and/or the patient is resident.
- Different types of alloy in the mouth can lead to galvanic reactions.
- Only use tioLogic® TWINFIT prosthetic components in combination with tioLogic® TWINFIT implants.

### 8. Warnings

Temporary restorations made of PEEK can be left in situ for a maximum of 180 days.

### 9. Storage and shelf life

Store in a dry place.

### 10. Information for single use products

Closure screws, gingiva formers, temporary abutments, impression posts and impression caps (closed impression-taking) are designed for single use only. Reconditioning of a closure screw, gingiva former, temporary abutment, impression post or impression cap that has been inserted previously (recycling) or reuse on patients is not permitted since it can then no longer be guaranteed that the article can be reprocessed safely or can function safely.