DE	Gebrauchsanweisung	· tioLogic® TWINFIT Verschlusschrauben, Gingivaformer, Abformpfosten, provisorische Aufbauten
EN	Instructions for use	· tioLogic® TWINFIT closure screws, gingiva formers, impression posts, temporary abutments
FR	Mode d'emploi	· Vis de couverture, conformateurs gingivaux, piliers d'empreinte, piliers provisoires tioLogic® TWINFIT
ES	Modo de empleo	· Tornillos de cierre, conformadores de encía, pilares de impresión, pilares provisionales tioLogic® TWINFIT
IT	Modalità d'uso	· Viti di chiusura, viti di guarigione, transfer d'impronta, abutment provvisori tioLogic® TWINFIT

1. Hersteller

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Dentaurum GmbH & Co. KG · Turnstr. 31 · 75228 Ispringen · Germany

2. General device description

Prosthetic abutments and components replace missing teeth, restore chewing function and thus improve the quality of life. The implants of the tioLogic® implant family are used to fixate corresponding abutment components, depending on the indication. The abutment components serve to fabricate prosthetic meso-constructions, which are then restored with a prosthetic superstructure.

Abutment components					
Series of abutments*	S	M		L	
Implant diameters in mm	3.3	3.7	4.2	4.8	5.5

*all of the abutment components listed in these Instructions for use are available in the S, M and L series of abutments.

tioLogic® TWINFIT closure screws are used to seal the freshly placed implant. The closure screws are screwed into the implant to prevent fluids from entering the implant and/or gingival tissue from growing into the implant during closed healing. The closure screws are supplied with platform connector geometry.

Closure screw			
Series of abutments	S	M	L

tioLogic® TWINFIT gingiva formers are intended for wound closure of the freshly placed implant while supporting gingiva management. The gingiva formers are screwed into the implant to prevent fluids from entering the implant and/or covering gingival tissue from growing into the implant during open healing. The gingiva formers are supplied with platform connector geometry.

Gingiva former			
Series of abutments	S	M	L
Heights in mm	2.0 / 4.0 / 6.0		

tioLogic® TWINFIT provisional abutments are used for time-limited long-term, non-functional, immediate restoration on implants and may be used when there is absolute primary stability and there is no recession of the implant site. The temporary abutments are supplied with platform connector geometry.

Temporary abutment			
Series of abutments	S	M	L
Length in mm	12.0		

Screws for tioLogic® TWINFIT temporary abutments are used to secure temporary abutments for time-limited long-term temporary, non-functional immediate restoration on (laboratory) implants. They are supplied with the respective temporary abutment.

Screws for tioLogic® TWINFIT temporary abutments			
Series of abutments	S	M	L
Screw length in mm	13.5		

tioLogic® TWINFIT impression components with a platform connector geometry are available for the open and closed impression-taking methods. The tioLogic® TWINFIT impression posts, open, are provided with the corresponding screw for impression posts. The tioLogic® TWINFIT impression post, closed, is supplied with a tioLogic® TWINFIT impression cap and a retaining screw.

Impression components			
Series of abutments	S	M	L
Length of impression post, open in mm	10.0 / 14.0		
Length of screw for impression post, open in mm	16.0 / 25.0		
Length of impression post, closed in mm	9.5		
Length of retaining screw in mm	14.0		

The tioLogic® TWINFIT implant system contains Implant system with specially coordinated abutments and components incl. instruments.

3. Intended purpose

Closure screws are intended to close the inserted implant or abutment during the healing period to prevent penetration of saliva.

Gingiva formers are intended to close the inserted implant during the healing period to prevent penetration of saliva and to shape the surrounding soft tissue for the subsequent prosthetic restoration.

Abutments are intended to connect the implants with the prosthetic restoration.

Screws are used to fix abutments on the implants.

Impression components are intended for three-dimensional recording and transfer of the intraoral position of the implants.

4. Intended users

The tioLogic® TWINFIT product range is supplied exclusively to qualified personnel such as surgeons, dentists and dental technicians. It should only be used by surgeons, dentists or dental technicians who are familiar with dental implantology, including diagnosis, preoperative planning, surgical techniques and prosthetic restorations.

5. Medical indications / contraindications

5.1 Indications

5.1.1 tioLogic® TWINFIT closure screws

- Implant closure with covered/closed healing

5.1.2 tioLogic® TWINFIT gingiva formers

- Implant closure with open healing
- Shaping of the gingiva / soft tissue during healing

5.1.3 tioLogic® TWINFIT temporary abutment

- Implant closure with open healing
- Shaping of the gingiva / soft tissue during healing
- Time-limited long-term provisional, non-functional immediate restoration
- Maximum time in situ up to 180 days

5.1.4 Screw for temporary abutment

- Use exclusively with the temporary abutment
- For final fitting, an unused AniTite screw is recommended!

5.1.5 tioLogic® TWINFIT impression components

- Implant-supported impression taking of single teeth and partially or fully edentulous jaws

5.2 Contraindications

5.2.1 tioLogic® TWINFIT closure screws

- Use as a temporary as well as a long-term restoration.
- Intended for the healing stage only.

5.2.2 tioLogic® TWINFIT gingiva formers

- Use as a temporary as well as a long-term restoration.
- Intended for the healing stage only.

5.2.3 tioLogic® TWINFIT temporary abutment

- Time in situ longer than 180 days
- Correction of angulation exceeding 20°, or exceeding 15° in the case of implants \varnothing 3.3 mm
- Primary bracing of abutments
- Single restoration with cantilever bridge unit
- Restorations with a length ratio to the length of the implant which exceeds 1:1.25

5.2.4 Screw for temporary abutment

- Use with abutments other than on tioLogic® ST and tioLogic® TWINFIT temporary abutments

5.2.5 tioLogic® TWINFIT impression components

- Use as restoration
- Use of the supplied screws on other prosthetic components

6. Patient target group

All skeletally mature patients undergoing implant treatment. Please consider the indications / contraindications.

7. Information on use, processing, fitting and insertion

Prior to use, the dental professional ensures that they have carefully read and understood the tioLogic® TWINFIT Instructions for use.

We strongly recommend that, prior to using Dentaaurum products, dental professionals attend a Dentaaurum training course to familiarize themselves with the Dentaaurum products.

Carry out thorough radiological, psychological and physical examinations before implantation and inform patients of the process in detail. Dentaaurum recommends full clinical, radiological, photographic and statistical documentation.

The closure screws, gingiva formers, impression posts, provisional abutments and components of the tioLogic® TWINFIT implant system can be documented, e.g. in the patient file or PatientPass, using the additional labels.

The implantologist should ensure the products cannot be aspirated during intra-oral use.

Precautions:

Before surgical implantation, it is important to carry out a careful patient anamnesis and, if required, a third-party anamnesis by the patient's doctor (GP/PCP) in order to clarify if 1) implantation may be difficult due to an anatomical situation, 2) there is the risk of a serious surgical problem or a general risk, 3) wound healing/osseointegration may be difficult, or 4) good hygiene and care of the implant, abutment and restoration may worsen. Do not fit dental implants if any of the above cases or a combination hereof is severe or uncontrollable.

Prosthetic restoration of the implant must not be performed until the soft tissue has healed without irritation. Furthermore, a radiographic check-up should be performed before the prosthetic restoration.

7.1 tioLogic® TWINFIT closure screw

After implant insertion and monitoring of the operation site, the wound can be closed by suturing. Ensure that the wound closure is saliva-proof and that there is good blood circulation.

To remove the closure screw, press the hex key 1.3 into the closure screw and pull it from the silicone lid. Insert the closure screw into the implant and tighten. Closure screws must fit flush on the implant to ensure that bone tissue and soft tissue growth cannot penetrate into the implant. The closure screws are inscribed with S, M or L corresponding to the series of abutments.

Closure screws are intended for single use only.

If a closure screw loses its sterility during the operation, it is possible to purchase closure screws individually as a non-sterile product.

IMPORTANT: Closure screws purchased separately should be cleaned and sterilized before they are inserted into the implant.

Tightening torque	
Closure screw	15 Ncm

7.2 tioLogic® TWINFIT gingiva formers

Gingiva formers, cylindrical or anatomically shaped, or 4Base abutments used directly – particularly gentle on the tissue – are available for the implantologist to ensure optimal gingiva management. Gingiva former with platform connector geometry are available for the S, M and L series of abutments.

The anatomical gingiva formers are designed to form a wide gingival contour. Depending on the type of prosthetic restoration, this can make it easier for the implantologist to fit the restoration.

The gingiva formers are selected according to the series of abutments, gingival height and insertion depth of the implant and are screwed into the implant using the AF 1.3 hex key. They are available for the series of abutments S, M or L and the different gingival heights (2.0 mm, 4.0 mm or 6.0 mm). Selection aids with conical and platform connector geometry for the abutment series S, M and L are available to help choose the gingiva former and 4Baseabutment on the model (gingival height and angulation).


Gingiva formers can also be used with open healing of the implant for preserving the soft tissue.

IMPORTANT: Clean and sterilize the gingiva formers and 4Base abutments before insertion in the implant.

If a temporary restoration is fitted, relieve the denture during gingiva forming. The impression should not be taken until the tissue is completely free of inflammation.

Tightening torque	
Gingiva former	15 Ncm

7.3 tioLogic® TWINFIT temporary abutment

It is possible to fit a time-limited long-term, non-functional immediate temporary restoration on implants if there is absolute primary stability and no recession of the implant site. In aesthetically relevant areas, the peri-implant structures are retained with a temporary abutment. An optimum impression can be taken after formation of the structures.

A temporary abutment with platform connector geometry is available for the S, M and L series of abutments. They are supplied non-sterile and are made from high-strength plastic (PEEK), which can be quickly and easily customized.

The temporary abutment can be faced directly with composite or fitted with a temporary crown or bridge. In both cases, secure the abutment intra-orally with the screw for the temporary abutment. Tighten the screw with the hex key AF 1.3. Mark the contours and adapt extra-orally. The restoration should only be shortened as far as the upper edge of the screw for the temporary abutment.

With a direct build-up of the facing, face the temporary abutment with composite extra-orally and then secure to the implant using the correct torque. Seal the screw opening with acrylic. With a crown restoration, fit the temporary abutment before you seal the screw aperture with wax and place the temporary restoration. The crown should only be retained with temporary cement.

Important:

The implant should **NOT** be subjected to mechanical stress

Temporary abutments can be in place for a maximum of 180 days

Tightening torque	
Screw for temporary abutment on the model / in the mouth	15 Ncm

7.4 tioLogic® TWINFIT impression components

The impression can be taken using either the open or closed technique. Relevant components with platform geometry for the abutment series S, M and L are available for both impression techniques.

IMPORTANT:

When positioning the impression post on the inserted implant, ensure it is seated without a gap. Please also note that the impression component should **NOT** come into contact with the individual impression tray.

In the case of removable dental restorations (restorations with 4Base abutments, ball abutments, tioLOC), the impression can also be taken with other special impression components over the respective primary abutments.

Silicone or polyether impression materials are recommended for impression-taking due to their high precision and elastic recovery.

Marking with interface M.

7.4.1 Open impression technique

Impression posts are available in lengths of 10.0 mm and 14.0 mm with the corresponding screws to cater for different occlusal spaces.

The impression posts are laser-marked with S, M or L on the retention surface and at the interface.

In order to make the abutment series more recognizable, the interface is additionally marked with dots:

- 1 dot is equivalent to abutment series S
- 2 dots are equivalent to abutment series M
- 3 dots are equivalent to abutment series L

After the impression has been taken, an individual tray is fabricated. This is strengthened and perforated in the region of the implants.

Groove on sure-grip screw

The closure screws, gingiva formers and the temporary restoration should be removed prior to taking the impression.

Push down the screw before fitting the impression post. This provides additional guidance when fitting the post. The inner connection is shorter with an open impression post to ensure a compression-free impression even with divergent axes. For exact implant-supported registration of the maxillomandibular relationship, the impression should be taken only on the basis of the platform connector geometry. Position the impression post according to the

S, M or L series of abutment on the implant until the rotational safety stop engages. A correct fit of the implant post on the implant shoulder is indicated when a visible mark on the screw is level with the upper edge of the impression post (screw should only be inserted and not tightened). If the rotational safety stop is not engaged, then the marking on the screw is not visible.


In addition to the marking, the red sure-grip screw also has a shortened thread which will only grip in the (laboratory) implant if the impression post has been inserted in the correct position into the connection.

If either of the two cases occurs, realign the impression post and check to ensure that it fits correctly (X-ray check for correct fit of the impression post). Then tighten the sure-grip screw using the AF 1.3 hex key or by hand.

When preparing the impression tray, ensure that there is no contact between the impression posts or screws and the tray at the perforations. The impression is to be taken with a material based on silicone or polyether. The impression posts are secured in the impression material by the retention. Ensure that the peri-implant region is accurately reproduced in the impression. Care must be exercised to ensure that the screws are still accessible after the impression has been taken. To remove the impression tray, loosen and retract the screws, and remove the impression tray. Send the impression tray together with the screws to the dental technician. Use the practice record card (REF 989-966-22) to provide the dental technician with all the relevant information. The respective gingiva formers and temporary restorations are refitted after impression-taking.

Tightening torque	
Screw for open impression post on the model/intra-orally	15 Ncm

7.4.2 Closed impression technique

For exact implant-supported registration of the maxillo-mandibular relationship, the impression should be taken only on the basis of the platform connector geometry.

Components for the closed impression technique include impression posts, retaining screws and impression caps. Impression posts and caps are laser-printed or marked with the S, M or L series of abutments.

After the impression has been taken, an individual tray is fabricated. It is designed to be large enough in the region of the implants so it does not come into contact with the impression posts, including the attached impression caps.

Remove the cover screws, the gingiva formers and the temporary restoration prior to taking the impression.

For easier insertion, the impression posts for closed-tray impressions can be mounted on positioning aids in the corresponding series of abutments (S, M, or L). Fit the impression post for closed impression-taking with the long tube into the implant. Then align the rotational safety stop so the PentaStop interlocks. The screw connection can be made occlusally through the positioning aid. The enclosed red retaining screw has a shortened thread which will only grip in the (laboratory) implant if the impression post has been inserted in the correct position into the connection. Align the impression post if necessary and check to ensure that it fits correctly (X-ray check of the proper fit of the impression post). Tighten the retaining screw with the hex key AF 1.3.

Fit the corresponding impression cap S, M or L according to the vertical retention grooves on the impression post until it perceptibly and audibly clicks into place.

IMPORTANT:

The impression caps are items for single use only and are NOT suitable for sterilization or multiple use. Multiple use results in transfer inaccuracies.

The orientation of the retention grooves ensures that they can be positioned without coming into contact with the adjacent teeth. The impression is to be taken with a material based on silicone or polyether. The impression caps are secured in the impression material by the retention. Ensure that the peri-implant region is accurately reproduced in the impression. Remove the impression tray after curing the impression compound. Send the impression posts with fixing screws together with the impression to the laboratory. Use the practice record card (REF 989-966-22) to provide the dental technician with all the relevant information. The respective gingiva formers are refitted after impression-taking.

Tightening torque	
Retaining screw on the model / in the mouth	15 Ncm

7.5 Additional compounds

The prosthetic screw (not included in scope of delivery) serves only to help the dental technician fix the abutment on the working cast temporarily.

A new AniTite screw should always be used for the final fitting and fixation of the abutment in the mouth at the prescribed torque.

Overview of additional components:

Torque ratchet
Hex key
Adapter, ISO shank, hexagon/ratchet

Caution:

Do not confuse the prosthetic screw with the AniTite screw. The prosthetic screw is colored silver and the AniTite screw is anodized gold.

8. Technical specification

Composition

Product	Composition
tioLogic® TWINFIT closure screw	3.7165
tioLogic® TWINFIT gingiva former	3.7165
tioLogic® TWINFIT temporary abutment	PEEK / 3.7165
tioLogic® TWINFIT impression post, open	3.7165 / VMQ
tioLogic® TWINFIT impression post, closed	3.7165 / POM
tioLogic® TWINFIT impression cap	POM
Positioning aid	1.4197

Information about the composition is included in the material list for our implantology products (REF 989-801-05), or see www.dentaurum.com/ifu. The SSCP is available at <https://ec.europa.eu/tools/eudamed> and www.dentaurum.com.

9. Adverse reactions and residual risks

Possible complications and risks during prosthetic treatment (temporary adverse reactions) are:

Pain, swellings, difficulties speaking, inflammation of the gums.

Possible complications and risks during prosthetic treatment (permanent adverse reactions) are:

fractured dentures, aesthetic problems.

10. Safety information and precautions

- Do not use the product if there is a known allergic reaction to one or more of the material components.
- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. Therefore, secure all articles used intra-orally against swallowing and/or aspiration.
- Different types of alloy in the oral cavity can lead to galvanic reactions.
- If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.
- Only use all tioLogic® TWINFIT prosthetic components in combination with tioLogic® TWINFIT implants! Use of third-party components can negatively impact the function and safety. In such cases, Dentaureum GmbH & Co. KG accepts no liability for replacement products or guarantees.
- Observe the indications for the individual products. Consider all indications that are not listed explicitly to be contraindications.
- Disinfection, cleaning and sterilizing should not be done in the original packaging.
- Products must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.
- Technical advice (oral and written) is based on the scientific and technical knowledge available at the time the product is put on the market. It does not release the dental professional from the responsibility of personally checking the suitability of the products for the intended indications and applications. Advice is only given as a non-binding recommendation and cannot be assumed to provide any form of assurance, liability or guarantee.

Severe incidents:

Report all severe incidents arising from the use of the product to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.

Severe incidents are identified as malfunctions and a drop in performance or deterioration of the properties of a medical device which resulted or may have resulted in a direct or indirect temporary or permanent deterioration in a patient's state of health, in the death of a patient, a dental professional or other persons or as posing a serious danger to public health.

11. Storage and shelf life

Store the tioLogic® TWINFIT implants, abutments and components always in a dry, dark place at room temperature in the original packaging. If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.

12. Disposal

Adhere strictly to the rules that apply for your region and the details outlined in the safety data sheets, if available. Products to be disposed of should be treated as surgical waste and be decontaminated in compliance with the relevant regulations for medical waste disposal.

13. Information for single use products

Prosthetic abutments and components are supplied non-sterile and are intended for single use on a single patient only. The product must be disinfected, cleaned and sterilized before first use.

Reconditioning of abutments and components (recycling) as well as reuse on patients are not permitted as safe reconditioning and/or functional safety cannot be guaranteed and may lead to a loss in mechanical, chemical and/or biological properties. Reuse may cause infection. Do not use implants after the expiry date stated on the label.

14. Information for multiple use products

The positioning aids are supplied non-sterile and are intended for multiple use. The product must be disinfected, cleaned and sterilized before being used initially and before being used again. Titanium scan abutments must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.

15. Information on disinfection / cleaning / sterilizing

The tioLogic® TWINFIT prosthetic abutments and components are supplied non-sterile. They are for single use on one patient only. They must be cleaned and disinfected before and after being used on a patient (e.g. transfer to dental laboratory). Dentaureum generally recommends that the product be sterilized in addition, as non-sterile components may cause infections. Only sterilized prosthetic components may be used if bleeding occurs.

16. Product overview and scope of delivery

Please refer to the tioLogic® TWINFIT product catalog for the tioLogic® TWINFIT product range.

17. Quality information

Dentaureum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaureum has no influence on the processing on site.

18. Explanation of symbols

Please refer to the label. See chapter "Additional information".

19. Information on MRT safety

The products within these Instructions for use are MRT safe.

20. Additional information

Additional information can be found at www.dentaureum.com/ifu

Instructions for use	REF
4Base abutments and components	989-365-93
Screw	989-804-98
tioLogic® TWINFIT Additional components	989-805-88
Torque ratchet	989-719-00
Preparation instructions for the implant system components	989-801-07
Explanation of symbols used on the label	989-313-00
MRI safety of Dentaureum products	989-804-43

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