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Abbreviated instructions for use

Dear customer,

Thank you for choosing a quality product from Dentaurum Implants.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit.

Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaurum.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurum Implants GmbH | Turnstr. 31 | 75228 Ispringen | Germany

2. General description of product, intended use and characteristics, technical description

For the prosthetic restoration of the tioLogic® TWINFIT implants, the tioLogic® TWINFIT implant system contains specifically coordinated instruments, abutments, and accessory components. Only original components of the tioLogic® TWINFIT implant system should be combined in accordance with the Instructions for use / user manuals.

The abutment is selected by the dentist/dental technician to suit the individual patient's needs.

The diameter of the implant and the gingival thickness determine the choice of abutment.

Selection aids with conical and platform connector geometry for the abutment series S, M and L are available to help choose the right abutment on the model (gingival height and angulation).

3. Area of application

The titanium abutments tioLogic® TWINFIT are prosthetic components for single tooth and bridge restorations on tioLogic® TWINFIT implants.

4. Medical indications/contraindications

Titanium abutments tioLogic® TWINFIT

Indications

- Individual crown and bridge work (cemented, individually screw-retained)
- Abutment on implant Ø 3.3 mm only in the area of the upper lateral incisors and lower incisors

Contraindications

- Primary bracing of abutments
- Correction of angulation exceeding 20°, or exceeding 15° in the case of implants Ø 3.3 mm
- Cast-on technique

5. Information on use, processing, initial implementation and handling

Titanium abutment tioLogic® TWINFIT

All titanium abutments tioLogic® TWINFIT are available in conical and platform connector geometry for the abutment series S, M and L.

The straight and angled (S 15°, M 20°, L 20°) titanium abutments can be slightly altered to adapt to the gingival line (straight 1.5 mm, 2.5 mm, 4.0 mm / angled 1.5 mm, 3.0 mm). These abutments have an exactly defined crown margin and an integrated rotational stop.

The angled titanium abutments are also available in different gingival heights (labial / palatal).

The universal titanium abutments are anatomically shaped and can be prepared as required for the individual case.

For digital manufacture of individual single tooth and bridge restorations, CAD/CAM titanium blocks (PreForms) are available for one-piece abutments and CAD/CAM titanium bases for hybrid abutments.

The prosthetic screw (not included in delivery) serves only to help the dental technician fix the abutment on the working cast. We recommend using a new AnoTite screw for the final fitting in the mouth.

Tightening torque

- Prosthetic screw, on model: manually
- AnoTite screw intra-orally: 30 Ncm

Attention: Do not confuse the prosthetic screw with the AnoTite screw. The prosthetic screw is colored silver and the AnoTite screw is anodized gold.

6. Composition

- Titanium abutments: titanium alloy (Ti6Al4V)

■ Prosthetic screws:

- AnoTite screw: titanium alloy (Ti6Al4V)
- Prosthetic screw: titanium alloy (Ti6Al4V)

■ Hex key:

- Hex key - ISO shank SW 1.3: 1.4197
- Hex key - ratchet SW 1.3: 1.4197

■ Hexagonal ball end key - ISO shank SW 1.3: 1.4305

7. Safety instructions

■ The product should not be used if there is a known allergic reaction to one or more of the material components.

■ Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. All articles used intra-orally should therefore be secured against swallowing and/or aspiration.

■ All serious incidents arising from the use of the product should be reported to the manufacturer and the competent authority in the country in which the user and/or the patient is resident.

■ Different types of alloy in the mouth can lead to galvanic reactions.

■ Only use tioLogic® TWINFIT prosthetic components in combination with tioLogic® TWINFIT implants.

8. Storage and shelf life

Store in a dry place.

9. Information for single use products

Titanium abutments are designed for single use only. Reconditioning of titanium abutments that have been inserted previously (recycling) or reuse on patients is not permitted since it can then no longer be guaranteed that the article can be reprocessed safely or can function safely.

10. Information for multiple use products

The hex key SW 1.3 and the hexagonal ball end key SW 1.3 are intended for multiple use. The product must be disinfected, cleaned and sterilized before reuse. Additional information can be found at www.dentaurum.com (Processing Instructions Instruments and Accessories REF 989-801-09).

11. Information on disinfection/cleaning/sterilizing

tioLogic® TWINFIT prosthetic components are delivered non-sterile. They are for single use on one patient only. They must be cleaned and disinfected before being used on a patient. Dentaurum generally recommends that the product be sterilized in addition. Only sterilized prosthetic components may be used if bleeding occurs. Additional information can be found at www.dentaurum.com (Processing Instructions Instruments and Accessories REF 989-801-09).

12. Product overview / scope of delivery, incl. REF

Please refer to the tioLogic® TWINFIT product catalog for the tioLogic® TWINFIT product range.

13. Quality information

Dentaurum Implants ensures faultless quality of its products. These recommendations are based upon Dentaurum's own experiences. The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as Dentaurum Implants has no influence on the processing on site.

14. Explanation of symbols

 Refer to the label. Additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00).

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