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Short instructions

Dear customer,

Thank you for choosing a quality product from Dentaurum Implants. It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative. As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaur.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurum Implants GmbH | Turnstr. 31 | 75228 Ispringen | Germany

2. General description of product, intended use and characteristics, technical description

The tioLogic® TWINFIT implant system contains specially coordinated instruments, abutments and accessories for the prosthetic restoration of tioLogic® TWINFIT implants. Only original components of the tioLogic® TWINFIT implant system should be combined in accordance with the instructions for use/ user manuals.

The abutment is selected by the dentist/dental technician to suit the individual patient’s needs.

The diameter of the implant and the gingival thickness determine the choice of abutment.

Selection aids with conical and platform connector geometry for the abutment series S, M and L are available to help choose the right abutment on the model (gingival height and angulation).

3. Area of application

The titanium abutments tioLogic® TWINFIT are prosthetic components for single tooth and bridge restorations on tioLogic® TWINFIT implants.

4. Medical indications/contraindications

Titanium abutments tioLogic® TWINFIT

Indications

- Individual crown and bridgework (cemented, individually screw-retained)
- The abutment can only be used on implants Ø 3.3 mm in the area of the lateral incisors and lower incisors

Contraindications

- Primary bracing of abutments
- Correction of angulation exceeding 20°, or exceeding 15° in the case of implants Ø 3.3 mm
- Cast-on technique

5. Information on use, processing, initial implementation and handling

Titanium abutment tioLogic® TWINFIT

All titanium abutments tioLogic® TWINFIT are available in conical and platform connector geometry for the abutment series S, M and L.

The straight and angled (S 15°, M 20°, L 20°) titanium abutments can be slightly altered to adapt to the gingival line (straight 1.5 mm, 2.5 mm, 4.0 mm/angled 1.5 mm, 3.0 mm). These abutments have an exactly defined crown margin and an integrated rotational stop.

The angled titanium abutments are also available in different gingival heights (labial/palatal).

The universal titanium abutments are anatomically shaped and can be prepared as required for the individual case.

For digital manufacture of individual single tooth and bridge restorations, CAD/CAM titanium blocks (PreForms) are available for one-piece abutments and CAD/CAM titanium bases for hybrid abutments.

tioLogic® CAD/CAM data records for 3shape™, dental wings and exocad® can be downloaded free of charge under www.dentaurum.com/digital. The data records were created and verified in collaboration with these manufacturers.

They include all the information needed for all types of restorations as a complete package, i.e. for hybrid abutments, one-piece abutments, bridges and bars.

After selecting the indication in the data record of the respective CAD software, the situation in question can be scanned with the help of scan abutments. The matching process and design are carried out according to the instructions of the software manufacturer and dental prosthetic rules.

Additional components

- Scan abutment tioLogic® TWINFIT

The prosthetic screw (not included in delivery) serves only to help the dental technician fix the abutment on the working cast. We recommend using a new AnoTite prosthetic screw for the final fitting in the mouth.

Tightening torque

- Prosthetic screw, on model: manually
- AnoTite screw, intra-orally: 30 Ncm

Attention: Do not confuse the prosthetic screw with the AnoTite screw. The prosthetic screw is colored silver and the AnoTite screw is anodized gold.

6. Composition

- Titanium abutments: titanium alloy (Ti6Al4V)
- CAD/CAM titanium base: titanium alloy (Ti6Al4V)
- CAD/CAM titanium block: titanium alloy (Ti6Al4V)
- Scan abutment: PEEK (polyether ether ketone)
- Prosthetic screws:
 - AnoTite screw: titanium alloy (Ti6Al4V)
 - Prosthetic screw: titanium alloy (Ti6Al4V)

7. Safety instructions

- The product should not be used if there is a known allergic reaction to one or more of the material components.

- Different types of alloy in the mouth can lead to galvanic reactions.

- Only use the tioLogic® TWINFIT prosthetic components in combination with tioLogic® TWINFIT implants!

8. Storage and shelf life

Store in a dry place.

9. Information for single use products

Titanium abutments are designed for single use only. Reconditioning of titanium abutments that have been inserted previously (recycling) or reuse on patients are not permitted.

10. Information on disinfection/cleaning/sterilizing

tioLogic® TWINFIT prosthetic components are delivered non-sterile. They are for single use on one patient only. They must be disinfected, cleaned and, where applicable, sterilized before and after being used on a patient.

11. Product overview/scope of delivery, incl. REF

Please refer to the tioLogic® TWINFIT product catalog for the tioLogic® TWINFIT product range.

12. Quality information

Dentaurum Implants ensures faultless quality of its products. These recommendations are based upon Dentaurum’s own experiences. The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as Dentaurum Implants has no influence on the processing on site.

13. Explanation of symbols

 Refer to the label. Additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00).

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