---

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Dentaurum GmbH & Co. KG · Turnstr. 31 · 75228 Ispringen · Germany

## 2. General device description

Prosthetic abutments and components replace missing teeth, restore chewing function and thus improve the quality of life. The implants of the tioLogic® implant family are used to fixate corresponding abutment components, depending on the indication. Prosthetic abutments and components serve to produce prosthetic meso-constructions, which are then restored with a prosthetic superstructure.

Abutment components						
Series of abutments*	S		M		L	
Implant diameters in mm	3.3	3.7	4.2	4.8	5.5	

\*all of the abutment components listed in these Instructions for use are available in the S, M and L series of abutments as well as with platform or connector geometry.

tioLogic® TWINFIT titanium frameworks are prefabricated blanks/components for machined production of customized, individual single-piece abutments. The titanium abutments are available in various basic designs, such as straight, angled, or anatomical shapes, as well as with different gingival heights. The titanium abutments are labeled according to the S, M, or L series of abutments, and are provided either with conical or platform connector geometry, along with a AnOTite screw for each.

Titanium abutments			
Conical, straight, GH in mm	1.5	2.5	4.0
Conical, straight, length in mm	8.2	9.5 / 12.0	11.2
Conical, angled, GH in mm	1.5	3.0	-
Conical, angled, length in mm	8.7	10.2	-
angulation	S 15°, M 20°, L 20°	S 15°, M 20°, L 20°	-
Platform, straight, GH in mm	1.0	2.5	4.0
Platform, straight, length in mm	8.0	9.5 / 12.0	11.0
Platform, angled, GH in mm	1.5	3.0	-
Platform, angled, length in mm	8.7	10.2	-
angulation	S 15°, M 20°, L 20°	S 15°, M 20°, L 20°	-

The tioLogic® TWINFIT implant system contains specially coordinated abutments and accessories, including instruments, for the insertion and prosthetic restoration of the implants.

## 3. Intended purpose

Abutments are intended to connect the implants with the prosthetic restoration.

## 4. Intended users

The tioLogic® TWINFIT product range is supplied exclusively to qualified personnel such as surgeons, dentists and dental technicians. It should only be used by qualified personnel who are familiar with dental implantology, including diagnosis, preoperative planning, surgical technique and prosthetic treatment.

## 5. Medical indications / contraindications

### 5.1 Indications

#### 5.1.1 tioLogic® TWINFIT titanium abutments

- For the production of individual crown and bridge techniques as a connecting piece between the implant and the prosthetic abutment (cemented, individually screw-retained)
- Adhesive bond with ceramic mesostructures
- Abutment on implant Ø 3.3 mm only for replacement in the area of the upper lateral incisors and lower incisors

### 5.2 Contraindications

- In patients with a known intolerance to individual components.
- The contact surface between the implant and abutment must not be blasted or otherwise mechanically processed
- Insufficient oral hygiene
- Bruxism

#### 5.2.1 tioLogic® TWINFIT titanium abutments

- Primary bracing of titanium blocks and titanium bases
- Correction of angulation exceeding 20°, or exceeding 15° in the case of implants Ø 3.3 mm
- Cast-on technique
- Direct laser welding
- Restorations with a length ratio to the length of the implant which exceeds 1:1.25
- Single crowns with a cantilever bridge unit
- Abutment on implant Ø 3.3 mm only for replacement in the area of the upper and lower posterior teeth

## 6. Patient target group

All skeletally mature patients undergoing implant treatment. Please consider the indications / contraindications.

## 7. Information on use, processing, fitting and insertion

Prior to use, the dental professional ensures that they have carefully read and understood the tioLogic® TWINFIT Instructions for use.

We strongly recommend that, prior to using Dentaurum products, dental professionals attend a Dentaurum training course to familiarize themselves with the Dentaurum products. Carry out thorough radiological, psychological and physical examinations before implantation and inform patients of the process in detail. Dentaurum recommends full clinical, radiological, photographic and statistical documentation.

The titanium abutments and components of the tioLogic® TWINFIT implant system can be documented, e.g. in the patient file or PatientPass, using the additional labels.

### Precautions:

Before surgical implantation, it is important to carry out a careful patient anamnesis and, if required, a third-party anamnesis by the patient's doctor (GP/PCP) in order to clarify if 1) implantation may be difficult due to an anatomical situation, 2) there is the risk of a serious surgical problem or a general risk, 3) wound healing/osseointegration may be difficult, or 4) good hygiene and care of the implant, abutment and restoration may worsen. Do not fit dental implants if any of the above cases or a combination hereof is severe or uncontrollable.

Prosthetic restoration of the implant must not be performed until the soft tissue has healed without irritation. Furthermore, a radiographic check-up should be performed before the prosthetic restoration.

### Important:

**The implantologist should ensure the products cannot be aspirated during intra-oral use.**

## 7.1 tioLogic® TWINFIT titanium abutments

The tioLogic® TWINFIT abutments with a conical and platform connector geometry are available and marked according to the S, M or L series of abutments.

The straight and angled (S 15°, M 20°, L 20°) titanium abutments can be slightly altered to adapt to the gingival line (straight 1.5 mm, 2.5 mm, 4.0 mm / angled 1.5 mm, 3.0 mm). These abutments have an exactly defined crown margin and an integrated rotational stop. The angled titanium abutments are also available in different gingival heights (labial / palatal).

The trimmable titanium abutments (GH 2.5 mm and length 12.0 mm) are anatomically shaped and can be prepared as required for the individual case.

All tioLogic® TWINFIT titanium abutments can be used in combination with cemented or screw-retained solutions, both for bridges and crowns in one-tooth restorations.

The occlusal space and the axial direction must be checked when the abutment is inserted. Both can be altered by marking where necessary.

Do not overheat the titanium abutment during preparation as this can lead to differences in surface hardening (alpha-case layer). This could complicate or handicap further processing. This should primarily be considered for trimmable titanium abutments, where it may be necessary to prepare or remove more material. Fix the mesostructure onto the model, cast, and finish.

A new AnoTite screw should be used for the final fitting in the mouth. Tighten the screws manually with hex screw head geometry with the hex key ratchet AF 1.3 or mechanically with the hex key ISO shank AF 1.3.

### Tip:

In order to achieve an accurate fit of the bridge later, fabricate the crown caps and connect them using the same fully burn-out resin and, once cured, splint them together using the same material. Next, fabricate the framework according to dental technical prerequisites (reduced final anatomical shape). The waxes from Dentaureum which burn without leaving a residue (StarWax range) are particularly suitable.

For easier decementation, a notch can be integrated into the metal framework design in the palatal / cervical region, which is not veneered (attachment for decementation tool, i.e. crown remover).

The crown is veneered according to the manufacturer's instructions (e.g. with ceraMotion®, Dentaureum GmbH & Co. KG). If the restoration framework is made from titanium, a suitable titanium ceramic should be used, such as ceraMotion® Ti from Dentaureum GmbH & Co. KG. Design the proximal areas so that the implant neck can be cleaned with an interdental brush. The metal alloy is cast and processed according to the relevant manufacturer's instructions.

### Important:

**It is essential that the framework has an absolutely passive fit on the titanium abutments.**

Tightening torque	
Prosthetic screw on the model	30 Ncm
AnoTite screw intra-orally	30 Ncm

To achieve optimum preloading of the AnoTite screw, the AnoTite screw must be retightened with the same torque approx. 5 minutes after insertion.

## 7.2 Additional components

The prosthetic screw (not included in scope of delivery) serves only to help the dental technician fix the abutment on the working cast temporarily.

A new AnoTite screw should always be used for the final fitting and fixation of the abutment in the mouth at the prescribed torque.

Overview of additional components:

Torque ratchet
Hex key
Adapter, ISO shank, hexagon/ratchet

### Caution:

Do not confuse the prosthetic screw with the AnoTite screw. The prosthetic screw is colored silver and the AnoTite screw is anodized gold.

## 8. Technical specification

### Composition

Product	Composition
tioLogic® TWINFIT titanium abutment	3.7165
AnoTite screw	3.7165
Prosthetic screw	3.7165

Information about the composition is included in the material list for our implantology products (REF 989-801-06), or see [www.dentaureum.com/ifu](http://www.dentaureum.com/ifu). The SSCP is available at <https://ec.europa.eu/tools/eudamed> and [www.dentaureum.com](http://www.dentaureum.com).

## 9. Adverse reactions and residual risks

Possible complications and risks during prosthetic treatment (temporary adverse reactions) are:

Pain, swellings, difficulties speaking, inflammation of the gums.

Possible complications and risks during prosthetic treatment (permanent adverse reactions) are:

fractured dentures, aesthetic problems.

## 10. Safety information and precautions

- Do not use the product if there is a known allergic reaction to one or more of the material components.
- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. Therefore, secure all articles used intra-orally against swallowing and/or aspiration.
- Different types of alloy in the oral cavity can lead to galvanic reactions.
- If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.
- Only use all tioLogic® TWINFIT prosthetic components in combination with tioLogic® TWINFIT implants! Use of third-party components can negatively impact the function and safety. In such cases, Dentaureum GmbH & Co. KG accepts no liability for replacement products or guarantees.
- Observe the indications for the individual products. Consider all indications that are not listed explicitly to be contraindications.
- Disinfection, cleaning and sterilizing should not be done in the original packaging.
- Products must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.
- Technical advice (oral and written) is based on the scientific and technical knowledge available at the time the product is put on the market. It does not release the dental professional from the responsibility of personally checking the suitability of the products for the intended indications and applications. Advice is only given as a non-binding recommendation and cannot be assumed to provide any form of assurance, liability or guarantee.

### Severe incidents:

Report all severe incidents arising from the use of the product to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.

Severe incidents are identified as malfunctions and a drop in performance or deterioration of the properties of a medical device which resulted or may have resulted in a direct or indirect temporary or permanent deterioration in a patient's state of health, in the death of a patient, a dental professional or other persons or as posing a serious danger to public health.

## 11. Storage and shelf life

Store the tioLogic® TWINFIT implants, abutments and components always in a dry, dark place at room temperature in the original packaging. If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.

## 12. Disposal

Adhere strictly to the rules that apply for your region and the details outlined in the safety data sheets, if available. Products to be disposed of should be treated as surgical waste and be decontaminated in compliance with the relevant regulations for medical waste disposal.

### 13. Information for single use products

Prosthetic abutments and components are supplied non-sterile and are intended for single use on one patient only. The product must be disinfected, cleaned and sterilized before first use.

Reconditioning of abutments and components (recycling) as well as reuse on patients are not permitted as safe reconditioning and/or functional safety cannot be guaranteed and may lead to a loss in mechanical, chemical and/or biological properties. Reuse may cause infection. Do not use implants after the expiry date stated on the label.

### 14. Information for multiple use products

The torque ratchet, hex keys, and ISO shank, hexagon/ratchet adapters are supplied non-sterile and are intended for multiple use. The product must be disinfected, cleaned and sterilized before being used initially and before being used again. Products must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.

### 15. Information on disinfection / cleaning / sterilizing

The tioLogic® TWINFIT prosthetic abutments and components are supplied non-sterile. They are for single use on one patient only, unless explicitly intended for multiple use. They must be cleaned and disinfected before and after being used on a patient (e.g. transfer to dental laboratory). Dentaurem generally recommends that the product be sterilized in addition, as the use of non-sterile components may cause infections. Only sterilized prosthetic components may be used if bleeding occurs.

### 16. Product overview and scope of delivery

Please refer to the tioLogic® TWINFIT product catalog for the tioLogic® TWINFIT product range.

### 17. Quality information

Dentaurem ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaurem has no influence on the processing on site.

### 18. Explanation of symbols

Please refer to the label. See chapter "Additional information".

### 19. Information on MRT safety

The products within these Instructions for use are MRT safe.

### 20. Additional information

Additional information can be found at [www.dentaurem.com/ifu](http://www.dentaurem.com/ifu)

Instructions for use	REF
Screw	989-804-98
tioLogic® TWINFIT Additional components	989-805-88
Torque ratchet	989-719-00
Preparation instructions for the implant system components	989-801-07
Explanation of symbols used on the label	989-313-00
MRI safety of Dentaurem products	989-804-43

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