DE

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Abbreviated instructions

Dear customer

Thank you for choosing a quality product from Dentaurum Implants. It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative. As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaurum.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurum Implants GmbH | Turnstr. 31 | 75228 Ispringen | Germany

2. General description of product, intended use and characteristics, technical description

The tiologic® TWINFIT implant system contains specially coordinated instruments, abutments and accessories for the prosthetic restoration of tiologic® TWINFIT implants. Only original components of the tiologic® TWINFIT implant system should be combined in accordance with the instructions for use / user manuals.

The abutment is selected by the dentist/dental technician to suit the individual patient’s needs.

The diameter of the implant and the gingival thickness determine the choice of abutment.

3. Area of application

The precious metal abutments tiologic® TWINFIT are prosthetic components for indicated single tooth, bridge and telescopic restorations on tiologic® TWINFIT implants.

4. Medical indications/contraindications

Precious metal abutments tiologic® TWINFIT

Indications

- Single tooth and bridge restorations and primary components for double crown restorations

Contraindications

- Single restoration with cantilever unit
- Restorations with a length ratio to the length of the implant which exceeds 1:1.25
- Restorations with an angle correction of more than 20° to the implant axis
- Primary bracing

5. Information on use, processing, initial implementation and handling

Precious metal abutment tiologic® TWINFIT

All precious metal abutments tiologic® TWINFIT are available in conical and platform connector geometry for the abutment series S, M and L.

The precious metal abutment consists of a cast-on base made from a precious metal alloy, an acrylic extension and an AnoTite screw. The abutment construction makes individualization easy and is extremely precise, due to the pre-fabricated inner joint. Precious metal abutments are processed in the laboratory according to conventional dental prosthetic rules.

Processing

In order to achieve a long term, stable, single restoration, an optimal proportion between crown and inserted implant length must be ensured. If this value is exceeded, then it is preferable to fabricate a bridge restoration on two or more abutments instead.

The precious metal abutment tiologic® TWINFIT is placed on the laboratory implant (pay attention to the rotational security) and fixed with the prosthetic screw L 9.0 mm (not included in delivery). The prosthetic screw is only used for processing the abutment in the laboratory. The acrylic extension is shortened and trimmed according to the amount of space available occlusally and anatomically. The precious metal abutment is removed, then the mesio construction is waxed-up. A precious metal alloy must be used to cast-on to the precious metal base and the instructions of the material manufacturer must be followed. If the cervical cast-on segment is to be bonded with ceramic material, it is important to ensure it has a wall thickness of at least 0.3 mm in order to avoid cracking within the ceramics. It must be free of wax, modeling acrylic and grease to prevent the melt from overflowing to the base surface of the abutment.

The mesio construction is fixed to the model and the crown is waxed and cast according to conventional PFM prerequisites. The crown is veneered according to

the manufacturer’s instructions for use (e.g. ceraMotion® Me, Dentaurum). The following procedure warrants an individual marginal design. The prosthetic screw (not included in delivery) serves only to help the dental technician fix the abutment on the working cast. We recommend using a new AnoTite prosthetic screw for the final fitting in the mouth.

Tightening torque

- Prosthetic screw, on model: manually
- AnoTite screw, intra-orally: 30 Ncm

Attention: Do not confuse the prosthetic screw with the AnoTite screw. The prosthetic screw is colored silver and the AnoTite screw is anodized gold.

6. Composition

- Precious metal abutment:
 - Cast-on base: precious metal alloy

Chemical composition (% by mass)	Au	60.0 %
	Pt	19.0 %
	Pd	20.0 %
	Ir	1.0 %
		Density
Physical & mechanical properties	Melting range	1400–1490 °C/ 2552–2714 °F
	Tensile strength	> 750 MPa
	Hardness	> 215 HV5
	0.2% yield strength	> 650 MPa
	Modulus of elasticity	136 GPa
	Elongation at rupture	>2.0 %
	CTE [25–500 °C/ 77–932 °F]	11.9 x 10 ⁻⁶ K ⁻¹
	CTE [25–500 °C/ 77–932 °F]	12.2 x 10 ⁻⁶ K ⁻¹
	Condition	15–75 % KV
	Precious metal abutment S	0.307 g
Precious metal abutment M	0.359 g	
Precious metal abutment L	0.482 g	

– Burn-out acrylic: polycarbonate

- Prosthetic screws:
 - AnoTite screw: titanium alloy (Ti6Al4V)
 - Prosthetic screw: titanium alloy (Ti6Al4V)

- Hex keys:
 - Hex key– ISO shank SW 1.3 (1.4197)
 - Hex key– ratchet SW 1.3 (1.4197)

7. Safety instructions

- The product should not be used if there is a known allergic reaction to one or more of the material components.

- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. All articles used intra-orally should therefore be secured against swallowing and/or aspiration.

- All serious incidents arising from the use of the product should be reported to the manufacturer and the competent authority in the country in which the user and/ or the patient is resident.

- Different types of alloy in the mouth can lead to galvanic reactions.

- Only use tiologic® TWINFIT prosthetic components in combination with tiologic® TWINFIT implants.

8. Storage and shelf life

Store in a dry place.

9. Information for single use products

Titanium abutments are designed for single use only. Reconditioning of precious metal abutments that have been inserted previously (recycling) or reuse on patients is not permitted since it can then no longer be guaranteed that the article can be reprocessed safely or can function safely.

10. Information for multiple use products

The hex keys SW 1.3 are intended for multiple use. They must be disinfected, cleaned and sterilized before being used again. Additional information can be found at www.dentaurum.com (Processing Instructions Instruments and Accessories REF 989-801-09).

11. Information on disinfection/ cleaning/sterilizing

tiologic® TWINFIT prosthetic components are delivered non-sterile. They are for single use on one patient only. They must be cleaned and disinfected before being used on a patient. Dentaurum generally recommends that the product be sterilized in addition. Only sterilized prosthetic components may be used if bleeding occurs.

Additional information can be found at www.dentaurum.com (Processing Instructions Instruments and Accessories REF 989-801-09).

12. Product overview/Scope of delivery, incl. REF

Please refer to the tiologic® TWINFIT product catalog for the tiologic® TWINFIT product range.

13. Quality information

Dentaurum Implants ensures faultless quality of its products. These recommendations are based upon Dentaurum’s own experiences. The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as Dentaurum Implants has no influence on the processing on site.

14. Explanatoin of symbols

Refer to the label. Additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00).

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