DE

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Abbreviated instructions

Dear customer

Thank you for choosing a quality product from Dentaurum Implants. It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit.

The Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative. As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaurum.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurum Implants GmbH | Turnstr. 31 | 75228 Ispringen | Germany

2. General description of product, intended use and characteristics, technical description

The tiologic® TWINFIT implant system contains specially coordinated instruments, abutments and accessories for the fabrication of tiologic® TWINFIT implants. Only original components of the tiologic® TWINFIT implant system should be combined in accordance with the instructions for use / user manuals. The abutment is selected by the dentist/dental technician to suit the individual patient's needs.

The diameter of the implant and the gingival thickness determine the choice of abutment.

Selection aids with conical and platform connector geometry for the abutment series S, M and L are available to help choose the right abutment on the model (gingival height and angulation).

3. Area of application

Ball and tioLOC abutments tiologic® TWINFIT are prosthetic components for implant retained, gingiva-supported restorations on tiologic® TWINFIT implants.

4. Medical indications/contraindications

4.1 Ball abutments tiologic® TWINFIT

Indications

- Implant retained, gingiva-supported overdenture
- Incorporation into an existing overdenture as temporary restoration or into a cast framework
- Prosthetic restoration on at least two implants

Contraindications

- Uneven number of implants per jaw
- Unfavorable positioning of implants which prevent a tangential rotation axis
- Restorations with mixed retention (tooth/implant)
- In combination with other retention elements
- Divergence of more than 20° between the implants
- Heavy mechanical strain with ø 3.3 mm implants

4.2 tioLOC abutments tiologic® TWINFIT

Indications

- Implant-retained, gingiva-supported coverdenture with low profile
- Incorporation into an existing overdenture
- Prosthetic restoration on two or four implant posts
- Special retainers (green and red) compensate for diverging implants of up to 40° (not applicable to tiologic® TWINFIT implants ø 3.3 mm).

Contraindications

- Uneven number of implants per jaw
- Unfavorable positioning of implants which prevent a tangential rotation axis
- Restorations with mixed retention (tooth/implant)
- In combination with other retention elements
- Divergence of more than 10° between the implants when using ø 3.3 mm tiologic® TWINFIT implants is contraindicated
- Intermaxillary space should be at least 4.0 mm

5. Information on use, processing, initial implementation and handling

5.1 Ball abutment tiologic® TWINFIT

Ball abutments with platform connector geometry are available for the S, M and L series of abutments. The gingival heights are 1.5 mm, 3.0 mm and 4.5 mm; referring to the distance between the upper edge of the implant and the lower edge of the socket for the PentaGrip insertion key tiologic® TWINFIT ISO shank of the ball abutment. This lower edge should be placed approx. 1.0 mm above the gingival line. To achieve optimum retention, all ball abutments should be positioned at the same level. The ball is 2.25 mm in diameter. The components of the ball abutment may not be modified. Only the matrix withdrawal force can be adjusted by activating the inner matrix (2 N – 12 N). The ball abutment is placed with the PentaGrip insertion key tiologic® TWINFIT ISO shank.

Tightening torque

- Ball abutment on a model: manually

- Ball abutment intra-orally: 35 Ncm

5.2 tioLOC abutment tiologic® TWINFIT

tioLOC abutments with platform connector geometry are available for the S, M and L series of abutments. The gingival heights are 1.5 mm, 3.0 mm and 4.5 mm, referring to the cylindrical area of the tioLOC abutment up to the lower edge of the retention element. Its upper edge of the tioLOC abutments should be placed approx. 1.0 mm above the gingival line. The head with retainer is identical for all tioLOC abutments. tioLOC abutment components must not be modified. The LOCATOR® matrix comprises a titanium body with inner retention insert. It is for polymerizing into an existing or new overdenture.

- Divergence of more than 10° between the implants when using ø 3.3 mm tiologic® TWINFIT implants is contraindicated

- Intermaxillary space should be at least 4.0 mm

The withdrawal force of the LOCATOR® matrix is regulated with an exchangeable retention insert.

The retention insert is available in the following versions:

transparent	0° – 10°	medium	2.270 g	
pink	0° – 10°	light	1.365 g	
blue	0° – 10°	very light	680 g	
green	0° – 20°	extended angulation	1.815 g	NOT for use with tiologic® TWINFIT implants ø 3.3 mm
red	0° – 20°	extended angulation	680 g	NOT for use with tiologic® TWINFIT implants ø 3.3 mm

The LOCATOR® matrix is supplied with a black retention ring inserted (polymerizing in), a white spacer ring (blocking out undercut) and the definitive retention rings in transparent, pink and blue. The retention inserts of LOCATOR® matrices can be exchanged with the tioLOC retention key. To do so, the pointed section of the tioLOC retention key is pressed into the lower side corner of the retention insert which is then pulled downwards and out.

The tioLOC abutment is inserted with the tioLOC insertion key.

Tightening torque

- tioLOC abutment on the model: manually

- tioLOC abutment intra-orally: 30 Ncm

6. Composition

- Ball abutment**: titanium alloy (Ti6Al4V)

- Laboratory implant: aluminium

- Matrix Dalbo®-Plus**: titanium alloy

- Lamellae retention insert**: Elitor®

- tioLOC abutment**: titanium alloy (Ti6Al4V)

- LOCATOR® matrix**: titanium alloy (Ti6Al4V)

- LOCATOR® retention insert**: DuPont™ Zytel® 101L NC10 Nylon (Polyamide 66)

- LOCATOR® laboratory implant**: titanium alloy (Ti6Al4V)

- tioLOC impression cap**: titanium alloy (Ti6Al4V), insert: DuPont™ Zytel® 101L NC10 Nylon (Polyamide 66)

- 4Base tioLOC cap**: (Ti6Al4V)

- PentaGrip insertion key – ISO shaft 4Base/ball abutment**: 1.4197

- Hex keys**:
 - Hex key – ISO shank SW 1.3: 1.4197
 - Hex key – ratchet SW 1.3: 1.4197

7. Safety instructions

- The product should not be used if there is a known allergic reaction to one or more of the material components.

- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. All articles used intra-orally should therefore be secured against swallowing and/or aspiration.

- All serious incidents arising from the use of the product should be reported to the manufacturer and the competent authority in the country in which the user and/ or the patient is resident.

- Different types of alloy in the mouth can lead to galvanic reactions.

- Only use tiologic® TWINFIT prosthetic components in combination with tiologic® TWINFIT implants.

8. Storage and shelf life

Store in a dry place.

9. Information for single use products

Ball and tioLOC abutments and matrices are designed for single use only. Reconditioning ball and tioLOC abutments that have been used previously (recycling) or reuse on patients is not permitted since it can then no longer be guaranteed that the article can be reprocessed safely or can function safely.

10. Information for multiple use products

The PentaGrip insertion key and the hex keys SW 1.3 are intended for multiple use. They must be disinfected, cleaned and sterilized before being used again. Additional information can be found at www.dentaurum.com (Processing Instructions Instrumnts and Accessories REF 989-801-09).

11. Information on disinfection/cleaning/sterilizing

tiologic® TWINFIT prosthetic components are delivered non-sterile. They are for single use on one patient only. They must be cleaned and disinfected before being used on a patient. Dentaurum generally recommends that the product be sterilized in addition. Only sterilized prosthetic components may be used if bleeding occurs.

Additional information can be found at www.dentaurum.com (Processing Instructions Instruments and Accessories REF 989-801-09).

12. Product overview /Scope of delivery, incl. REF

Please refer to the tiologic® TWINFIT product catalog for the tiologic® TWINFIT product range.

13. Quality information

Dentaurum Implants ensures faultless quality of its products. These recommendations are based upon Dentaurum's own experiences. The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as Dentaurum Implants has no influence on the processing on site.

14. Explanation of symbols

 Refer to the label. Additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00).

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Dentaurum Implants GmbH | Turnstr. 31 | 75228 Ispringen | Allemagne