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Dentaurum GmbH & Co. KG · Turnstr. 31 · 75228 Ispringen · Germany

2. General device description

Prosthetic abutments and components replace missing teeth, restore chewing function and thus improve the quality of life. The implants of the tioLogic® implant family are used to fixate corresponding abutment components, depending on the indication. Prosthetic abutments and components serve to produce prosthetic meso-constructions, which are then restored with a prosthetic superstructure.

Abutment components						
Series of abutments*	S		M		L	
Implant diameters in mm	3.3	3.7	4.2	4.8	5.5	

*all of the abutment components listed in these Instructions for use are available in the S, M and L series of abutments.

The ball / tioLOC abutments are only available with a platform connector geometry.

tioLogic® TWINFIT ball and tioLOC abutments are prefabricated components for fabricating implant-retained tissue-borne prosthetic restorations. The ball and tioLOC abutments are available in different gingival heights to suit the varying gingival conditions of patients. The ball and tioLOC abutments are marked according to the S, M or L series of abutments. They are supplied exclusively with a platform connector geometry.

Ball abutment						
Series of abutments	S		M		L	
GH in mm	1.5 / 3.0 / 4.5					
Ball diameter in mm	2.25					

tioLOC abutment						
Series of abutments	S		M		L	
GH in mm	1.5 / 3.0 / 4.5					

Ball abutment laboratory implants are used to create working casts. They accurately reflect the position of the tioLogic® TWINFIT ball abutments. Processing and use take place exclusively in the (practice) laboratory or in the manufacturing center and not on the patient. This is not a medical device.

tioLogic® TWINFIT 4Base tioLOC cap, titanium, is intended for use in the fabrication or supplementation of an implant-retained, tissue-borne restoration with tioLOC abutments. It is possible to place a 20° or 30° angulated 4Base abutment on the implant if the implants diverge strongly. The 4Base tioLOC cap titanium on the angled 4Base abutment enables integration into the LOCATOR® work. Existing coverdentures can be modified or remade. The 4Base tioLOC cap, titanium is marked with a "4" for 4 Base.

The tioLogic® TWINFIT implant system contains Implant system with specially coordinated abutments and components incl. instruments.

3. Intended purpose

Abutments are intended to connect the implants with the prosthetic restoration.

Ball abutment laboratory implants are designed to represent the position of the ball abutment (axial inclination and height) in a working cast.

4. Intended users

The tioLogic® TWINFIT product range is supplied exclusively to qualified personnel such as surgeons, dentists and dental technicians. It should only be used by qualified personnel who are familiar with dental implantology, including diagnosis, preoperative planning, surgical technique and prosthetic treatment.

5. Medical indications / contraindications

5.1 Indications

5.1.1 tioLogic® TWINFIT ball abutments

- Resilient anchoring of implant-retained, tissue-borne coverdentures on two implants to ensure the tangential axis of rotation
- Anchoring of implant-retained full dentures on two implants
- Incorporation in an existing coverdenture as an interim denture or with partial denture framework possible
- Prosthetic restoration on at least two implants

5.1.2 Ball abutment laboratory implant

- Precise reproduction of positions for cast models (dental stone)

5.1.3 tioLogic® TWINFIT tioLOC abutments

- Anchoring of implant-retained, tissue-borne coverdentures with very
- Low profile height
- Incorporation into an existing coverdenture is possible
- Prosthetic restoration on two or four implant posts

5.1.4 tioLogic® TWINFIT 4Base tioLOC cap, titanium

- For use on straight or angled 4Base abutments
- Integration into implant-retained, tissue-borne dental restorations and combination with tioLOC abutments and matrices

5.2 Contraindications

- In patients with a known intolerance to individual components.
- The contact surface between the implant and abutment must not be blasted or otherwise mechanically processed
- Insufficient oral hygiene
- Bruxism

5.2.1 tioLogic® TWINFIT ball abutments

- Uneven number of implants per jaw
- Unfavorable positioning of implants which prevent a tangential rotation axis
- Restorations with mixed retention (tooth/implant)
- Combination with other retention elements
- Divergence of more than 20° between the implants
- Heavy mechanical strain with \varnothing 3.3 mm implants

5.2.2 Ball abutment laboratory implant

- Use of laboratory implants with defective or damaged connector geometry
- Use on patients

5.2.3 tioLogic® TWINFIT tioLOC abutments

- Uneven number of implants per jaw
- Unfavorable positioning of implants which prevent a tangential rotation axis
- Restorations with mixed retention (tooth/implant)
- Combination with other retention elements
- Implant divergences of more than 10° using tioLogic® TWINFIT implants \varnothing 3.3 mm are contraindicated
- Intermaxillary space should be at least 4.0 mm

5.2.4 tioLogic® TWINFIT 4Base tioLOC cap, titanium

- Use as single tooth restoration
- Restorations with a length ratio to the length of the implant which exceeds 1:1.25
- Use of other abutments than on 4Base abutments
- Combination of other restoration variants than with tioLOC

6. Patient target group

All skeletally mature patients undergoing implant treatment. Please consider the indications / contraindications.

7. Information on use, processing, fitting and insertion

Prior to use, the dental professional ensures that they have carefully read and understood the tioLogic® TWINFIT Instructions for use.

We strongly recommend that, prior to using Dentaurum products, dental professionals attend a Dentaurum training course to familiarize themselves with the Dentaurum products.

Carry out thorough radiological, psychological and physical examinations before implantation and inform patients of the process in detail. Dentaurum recommends full clinical, radiological, photographic and statistical documentation.

The ball and tioLOC abutments and components of the tioLogic® TWINFIT implant system can be documented, e.g. in the patient file or PatientPass, using the additional labels.

Precautions:


Before surgical implantation, it is important to carry out a careful patient anamnesis and, if required, a third-party anamnesis by the patient's doctor (GP/PCP) in order to clarify if 1) implantation may be difficult due to an anatomical situation, 2) there is the risk of a serious surgical problem or a general risk, 3) wound healing/osseointegration may be difficult, or 4) good hygiene and care of the implant, abutment and restoration may worsen. Do not fit dental implants if any of the above cases or a combination hereof is severe or uncontrollable.

Prosthetic restoration of the implant must not be performed until the soft tissue has healed without irritation. Furthermore, a radiographic check-up should be performed before the prosthetic restoration.

Important:

The implantologist should ensure the products cannot be aspirated during intra-oral use.

7.1 tioLogic® TWINFIT ball abutments


Ball anchors function in such a manner that it is advisable to support the prosthetic restoration on at least two implants. The implants must not diverge by more than 20°. Care should be taken to avoid an excessive mechanical loading when using ball head abutments together with tioLogic® TWINFIT \varnothing 3.3 mm implants. The tioLogic® TWINFIT ball abutments are available with platform connector geometry in gingival heights of 1.5 mm, 3.0 mm and 4.5 mm for S, M and L abutments. The gingival height is the distance

between the uppermost edge of the implant and the lowermost edge of the cylindrical area of the ball abutments. This upper edge should be approx. 1.0 mm above the gingival line. To achieve optimum retention, the ball abutments should each be positioned at the same level. The ball head is 2.25 mm in diameter. No components of the ball abutment may be modified. Only the matrix withdrawal force can be adjusted by activating the inner matrix. The ball abutment is placed with the tioLogic® TWINFIT PentaGrip insertion key ISO shank.

Working procedures

- In the direct technique, process the matrix into an existing denture directly in the patient's mouth without fabricating a model. Neither a tioLogic® TWINFIT impression post nor laboratory implant are required. The tioLogic® TWINFIT ball abutments are attached to the tioLogic® TWINFIT implants intra-orally. During the clinical procedure, protect the undercuts with the help of a rubber dam. It is put over the matrix to prevent cold-curing resin filling the undercuts in the mouth which would prevent the denture being removed after the resin has cured. Observe the instructions for use of the resin.
- The indirect technique requires an impression and a model. The matrix is processed into the denture in the laboratory.

With the indirect version, the impression can also be taken using the closed technique directly over the ball abutments fixed in the mouth. To do this, the correct tioLogic® TWINFIT ball abutments for the gingival height and S, M or L abutment should be selected, the tioLogic® TWINFIT closure screws or tioLogic® TWINFIT gingiva formers removed and the tioLogic® TWINFIT ball abutment fixated in the implant. Apply impression material around all ball head abutments, fill the closed impression tray with

impression material and take the impression. Remove the tray after the impression material is cured. Ensure that the impression has captured the ball abutments exactly. The ball abutment laboratory implant is then repositioned in the impression in the laboratory. The flat surfaces beneath the ball ensure that the implant axes are transferred precisely. Then pour dental stone into the impression. Cast the dental arch as usual and create the base. The laboratory implant must fit precisely on the model. A ball abutment laboratory implant is available for all series of abutments as all balls have a diameter of 2.25 mm.

Processing and integration of the following matrix/inner matrix into the denture should be carried out using the corresponding components and according to the manufacturer's instructions.

REF 386-630-00  Dalbo®-PLUS matrix TE basic complete

REF 386-631-00  Dalbo®-PLUS Lamella retention insert E

REF 387-663-00  Dalbo®-PLUS screwdriver / activator

Recall:

Dentures and their retention units must be monitored at six-monthly intervals. Take the following points into account:

- Eliminate unfavorable movements of the denture (optimize relining of the denture, activate or replace the matrices)
- Check the fit of the ball abutments on the implant (tighten if necessary)
- Oral hygiene (remove plaque and dental calculus and, if necessary, re-instruct the patient on cleaning implants)

Tightening torque

Ball abutment on the model: / in the mouth	35 Ncm
--	--------

7.2 tioLogic® TWINFIT tioLOC abutments

The tioLOC-technique involves an implant-retained, tissue-borne restoration with a very low profile height. Existing coverdentures can be modified or remade. The tioLOC abutment can be used for tissue-borne coverdentures with two to four implant abutments. The intermaxillary space should be at least 4.0 mm. The tioLOC components must not be modified. Only the tightening forces of the matrices can be adjusted as required using the exchangeable retention inserts.

The tioLogic® TWINFIT tioLOC abutments are available with a platform connector geometry in gingival heights of 1.5 mm, 3.0 mm and 4.5 mm for S, M and L series of abutments. The gingival height refers to the cylindrical section of the tioLOC abutment. Its upper edge should be placed approx. 1.0 mm above the gingival line. To achieve optimum retention, the tioLOC abutments should be positioned at the same level. The head with retainer element is identical for all tioLOC abutments. The tioLOC abutment is inserted with the insertion key tioLOC abutment. All tioLOC abutment components may not be modified.

Working procedures

The direct technique involves processing the matrix into an existing denture directly in the patient's mouth

without creating a model. Neither an impression post nor a laboratory implant are required. The tioLogic® TWINFIT tioLOC abutments are attached intra-orally to the corresponding tioLogic® TWINFIT implants. During the clinical procedure, protect the undercuts with the help of a rubber dam. Place this over the matrix. This prevents the undercuts in the mouth from filling with cold-curing polymer and the denture from becoming impossible to remove after curing. Instructions for the resin should be adhered to strictly.

The indirect technique requires an impression and a model. The matrix is processed into a denture in the laboratory.

Processing and integration of the following matrices/caps into the denture should be carried out using the corresponding components and according to the manufacturer's instructions.

REF 346-512-00		Impression caps, 2 pieces	∅ 4.0 mm	L 10.0 mm
REF 346-514-00		Laboratory analogs 2 pieces	∅ 4.0 mm	L 10.0 mm
REF 346-516-00		Set of retention caps: 1 stainless steel case for caps, 1 protective plate, 1 laboratory cap, 3 retention caps (with internal pin)		
REF 346-518-00		Black laboratory caps (with pin), 4 pieces	GH 1.0 mm	
REF 346-520-00		Yellow laboratory caps (with pin), 4 pieces		700 g 0° - 20°
REF 346-522-00		Pink laboratory caps (with pin), 4 pieces		900 g 0° - 20°
REF 346-524-00		Transparent laboratory caps (with pin), 4 pieces		1.500 g 0° - 20°
REF 346-526-00		Violet laboratory caps (with pin), 4 pieces		2.700 g 0° - 20°
REF 346-530-00		Protective plate 10 pieces		
REF 347-745-00		Insertion/removal tool for elastic caps		

Recall:

Dentures and their retention units must be monitored at six-monthly intervals. Take the following points into account:

- Eliminate unfavorable movements of the denture (optimize relining of the denture, activate or replace the matrices)
- Check the fit of the ball abutments on the implant (tighten if necessary)
- Oral hygiene (remove plaque and dental calculus and, if necessary, re-instruct the patient on cleaning implants)

Tightening torque	
tioloc Abutment on the model / in the mouth	30 Ncm

7.3 tioLogic® TWINFIT 4Base tioloc cap, titanium

 4Base tioloc cap, titanium, is intended for use in the fabrication or supplementation of an implant-retained, tissue-borne restoration with tioloc abutments. They fit all three 4Base abutments (S, M and L). The seating surface of all 4Base abutments is identical (∅ 4.6 mm). If the implants diverge considerably, it is possible to place a 20° or 30° angled 4Base abutment on the implant. Then place the 4Base tioloc cap titanium on the 4Base abutment and screw it tight using the tioloc insertion key. The 4Base tioloc cap titanium on the angled 4Base abutment enables integration into the LOCATOR® work. Existing coverdentures can be modified or remade. The 4Base tioloc cap, titanium is marked with a 4 for 4Base. The tioloc abutment is inserted with the insertion key tioloc abutment.

Tightening torque	
4Base tioloc cap, titanium on the model / intra-orally	25 Ncm

7.4 Additional components

The prosthetic screw (not included in scope of delivery) serves only to help the dental technician fix the abutment on the working cast temporarily.

A new AnoTite screw should always be used for the final fitting and fixation of the abutment in the mouth at the prescribed torque.

Overview of additional components:

Torque ratchet
tiologic® TWINFIT PentaGrip insertion key, ISO shank
tioloc Insertion key, ratchet
Hex key
Adapter, ISO shank, hexagon/ratchet
Prosthetic screw on the model

Caution:

Do not confuse the prosthetic screw with the AnoTite screw. The prosthetic screw is colored silver and the AnoTite screw is anodized gold.

8. Technical specification**Composition**

Product	Composition
tiologic® TWINFIT ball abutment	3.7165
Ball abutment laboratory implant	Al alloy
tiologic® TWINFIT tioloc abutment	3.7165
tiologic® TWINFIT 4Base tioloc cap, titanium	3.7165
AnoTite screw	3.7165

Information about the composition is included in the material list for our implantology products (REF 989-801-05), or see www.dentaurum.com/ifu. The SSCP is available at <https://ec.europa.eu/tools/eudamed> and www.dentaurum.com.

9. Adverse reactions and residual risks

Possible complications and risks during prosthetic treatment (temporary adverse reactions) are:

Pain, swellings, difficulties speaking, inflammation of the gums.

Possible complications and risks during prosthetic treatment (permanent adverse reactions) are:

fractured dentures, aesthetic problems.

10. Safety information and precautions

- Do not use the product if there is a known allergic reaction to one or more of the material components.
- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. Therefore, secure all articles used intra-orally against swallowing and/or aspiration.
- Different types of alloy in the oral cavity can lead to galvanic reactions.
- If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.
- Only use all tioLogic® TWINFIT prosthetic components in combination with tioLogic® TWINFIT implants! Use of third-party components can negatively impact the function and safety. In such cases, Dentaurum GmbH & Co. KG accepts no liability for replacement products or guarantees.
- Observe the indications for the individual products. Consider all indications that are not listed explicitly to be contraindications.
- Disinfection, cleaning and sterilizing should not be done in the original packaging.
- Products must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.
- Technical advice (oral and written) is based on the scientific and technical knowledge available at the time the product is put on the market. It does not release the dental professional from the responsibility of personally checking the suitability of the products for the intended indications and applications. Advice is only given as a non-binding recommendation and cannot be assumed to provide any form of assurance, liability or guarantee.

Severe incidents:

Report all severe incidents arising from the use of the product to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.

Severe incidents are identified as malfunctions and a drop in performance or deterioration of the properties of a medical device which resulted or may have resulted in a direct or indirect temporary or permanent deterioration in a patient's state of health, in the death of a patient, a dental professional or other persons or as posing a serious danger to public health.

11. Storage and shelf life

Store the tioLogic® TWINFIT implants, abutments and components always in a dry, dark place at room temperature in the original packaging. If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.

12. Disposal

Adhere strictly to the rules that apply for your region and the details outlined in the safety data sheets, if available. Products to be disposed of should be treated as surgical waste and be decontaminated in compliance with the relevant regulations for medical waste disposal.

13. Information for single use products

Prosthetic abutments and components are supplied non-sterile and are intended for single use on a single patient only. The product must be disinfected, cleaned and sterilized before first use.

Reconditioning of abutments and components (recycling) as well as reuse on patients are not permitted as safe reconditioning and/or functional safety cannot be guaranteed and may lead to a loss in mechanical, chemical and/or biological properties. Reuse may cause infection. Do not use implants after the expiry date stated on the label.

14. Information for multiple use products

The torque ratchet, tioLogic® TWINFIT PentaGrip insertion key, ISO shank, tioLOC insertion key, ratchet, hex key and the ISO shank-hex/ratchet adapter are supplied non-sterile and are intended for multiple use. The product must be disinfected, cleaned and sterilized before being used initially and before being used again. Titanium scan abutments must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.

15. Information on disinfection / cleaning / sterilizing

The tioLogic® TWINFIT prosthetic abutments and components are supplied non-sterile. They are for single use on one patient only, unless explicitly intended for multiple use. They must be cleaned and disinfected before and after being used on a patient (e.g. transfer to dental laboratory). Dentaurem generally recommends that the product be sterilized in addition, as non-sterile components may cause infections. Only sterilized prosthetic components may be used if bleeding occurs

16. Product overview and scope of delivery

Please refer to the tioLogic® TWINFIT product catalog for the tioLogic® TWINFIT product range.

17. Quality information

Dentaurem ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaurem has no influence on the processing on site.

18. Explanation of symbols

Please refer to the label. See chapter "Additional information".

19. Information on MRT safety

The products within these Instructions for use are MRT safe.

20. Additional information

Additional information can be found at www.dentaurem.com/ifu

Instructions for use	REF
tioLogic® TWINFIT Additional components	989-805-88
Torque ratchet	989-719-00
4Base abutments and components	989-365-93
Preparation instructions for the implant system components	989-801-07
Explanation of symbols used on the label	989-313-00
MRI safety of Dentaurem products	989-804-43

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