DE Kurzanweisung

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Abbreviated Instructions for use

Dear customer

Thank you for choosing a quality product from Dentaurum Implants.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit.

Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaurum.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurum Implants GmbH | Turnstr. 31 | 75228 Ispringen | Germany

2. General description of product, intended use and characteristics, technical description

tiologic® TWINFIT implants are designed for insertion in the endosteal region of the maxilla or mandible. Depending on the indication, appropriate transgingival abutments are secured on the implants and fitted with a prosthetic superstructure.

The tiologic® TWINFIT implant system has specially coordinated instruments, abutments, and accessories available for implant insertion and prosthetic restoration. Only original components of the tiologic® TWINFIT implant system should be used in combination in accordance with the Instructions for use / user manuals.

3. Area of application

tiologic® TWINFIT implants can be used in both the mandible and maxilla for surgical immediate implantation, delayed immediate implantation and delayed implantation using either the one-stage or two-stage technique.

4. Medical indications/contraindications

Indications

Indications for implant placement are small and large bounded saddles (one-tooth restorations, increasing the number of abutments), a shortened dentition and an edentulous jaw in the maxilla or mandible. Before deciding to use an implant, the risks and advantages / disadvantages of dental implants and alternative methods of treatment should be considered.

In all cases, the diameter and length of the tioLogic® TWINFIT implant should be in proportion to the prosthetic restoration. If practical with the individual oral situation, implants with a minimum diameter of 4.2 mm should generally be used for restorations that subject the implant and superstructure to high mechanical loading. The tiologic® TWINFIT implant type S ø 3.3 mm is available for patients with narrow alveolar ridges. Due to the smaller diameter and lower load capacity (compared to the tiologic® TWINFIT implant type M ø 4.2 mm), these implants have a limited range of indications. In fully edentulous cases, four or more tiologic® TWINFIT implants type S ø 3.3 mm must be inserted with a splinted bar restoration. They must be combined with tiologic® TWINFIT implants ø 4.2 mm, ø 4.8 mm or ø 5.5 mm and splinted in partially edentulous cases with implant supported restorations. For restorations with single crowns, tiologic® TWINFIT implant type S ø 3.3 mm should only be used for the lower incisors or the upper lateral incisors and only with a minimum length of 11.0 mm. Restorations with single crowns on tiologic® TWINFIT implants ø 3.7 mm, ø 4.2 mm, ø 4.8 mm or ø 5.5 mm require a minimum implant length of 9.0 mm.

Contraindications

General contraindications:

General contraindications for dental surgery procedures apply. These include: reduced immunodeficiency / steroid treatment / blood coagulation disorders / uncontrolled endocrine diseases / rheumatic disorders / bone system diseases / cirrhosis of the liver / drug, alcohol or tobacco abuse / depression, psychopathic disorders / poor patient compliance / chronic underlying inflammatory diseases / incomplete physical growth of patient.

Local / personal contraindications:

Osteomyelitis / radiotherapy in the head region / recurring mucosal diseases / temporomandibular joint dysfunctions / parafunctions / lack of vertical or horizontal bone availability, jaw defects, inadequate bone quality / poor oral hygiene.

5. Information on use, processing, initial implementation and handling

The tiologic® TWINFIT product range is supplied exclusively to doctors, dentists and dental technicians. It should only be used by doctors, dentists or dental technicians who are familiar with dental implantological procedures, including diagnosis, preoperative planning, surgical technique and prosthetic treatment.

Before use, operators should ensure that they have carefully read and understood the full tiologic® TWINFIT Instructions for use / manuals. These abbreviated Instructions for use do not replace the Instructions for use / manuals for the tiologic® implant system. As the instructions and manuals cannot provide all information for immediate use, we strongly recommend that, before using the system, operators attend a tiologic® TWINFIT system training course offered by Dentaurum Implants to learn the correct techniques.

■ Please refer to the **Product Catalog** and the **Surgery Manuals** for information on precautions and the selection of components for the clinical procedure.

■ Refer to the **Product Catalog** and the **Prosthetic Manuals** for information on precautions and the selection of components for the prosthetic procedure.

Before using this product, the operator must examine the patient thoroughly and explain the procedure in detail. Dentaurum Implants recommends full clinical, radiological, photographic and statistical documentation.

The tiologic® TWINFIT implant system components can be documented, for example in the patient file or PatientPass, using the additional labels.

The operator should ensure the products cannot be aspirated during intra-oral use.

6. Composition

■ Titanium Grade 4 (pure titanium) DIN EN ISO 5832-2

Chemical composition (% by mass)	O	0.4 % max.
	Fe	0.5 % max.
	C	0.1 % max.
	N	0.05 % max.
	H	0.012 % max.
	Ti	Residue
Physical and mechanical properties	0.2% yield strength	520 MPa min.
	Tensile strength	680 MPa min.
	Elongation at rupture	10 % min.

■ Titanium Grade 5 (titanium alloy) DIN EN ISO 5832-3

Chemical composition (% by mass)	Al	5.5 % – 6.75 %
	V	3.5 % – 4.5 %
	Fe	0.3 % max.
	C	0.08 % max.
	N	0.05 % max.
	H	0.015 % max.
	O	0.2 % max.
	Ti	Residue
Physical and mechanical properties	0.2% yield strength	780 MPa
	Tensile strength	860 MPa
	Elongation at rupture	10 % min.

7. Side effects

tiologic® TWINFIT implant types are not known to have any side effects or produce any interaction.

It cannot, however, be ruled out that in rare cases allergies to components used in the materials of the tiologic® TWINFIT implant system may occur or that there may be electrochemically-induced discomfort.

8. Safety instructions

■ The product should not be used if there is a known allergic reaction to one or more of the material components.

■ Different types of alloy in the mouth can lead to galvanic reactions.

■ Only use the tiologic® TWINFIT prosthetic components in combination with tiologic® TWINFIT implants.

9. Warnings

Thorough placement of dental implants has a high rate of success and implants have a long durability, successful treatment cannot be guaranteed. The operator should note and document any problematic cases and inform the manufacturer.

An inadequate number of implants, implants with insufficient length or diameter, unfavorable positioning of the implants or a statically-poor prosthetic restoration can cause fatigue fracture in implants, abutments, prosthetic restorations and screws under bio-mechanical loading. Placement of the implants and fabrication of the prosthetic restoration should take into account the individual oral situation of a patient to avoid overloading the components.

Using tiologic® implant system components in combinations contrary to the Instructions for use / manuals can also cause mechanical failure, damage to the tissue or unsatisfactory aesthetic results.

10. Storage and shelf life

tiologic® TWINFIT implants should always be stored in a dry, dark place at room temperature in the original packaging. The implants should not be used after the expiry date (see label).

11. Notes on sterile packaging

All tiologic® TWINFIT implants are supplied sterile with the appropriate closure screws and with drill stops of appropriate lengths and widths in gamma-sterilized double packaging (foil and blister packaging) that protects the content against contamination. The packaging should not be removed until shortly before insertion.

Dentaurum Implants guarantees the sterility of the tiologic® TWINFIT implants until the expiration date on the packaging, provided that the original packaging is undamaged. No guarantee can be given for the sterility after the expiration date and the tiologic® TWINFIT implant may no longer be used on patients after this date. If the sterile packaging is damaged, the product may not be used or re-sterilized. **Packaging that has been opened may not be re-sterilized.** A tiologic® TWINFIT implant that has been removed from the packaging but has not been inserted, may also no longer be used or re-sterilized.

All other components of the tiologic® TWINFIT implant system are supplied non-sterile.

12. Information for single use products

tiologic® TWINFIT implants are designed for single use only. Reconditioning of a tiologic® TWINFIT implant that has already been used (recycling) and its reuse on a patient are not permitted.

13. Information for multiple-use products

Rotating instruments and many accessories are intended for multiple use. The product must be disinfected, cleaned and sterilized before being used again. This applies in particular to components that can be disassembled, e. g. the torque wrench. Rotary instruments – used with proper care and provided that they are not damaged or contaminated – can be reused in dense bone 30 to 40 times; any further reuse or the use of damaged and / or contaminated instruments should be avoided. The operator is responsible for ensuring the instruments are in good condition. No liability is accepted if these instructions are disregarded.

14. Information on disinfection/cleaning/sterilizing

Please refer to the Surgery Manual tiologic® TWINFIT at www.dentaurum.com for information on disinfection.

15. Notes on compliance / Conformity declaration

Development, clinical testing, production and quality control of the tiologic® TWINFIT product range are executed in accordance with the Medical Device Directive 93/42/EEC.

Sections 9 and 10 of our General Terms of Delivery and Payment apply with regard to warranty or liability – unless otherwise provided in these abbreviated instructions.

Warranty and liability are rendered void in particular if the products are not used by the operator or a third party in accordance with the Instructions for use, this also applies if the tiologic® TWINFIT product range is used in combination with products from other manufacturers which have not been specifically recommended for use by Dentaurum Implants. Dentaurum Implants has no control over processing and application of the product. These are the sole responsibility of the user.

Technical advice (oral and written) is based on the scientific and technical knowledge available when the product is put on the market. It does not release the user from the responsibility of personally checking the suitability of the products for the intended indication and application. Advice is only given for a non-binding recommendation, which cannot be assumed to provide any form of assurance or guarantee.

All products are subject to continuous development based on current scientific knowledge and we reserve the right to make changes in the construction, design or material of the products.

16. Product overview / scope of delivery, incl. REF

Please refer to the tiologic® TWINFIT product catalog for the product range for tiologic® TWINFIT. **Note: Not all components are available in every country.**

17. Quality information

Dentaurum Implants ensures faultless quality of its products. The recommendations in these instructions are based upon Dentaurum Implant’s own experiences. The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as Dentaurum Implants has no influence on the processing on site.

18. Explanation of symbols

Refer to the label. Additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00).

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