DE **Gebrauchsanweisung** · tiologic[®] TWINFIT Implantate

EN **Instructions for use** · tiologic[®] TWINFIT implants

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Dentaurum GmbH & Co. KG · Turnstr. 31 · 75228 Ispringen · Germany

2. General device description

tioLogic® TWINFIT dental implants are designed for insertion in the endosteal region of the maxilla or mandible. They are used for anchoring functional and aesthetic oral rehabilitations in partially edentulous and fully edentulous patients.

tioLogic® TWINFIT implant data					
Implant diameters in mm	3.3	3.7	4.2	4.8	5.5
Color coding	Brown	Yellow	Red	Blue	Green
Implant lengths in mm	9/11/13/15	7/9/11/13/15			9/11/13/15

Depending on the indication, appropriate transgingival abutments are secured on the tioLogic® TWINFIT implants and form the connection between the implant and a prosthetic superstructure. Additional information on the abutments available at www.dentaurum.com/ifu

The implants can be used for surgical immediate implantation, delayed immediate implantation, and later implantation, either one- or two-phase. tioLogic® TWINFIT implants have a cylindrical / conical external geometry with a blasted / etched surface. The surface of the implant shoulder is machined. The internal geometry of the tioLogic® TWINFIT implant allows both a platform and a conical connection and has 5 grooves to facilitate the positioning and rotational security of prosthetic abutment components.

The tioLogic® TWINFIT implant system contains specially coordinated abutments and accessories, including instruments, for the insertion and prosthetic restoration of the implants.

3. Intended purpose

Implants are designed for insertion in the endosteal region of the maxilla or mandible. They are used for anchoring functional and aesthetic oral rehabilitations.

4. Intended users

The tioLogic® TWINFIT product range is supplied exclusively to qualified personnel such as surgeons, dentists and dental technicians. It should only be used by qualified personnel who are familiar with dental implantology, including diagnosis, preoperative planning, surgical technique and prosthetic treatment.

5. Medical indications / contraindications

5.1 Indications

Indications for implant placement are small and large-bounded saddles (one-tooth restorations, increasing the number of abutments), a shortened dental arch and an edentulous jaw in the maxilla or mandible. Consider the risks and advantages / disadvantages of dental implants and alternative methods of treatment before deciding on the indication for an implantation.

Whenever implants are used, the implant diameter and length of the tioLogic® TWINFIT implant types should be in proportion to the prosthetic restoration. If practical with the individual oral situation, generally use implants with a minimum diameter of 4.2 mm for restorations with a high mechanical loading on implants and superstructures.

The tioLogic® TWINFIT implant types S \varnothing 3.3 mm are available for patients with indications with a narrow alveolar ridge. Due to the smaller diameter and lower mechanical load capacity (compared to the tioLogic® TWINFIT implant type M \varnothing 4.2 mm), these implants have a limited range of indications.

In fully edentulous cases, four or more tioLogic® TWINFIT implants with \varnothing 3.3 mm must be inserted with a splinted bar construction.

In partially edentulous cases, combine tioLogic® TWINFIT implants \varnothing 3.3 mm for implant-supported restorations with tioLogic® TWINFIT implants \varnothing 4.2 mm or \varnothing 4.8 mm or \varnothing 5.5 mm, and the prosthetic finishing should be fixed and splinted.

In single crown restorations, use tioLogic® TWINFIT implants \varnothing 3.3 mm only for the lower incisors or the upper lateral incisors and only with a minimum 11.0 mm implant length. Single restorations on tioLogic® TWINFIT \varnothing 3.7 mm, \varnothing 4.2 mm, \varnothing 4.8 mm or \varnothing 5.5 mm implant types require a minimum 9.0 mm implant length.

The healing period for tioLogic® TWINFIT implants \varnothing 3.3 mm is a minimum of 12 weeks.

Only use tioLogic® TWINFIT implants with a length of 7.0 mm if there is no space available for longer implants.

5.2 Contraindications

General contraindications:

Insufficient bone or tissue available and/or inadequate bone quality, local root remnants, bone disease and wound healing disorders, local infection of implant site, serious therapy-resistant functional disorders, uncontrolled diabetes mellitus, long-term immunosuppressive therapy, connective tissue disease, blood disorders (e.g. leukemia, hemophilia), intra-oral infection or malignancies, uncontrolled parafunctional habits, non-treatable occlusal or articulation conditions, serious psychological illnesses, xerostomia and titanium allergies.

Local / personal contraindications:

Incomplete growth of jaw, unfavorable anatomical bone situations, previously irradiated bone, temporomandibular joint disorders and treatable pathologic diseases to the jaw, severe systemic diseases, reduced immune defense and leukocyte disfunctions leading to increased risk of infection, endocrine disorders, medical anticoagulation / hemorrhagic diathesis, arteriosclerosis and stroke, hypertension, heart attack, illnesses with periodic use of steroids, hepatitis, diabetes mellitus and pregnancy, metabolic bone disorders, all illnesses affecting bone regeneration or blood microcirculation, rheumatic diseases and drug/alcohol/tobacco addiction, lack of patient compliance, inadequate oral hygiene, periodontitis, bruxism, parafunctional habits and oral mucosal abnormalities.

It is prohibited to rework the contact surface of implant abutments and the screw aperture. It is otherwise not possible to guarantee a safe connection.

In single crown restorations, use tioLogic® TWINFIT implants \varnothing 3.3 mm only for the lower incisors or the upper lateral incisors. A single restoration of the upper central incisors, canines, premolars and molars with tioLogic® TWINFIT implants \varnothing 3.3 mm is contraindicated. Double crown constructions and implant divergences of more than 10° are also prohibited for this implant diameter with tioLOC abutments. Immediate loading of single crown restorations is not recommended for implants with a length of 7.0 mm.

6. Patient target group

All skeletally mature patients undergoing implant treatment. Please consider the indications / contraindications.

7. Information on use, processing, fitting and insertion

Prior to use, the dental professional ensures that they have carefully read and understood the tioLogic® TWINFIT Instructions for use.

We strongly recommend that, prior to using Dentaurum products, dental professionals attend a Dentaurum training course to familiarize themselves with the Dentaurum products. Carry out thorough radiological, psychological and physical examinations before implantation and inform patients of the process in detail. Dentaurum recommends full clinical, radiological, photographic and statistical documentation.

There are self-adhesive labels in the outer packaging with REF and LOT numbers for documentation in the PatientPass and the surgical protocol.

The implantologist should ensure the products cannot be aspirated during intra-oral use.

Precautions:

Before surgical implantation, it is important to carry out a careful patient anamnesis and, if required, a third-party anamnesis by the patient's doctor (GP/PCP) in order to clarify if 1) implantation may be difficult due to an anatomical situation, 2) there is the risk of a serious surgical problem or a general risk, 3) wound healing/osseointegration may be difficult, or 4) good hygiene and care of the implant, abutment and restoration may worsen. No dental implants should be fitted if any of the above cases or a combination hereof is severe or uncontrollable.

7.1 Preparation for surgery

Ensure careful and validated hygiene measures throughout the entire operation and document these. Prepare the treatment room, instrumentarium and patient under these aspects. Check all surgical instruments required for the operation to ensure they are complete, functional and sterile. We recommend having several implants and surgical instruments available as a precaution.

All ADVANCED drills and thread taps are supplied non-sterile. Please sterilize them before use. They should be thoroughly cleaned, disinfected, conditioned and sterilized before being used for the first time (factory new) and immediately after each use. Rotary instruments should always be checked to ensure that they are sharp, in good condition and the markings are legible, as they have a limited service life. Instruments can become blunt as a result of use and cleaning. Only instruments that are sharp and free from defects should be used. Rotary instruments can be reused in hard bone 15 to 20 times if used with proper care and provided that they are not damaged or contaminated. Avoid any further reuse or the use of damaged and / or contaminated instruments. The dental professional bears responsibility. No liability is accepted if these instructions are disregarded.

7.2 Notes on operation planning

Implants need a distance of at least 1.5 mm to the neighboring teeth. The distance between implants must be at least 3.0 mm. Select the implant diameter and length so that adequate bone around the implant is available (min. 1.0 mm).

For technical reasons, the depth drill and the stepped countersink are 1.0 mm longer (apically) than the given preparation length. This should be taken into account during diagnostics, planning and preparation. The distance to the mandibular nerve and inferior alveolar nerve should be at least 1.5 mm!

7.3 General information on implant site preparation

Careful preparation of the hard and soft tissue in the bone is important for optimum healing of the implant site. Avoid thermal or mechanical trauma at all costs. To this end, keep the temperature as low as possible when preparing the site and do not exceed the maximum implant insertion torque (max. 40 Ncm).

For this reason, a low-speed preparation is recommended for certain drills. Generally speaking, decrease the torque as the drill diameter increases.

Article		Low-speed preparation	Maximum torque
Marking drill		max. 800 rpm	max. 800 rpm
Depth drill		max. 800 rpm	max. 800 rpm
Surface cutter		50 rpm– 80 rpm	max. 500 rpm
Stepped countersink		50 rpm– 80 rpm	max. 500 rpm
Expanders		50 rpm– 80 rpm	max. 500 rpm
Thread tap		max. 10 rpm or manually	max. 10 rpm or manually


Drilling protocol - tiologic® TWINFIT

- Adapt preparation protocol according to the indication and situation of each patient if required.
- The initial situation is the preparation protocol
- Adapt if required

- = max. Implanttiefe
- = min. depth marking 7
- = optional use
- = refer to Instructions for use

View surgical tray

Step	Article	Implant diameters																																																	
		Ø 3.3	Ø 3.7	Ø 4.2	Ø 4.8	Ø 5.5																																													
1	Marking drill																																																		
2	Depth drill																																																		
3	Surface cutter																																																		
4 a	Preceding stepped countersink																																																		
4	Stepped countersink																																																		
Preparation according to bone quality		<table border="1"> <thead> <tr> <th colspan="3">Preparation according to bone quality</th> <th colspan="3">Preparation according to bone quality</th> <th colspan="3">Preparation according to bone quality</th> <th colspan="3">Preparation according to bone quality</th> <th colspan="3">Preparation according to bone quality</th> </tr> <tr> <th>D4-D3</th> <th>D3-D2</th> <th>D2-D1</th> <th>D4-D3</th> <th>D3-D2</th> <th>D2-D1</th> <th>D4-D3</th> <th>D3-D2</th> <th>D2-D1</th> <th>D4-D3</th> <th>D3-D2</th> <th>D2-D1</th> <th>D4-D3</th> <th>D3-D2</th> <th>D2-D1</th> </tr> <tr> <td>soft</td> <td>medium</td> <td>hard</td> <td>soft</td> <td>medium</td> <td>hard</td> <td>soft</td> <td>medium</td> <td>hard</td> <td>soft</td> <td>medium</td> <td>hard</td> <td>soft</td> <td>medium</td> <td>hard</td> </tr> </thead> </table>					Preparation according to bone quality			Preparation according to bone quality			Preparation according to bone quality			Preparation according to bone quality			Preparation according to bone quality			D4-D3	D3-D2	D2-D1	D4-D3	D3-D2	D2-D1	D4-D3	D3-D2	D2-D1	D4-D3	D3-D2	D2-D1	D4-D3	D3-D2	D2-D1	soft	medium	hard	soft	medium	hard	soft	medium	hard	soft	medium	hard	soft	medium	hard
Preparation according to bone quality			Preparation according to bone quality			Preparation according to bone quality			Preparation according to bone quality			Preparation according to bone quality																																							
D4-D3	D3-D2	D2-D1	D4-D3	D3-D2	D2-D1	D4-D3	D3-D2	D2-D1	D4-D3	D3-D2	D2-D1	D4-D3	D3-D2	D2-D1																																					
soft	medium	hard	soft	medium	hard	soft	medium	hard	soft	medium	hard	soft	medium	hard																																					
5	Expanders																																																		
6	Thread tap																																																		
		Ø 3.3	Ø 3.7	Ø 4.2	Ø 4.8	Ø 5.5																																													

7.4 Implant site preparation

Select the rotary instruments and their drilling sequence and depth for the preparation of the implant site depending on the bone quality. A drilling protocol for different bone qualities (soft, medium, hard) is available for the implantologist. Determination of the bone quality rests on the implantologist. Use the drills with intermittent drilling and ensure sufficient cooling using a sterile, cooled physiological saline solution (5 °C / 41 °F) as the diameter increases.

Carry out implant site preparation using the ADVANCED for tioLogic® drilling system, bearing in mind the drilling protocol and the color marking on the trays, as follows:

- **Drillpoint marking:**
creates an incision and marks the insertion point.
- **Option:**
slight smoothing of thin crestal bone conditions with a round bur.
- **Depth drilling ADVANCED:**
Determine the depth and direction of the implant site. The depth drill (∅ 2.0 mm) does not have an integrated depth stop. Depth markings (7.0 mm, 9.0 mm, 11.0 mm, 13.0 mm and 15.0 mm) on the cutting edges of the depth drill indicate when the planned implant length has been reached.
- **Paralleling pin:**
Monitors drilling depth and direction. Can also serve as orientation using the thread of an implant already in place. Secure the paralleling pin with a cord.
- **Surface cutting ADVANCED:**
Marking of bone based on diameter to monitor crestal bone availability around the implant. By using the surface cutter, the insertion depth may be larger than planned.
→ re-examination and adaption of implant length if required
- **Stepped countersinking ADVANCED:**
Gradual expansion of the bone according to the implant contour with increasing diameter until the planned diameter for the implant site has been reached. The stepped countersinks do not have a depth stop. Depth markings (7.0 mm, 9.0 mm, 11.0 mm, 13.0 mm and 15.0 mm) on the cutting edges of the stepped countersink indicate when the planned implant length has been reached. Remove the bone chip that has been collected in the drill. The bone chip may be used as autologous material, and external cooling is improved.
- **Expanding ADVANCED:**
Penetration and stress reduction of the cortical bone and bone expansion according to the implant contour. The expanders should be used with a minimum insertion depth of 7.0 mm and a maximum insertion depth the length of the planned implant. The expanders do not have a depth stop. Depth markings (7.0 mm, 9.0 mm, 11.0 mm, 13.0 mm and 15.0 mm) on the cutting edges of the expander indicate when the planned implant length has been reached. Remove the bone chip that has been collected in the drill. The bone chip may be used as autologous material, and external cooling is improved.

IMPORTANT:

The surface cutter, stepped countersink and expander should be used in low-speed to keep the temperature as low as possible to avoid thermal trauma.

- **Thread cutting ADVANCED:**
It may be necessary to use a thread tap to reduce the implant insertion torque in hard bone. The thread taps do not have a depth stop. Depth markings (7.0 mm, 9.0 mm, 11.0 mm, 13.0 mm and 15.0 mm) on the cutting edges of the thread tap indicate when the planned implant length has been reached. The thread tap can be used either with a machine (max. 10 rpm) or manually with the sure-grip wheel or the torque ratchet.

IMPORTANT:

All drills, except for the marking drill, the surface cutter and the thread tap, can be used with depth-stop sleeves. Availability of a variable drill depth which is dependent on the bone quality is limited when using depth-stop sleeves.

7.5 tioLogic® TWINFIT implant packaging:

The implant is optimally stored and protected in a triangular-shaped holder in the gamma-sterilized blister packaging. The lid of the triangular-shaped holder is the same color as for the coding-coding of the implant diameters. The tioLogic® TWINFIT implant is removed from the holder using the tioLogic® TWINFIT insertion key. This ensures that the implant is not touched before it is inserted. The blister packaging is also protected by outer packaging.

The label on both the outer packaging and the blister packaging states the article number, the name, length and diameter of the implant, the date of manufacture, the sterility expiry date and the LOT number. There are self-adhesive labels in the outer packaging with REF and LOT numbers for documentation in the PatientPass and the surgical protocol.

tioLogic® TWINFIT implants should only be stored in the original packaging at room temperature in a dark, dry place. The implant should not be used after the expiry date (see label). The contents only remain sterile as long as the packaging is undamaged. Do not use the product if the sterile packaging is damaged.

Opening the sterile packaging and removal of implant:

Transfer the blister packaging with the implant, closure screw and depth stops into the sterile area. Peel back the cover of the blister packaging and remove the sterile inner container. Remove the silicone seal with the closure screw.

Note: Do not tilt as the implant lies loose in the inner container.

Grasp the implant with the handpiece tioLogic® TWINFIT insertion key. The insertion key snaps into place on the implant when turned slightly. The insertion key must engage fully in the rotational security. The implant is otherwise not secure and the inner connection of the implant may be damaged during insertion. The insertion key is designed to ensure contact-free insertion, regardless of the indication.

7.6 Implant insertion:

Following removal from the packaging with the tioLogic® TWINFIT insertion key, the implant can be inserted into the prepared implant site either manually or using the handpiece. During insertion, torque must not exceed 40 Ncm. The motor speed during handpiece insertion should not exceed 10 rpm. Excessive torque or motor speed can damage the implant site.

The five markings on the insertion key correspond to the five antirotational mechanisms on the implant and allow alignment of these with regard to the subsequent prosthetic restoration. To position the implant at bone level, the implant is inserted into the bone as far as the lower edge of the polished cervical area, i.e. its final position is slightly transcrestal (0.3 mm). If the implant should be placed subcrestally, adapt the drilling protocol and the drill depth accordingly. Remove the insertion key after implant insertion from the implant (removal should be done carefully in cases of reduced horizontal bone).

Ensure that epithelial tissue does not enter the implant site during implant insertion. Clean the implant internal thread and seal the implant.

7.7 Temporary closure:

To remove the closure screw, press the hex key 1.3 into the closure screw and pull it from the silicone lid. Insert the closure screw into the implant and tighten. Closure screws must fit flush on the implant to ensure that bone tissue cannot penetrate into the implant. The closure screws are inscribed with S, M or L corresponding to the series of abutments. Closure screws are intended for single use only.

If open healing is planned or indicated, insert the relevant gingiva former (S, M or L) into the implant instead of a closure screw.

Tightening torque	
Closure screw	15 Ncm
Gingiva former	15 Ncm

7.8 Wound closure:

After checking the operation site, close the wound by suturing. Ensure that the wound closure is saliva-proof and that there is good blood circulation. For open healing, ensure that the tissue is sutured close to the gingiva former. After the implant insertion is completed, x-ray the site to monitor the fit and position of the implant.

Important information and implant care by the patient following implantation:

Cool externally and avoid mechanical impact on the operation site directly after the operation. A basic requirement for long-term success of the dental implant is the patient's meticulous oral hygiene and regular check-ups with the dentist or implantologist. Implant data and appointments are documented at the dental office in the patient's file and the PatientPass.

- Temporary denture (non-implant-borne):
- Do not fit a temporary prosthetic restoration until at least 14 days after implant insertion. Always ensure that there is no mechanical load on the fitted implant. The restoration should be relieved over the implants and fitted with soft lining. If there are residual teeth, a temporary prosthetic restoration is generally fabricated on the abutment teeth prior to implant placement, or an existing denture is converted.

7.9 Immediate restoration (temporary abutment):

It is possible to fit a long-term, non-functional immediate temporary restoration on implants if there is absolute primary stability and no recession of the implant site. In aesthetically relevant areas, the peri-implant structures are retained with a temporary abutment. An optimum impression can be taken after formation of the structures.

Temporary abutments with platform connector geometry are available for the S, M and L series of abutments. They are supplied non-sterile and are made from high-strength plastic (PEEK), which can be quickly and easily customized. Temporary abutments can be in place for a maximum of 180 days. Either face the temporary abutment directly with acrylic or fit with a temporary crown / bridge. In both cases, secure the abutment intra-orally with the screw for the temporary abutment.

Tightening torque	
temporary abutment on the model / in the mouth	15 Ncm

7.10 Healing stage:

The healing stage in the mandible is normally 3 months and in the maxilla 6 months. It can be shorter or longer depending on the bone quality, the surgical procedure used and the patient's anatomy.

If examinations after the healing stage indicate osseointegration of the implant, begin with the prosthetic restoration.

Include the patient in a regular recall program after the restoration has been fitted to ensure the long-term success of the implant restoration. Give patients instructions on the appropriate hygiene regime for the implants and prosthetic restoration.

7.11 Implant exposure and gingiva management:

Expose the implant after the healing stage. The patient should be prepared in the same way as for other surgical procedures. The patient is given a local anaesthetic. The implants can be exposed using different techniques and aids, e.g. by using a scalpel or laser. If a scalpel or a laser is used, the peri-implant tissue (attached gingiva) is conserved and optimum aesthetic results are achieved (gingiva management).

7.12 Gingiva forming:

Gingiva formers, anatomically shaped or cylindrical, or 4Base abutments used directly – particularly gentle on the tissue – are available for the implantologist to ensure optimal gingiva management. The anatomical gingiva formers are designed to form a wide gingival contour. Depending on the type of prosthetic restoration, this can make it easier for the implantologist to fit the restoration. Select the gingiva formers according to the series of abutments, gingival height and insertion depth of the implant. They are available with a platform connector geometry for the series of abutments S, M or L in different gingival heights (laser-marked). Gingiva formers can also be used with open healing of the implant for specific indications and for preserving the soft tissue.

Important: Clean and sterilize the gingiva formers and 4Base abutments before insertion in the implant.

If a temporary restoration is fitted, relieve the denture during gingiva forming. The impression should not be taken until the tissue is completely free of inflammation.

Tightening torque	
Gingiva former	15 Ncm

8. Technical specification

Composition

Chemical composition of the implants (% by mass)										
tioLogic® TWINFIT Ø in mm	Material Number (DIN, AISI)	Material Description	Al	V	Fe	C	N	H	O	Ti
4.8/ 5.5	3.7065 (DIN EN ISO 5832-2)	Pure titanium (Ti)	–	–	≤ 0.5	≤ 0.08	≤ 0.05	≤ 0.0125	≤ 0.4	Resi- due
3.3/ 3.7/ 4.2	3.7165 (DIN EN ISO 5832-3)	Titanium alloy (Ti6Al4V)	5.5 – 6.75	3.5 – 4.5	≤ 0.3	≤ 0.08	≤ 0.05	≤ 0.015	≤ 0.2	Resi- due

Product	Composition
tioLogic® TWINFIT closure screw	3.7165
tioLogic® TWINFIT gingiva former	3.7165
tioLogic® TWINFIT temporary abutment	PEEK / 3.7165

Information about the composition is included in the material list for our implantology products (REF 989-801-05), or see www.dentaurum.com/ifu. The SSCP is available at <https://ec.europa.eu/tools/eudamed>.

9. Adverse reactions and residual risks

Possible complications and risks during surgical treatment (temporary adverse reactions) are:

Pain, swellings, difficulties speaking, inflammation of the gums.

Possible complications and risks during surgical treatment (permanent adverse reactions) are:

Chronic pain in connection with the implant, permanent paresthesia, dysesthesia, nerve damage, exfoliation, hyperplasia, local or systemic infections, oroantral or oronasal fistulas, alveolar bone loss in maxilla/mandible, irreversible damage to adjacent teeth, fractures to the implant, jaw, bone or restoration, aesthetic problems.

10. Safety information and precautions

- Do not use the product if there is a known allergic reaction to one or more of the material components.
- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. Therefore, secure all articles used intra-orally against swallowing and/or aspiration.
- Different types of alloy in the oral cavity can lead to galvanic reactions.
- Do not use the product if the sterile packaging is damaged.
- Only use tioLogic® TWINFIT prosthetic components in combination with tioLogic® TWINFIT implants. Use of third-party components can negatively impact the function and safety. In such cases, Dentaurum GmbH & Co. KG accepts no liability for replacement products or guarantees.
- Use drills, thread taps and components only with compatible implant systems and diameters. Using these tools with other implant systems and diameters can cause damage to the bone and tissue and can produce aesthetically unsatisfactory results.
- Observe the indications for the individual products. Consider all indications that are not listed explicitly to be contraindications.
- Disinfection, cleaning and sterilizing should not be done in the original packaging.
- Products must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.
- Technical advice (oral and written) is based on the scientific and technical knowledge available at the time the product is put on the market. It does not release the dental professional from the responsibility of personally checking the suitability of the products for the intended indications and applications. Advice is only given as a non-binding recommendation and cannot be assumed to provide any form of assurance, liability or guarantee.
- All products are subject to continuous development based on current scientific knowledge, and we reserve the right to make changes in the construction, design or material of the products.
- Under certain circumstances such as, for example, non-observance of the indications or the surgical procedure (Instructions for use), it may be necessary to remove an implant. The dental professional must assess the degree of osseointegration at the planned time for removal.
- Possible reasons for the removal of non-osseointegrated implants:
- Surgical complications, insufficient osseointegration, peri-implantitis, suboptimal position of the implant, infection of the surrounding tissue, patient's request.

Severe incidents:

Report all severe incidents arising from the use of the product to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.

Severe incidents are identified as a drop in performance or deterioration of the properties of a medical device which resulted or may have resulted in a direct or indirect temporary or permanent deterioration in a patient's state of health, in the death of a patient, a dental professional or other persons or as posing a serious danger to public health.

11. Storage and shelf life

Store the tioLogic® TWINFIT implants, abutments and components always in a dry, dark place at room temperature in the original packaging. Do not use the product if the sterile packaging is damaged. Do not use implants after the expiry date stated on the label.

12. Disposal

Adhere strictly to the rules that apply for your region and the details outlined in the safety data sheets, if available. Products to be disposed of should be treated as surgical waste and be decontaminated in compliance with the relevant regulations for medical waste disposal.

13. Information for single use products

tioLogic® TWINFIT implants are gamma-sterilized on delivery. Prosthetic abutment components are delivered non-sterile. All implants and prosthetic abutment components are for single use on one patient only.

Reconditioning of a tioLogic® TWINFIT implant and abutment components (recycling) or reuse on patients are not permitted as safe reconditioning and/or functional safety cannot be guaranteed and may lead to a loss in mechanical, chemical and/or biological properties. Reuse may cause infection. Do not use implants after the expiry date stated on the label.

14. Information for multiple use products

All drills, thread taps and instruments for implantation are delivered non-sterile and are intended for multiple use. The product must be disinfected, cleaned and sterilized before being used initially and before being used again. Drills, instruments and components must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.

15. Information on disinfection / cleaning / sterilizing

The tioLogic® TWINFIT prosthetic abutment components are supplied non-sterile. They are for single use on one patient only. They must be cleaned and disinfected before and after being used on a patient (e.g. transfer to dental laboratory). Dentaaurum generally recommends that the product be sterilized in addition, as non-sterile components may cause infections. Only sterilized prosthetic components may be used if bleeding occurs.

16. Product overview and scope of delivery

Please refer to the tioLogic® TWINFIT product catalog for the tioLogic® TWINFIT product range.

17. Quality information

Dentaaurum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaaurum has no influence on the processing on site.

18. Explanation of symbols

Please refer to the label. See chapter "Additional information".

19. Information on MRT safety

The products within these Instructions for use are MRT safe.

20. Additional information

Additional information can be found at www.dentaaurum.com/ifu

Instructions for use	REF
tioLogic® TWINFIT closure screws, gingiva formers, impression posts, temporary abutments	989-365-90
Torque ratchet	989-719-00
Preparation instructions for the implant system components	989-801-07
Explanation of symbols used on the label	989-313-00
MRI safety of Dentaaurum products	989-804-43

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