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Placement of the implants and fabrication of the prosthetic restoration should take into account the particular oral situation to avoid overloading the components. Using any components in combination contrary to the instructions for use/manual or that are not described as being specifically combinable can cause mechanical failure, damage to the tissue or unsatisfactory aesthetic results.

CITO mini® implants are not known to have any side effects or produce any interaction. It cannot, however, be ruled out that in rare cases allergies to components used in the materials of the CITO mini® implant system may occur or that there may be electrochemically induced discomfort.

### 6. Further information

Though placement of dental implants has a high rate of success and implants have a long durability, successful treatment cannot be guaranteed. The operator should note and document any problem cases and inform the manufacturer.

An inadequate number of implants, implants with insufficient length or diameter, unfavorable positioning of the implants or a statically poor prosthetic restoration can cause premature implant loss and fatigue fractures in implants and prosthetic components under biomechanical loading.

Placement of the implants and fabrication of the prosthetic restoration should take into account the particular oral situation to avoid overloading the components. Using any components in combination contrary to the instructions for use/manual or that are not described as being specifically combinable can cause mechanical failure, damage to the tissue or unsatisfactory aesthetic results.

CITO mini® implants are not known to have any side effects or produce any interaction. It cannot, however, be ruled out that in rare cases allergies to components used in the materials of the CITO mini® implant system may occur or that there may be electrochemically induced discomfort.

All products are subject to continuous development based on current scientific knowledge and we reserve the right to make changes in the construction, design or material of the products.

### 7. Storage

CITO mini® implants should always be stored in a dry, dark place at room temperature in the original packaging. The implants should not be used after the sterility expiry date (see label).

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Dentaurum Implants GmbH | Turnstr. 31 | 75228 Ispringen | Germany

### 2. Brief description

CITO mini® implants are designed for endosteal insertion in the maxilla or mandible. The one-piece implants are provided with a Ø 1.8 mm ball, on top of which the female parts of the connection are fixated, once they are polymerized into the prosthesis of the patient. The CITO mini® implant system contains specially coordinated instruments, abutments and accessories for placement of the implants and fabricating the prosthetic restoration. Only original components of the CITO mini® implant system should be combined in accordance with the instructions for use / user manual.

### 3. Indications

Permanent fixation of partial and total prostheses.

### 4. Contraindications

#### General contraindications:

General contraindications for dental surgery procedures apply. These include: reduced immunodeficiency/steroid treatment/blood coagulation disorders/uncontrolled endocrine diseases/rheumatic disorders/bone system diseases/cirrhosis of the liver/drug, alcohol or tobacco abuse/depression, psychopathic disorders/poor patient compliance/chronic inflammatory underlying diseases.

#### Local / personal contraindications:

Osteomyelitis/radiotherapy in the head region/recurrent mucosal diseases/temporo-mandibular joint dysfunctions/parafunctions/lack of vertical or horizontal bone availability, jaw defects, inadequate bone quality/poor oral hygiene.

### 5. Further information

Though placement of dental implants has a high rate of success and implants have a long durability, successful treatment cannot be guaranteed. The operator should note and document any problem cases and inform the manufacturer.

An inadequate number of implants, implants with insufficient length or diameter, unfavorable positioning of the implants or a statically poor prosthetic restoration can cause premature implant loss and fatigue fractures in implants and prosthetic components under biomechanical loading.

Placement of the implants and fabrication of the prosthetic restoration should take into account the particular oral situation to avoid overloading the components. Using any components in combination contrary to the instructions for use/manual or that are not described as being specifically combinable can cause mechanical failure, damage to the tissue or unsatisfactory aesthetic results.

CITO mini® implants are not known to have any side effects or produce any interaction. It cannot, however, be ruled out that in rare cases allergies to components used in the materials of the CITO mini® implant system may occur or that there may be electrochemically induced discomfort.

### 6. Packaging and sterility

All CITO mini® implants are supplied singly with the appropriate closure screws in gamma-sterilized double packaging (foil and blister packaging) that protects the content against contamination. Outer packaging provides added protection to the double packaging. The product is only sterile if the packaging is not damaged. If the double packaging is visibly damaged, the product must not be used. The CITO mini® implants are single-use products. All other components of the CITO mini® implant system are supplied non-sterile. The surgical components should be thoroughly disinfected, cleaned and sterilized before use. This applies in particular to components that can be disassembled, e. g. the torque wrench.

### 10. Explanation of symbols

Refer to the label. Additional information can be found at [www.dentaurum.de](http://www.dentaurum.de) (Explanation of symbols REF 989-313-00).

Date of Information: 05/15

Subject to modifications

### 7. Storage

CITO mini® implants should always be stored in a dry, dark place at room temperature in the original packaging. The implants should not be used after the sterility expiry date (see label).

### 8. Application / Availability / Precautions / Documentation

The CITO mini® product range is supplied exclusively to doctors, dentists and dental technicians. It should only be used by doctors, dentists or dental technicians who are familiar with dental implantological procedures, including diagnosis, preoperative planning, surgical techniques and prosthetic treatment.

Before use, operators should ensure that they have carefully read and understood the full CITO mini® instructions for use/manuals. This quick reference guide does not replace the instructions for use/manuals for the CITO mini® implant system. As the instructions and manuals cannot provide all information for immediate use, we strongly recommend that, before using the system, operators attend a CITO mini® training course offered by Dentaurum Implants to learn the correct techniques.

When choosing the number of implants to be inserted, it is important to consider the bone quality and whether the implants will be loaded immediately. To fixate total prostheses, Dentaurum Implants recommends at least 4 implants in the mandible and 6 implants in the maxilla.

When positioning the implants, adhere to the minimum distance between the implants. Consider the size of the female part of the connection.

Immediate implant loading is only possible, if a torque of 35 Ncm was reached during insertion. The torque of 45 Ncm should not be exceeded during insertion.

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Development, clinical testing, production and quality control of the CITO mini® product range are completed in accordance with the Medical Device Directive 93/42/EEC.

Sections 9 and 10 of our General Terms of Delivery and Payment apply with regard to warranty or liability – unless otherwise provided in this quick reference guide.

Warranty and liability are rendered void in particular if the products are not used by the operator or a third party in accordance with the instructions for use; this also applies if the CITO mini® product range is used in combination with products of other manufacturers, which have not been specifically recommended for use by Dentaurum Implants.

Dentaurum Implants has no control over processing and application of the product, which are the sole responsibility of the user. Technical advice (oral and written) is based on the scientific and technical knowledge available when the product is put on the market. It does not release the user from the responsibility of personally checking the suitability of the products for the intended indication and application. Advice is only given as a non-binding recommendation, which cannot be assumed to provide any form of assurance or guarantee.

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