CITO mini®

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Dentaurum GmbH & Co. KG · Turnstr. 31 · 75228 Ispringen · Germany

2. General device description

CITO mini® dental implants are designed for insertion in the endosteal region of the maxilla or mandible. They are used for anchoring functional and aesthetic oral rehabilitations in partially edentulous and fully edentulous patients.

| CITO mini® implant data | | | |
|-------------------------|-------|-------|-------|
| Implant diameters in mm | 1.8 | 2.2 | 2.5 |
| Implant lengths in mm | 11/13 | 11/13 | 11/13 |

The one-piece implants are fitted with a ball \varnothing 1.8 mm, on top of which matrices of the system are fixated once they have been polymerized into the patient's dental prosthesis. Depending on the indication, there are matrices with different withdrawal forces and freedom of movement for divergencies between the implants.

The implants can be used for surgical immediate implantation, delayed immediate implantation, and later implantation, either one- or two-phase. CITO mini® implants have a cylindrical / conical external geometry with a blasted / etched surface. The surface of the implant shoulder is machined.

The CITO mini® implant system contains specifically coordinated abutments and components, including instruments, for the insertion and prosthetic restoration of the implants.

3. Intended purpose

Implants are designed for insertion in the endosteal region of the maxilla or mandible. They are used for anchoring functional and aesthetic oral rehabilitations.

4. Intended users

The CITO mini® product range is supplied exclusively to qualified personnel such as surgeons, dentists and dental technicians. The products should only be used by qualified personnel who are familiar with dental implantology, including diagnosis, preoperative planning, surgical technique and prosthetic treatment.

5. Medical indications / contraindications

5.1 Indications

CITO mini® implant types are used to stabilize prostheses in the edentulous or partially edentulous maxilla and mandible. Consider the risks and advantages / disadvantages of dental implants and alternative methods of treatment before deciding on the indication for an implantation.

Whenever implants are used, the implant diameter and length of the CITO mini® implant types should be in proportion to the prosthetic restoration.

When choosing the number of implants to be inserted, consider the bone quality and whether the implants will be loaded immediately.

To fixate full dentures, Dentaurum recommends at least four implants in the mandible and six implants in the maxilla. For restorations in the maxilla, use CITO mini® implants \varnothing 2.2 mm or \varnothing 2.5 mm.

When positioning the implants, pay attention to keeping an adequate distance between the implants and to ensuring that the static force on the bone is optimally applied. The size of the matrix should also be considered.


Immediate implant loading is only possible if a torque of 35 Ncm was reached during insertion. Torque should not exceed 40 Ncm during insertion.

Planning scheme CITO mini® in the maxilla

- strategic pillar
- optional pillar
- existing tooth

The number of pillars required is based on the consensus paper from the German implantological societies BDO, DGMKG, DGI, DGZI in 2009 (mandible = 4 pillars / maxilla = 6 pillars for removable dentures)^{1,2}.

The number and the position of the optional pillars should be considered separately in each case and chosen at the discretion of the implantologist.

Literature: ¹ Mundt T, Lucas C, Biffar R, Heinemann F. Stabilisierung von Teilprothesen mit Mini-Implantaten – zwei Fallberichte [Stabilization of partial dentures with mini implants - two case studies]. Dtsch Zahnärztl Z 2015;70:416-424.
² Mundt T, Al Jaghisi A, Schwahn B, Hilgert J, Lucas C, Biffar R, Schwahn C, Heinemann F. Immediate versus delayed loading of strategic mini dental implants for the stabilization of partial removable dental prostheses: a patient cluster randomized, parallel-group 3-year trial. BMC Oral Health. 2016 Jul 30;17(1):30. doi: 10.1186/s12903-016-0259-z

Planning scheme CITO mini® in the mandible

- strategic pillar
- optional pillar
- existing tooth

The number of pillars required is based on the consensus paper from the German implantological societies BDO, DGMKG, DGI, DGZI in 2009 (mandible = 4 pillars / maxilla = 6 pillars for removable dentures) ^{1,2}.

The number and the position of the optional pillars should be considered separately in each case and chosen at the discretion of the implantologist.

| Class 0 | Class 1 | Class 2 | Class 3 | Class 4 | Class 5 |
|--|---|---|--|---|---|
| <ul style="list-style-type: none"> edentulous | <ul style="list-style-type: none"> only anterior teeth | <ul style="list-style-type: none"> no canine 1 tooth from pos. 4 possibly anterior teeth | <ul style="list-style-type: none"> no canine min. 2 teeth from pos. 4 possibly anterior teeth | <ul style="list-style-type: none"> canine possibly anterior teeth | <ul style="list-style-type: none"> canine min. 1 tooth from pos. 4 possibly anterior teeth |
| | | | | | |

Literature: ¹ Mundt T, Lucas C, Biffar R, Heinemann F. Stabilisierung von Teilprothesen mit Mini-Implantaten – zwei Fallberichte [Stabilization of partial dentures with mini implants - two case studies]. Dtsch Zahnärztl Z 2015;70(4):16-424.
² Mundt T, Al Jaghni A, Schwahn B, Hilgert J, Lucas C, Biffar R, Schwahn C, Heinemann F. Immediate versus delayed loading of strategic mini dental implants for the stabilization of partial removable dental prostheses: a patient cluster randomized, parallel-group 3-year trial. BMC Oral Health. 2016 Jul 30;17(1):30. doi: 10.1186/s12903-016-0259-z

Planning scheme CITO mini® (REF 989-502-45)

5.2 Contraindications

General contraindications:

Insufficient bone or tissue available and/or inadequate bone quality, local root remnants, bone disease and wound healing disorders, local infection of implant site, serious therapy-resistant functional disorders, uncontrolled diabetes mellitus, long-term immunosuppressive therapy, connective tissue disease, blood disorders (e.g. leukemia, hemophilia), intra-oral infection or malignancies, uncontrolled parafunctional habits, non-treatable occlusal or articulation conditions, serious psychological illnesses, xerostomia and titanium allergies.

Local / personal contraindications:

Incomplete growth of jaw, unfavorable anatomical bone situations, previously irradiated bone, temporomandibular joint disorders and treatable pathologic diseases to the jaw, severe systemic diseases, reduced immune defense and leukocyte disfunctions leading to increased risk of infection, endocrine disorders, medical anticoagulation / hemorrhagic diathesis, arteriosclerosis and stroke, hypertension, heart attack, illnesses with periodic use of steroids, hepatitis, diabetes mellitus and pregnancy, metabolic bone disorders, all illnesses affecting bone regeneration or blood microcirculation, rheumatic diseases and drug/alcohol/tobacco addiction, lack of patient compliance, inadequate oral hygiene, periodontitis, bruxism, parafunctional habits and oral mucosal abnormalities.

An inadequate number of implants, implants with insufficient length or diameter, unfavorable positioning of the implants or a statically poor prosthetic restoration can cause premature implant loss and fatigue fractures in implants and abutments under biomechanical loading.

NOTE: The insertion of CITO mini® implants in soft bone is not indicated.

6. Patient target group

All skeletally mature patients undergoing implant treatment. Please consider the indications / contraindications.

7. Information on use, processing, fitting and insertion

Prior to use, the dental professional ensures that they have carefully read and understood the CITO mini® Instructions for use.

We strongly recommend that, prior to using Dentaurem products, dental professionals attend a Dentaurem training course to familiarize themselves with the Dentaurem products.

Carry out thorough radiological, psychological and physical examinations before implantation and inform patients of the process in detail. Dentaurem recommends full clinical, radiological, photographic and statistical documentation.

There are self-adhesive labels in the outer packaging with REF and LOT numbers for documentation in the PatientPass and the surgical protocol.

The implantologist should ensure the products cannot be aspirated during intra-oral use.

Precautions:

Before surgical implantation, it is important to carry out a careful patient anamnesis and, if required, a third-party anamnesis by the patient's doctor (GP/PCP) in order to clarify if 1) implantation may be difficult due to an anatomical situation, 2) there is the risk of a serious surgical problem or a general risk, 3) wound healing/osseointegration may be difficult, or 4) good hygiene and care of the implant, abutment and restoration may worsen. No dental implants should be fitted if any of the above cases or a combination hereof is severe or uncontrollable.

7.1 Preparation for surgery

Ensure that hygiene measures are carefully followed, documented and validated throughout the surgical procedure. Prepare the treatment room, instrumentarium and patient under these aspects. Check all surgical instruments required for the operation to ensure they are complete, functional and sterile. We recommend having several implants and surgical instruments available as a precaution.

All rotary CITO mini® drills and marking drills are supplied non-sterile and must be sterilized before use. They should be thoroughly cleaned, disinfected, conditioned and sterilized before being used for the first time (factory new) and immediately after each use. Rotary instruments should always be checked to ensure that they are sharp, in good condition and the markings are legible, as they have a limited service life. Instruments can become blunt as a result of use and cleaning. Only instruments that are sharp and free from defects should be used. Rotary instruments can be reused in hard bone 15 to 20 times if used with proper care and provided that they are not damaged or contaminated. Avoid any further reuse or the use of damaged and / or contaminated instruments. The dental professional bears responsibility. No liability is accepted if these instructions are disregarded.

7.2 Notes on operation planning

Implants need a distance of at least 1.5 mm to the neighboring teeth. The distance between the implants must be at least 3.0 mm. Select the implant diameter and length so that there is adequate bone around the implant (min. 1.0 mm).

For technical reasons, the depth drill and the stepped countersink are 0.6 mm longer (apically) than the given preparation length. This should be taken into account during diagnostics, planning and preparation. The distance to the mandibular nerve and inferior alveolar nerve should be at least 1.5 mm!

7.3 General information on implant site preparation

Careful preparation of the hard and soft tissue in the bone is important for optimum healing of the implant site. Avoid thermal or mechanical trauma at all costs. To this end, keep the temperature as low as possible when preparing the site and do not exceed the maximum implant insertion torque (max. 40 Ncm).

| Article | | Maximum torque |
|------------------------------------|---|----------------|
| Marking drill |  | max. 800 rpm |
| CITO mini® stepped drill, ø 1.8 mm |  | max. 500 rpm |
| CITO mini® stepped drill, ø 2.2 mm |  | max. 500 rpm |
| CITO mini® stepped drill, ø 2.5 mm |  | max. 500 rpm |

7.4 Implant site preparation

Select the rotary instruments and their drilling sequence and depth for the preparation of the implant site depending on the bone quality. A drilling protocol for different bone qualities (soft, medium, hard) is available for the implantologist. Determination of the bone quality rests on the implantologist. Use the drills with intermittent drilling and ensure sufficient cooling using a sterile, cooled physiological saline solution (5 °C / 41 °F) as the diameter increases.

Carry out implant site preparation using the CITO mini® drilling system, bearing in mind the drilling protocol, as follows.

- **Drillpoint marking:** creates an incision and marks the insertion point.
- **CITO mini® stepped drill:**

For technical reasons, the stepped drill is 0.6 mm longer than the given preparation length. Take this into account during diagnosis and preparation. The stepped drill has no integrated depth stop. Depth markings (6.0 mm, 11.0 mm, 13.0 mm and 15.0 mm) on the cutting edges of the depth drill indicate when the planned implant length has been reached. Drill with external cooling using a sterile, cooled physiological saline solution (5 °C / 41 °F). Drilling should be intermittent without applying pressure to ensure that the tip of the drill can cool.


IMPORTANT:
Should the bone be very hard, it may be necessary to use a stepped drill one size larger to ream the cortical bone.

Drilling protocol

Drilling protocol with CITO mini® drills. Taking different bone qualities into account.

The drilling protocol may need to be adapted depending on the indication and the individual situation of the patient. The insertion of CITO mini® implants in soft bone is not indicated.

CITO mini®

| | Soft bone quality | Medium bone quality | Hard bone quality |
|---|--|---|---|
| |  |  |  |
| | ø 1.8 ³ ø 2.2 ø 2.5 | ø 1.8 ³ ø 2.2 ø 2.5 | ø 1.8 ³ ø 2.2 ø 2.5 |
| Marking drill ² | | X X X | X X X |
| Stepped drill for ø 1.8 mm ¹ | | min. 6.0 mm | X ⁴ |
| Stepped drill for ø 2.2 mm ¹ | | min. 6.0 mm | X ⁴ |
| Stepped drill for ø 2.5 mm ¹ | | min. 6.0 mm | X ⁴ |

The insertion of CITO mini® implants in soft bone is not indicated.

¹ The insertion depth / length of the stepped drill depends on the desired primary stability and the bone quality. The depth markings / scaling help to avoid damaging neighboring structures (e.g. nerves) due to perforations that are too deep.
² Optional use.
³ Not intended for restorations in the maxilla.
⁴ Should the bone be very hard, it may be necessary to use a stepped drill one size larger to ream the cortical bone.

7.5 CITO mini® implant packaging:

The implant is optimally stored and protected in a triangular-shaped holder in the gamma-sterilized blister packaging. The CITO mini® implant is removed from the holder using the CITO mini® PentaGrip insertion key. This ensures that the implant is not touched before it is inserted. The blister packaging is also protected by outer packaging.

The label on both the outer packaging and the blister packaging states the article number, the name, length and diameter of the implant, the date of manufacture, the sterility expiry date and the LOT number. There are self-adhesive labels in the outer packaging with REF and LOT numbers for documentation in the PatientPass and the surgical protocol.

CITO mini® implants should only be stored in the original packaging at room temperature in a dark, dry place. Do not use the implant after the expiry date (see label). The contents only remain sterile as long as the packaging is undamaged. Do not use the product if the sterile packaging is damaged.

Opening the sterile packaging and removal of implant:

Transfer the blister packaging with the implant into the sterile area. Peel back the cover of the blister packaging and remove the sterile inner container. Remove the silicone seal.

Note: Do not tilt as the implant lies loose in the inner container.

Grasp the implant with the handpiece CITO mini® PentaGrip insertion key. The insertion key snaps into place on the implant when turned slightly. The insertion key must engage fully on the antirotational mechanism. The implant is otherwise not secure and may be damaged during insertion. The CITO mini® PentaGrip insertion key is designed to ensure contact-free insertion, regardless of the indication.

7.6 Implant insertion:

Once grasped by the PentaGrip insertion key, the CITO mini® implant can be inserted into the prepared implant site. The PentaGrip insertion key can be extended using a drill extension and can be used manually by attaching the adapter. A torque of 40 Ncm should not be exceeded with any insertion procedure.

It is necessary to achieve a torque of min. 35 Ncm for immediate loading of the implant. Maintain a maximum motor speed of 10 rpm during handpiece insertion. Excessive torque or motor speed can damage the implant site.

For manual insertion, an adapter ISO shank hexagon/ratchet is available for use with the torque ratchet or sure-grip wheel. The adapter is placed on the PentaGrip insertion key. The CITO mini® implant is inserted up to the planned and prepared position in the bone. Ensure that epithelial tissue does not enter the implant site during implant insertion. If the implant is difficult to insert, the implant site should be flushed and prepared again. Remove the insertion key after implant insertion from the rotational security of the CITO mini® implant (removal should be done carefully in cases of reduced horizontal bone).

Important information and implant care by the patient following implantation:

Cool externally and avoid mechanical impact on the operation site directly after the operation. A basic requirement for long-term success of the dental implant is the patient's meticulous oral hygiene and regular check-ups with the dentist or implantologist. Implant data and appointments are documented at the dental office in the patient's file and the PatientPass.

7.7 Temporary denture (non-implant-borne):

Do not fit a temporary prosthetic restoration until at least 14 days after implant insertion. Always ensure that there is no mechanical load on the fitted implant. The restoration should be relieved over the implants and fitted with soft lining. If there are residual teeth, a temporary prosthetic restoration is generally fabricated on the abutment teeth prior to implant placement, or an existing denture is converted.

7.8 Healing stage:

The healing stage in the mandible is normally 3 months and in the maxilla 6 months. It can be shorter or longer depending on the bone quality, the surgical procedure used and the patient's anatomy.

If examinations after the healing stage indicate osseointegration of the implant, begin with the prosthetic restoration.

Patients should be entered into a regular recall program after the restoration has been fitted to ensure the long-term success of the implantological restoration. Patients should be given instructions on the appropriate hygiene regime for the implants and prosthetic restoration.

8. Technical specification

Composition

| Chemical composition of the implants (% by mass) | | | | | | | | | | |
|--|----------------------------------|-----------------------------|------------------|-----------------|---------------|----------------|----------------|-----------------|---------------|--------------|
| CITO mini® ø in mm | Material Number (DIN, AISI) | Material Description | Al | V | Fe | C | N | H | O | Ti |
| 1.8/ 2.2/ 2.5 | 3.7165 (DIN EN ISO 5832-3) | Titanium alloy (Ti6Al4V) | 5.5 – 6.75 | 3.5 – 4.5 | ≤ – 0.3 | ≤ – 0.08 | ≤ – 0.05 | ≤ – 0.015 | ≤ – 0.2 | Resi- due |

Information about the composition is included in the material list for our implantology products (REF 989-801-05), or see www.dentaurum.com/ifu. The SSCP is available at <https://ec.europa.eu/tools/eudamed> and www.dentaurum.com.

9. Adverse reactions and residual risks

Possible complications and risks during surgical treatment (temporary adverse reactions) are:

Pain, swellings, difficulties speaking, inflammation of the gums.

Possible complications and risks during surgical treatment (permanent adverse reactions) are:

Chronical pain in connection with the implant, permanent paresthesia, dysesthesia, nerve damage, exfoliation, hyperplasia, local or systemic infections, oroantral or oronasal fistulas, alveolar bone loss in maxilla/mandible, irreversible damage to adjacent teeth, fractures to the implant, jaw, bone or restoration, aesthetical problems.

10. Safety information and precautions

- Do not use the product if there is a known allergic reaction to one or more of the material components.
- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. Therefore, secure all articles used intra-orally against swallowing and/or aspiration.
- Different types of alloy in the oral cavity can lead to galvanic reactions.
- Do not use the product if the sterile packaging is damaged.
- Only use CITO mini® prosthetic components in combination with CITO mini® implants. Use of third-party components can negatively impact the function and safety. In such cases, Dentaurum GmbH & Co. KG accepts no liability for replacement products or guarantees.
- Use drills, thread taps and components only with compatible implant systems and diameters. Using these tools with other implant systems and diameters can cause damage to the bone and tissue and can produce aesthetically unsatisfactory results.
- Observe the indications for the individual products. Consider all indications that are not listed explicitly to be contraindications.
- Disinfection, cleaning and sterilizing should not be done in the original packaging.
- Products must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.
- Technical advice (oral and written) is based on the scientific and technical knowledge available at the time the product is put on the market. It does not release the dental professional from the responsibility of personally checking the suitability of the products for the intended indications and applications. Advice is only given as a non-binding recommendation and cannot be assumed to provide any form of assurance, liability or guarantee.
- All products are subject to continuous development based on current scientific knowledge, and we reserve the right to make changes in the construction, design or material of the products.
- Under certain circumstances such as, for example, non-observance of the indications or the surgical procedure (Instructions for use), it may be necessary to remove an implant. The dental professional must assess the degree of osseointegration at the planned time for removal.
- Possible reasons for the removal of non-osseointegrated implants:
- Surgical complications, insufficient osseointegration, peri-implantitis, suboptimal position of the implant, infection of the surrounding tissue, patient's request.

Severe incidents:

Report all severe incidents arising from the use of the product to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.

Severe incidents are identified as a drop in performance or deterioration of the properties of a medical device which resulted or may have resulted in a direct or indirect temporary or permanent deterioration in a patient's state of health, in the death of a patient, a dental professional or other persons or as posing a serious danger to public health.

11. Storage and shelf life

Store the CITO mini® implants, abutments and components always in a dry, dark place at room temperature in the original packaging. If the packaging shows serious damage, the product must be checked for integrity and cleanliness before use, and if necessary must be disposed of. Do not use implants after the expiry date stated on the label.

12. Disposal

Adhere strictly to the rules that apply for your region and the details outlined in the safety data sheets, if available. Products to be disposed of should be treated as surgical waste and be decontaminated in compliance with the relevant regulations for medical waste disposal.

13. Information for single use products

CITO mini® implants are gamma-sterilized on delivery. Prosthetic abutment components are delivered non-sterile. All implants and prosthetic abutment components are for single use on one patient only.

Reconditioning of a CITO mini® implant and abutment components (recycling) or reuse on patients are not permitted as safe reconditioning and/or functional safety cannot be guaranteed and may lead to a loss in mechanical, chemical and/or biological properties. Reuse may cause infection. Do not use implants after the expiry date stated on the label.

14. Information for multiple use products

All CITO mini® drills and instruments for implantation are delivered non-sterile and are intended for multiple use. The product must be disinfected, cleaned and sterilized before being used initially and before being used again. Drills, instruments and components must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.

15. Information on disinfection / cleaning / sterilizing

The CITO mini® prosthetic abutment components are supplied non-sterile. They are for single use on one patient only. They must be cleaned and disinfected before and after being used on a patient (e.g. transfer to dental laboratory). Dentaureum generally recommends that the product be sterilized in addition, as non-sterile components may cause infections. Only sterilized prosthetic components may be used if bleeding occurs.

16. Product overview and scope of delivery

Refer to the CITO mini® order form for the CITO mini® product range.

17. Quality information

Dentaureum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaureum has no influence on the processing on site.

18. Explanation of symbols

Please refer to the label. See chapter "Additional information".

19. Information on MRT safety

The products within these Instructions for use are MRT safe.

20. Additional information

Additional information can be found at www.dentaureum.com/ifu

| Instructions for use | REF |
|--|------------|
| CITO mini® matrices and components | 989-804-97 |
| Torque ratchet | 989-719-00 |
| Preparation instructions for the implant system components | 989-801-07 |
| Explanation of symbols used on the label | 989-313-00 |
| MRI safety of Dentaureum products | 989-804-43 |

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