remanium® **R**

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**Dear customer,**

Thank you for choosing a quality product from Dentaureum.

It is essential to read these Processing Instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit.

It is not possible to describe every eventuality and possible application in these instructions. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at [www.dentaureum.com](http://www.dentaureum.com), even though you may frequently use the same product.

**Manufacturer**

Dentaureum GmbH & Co. KG | Turnstr. 31 | 75228 Ispringen | Germany

**Quality**

Dentaureum ensures faultless quality of its products. These recommendations are based upon Dentaureum's own experiences. The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as we, Dentaureum, have no influence on the processing on site.

# Contents

<b>1. General instructions for processing</b> .....	<b>14</b>
<b>2. Wax up</b> .....	<b>14</b>
<b>3. Spruing</b> .....	<b>15</b>
<b>4. Calculation of metal quantity</b> .....	<b>15</b>
<b>5. Investing</b> .....	<b>15</b>
<b>6. Wax burn-out and preheating</b> .....	<b>15</b>
<b>7. Melting and casting</b> .....	<b>16</b>
7.1. Crucible .....	16
7.2. Reuse of casting buttons .....	16
7.3. Cast-on .....	16
7.4. Casting equipment .....	16
<b>8. High-frequency melting with vacuum pressure casting units</b> .....	<b>16</b>
8.1. Premelting .....	16
8.2. Melting .....	16
<b>9. High-frequency melting with centrifugal casting units</b> .....	<b>17</b>
<b>10. Melting with torch</b> .....	<b>17</b>
10.1. Observing the melt during torch melting .....	17
<b>11. Overheating</b> .....	<b>17</b>
<b>12. Cooling of the casting ring</b> .....	<b>18</b>
<b>13. Devesting and sandblasting</b> .....	<b>18</b>
<b>14. Finishing</b> .....	<b>18</b>
<b>15. Finishing and polishing steps</b> .....	<b>19</b>
<b>16. Preparation of surfaces for ceramic veneering</b> .....	<b>19</b>
<b>17. Oxide firing</b> .....	<b>20</b>
<b>18. Ceramic veneering</b> .....	<b>20</b>
<b>19. Soldering</b> .....	<b>20</b>
19.1. Soldering after firing .....	21
<b>20. Laser welding</b> .....	<b>21</b>
<b>21. Cleaning</b> .....	<b>21</b>
<b>22. Contraindications and adverse reactions</b> .....	<b>21</b>

## 1. General instructions for processing

These Processing Instructions deal with the main steps when processing non-precious remanium® crown and bridge alloys.


Individual additional information regarding our non-precious remanium® crown and bridge alloys which may differ from these instructions can be found in the information supplied in each individual alloy package or in the internet under [www.dentaurum.de](http://www.dentaurum.de).

Our dental technical team is available to answer any questions that may arise. Dental Technical Hotline Tel No: + 49 72 31 / 803-410.

Please refer to the information at the end of these instructions for information on contraindications and adverse reactions.

## 2. Wax-up

The wax pattern is fabricated in the usual manner based on standard dental technical rules. Use only waxes or modeling materials that burn without residue, e.g. StarWax by Dentaurum.

Before waxing-up, the dies should be coated with a die spacer. The shape of the metal copings should correspond in a reduced scale to the final shape of the finished crowns. Missing portions of teeth should be compensated as in fig. ① to ⑧.

In order to ensure that the metal fills the investment mould completely, the copings must be waxed to a minimum thickness of 0.4 mm.

A uniform thickness of ceramic ensures stress-free adhesion to metal, see fig. ① to ⑧.


Apply ceramic to the framework in an equal thickness. The metal on the lingual side should finish below the incisal area of the crown in order to prevent chipping. Areas to be soldered should be waxed up spaciouly, see ⑤, ⑦ and ⑧.

Rounding the metal / ceramic interfaces, see ① to ⑦, will provide neat junctures and aesthetic shading.

Bridge units must be designed with either a collar or at least inter-dental reinforcement (thermal dissipation).


The outstanding properties of non-precious remanium® crown and bridge alloys allow very thin connections to be used between the crowns and bridge units. The teeth can later be shaped more individually and a blocked appearance can be avoided.

### 3. Spruing

#### For bridges: Wax pattern with sprue bars

Casting button sprue	ø 3.5 – 4 mm
Runner bar	ø 4 – 5 mm
Connection to the crowns	ø 2.5 – 3 mm
	Length 2.5 – 3 mm

#### Use the direct method for single crowns and small bridges

Sprue	ø 3 mm
	Length 15 – 20 mm
Voluminous crowns	ø 3.5 – 4 mm
	Length 15 – 20 mm

### 4. Calculation of metal quantity

General rule: Wax weight of pattern including sprues multiplied by the density of the alloy (8.2 – 8.6 g / cm<sup>3</sup>) plus 1 to 2 cylinder (approx. 10 g) for the casting button.

### 5. Investing

Use special investment material for non-precious alloys e.g. rema<sup>®</sup> CC, Castorit<sup>®</sup> all speed, Castorit<sup>®</sup>-super C, Trivest or rema<sup>®</sup> TT (Dentaurum) with a high overall expansion.

A casting ring system using corresponding ring liners is recommended in order to have exact control over the final fit. "Kera-Vlies<sup>®</sup>" in 1 or 2 mm thickness is especially recommended.

Silicone rings can also be used for investing. Please refer to the recommendations made for the investment material in question.

### 6. Wax burn-out and preheating

For specific information concerning conventional and speed preheating of the investment material please see the relevant manufacturer's Instructions for use.

In general the final temperature for all remanium<sup>®</sup> alloys should be 950 °C / 1742 °F.

This standard temperature may, however, vary according to the type of casting system used. The final temperature should be held for between 30 and 60 minutes according to the size of the casting ring.

## 7. Melting and casting

### 7.1. Crucible

Use a ceramic-based crucible (magnesium oxide, silicon oxide or aluminum oxide) that is suitable for non-precious alloys.

Do not use graphite crucibles.

Preheat empty crucible in the furnace. Use one crucible per alloy and clean after each casting.

Do not use special flux powder!

For larger quantities of metal, a preliminary melting phase (casting delay phase) is recommended.

### 7.2. Reuse of casting buttons

When the alloy is melted, some of the elements which later form the adhesive oxide are lost. If the alloy (e. g. in the form of a casting button) is re-melted several times, these elements may be diminished until they are completely destroyed. If the metal framework is later to be coated with ceramic, then only new material should be used for casting. This is the only way to ensure clear batch traceability.

### 7.3. Cast-on

Metals used for cast-on and high platinum content attachments cannot be used in combination with non-precious remanium® crown and bridge alloys.

### 7.4. Casting equipment

Melting and casting techniques suitable for casting non-precious remanium® crown and bridge alloys include high-frequency melting with vacuum pressure or centrifugal casting units, arc melting with vacuum pressure casting, or torch melting with propane/oxygen or acetylene/oxygen or other gas mixtures capable of reaching the casting temperature.

## 8. High-frequency melting with vacuum pressure casting units

Set vacuum to 250 – 450 mbar (recommended pressure: 450 mbar).

### 8.1. Premelting

Premelt the ingots until they are red hot (no visible rounding at the corners of the ingots).

### 8.2. Melting

Melt until the oxide layer of the molten metal breaks open, then cast immediately.

Observe the melting process and initiate the casting process manually. Oxides formed during melting (see melting illustrations in the operating instructions for the pressure casting machine) remain as slag in the crucible after casting and can be easily removed. If the correct melting and casting technique is used, there is only a minimum residue in the crucible and the molten metal is not overheated.

Automatic casting with a pyrometer control is unsuitable, as the pyrometer cannot always accurately determine the right time for casting.

## 9. High-frequency melting with centrifugal casting units


Wait until the last ingot has collapsed and the oxide layer of the molten metal has broken, then cast immediately.

We recommend premelting the metal if a large amount is to be cast.

## 10. Melting with torch

When melting, the remanium® crown and bridge alloys show a clean melt without sparks.

**Caution! Use torch with shower head nozzle!**

### a) Propane / oxygen melting

Select maximum torch setting as stated by manufacturer.  
Oxygen 2 – 3 bar.

### b) Acetylene / oxygen melting

Follow the instructions of the torch manufacturer.  
Acetylene: 0.7 bar / oxygen 3 bar.

The blue core of the flame at the torch head should be 4 – 5 mm in length. Preheat the empty crucible. Melt the metal in the crucible with a loose circular movement of the flame. The distance at which the torch nozzle should be held depends upon the type of nozzle and the gas mixture used.

### 10.1. Observing the melt during torch melting

When melting with a torch, an oxide layer forms when the last piece of metal has collapsed. With a loose circular movement of the flame, continue the melting process until the metal under the oxide skin is seen to move under the pressure of the flame. Begin the casting process **before** the oxide layer cracks open.

## 11. Overheating

**As with all melting and casting techniques, overheating the metal must be avoided.**

Overheating results in blowholes, micro-porosities and grain coarsening. These casting imperfections are frequently reasons for cracks in the veneering ceramic.

## 12. Cooling of the casting ring

Remove the casting ring carefully and allow to cool.

**Caution! Danger of fire. Always place hot casting rings onto heat-resistant surfaces.**

Ensure that a safe distance from flammable objects is observed.

## 13. Devesting and sandblasting

### **Avoid dust inhalation!**

Please observe the recommendations made on safety precautions in the Instructions for use enclosed with the product. Moisten the casting ring before devesting in order to reduce the amount of dust produced.

Do not strike the casting button with a hammer.

Sandblast using pure aluminum oxide, grain size approx. 125 µm.

## 14. Finishing

### **General information:**

**Please observe the recommendations made on safety precautions to be taken when processing the alloy mechanically. These can be found in the Instructions for use enclosed with the product.**

The higher physical characteristics of the non-precious alloys require grinding and polishing tools different to those normally used for processing precious metal alloys.

The grinding tools selected should only be used for one alloy.

Observe the minimal coping thickness of 0.2 – 0.3 mm.

Due to the high physical properties of the non-precious alloys, both the wall thickness of the crown and the thickness of the bridge connecting points can be reduced by approx. 30 % in comparison to average precious metal alloys.

Always grind in one direction without applying too much pressure. Avoid overlapping of the metal. Remove irregularities such as porosities and inclusions.

## 15. Finishing and polishing steps

Steps	Tools	REF	Remarks
<b>Removal of casting sprues</b>	Supercut STM separating disc	<b>130-111-00</b>	High cutting efficiency and long life span.
	Supercut STM separating disc	<b>130-113-00</b>	High cutting efficiency and long life span.
<b>Grinding of sprue attachments</b>	Grinding wheel	<b>131-322-50</b>	Highly suitable for grinding the stump of the sprue and for coarse grinding.
<b>Finishing</b>	Supercut STM separating disc	<b>130-112-00</b>	Rough finishing, e. g. of interdental spaces.
	Hard metal bur, with cross cut	<b>123-582-00</b>	General finishing with coarse tungsten carbide burs. Fine grinding can be done with fine tungsten carbide burs.
		<b>123-584-00</b>	
<b>123-585-00</b>			
<b>123-601-00</b>			
Aloxine-Point	<b>135-852-00</b> <b>135-853-00</b>	Pure aluminum oxide abrasives for smooth, even transitions, especially metal occlusal surfaces.	
<b>Buffing</b>	Silichrom polisher	<b>138-640-00</b> <b>138-645-00</b>	Abrasive performance, rough.
	Grey rubber polishers	<b>138-102-00</b> <b>138-302-00</b>	Abrasive performance, medium.
	Green rubber polishers	<b>138-101-00</b> <b>138-301-00</b>	Abrasive performance, fine.
<b>Polishing</b>	Polishing brushes	<b>141-800-00</b>	Universal polish with hand tool.
	Tiger Brillant polishing paste	<b>190-350-00</b>	Fast acting high shine polishing paste.
	Tiger Starshine Universal finishing paste	<b>190-301-00</b>	The ideal paste for final finishing.
<b>Polishing of crown's inner surfaces</b>	AL2 O3 -50µm + polishing beads	<b>128-017-00</b> <b>128-211-00</b>	For use with microblaster. Caution! Cover ceramic edges with wax.

## 16. Preparation of surfaces for ceramic veneering

Grind the metal surfaces with hard metal burs to ensure a smooth transition. Blast the surfaces with microblaster polishing beads using pure aluminum oxide and medium grain size (125 µm) at a low pressure of 2 – 3 bar. Clean with ultrasonic cleaner in distilled water.

## 17. Oxide firing

Oxide firing is not required with remanium® bonding alloys. It is advisable to fire an oxide bake in order to visually check the framework condition (5 minutes without vacuum at the firing temperature of the opaque, unless otherwise indicated by the ceramic manufacturer). Following the bake, the oxide should be removed by blasting with single-use aluminum oxide abrasive, grain size 125 µm and low pressure of 2 – 3 bar. Afterwards the surface must be cleaned again.

## 18. Ceramic veneering

remanium® ceramic bonding alloys can be veneered using conventional ceramics such as ceraMotion® Me or ceraMotion® LFC (Dentaurum) as long as the CTE is suitable for the alloy (see alloy Instructions for use).

The coefficient of thermal expansion (CTE) of many ceramic varieties changes, depending on the amount of times the ceramic is fired and the length of time within the firing furnace.

The longer and the more often the ceramic is fired, the higher the CTE. The metal alloy CTE remains constant; the adhesion takes place in the time required for the temperature to cool down from the firing temperature to the standby temperature.

The opaque material can be applied to the sandblasted and cleaned metal framework according to the ceramic manufacturer's recommendations.

It is important to dry the opaque thoroughly. Always ensure the ceramic furnace is calibrated.

The ceramic is then further fired according to the Instructions for use.

If no other information is given by the ceramic manufacturer, the object is cooled as recommended in the alloy Instructions for use.

**Important: After each firing, brush the framework under running water and allow to dry.**

## 19. Soldering

Try to avoid soldered joints if possible.

If however, it is necessary to solder, please ensure a suitable solder with the correct composition and melting temperature which is appropriate for the alloy is used.

Coat the solder objects with an adequate amount of flux!

**Caution! Avoid letting the flux dry out!**

Heat the framework until it is red hot. Position the parts that have been dipped in the flux and heat to flowing.

The following are suitable solders and flux for remanium® alloys:

Solder	CoCrMo-Sold 1	approx. 3 g	<b>REF 102-306-00</b>
Flux	rema® Flux 1	25 g	<b>REF 102-304-00</b>

**Caution: Soldered objects should cool slowly after the ceramic has been fired.**

### 19.1. Soldering after firing

Soldering after firing is not recommended as it reduces corrosive resistance and the ability of precious metal solders to fuse with non-precious alloys is low. Recommended alternatives are joining techniques such as laser welding, TIG welding and adhesive joining.

## 20. Laser welding

The laser welding technique produces a solder-free, mechanically strong and corrosion-resistant joint.

It is important to observe the framework design, the surface structures, the welding sequence and the welding parameters of each individual laser machine. Suitable welding wires are available as filler materials for all remanium® alloys.

CoCr welding wire:

ø 0.25 mm	Roll 2 m long	<b>REF 528-215-10</b>
ø 0.35 mm	Roll 2 m long	<b>REF 528-210-10</b>
ø 0.5 mm	Roll 2 m long	<b>REF 528-200-10</b>

## 21. Cleaning

Crowns and bridges made in remanium® crown and bridge alloys can be cleaned, after finishing, with the ultrasonic solution P. Refer to the cleaning recommendations on the labels.

Recommended ultrasonic solution after polishing:

Ultrasonic cleaning solution	1000ml	<b>REF 094-610-00</b>
------------------------------	--------	-----------------------

## 22. Contraindications and adverse reactions

Signs of intolerance to non-precious remanium® crown and bridge alloys are extremely rare if the manufacturer's Instructions for use are adhered to.

If the patient has a proven allergy against any component within the alloy, this alloy must not be used for safety reasons.

There have been individual reportings of local irritations which were electrochemically induced.

If various alloys have been used, it is possible that galvanic effects may occur.

Please observe the information contained in the enclosed Instructions for use.

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