Gebrauchsanweisung | Instructions for use |
Mode d'emploi | Modo de empleo | Modalità d'uso

CE 0483

remanium® 

Edelmetallfreie Legierungen für die Modellgusstechnik
Non-precious metal alloys for the cast partial denture technique
Alliages exempts de métaux précieux pour la coulée des châssis métalliques
Aleaciones exentas de metales preciosos para técnica de esqueléticos
Leghe non preziose per scheletrati

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Dear Customer,

Thank you for having chosen a quality product from Dentaureum.

In order to use this product at its best for you and your patients, it is important to study and follow these directions for use carefully.

The written instructions cannot cover all eventualities during operation. For this reason our Hotline is available to answer any other questions and ideas that may arise.

Due to a constant development we recommend, even when you use the same products frequently, to study the enclosed latest issue of the instructions for use or refer to our website at www.dentaureum.com.

Manufacturer

Dentaureum GmbH & Co. KG | Turnstr. 31 | 75228 Ispringen | Germany

With reference to our quality standards

Dentaureum ensures a faultless quality of the products manufactured by us. These recommendations are based upon our own experiences. The user himself is responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence on the processing on site.

Contents

1. General instructions for use	14
2. Waxing-up	14
3. Spruing	14-15
4. Investing	14
5. Casting ring systems	16
6. Wax burn-out and preheating	16
7. Melting and casting	16
7.1. Melting crucible	16
7.2. Metal quantity	16
7.3. Reuse of casting cones	16
7.4. Cast-on	17
7.5. Casting equipment	17
8. Induction melting with vacuum pressure casting machines	17
8.1. Pre-melt	17
8.2. Melt	17
9. Induction melting with centrifugal casting machines	17
10. Melting with torch	18
10.1. Observing the melt during torch casting	18
11. Overheating	18
12. Cooling of the casting ring	18
13. Devesting and sandblasting	19
14. Finishing	19
14.1. Processing and polishing steps	19
14.2. Electrolytic polishing	20
15. Laser welding	20
16. Soldering	20
17. Cleaning	21
18. Contraindications and side effects	21

1. General instructions for use

These are the instructions for use when working with non-precious remanium® cast partial denture alloys.

Any deviations to the norm or specialized additional information about our non-precious remanium® cast partial denture alloys can be found in the specific alloy supplement supplied in each individual alloy package or in the internet under www.dentaurum.com.

Our dental technical team is available to answer any questions that may arise. Dental Technical Hotline Tel No: +49 72 31/803-410.

Information on contraindications and side effects can be found at the end of these instructions for use.

2. Waxing-up

The outstanding physical qualities of the non-precious remanium® cast partial denture alloys enable the framework thickness to be kept to a minimum. A thickness of 0.4 mm is sufficient for full metal dentures. If the transversal connections on partial dentures are slender in design then the main body structure must be stronger. A cross section diameter of 2 x 4 mm is recommended for sublingual bars. Waxes and acrylics which burn out without leaving a residue are essential for high quality casting results.

3. Spruing

Round profiles \varnothing 3.5 – 4.0 mm. Length 20 mm minimum. For voluminous wax ups \varnothing 4.0 mm. For upper jaw 2 – 4 sprues, for lower jaw 2 – 3 sprues.

The sprues should be attached to the pattern in such a way as to determine the flow direction of the molten metal. It is possible to cast from above as well as through the model itself. The outline on page 15 shows recommendations of sprue diameter and position depending on different types of design and construction.

4. Investing

Always use phosphate bonded investment materials such as rema® Exakt, rema® Exakt F or rema® dynamic S with an ideally suited expansion. rema® dynamic S is an investment material recommended for quick casting methods with which it is possible to immediately place the ring into the furnace at final temperature without risking a bad fit.

Spruing

3.5 mm

3.5 mm

3.5 mm

3.5 mm

3.5 mm

3.5 mm

3.5 mm

3.5 mm

3.5 mm

3.5 mm

3.5 mm

3.5 mm

3.5 mm

3.5 mm

4 mm

4 mm

2.5 mm / 4 mm

4 mm

3.5 mm

3.5 mm

4 mm


2.5 mm / 4 mm


2.5 mm / 4 mm

2.5 mm / 4 mm

3.5 mm

3.5 mm

Length min. = 20 mm

Length max. = 40 mm

5. Casting ring systems

It is recommended to use reusable casting rings when investing: such as rema® Form quick investment system or Neo-Star investment system. The model base must always be surrounded by investment material. This prevents the model from being separated from the ring during the casting process.

6. Wax burn-out and preheating

No matter which investment material is used it is always imperative to follow the relevant manufacturer's instructions for use. The ring temperature of investment should be kept as low as possible. The temperatures range between 850 °C/1562 °F and 1050 °C/1922 °F according to the type of investment, casting machine and pattern design used.

A lower temperature provides a more homogenous metal structure and reduces oxidation. When deciding on a preheating temperature it is also important to take into account the pull-off force of the casting machine, the thickness of the wax pattern and the investment material's gas permeability.

7. Melting and casting

7.1. Melting crucible

Use a ceramic based crucible (i.e. magnesium oxide, silicon oxide or aluminium oxide).

Do not use carbon crucibles.

Preheat empty melting crucible in the furnace. Use one crucible per alloy and clean after each casting.

Do not use special flux powder!

For larger quantities of metal, a preliminary melting phase (casting delay phase) is recommended.

7.2. Metal quantity

In general, one metal casting requires 4 – 5 casting ingots at 6 g each.

7.3. Reuse of casting cones

When the casting cone is to be re-used one more time, at least $\frac{2}{3}$ of new material must be added. If the metal cone from the second casting is then to be used again the metal structure may change to such a degree that it no longer has any resemblance to the alloy's original constituents and physical data.

7.4. Cast-on

Metals used for cast-on and high platinum content attachments cannot be used in combination with non-precious remanium® cast partial denture alloys.

7.5. Casting equipment

Melting and casting techniques suitable for casting non-precious remanium® cast partial denture alloys are those which are capable of reaching the required casting temperatures. These include induction casting machines with vacuum-pressure or centrifuge, arc melting with vacuum pressure casting, or torch melting with Propane/Oxygen or Acetylene/Oxygen or other gas mixtures strong enough to reach the casting temperature.

8. Induction melting with vacuum pressure casting machines

Set vacuum to 250 – 450 mbar (recommended pressure: 450 mbar).

8.1. Premelting

Premelt the ingots until they are red hot (no visible rounding at the corners of the ingots).

8.2. Melting

Depending on the type of alloy, melt the metal until the skin which forms over the molten metal splits **then cast immediately** or wait for 1 – 2 seconds after the oxide skin layer has split open and then cast (see the relevant alloy information sheet).

Observe the melting procedure and initiate the casting procedure manually. Oxides formed during melting (cf. melting illustrations in the operating instructions for the pressure casting machine) remain as slag in the crucible after casting and are easily removed. If the correct melting and casting technique is used, there is only a minimum residue in the crucible and the molten metal is not overheated.

Automatic casting with a pyrometer control is unsuitable, as the pyrometer cannot always accurately determine the right time for casting.


9. Induction melting with centrifugal casting machines

After the last metal ingot has melted, wait until the oxide skin layer splits open (shiny molten state) **then cast immediately** or wait for 1 – 2 seconds after the oxide skin layer has split open and then cast – (see the relevant alloy information sheet).

We recommend premelting the metal if a large amount is to be cast.

10. Melting with torch

When melting, the remanium® cast partial denture alloys show a clean molten mass without sparks. **Caution! Use torch with shower head nozzle!**

a) Propane / Oxygen Melting

Select maximum burner setting as stated by manufacturer. Oxygen 2 – 3 bar, propane 1 bar.

b) Acetylene / Oxygen Melting

Follow instructions of burner manufacturer. Acetylene: 0.7 bar, oxygen 3 bar.

The blue core of the flame at the burner head should be 4 – 5 mm in length. Pre-heat the empty crucible. Melt the metal in the crucible with a loose circular movement of the flame. The distance at which the torch nozzle should be held depends upon the type of nozzle and the gas mixture used.

10.1. Observing the melt during torch casting

When melting with a torch, an oxide layer forms when the last metal piece has collapsed. With a loose circular motion of the flame, continue the melting process until the metal under the oxide skin is seen to move under the pressure of the flame. Begin the casting process **before** the oxide skin cracks open.

11. Overheating

As with all melting and casting techniques, **overheating the metal must be avoided**. Overheating results in contraction holes, micro porosities and grain coarsening. These casting imperfections are frequently reasons for cracking of a cast partial.

12. Cooling of the casting ring

Remove the casting ring carefully and allow to bench cool. In order to achieve the best mechanical qualities and a stress free material allow the metal to cool down slowly to room temperature. After at least 30 minutes of bench cooling, the casting ring may if necessary be carefully quenched in water.

Caution! Danger of fire! Always place hot casting rings onto heat resistant surfaces.

Ensure of a safe distance from flammable objects.

13. Devesting and sandblasting

Avoid dust inhalation!

Moisten the casting ring before devesting in order to reduce the amount of dust produced!

Do not strike the casting button with a hammer!

Sandblast using pure aluminium oxide, grain size approx. 125 µm to 250 µm.

Wear a protection mask!

14. Finishing

14.1. Processing and polishing steps

Steps	Tools	REF	Remarks
Removal of casting sprues	Supercut ST separating disc	130-102-00	High cutting power and durability for high speed grinders.
	Supercut STM cutter disc	130-111-00	High cutting power and durability for handpieces.
Grinding of sprue attachments	Grinding wheels	131-322-50	Ideal for sprue stumps and pre-grinding.
Finishing – rough	crosscut carbide burs	123-581-00 123-583-00 123-584-00 123-585-00	General use in inaccessible areas.
Finishing – medium	Diasint burs	137-826-00 137-827-00 137-831-00	Sintered diamond burs. No splintering due to metallic bond.
Finishing – fine	Grinding points – pink	132-802-00 132-803-00	Ceramic bonded corundum grinding points.
Electrolytic polishing	Electropolishing solution	128-301-00	Remember to protect and cover attachments.
Buffing	Silichrom polisher, cylinders	138-640-00	Rough effect.
	Silichrom polisher, discs	138-645-00	
	Rubber polisher, green, discs	138-101-00	Fine effect.
	Rubber polisher, green, cylinders	138-301-00	
Polishing	Polishing brushes	141-712-00	Triple row brushes for polishing lathes.
	Polishing paste Tiger brillant	190-350-00	Quick, high polishing paste.
	Tiger Starshine Universal Finishing Paste	190-301-00	High polishing Finishing Paste.

General information:

Caution: Always use suction units when grinding, cutting and polishing!
Wear protective glasses and a protection mask!

Dentaurum's remanium® cast partial denture alloys can be finished easily and safely with special finishing and polishing instruments. Suitable instruments and materials are available in the Dentaurum range of equipment.

Successful fitting and finishing of Dentaurum's remanium® cast partial denture alloys depend greatly upon the accuracy of the wax pattern and casting.

Use the specific instruments for only one alloy.

14.2. Electrolytic polishing

The duration of the polishing times depend upon the type of polishing unit and electrolyte used.

General duration: Polishing time without movement in the solution approx. 2 x 5 minutes (after half time change the position of the object in the solution) at 4 Amps, with movement approx. 10 minutes at 4 Amps.

Polishing baths should be warmed in order to increase their polishing ability 30 °C/86 °F to 40 °C/104 °F. Hot baths of > 50 °C/122 °F) become damaged and diminish the polishing effect.

15. Laser welding

All remanium® cast partial denture alloys contain next to no carbon and are therefore suitable for laser welding. remanium® GM 900, is entirely free of carbon and consequently especially suited.

With the laser welding technique it is possible to achieve solder free, corrosion resistant and mechanically strong joints. It is important to observe the geometry of the object, the surfaces, the welding sequence and the individual welding parameters recommended for the laser welding machine. Suitable laser welding wires in the same carbon free material are available for all remanium® metal alloys.

CoCr laser welding wire:

ø 0,25 mm	2 m coil	REF 528-215-10
ø 0,35 mm	2 m coil	REF 528-210-10
ø 0,5 mm	2 m coil	REF 528-200-10

16. Soldering

Try to avoid soldering if possible in order to prevent a material mixture. If however, it is necessary to solder, please ensure a suitable solder with the correct composition and melting temperature which is appropriate for the alloy used.

Paint the area to be soldered with sufficient flux!

Caution! Do not allow the flux to dry!

Heat the framework until it is red hot. Set the piece of solder covered in flux in place and heat it directly until it melts and flows!

Suitable solders and flux for remanium® alloys:

Solder	CoCrMo-Sold 1	approx. 3 g	REF 102-306-00
Flux	rema® Flux 1	25 g	REF 102-304-00

Alternative:

Cobalt-Chrome Solder sticks incl. flux	6 rods	REF 142-100-50
--	--------	-----------------------

17. Cleaning

remanium® cast partial dentures can be polished using conventional polishing materials suitable for prosthetics. Observe the polishing directions on the label. It is recommended not to use household polishing materials as these could cause the alloy to discolour.

Recommended ultrasonic cleaning after polishing:

Ultrasonic cleaning liquid	1000 ml	REF 094-610-00
----------------------------	---------	-----------------------

18. Contraindications and side effects

Signs of intolerance to non-precious remanium® cast partial denture alloys after correct processing are rare.

If the patient has a proven allergy against any component within the alloy, for safety reasons this alloy must not be used.

Some patient cases have shown traces of local irritation caused by electrochemical reactions.

If various alloys have been used, it is possible that galvanic effects may occur.

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