Orthocryl® LC

CE 0483

- Gebrauchsanweisung · **Lichthärtender Kunststoff**
Instructions for use · **Light-curing acrylic**
Mode d'emploi · **Résine photopolymérisable**
Modo de uso · **Acrílico fotopolimerizable**
Modalità d'uso · **Resina fotopolimerizzabile**

D
DENTAURUM
1886

1. Hersteller

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Notizen

This image shows a full page of primary-ruled paper. It features multiple sets of horizontal dashed lines spaced evenly down the page, providing a guide for handwriting practice. The paper is otherwise blank, with no margins or additional markings.

1. Manufacturer

Dentaurum GmbH & Co. KG
Turnstraße 31 · 75228 Ispringen · Germany

2. Product description

Orthocryl® LC is a light-curing acrylic that polymerizes under light in a specific wavelength range. Orthocryl® LC does not contain methyl methacrylate or dibenzoyl peroxide!

3. Intended purpose

Orthodontic acrylics are used for the fabrication of orthodontic appliances.

4. Intended users

All products are intended for dental professionals only.

5. Indications

Misalignment of teeth and jaws that are to be treated with orthodontic appliances.

6. Contraindications

For users / patients who are allergic to the material components in light-curing acrylics.
See 13. Composition.

7. Patient target group

Appliances fabricated by the specialist user are used for the treatment of children, adolescents and adults.

8. Processing steps

To avoid direct contact with Orthocryl® LC, it is recommended to wear face protection / safety glasses (REF 162-920-00) as well as special protective gloves i.e., Dermatril® disposable protective gloves made of special nitrile (REF 162-907-00, 162-908-00 and 162-909-00), also refer to the safety data sheet.


- Processing of Orthocryl® LC requires an injector, the cartridge with the polymer and the cartridge tip (Fig. 1).
- Remove the small blue protective cap on the Orthocryl® LC cartridge by turning it anti-clockwise (Fig. 2).
- Remove the large blue protective cap on the Orthocryl® LC cartridge by pushing up the protruding flap (Fig. 3).
- If necessary, shorten the interchangeable cartridge tip (REF 160-421-00) to the length required and rotate onto the opening of the cartridge as far as it will go (Fig. 4).
- The button on the upper side of the injector releases the forward and back feed of the injector plunger. Hold the button down and bring the plunger into the initial position.
- Put the cartridge onto the injector and lock with a 90° turn (Fig. 5).
- Push the plunger forward whilst holding the button down until resistance is felt.
- Apply the emerging Orthocryl® LC onto wires, screws and model with care.

8.1 Separate stone model/plastic model

Both stone and plastic models **must** be given a layer of separating medium before Orthocryl® LC is applied.

When working with stone models, we recommend using a special separating medium (REF 162-800-00) in order to attain a smooth, silky basal surface on the appliance and to avoid discolorations and plaster inclusions. It is not necessary to ventilate (soak) the model beforehand.

Plastic models **must** be given a layer of separating medium before Orthocryl® LC is applied in order to prevent the plastics from sticking. We recommend the 3D separating medium (REF 162-801-00), not available in the USA.

8.1.1 Application of the separating medium (REF 162-800-00) on stone models

- Apply the separating medium from the dosing bottle (REF 162-100-00) directly on the model and distribute evenly using a brush.
- Stand the model on its back (tuber plane) so that excess separating medium can drain away.
- After approx. 30 seconds, carefully blow away separating medium residues or droplets with compressed air and, if necessary, remove residues from the screw and wire retentions with a dry brush.
- Processing of the acrylic can start immediately afterwards!

8.1.2 Application of the 3D separating medium (REF 162-801-00) on plastic models

For the fabrication of braces, we recommend high-quality full plastic models produced using the DLP process (Digital Light Processing), e.g., prime4me® base model (REF 430-110-10) or the FFF process (Fused Filament Fabrication) with FILAPRO® (440-000-00) and the filaments FILAPRO® model gray (440-200-00), model white (440-201-00) and model cream (440-202-00). Besides being highly accurate and therefore ensuring very precise fitting, these models allow for the usual process in the laboratory.

Due to the large variety of materials for model printing, it is essential to test the suitability of the 3D separating medium before use.

- Apply the 3D separating medium evenly with a separating brush over the entire model before the expansion screw(s) and wire elements are finally fixed in place and allow the separating medium to dry thoroughly. To this end, apply the separating medium from the dosing bottle (REF 162-100-00) directly on the model and apply evenly using a brush.
- Stand the model on its back (tuber plane) so that excess separating medium can drain away if necessary.
- Allow the separating layer to dry thoroughly. The separating medium is completely dry when it no longer feels sticky.
- Repeat the process in order to achieve an even separating layer everywhere.
- In order to accelerate drying of the separating medium, careful blowing with compressed air or a household hair dryer can be helpful. For faster drying of the separating medium, the model can also be placed in a light-curing unit for a few minutes.

*The time required in the light-curing unit may vary depending on the size of the workpiece and the amount of separating medium. The drying time also depends on the type of device – with the Solidilite LED unit from Shofu, up to 3 minutes should be sufficient.

- If required, block out undercuts in the tooth/jaw area and fix the expansion screw(s) and wire elements as usual on the dried isolated plastic model.
- Processing of the acrylic can start immediately afterwards!
- For plastic models with the FILAPRO®, wait until the model has cooled to room temperature before removing the appliance.

9. Processing Orthocryl® LC

Example: expansion plate

Orthocryl® LC is generally applied to the model that is dry and isolated but not hydrated. Mark the position of the expansion screws on the model and drill a hole or a groove with a conical bur (e.g., REF 123-598-00) for the screw holder. Fill the hole on the model with heated wax (e.g., REF 119-150-00) and insert the expansion screw with the lower, basal part of the holder. Remove the expansion screw after the wax has hardened.

Fix bent wire elements to the model with sticky wax (REF 122-300-00).

Isolate the model as described under Section 8.1.

Processing of the acrylic can start immediately afterwards!

Completely surround the removed expansion screw with Orthocryl® LC (Fig. 6). Please observe that the expansion screw is completely embedded in acrylic without any air bubbles. Then reinsert the expansion screw into the wax-lined hole. The wax in the hole ensures that the expansion screw settles in the same position as determined before. With wire retentions, apply Orthocryl® LC to the space between wire and stone model first and then embed the complete wire (Fig. 7).

Fig. 6

Fig. 7


Fig. 8



Fig. 9

Fig. 10

Next, apply Orthocryl® LC until the desired form and thickness of the appliance (maxilla approx. 3 mm / mandible approx. 4 mm – 5 mm) is achieved (Fig. 8). Should this process require more than a single step, dip the tip of the cartridge slightly into the material already applied before each new step. This prevents air inclusions between the layers.

Eliminating air bubbles

If air bubbles form during application, these can be opened with a pointed instrument, e.g., a probe, and filled. Alternatively, one can also inject additional Orthocryl® LC from underneath the air bubble so that it moves towards the surface and disintegrates.

Orthocryl® Disco glitter

Spray a small amount of Orthocryl® Disco glitter on the first layer of Orthocryl® LC directly from the spray bottle. Then, cover it with another layer of Orthocryl® LC.

Application of decals

The decals – for visually enhancing appliances – work particularly well with highly transparent acrylics. The decals and the acrylic do not bond together chemically so it is important to embed them completely in the acrylic! Avoid abrading the decals during processing. Place the cut decal on the first Orthocryl® LC layer, at about the height of the wire retentions. Then, cover it with another layer of Orthocryl® LC.

9.1. Processing of the Orthocryl® LC bite paste

Example: activator

After the maxillary and mandibular models have been placed in the fixator with the aid of the construction bite*, separate the two halves of the fixator again and prepare the models as described in Section 9.

Orthocryl® LC is then applied to the palatal side of the maxilla and to the lingual side of the mandible. The Orthocryl® LC bite paste is then taken from the jar with a spatula or a wax knife and applied onto the occlusal or the incisal part of the mandibular model. To ensure that the material enters all areas of the structure (bite ramp) and to avoid air bubbles, make sure to apply an excess amount of bite paste.

After applying the bite ramps, reassemble the two halves of the fixator, remove any excess material, and place the closed fixator in the light-curing device.

*When placed in a space-saving manner (close to the round telescopic bars) into the fixator 072-004-00, average-sized models fit comfortably into the Solidilite LED device from Shofu.

10. Polymerization

The specifications on polymerization times refer to our recommended light-curing unit Solidilite LED from SHOFU Dental GmbH. Please observe the instructions for use for this unit.

Polymerization of Orthocryl® LC in the light-curing unit takes place at a wave length of 400 nm – 550 nm*. The time required in the light-curing unit may vary depending on the size of the workpiece and the amount of separating medium. The type of device also plays a role. The polymerization time is 3 minutes for single jaws and 6 minutes for bimaxillary devices. Then carefully remove the appliance from the model (Fig. 10) and expose again for 1 minute or 3 minutes – with the basal side facing upwards. In the case of bimaxillary units, first cure the maxilla and then the mandible. The two pieces are connected using bite paste and renewed curing.

*If a halogen lamp is to be used (e.g., Solidilite V), use a light output of at least 4x 150 watts = 600 watts. Take into account the age of the lamp and the Instructions for use. We recommend switching to curing devices with LED lamps.

Caution!

The oxygen inhibition layer (smear layer) present on the fully polymerized appliance must be removed, see Section 11.

11. Removing the inhibition layer

Completely remove the oxygen inhibition layer (smear layer). To do this, apply Orthocryl® LC cleaning liquid (REF 160-430-00) on a cloth and remove the smear layer.

Caution!

To exclude any risk to patients, it is essential to ensure that the oxygen inhibition layer (smear layer) is removed completely (see safety data sheet). Direct application of Orthocryl® LC and the bite paste in the mouth is not permitted!

12. Finishing and polishing

When finishing Orthocryl® LC appliances, proceed as with Orthocryl® cold-curing polymer. However, do not use cross-cut hard metal burs. Hard metal burs with transversal toothing or burs for soft acrylics are better suited. Employ extraction during finishing.

Pre-polish the appliance with silicone polishers or sand paper, and then with pumice powder and polishing brushes. A high gloss finish is achieved using Edelweiß polishing material (REF 190-100-00 or REF 190-200-00) and a linen polishing pad. Wear the prescribed protective equipment (face protection / eye protection).

Note:

Due to stabilizing fillers, Orthocryl® LC appliances exhibit a slight milky haze.

13. Composition

Mixture of various monomers and oligomers. Contains fillers and dyes. Refer also to the safety data sheets which are delivered automatically to customers with every first order and which are available at any time on request!

14. Information on dangers and safety

Causes irritation to the skin. Can cause allergic reactions. Causes serious irritation to the eyes. Avoid breathing dust/fume/gas/mist/vapors/spray. Do not wear contaminated work clothes outside the workplace. Wear protective gloves/protective clothing/eye protection/face protection/hearing protection. If in contact with skin: Wash with plenty of water. If skin irritation or rash occurs: Get medical advice/attention. If in contact with eyes: Rinse carefully with water for several minutes. Remove any contact lenses if possible. Continue rinsing. In case of persistent eye irritation: Get medical advice/attention. If you feel unwell, call the POISON INFORMATION CENTER or a doctor.

For curing Orthocryl® LC, use light-curing units that emit light at a wavelength of 400 nm – 550 nm.

In order to avoid direct skin contact with Orthocryl® LC during processing, we recommend wearing protective gloves (e.g., Dermatril® disposable gloves made of nitrile, REF 162-907-00, REF 162-908-00, REF 162-909-00), also see safety data sheet.

Direct application of Orthocryl® LC in the mouth is not permitted!

The product should not be used if there is a known allergic reaction to one or more of the material components. Should the dental professional and/or the patient become aware of serious problems arising from the use of the product, inform the manufacturer and the competent authority in the country in which the dental professional and/or the patient is resident. For further information, please refer to the safety data sheet.

15. Storage and shelf life

Since Orthocryl® LC is a light-curing acrylic, close cartridges after use and store away from light.

The shelf life in the closed original package is 48 months from the date of manufacture.

Protect from contamination and sunlight as they may cause premature polymerization.

16. Cleaning and caring for appliances

Appliances should be cleaned regularly and stored in fresh water when not in use. This procedure prevents gap corrosion or brown discoloration of the wire elements.

For further tips on care, see also "Oral hygiene during orthodontic treatment" (REF 989-328-10), which is available for download at www.dentaurum.com.

17. Disposal

When disposing of the product, please adhere to the national regulations and the details outlined in the safety data sheets. Avoid disposing of loose chips and dust particles.

18. Product range

Please refer to the Orthodontics or Dental Prosthetics Catalog for information on the Orthocryl® product range.

19. General information

The SSCP is available at <https://ec.europa.eu/tools/eudamed> and www.dentaurum.com.

For more information on our products, please visit the internet at www.dentaurum.com. You will find these Instructions for use and information on the current status at www.dentaurum.com/ifu.

20. Explanation of symbols used on the label

Observe the label. Additional information can be found on the internet at www.dentaurum.com/ifu
(Explanation of symbols used on the label REF 989-313-00).

Notes

This image shows a full page of white paper with horizontal dashed lines, typical of primary school handwriting practice paper. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

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