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Dentaurum GmbH & Co. KG · Turnstr. 31 · 75228 Ispringen · Germany

2. General device description

Abutment components					
Series of abutments*	S	M		L	
Implant diameters in mm	3.3	3.7	4.2	4.8	5.5

*all of the screw removal sets listed in these Instructions for use are available for the S, M and L series of abutments.

The tioLogic® ST and the tioLogic® TWINFIT screw removal sets correspond to the inner geometry of the respective implant system and serve to safely remove the fragments of a broken screw from the implant inner thread of tioLogic® ST or tioLogic® TWINFIT implants. The implant remains intact and can be prosthetically restored again. The screw removal sets are intended for single use only.

Implant system	Implant diameter	Product name	REF
tioLogic® ST	3.3 mm	tioLogic® ST screw removal set S	387-606-00
	3.7 mm	tioLogic® ST screw removal set M	387-607-00
	4.2 mm		
	4.8 mm	tioLogic® ST screw removal set L	387-608-00
	5.5 mm		
tioLogic® TWINFIT	3.3 mm	tioLogic® TWINFIT screw removal set S	347-606-00
	3.7 mm	tioLogic® TWINFIT screw removal set M	347-607-00
	4.2 mm		
	4.8 mm	tioLogic® TWINFIT screw removal set L	347-608-00
	5.5 mm		

3. Intended purpose

tioLogic® ST and tioLogic® TWINFIT screw removal sets are intended for removal of fragments of broken screws from the implant thread of tioLogic® ST or tioLogic® TWINFIT implants.

4. Intended users

The tioLogic® ST and tioLogic® TWINFIT product ranges are supplied exclusively to qualified personnel such as surgeons, dentists and dental technicians. They should only be used by qualified personnel who are familiar with dental implantology, including diagnosis, preoperative planning, surgical technique and prosthetic treatment.

5. Medical indications / contraindications

5.1 Indications

5.1.1 tioLogic® ST and tioLogic® TWINFIT screw removal sets

- Removal of screw fragments from tioLogic® ST implants with the tioLogic® ST screw removal set
- Removal of screw fragments from tioLogic® TWINFIT implants with the tioLogic® TWINFIT screw removal set
- Use on osseointegrated implants

5.2 Contraindications

5.2.1 tioLogic® ST and tioLogic® TWINFIT screw removal sets

- Use of tioLogic® ST screw removal sets with tioLogic® TWINFIT implants
- Use of tioLogic® TWINFIT screw removal sets with tioLogic® ST implants
- Disassembly of the screw removal set
- Multiple use of the screw removal set
- All indications not listed are contraindicated
- Use on non-osseointegrated implants

6. Patient target group

All skeletally mature patients undergoing implant treatment. Please consider the indications / contraindications.

7. Information on use, processing, fitting and insertion

Prior to use, dental professionals should ensure that they have carefully read and understood the full tioLogic® ST and tioLogic® TWINFIT Instructions for use.

We strongly recommend that, prior to using Dentaurum products, dental professionals attend a Dentaurum training course to familiarize themselves with the Dentaurum products.

Carry out thorough radiological, psychological and physical examinations before implantation and inform patients of the process in detail. Dentaurum recommends full clinical, radiological, photographic and statistical documentation.

Precautions:

Before surgical implantation, it is important to carry out a careful patient anamnesis and, if required, a third-party anamnesis by the patient's doctor (GP/PCP) in order to clarify if 1) implantation may be difficult due to an anatomical situation, 2) there is the risk of a serious surgical problem or a general risk, 3) wound healing/osseointegration may be difficult, or 4) good hygiene and care of the implant, abutment and restoration may worsen. Do not fit dental implants if any of the above cases or a combination hereof is severe or uncontrollable.

Prosthetic restoration of the implant must not be performed until the soft tissue has healed without irritation. Furthermore, a radiographic check-up should be performed before the prosthetic restoration.

Important:

The implantologist should ensure the products cannot be aspirated during intra-oral use.

7.1 tioLogic® ST and tioLogic® TWINFIT screw removal sets


Screw removal set tioLogic® ST implants

Screw removal set tioLogic® TWINFIT implants

All tioLogic® ST and tioLogic® TWINFIT screw removal sets are available for the S, M and L series of abutments. They each include a pre-installed drill. On the one hand this is the centering tool ①, on the other the left-cutting drill ②. The screw removal sets are intended for single use only.

Step 1: Insert the screw removal set with the centering tool (marked ① on the instrument head) into the antirotational mechanism of the implant so that it “locks” into place. The correct position in the antirotational mechanism ensures that the centering tool is guided properly so that it will not damage the inner connection of the implant. Move the centering tool (manually) until it makes contact and turn **anti-clockwise**.

Step 2: Insert the screw removal set with the flat, thin left-cutting drill ② (marked on the instrument head) into the antirotational mechanism of the implant so that it “locks” into place. The correct position in the antirotational mechanism ensures that the drill is guided properly so that it will not damage the inner connection of the implant. Move the left-cutting drill (manually) until it makes contact and turn **anti-clockwise**.

8. Technical specification

Composition

Product	Composition
Screw removal set	1.4305 / 3.7165 / 1.3343

Information about the composition is included in the material list for our implantology products (REF 989-801-05), or see www.dentaurum.com/ifu. The SSCP is available at <https://ec.europa.eu/tools/eudamed> and www.dentaurum.com.

9. Adverse reactions and residual risks

Possible complications and risks during prosthetic treatment (temporary adverse reactions) are:

Pain, swellings, difficulties speaking, inflammation of the gums.

Possible complications and risks during prosthetic treatment (permanent adverse reactions) are:

fractured dentures, aesthetic problems.

10. Safety information and precautions

- Do not use the product if there is a known allergic reaction to one or more of the material components.
- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. Therefore, secure all articles used intra-orally against swallowing and/or aspiration.
- Different types of alloy in the oral cavity can lead to galvanic reactions.
- If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.
- Only use all tioLogic® ST prosthetic components in combination with tioLogic® ST implants! Use of third-party components can negatively impact the function and safety. In such cases, Dentaureum GmbH & Co. KG accepts no liability for replacement products or guarantees.
- Only use all tioLogic® TWINFIT prosthetic components in combination with tioLogic® TWINFIT implants! Use of third-party components can negatively impact the function and safety. In such cases, Dentaureum GmbH & Co. KG accepts no liability for replacement products or guarantees.
- Observe the indications for the individual products. Consider all indications that are not listed explicitly to be contraindications.
- Disinfection should not be done in the original packaging.
- Products must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.
- Technical advice (oral and written) is based on the scientific and technical knowledge available at the time the product is put on the market. It does not release the dental professional from the responsibility of personally checking the suitability of the products for the intended indications and applications. Advice is only given as a non-binding recommendation and cannot be assumed to provide any form of assurance, liability or guarantee.

Severe incidents:

Report all severe incidents arising from the use of the product to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.

Severe incidents are identified as malfunctions and a drop in performance or deterioration of the properties of a medical device which resulted or may have resulted in a direct or indirect temporary or permanent deterioration in a patient's state of health, in the death of a patient, a dental professional or other persons or as posing a serious danger to public health.

11. Storage and shelf life

Store the tioLogic® TWINFIT implants, abutments and components always in a dry, dark place at room temperature in the original packaging. If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.

12. Disposal

Adhere strictly to the rules that apply for your region and the details outlined in the safety data sheets, if available. Products to be disposed of should be treated as surgical waste and be decontaminated in compliance with the relevant regulations for medical waste disposal.

13. Information for single use products

tioLogic® ST and tioLogic® TWINFIT screw removal sets are supplied non-sterile and are intended for single use on a single patient only. The product must be disinfected before first use. The screw removal sets may not be sterilized.

Reconditioning of abutments and components (recycling) as well as reuse on patients are not permitted as safe reconditioning and/or functional safety cannot be guaranteed and may lead to a loss in mechanical, chemical and/or biological properties. Reuse may cause infection. Do not use implants after the expiry date stated on the label.

14. Information on disinfection / cleaning / sterilizing

The tioLogic® ST and tioLogic® TWINFIT screw removal sets are supplied non-sterile. They are for single use on one patient only. They must be disinfected before being used on a patient. They may not be sterilized.

15. Product overview and scope of delivery

Please refer to the tioLogic® ST product catalog for the tioLogic® ST product range.

Please refer to the tioLogic® TWINFIT product catalog for the tioLogic® TWINFIT product range.

16. Quality information

Dentaureum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaureum has no influence on the processing on site.

17. Explanation of symbols

Please refer to the label. See chapter "Additional information".

18. Additional information

Additional information can be found at www.dentaureum.com/ifu

Instructions for use	REF
Preparation instructions for the implant system components	989-801-07
Explanation of symbols used on the label	989-313-00

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