
tomas®-pin SD tomas®-pin EP

DE **Gebrauchsanweisung** · tomas®-pin SD, tomas®-pin EP

EN **Instructions for use** · tomas®-pin SD, tomas®-pin EP

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

Dentaurum GmbH & Co. KG · Turnstr. 31 · 75228 Ispringen · Germany

2. General description

A temporary skeletal anchorage option is created for orthodontic treatment with the aid of an endosteally anchored mini implant (tomas®-pin). The tomas®-pin SD and the tomas®-pin EP have a self-drilling (SD = self-drilling) thread (Ø 1.6 mm) in different lengths (see point 10). The gingival collar has a height of 2.0 mm and a maximum Ø of 2.8 mm. The head of the tomas®-pin SD (Fig. 1) has a cross slot (slot width: 0.56 mm / 22; slot depth: max. 1.15 mm). The tomas®-pin EP (Fig. 2) has a mushroom-shaped head. Different orthodontic appliances can be connected to the head of the tomas®-pin, according to the indication, to achieve or support the required tooth movement.

tomas® is a system of coordinated components for the insertion and coupling of the mini implants. The fracture torque for the tomas®-pin is 35 Ncm.

2.1 Materials used:

tomas®-pin implants	Titanium alloy according to ISO 5832-3
tomas®-abutments	Stainless steel 1.4404 and 1.4301
tomas®-abutments with wire element	Stainless steel 1.4404 and 1.4301, and cobalt-nickel-chromium alloy according to ISO 5832-8

The material composition is included in the orthodontics materials list; please refer to the catalog or www.dentaurum.com. The SSCP is available at <https://ec.europa.eu/tools/eudamed> and www.dentaurum.com.

3. Intended purpose

The implantable products of the tomas®-pin system have the following intended purpose:

- Orthodontic endosteal anchorage elements are inserted into the jawbone to serve as a temporary anchor.

3.1 Indications and contraindications

The tomas®-pin serves as a temporary orthodontic anchor for the skeletal support of orthodontic appliances. It can also be used as part of pre-prosthetic treatment.

3.2 Intra-oral contraindications:

- Unfavorable anatomical jaw positions (e.g., reduced bone availability), inadequate bone quality, pathological findings in bone (e.g., cysts),
- Pathological and/or recurring findings in the mucosa (e.g., lichen),
- Inadequate compliance (e.g., oral hygiene, psychoses, drug, alcohol or tobacco abuse),
- Long-standing therapy-resistant functional disorders,
- Untreated periodontitis.

3.3 Temporary contraindications:

- Temporary medication intake (e.g., anticoagulants, immune suppressants, steroids), acute inflammatory diseases (e.g., maxillary sinusitis), pregnancy.

3.4 General medical contraindications

- Reduced immune defense, poorly managed diabetes, treatment with antiresorptive substances which may cause medication-related osteonecrosis of the jaw (MRONJ),
- Diseases relating to the bone structure,
- Radiotherapy in the head region, active tumors,
- Severe internal diseases (e.g., cirrhosis of the liver),
- Blood coagulation disorders,
- Rheumatic disorders,
- Titanium or nickel allergy,

Adverse reactions and residual risks

- Pain,
- Swelling,
- Excessive bleeding,
- Wound healing disorders,
- Broken drill tip,
- Sensitization to one of the materials,
- Limitation of speech and chewing function,
- Soft tissue or bone infections,
- Nerve damage (paresthesia, dysesthesia, hypo- or hyperesthesia),
- Damaged adjacent teeth,
- Implant fracture or loss,
- Swallowing or aspiration.

4. Precautions

The tomas®-pin should only be inserted by orthodontists, dentists, oral surgeons or maxillo-facial surgeons. Before use, dental professionals should ensure that they have carefully read and observed the Instructions for use supplied.

We recommend dental professionals complete a relevant training course prior to using the tomas® system.

When using tomas®, only use the appropriate original components.

All patients must be examined thoroughly and informed adequately of the procedure, risks and possible complications before treatment with the tomas®-pin!

The product may only be used if it has been correctly prepared. There must be sufficient space at the insertion site for the respective instruments.

Sufficient distance between the implant and tooth roots, vessels and nerves must be ensured. Spatial conditions should be checked in advance using imaging diagnostics.

5. Form of presentation

The tomas®-pin SD and tomas®-pin EP are delivered sterile or non-sterile.

All other tomas® components are supplied non-sterile and must be reprocessed and sterilized before use on the patient.

5.1 Sterile pins – tomas®-pin SD and tomas®-pin EP

These tomas®-pin versions are delivered sterile, in a glass tube wrapped in blister foil packaging. The packaging should not be removed until shortly before insertion.

Dentaurum guarantees the sterility of the tomas®-pin until the expiry date on the packaging, provided that the original packaging is undamaged. Once the date has been exceeded, the sterility of the products can no longer be guaranteed and the tomas®-pin must then no longer be used on a patient.

If the sterile packaging is damaged, do not re-sterilize. Packaging that has been opened may not be re-sterilized! A tomas®-pin that has been removed from the packaging but has not been inserted, may also no longer be used or re-sterilized.

5.2 Non-sterile pins – tomas®-pin SD N and tomas®-pin EP N

The non-sterile tomas®-pins must be cleaned, disinfected and sterilized prior to clinical use according to the Dentaurum reprocessing instructions (www.dentaurum.com/ifu).

6. Insertion of the tomas®-pin SD and tomas®-pin EP

6.1 General information

Adhere to all essential hygienic measures for invasive surgery during insertion of the tomas®-pin.

When using the tomas®-pin as an anchorage unit, ensure that the patient's head and surrounding soft tissue are not subjected to detrimental mechanical influences (e.g., movement of the mucosa, interference from bands, habits).

6.2 Procedure

Selection of the insertion site. Subsequent loading should, as far as possible, be at an angle of 90° or at right angles to the long axis of the tomas®-pin.

- Rinse with disinfecting mouth rinse,
- local anesthesia,
- **pilot drilling in the case of thick compact bone,**
- **when using a 14 mm tomas®-pin, the compact bone must be perforated prior to insertion. When using a 12 mm tomas®-pin, this is recommended,**
- the maximum motor speed for pilot drilling is 600 rpm; ensure that adequate cooling is provided using a sterile, physiological saline solution,
- pick up the tomas®-pin with the insertion instrument,
- screw the tomas®-pin into the bone,
- the tomas®-pin can be loaded immediately,
- coupling the tomas®-pin with the orthodontic appliance.

6.2.1 Picking up the tomas®-pin (sterile version)

- Remove the glass tube from the sterile packaging.
- When the glass tube is opened, the tomas®-pin sits in the metal holder of the cap and is secured in position by a silicone plug against accidental falling out. Hold the silicone plug upwards and carefully remove it.
- Insert the insertion instrument into the metal holder and press until the head of the tomas®-pin perceptibly engages in the retention of the instrument. It must sit firmly in the insertion instrument.
- The tomas®-pin is now removed directly from the metal holder using the insertion instrument and inserted at the intended site.

6.2.2 Picking up the tomas®-pin (non-sterile version)

Place the insertion instrument along the longitudinal axis on the head of the tomas®-pin and press until it perceptibly engages in the retention of the instrument. It must sit firmly in the insertion instrument.

6.2.3 Insertion in the bone

Insertion in the bone can be either manually or using a handpiece. The tomas®-screwdriver or the tomas®-applicator with tomas®-wheel can be used as a manual insertion instrument. The tomas®-torque ratchet is available as an additional aid for the manual procedure with a torque limited to max. 15 Ncm. We recommend a torque of 10 Ncm.

The tomas®-driver is available for insertion using a handpiece. The torque should be set at 15 Ncm max. and the motor speed at 25 rpm max. The tomas®-pin is inserted such that the collar of the implant rests against the gingiva.

7. Removal of the tomas®-pin

- Local anesthesia (optional),
- Remove the coupling elements,
- Unscrew the tomas®-pin manually with one of the insertion instruments mentioned in Chapter 6.2.3.

8. Composition of the tomas®-pin

The tomas®-pin is made of a nickel-free titanium-aluminum-vanadium alloy according to ISO 5832-3.

9. Storage and shelf life

The tomas®-pin should always be stored in a dry, dark place at room temperature in the original packaging.

The tomas®-pin should not be used after the expiry date (see label).

10. Products for single use

tomas®-pins are intended for single use only. Reconditioning of a tomas®-pin that has already been used (recycling) and its reuse on a patient are not permitted.

11. Further information

- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. Therefore, secure all articles used intra-orally against swallowing and/or aspiration.
- Use of third-party components can negatively impact the function and safety. In such cases, Dentaaurum GmbH & Co. KG accepts no liability for replacement products or guarantees.
- Instruments must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.
- Should the dental professional and/or the patient become aware of serious problems arising from the use of the product, inform the manufacturer and the competent authority in the country in which the dental professional and/or the patient is resident.

12. Disposal


When disposing of the product, please adhere to the national regulations.

13. Quality information

Dentaaurum ensures faultless quality of its products. The dental professional is solely responsible for the further processing of the products. Liability for failures cannot be taken, as Dentaaurum has no influence on the further processing on site.


14. Explanation of symbols

Refer to the label. Additional information can be found on the internet at www.dentaaurum.com (Explanation of symbols used on the label REF 989-313-00).

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