ES

Estimado cliente

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EN

Dear customer

Thank you for choosing a quality product from Dentaurum.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit.

Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the Internet at www.dentaurum.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurum GmbH & Co. KG
Turnstr. 31 | 75228 Ispringen | Germany

2. Notes on application

Herbst IV is a further development of the standard Herbst I hinge. It has a ball and socket hinge that allows the patient greater freedom of movement, increasing comfort and tolerance during treatment. Retention using lock washers (C-clips) facilitates fitting by the orthodontist and replaces the relatively awkward insertion of screws.

3. Fabrication of the Herbst IV appliance

Pour les impressions, which include the waxed out bands fitted on the first premolars and molars, with Class IV dental stone, e.g. Rapidur® (REF 164-450-00). Mount the models in the correct position in the orthodontic FKO-split-post fixator (REF 072-004-00) with plaster using the construction bite recorded by the operator. Bend remaloy® wire, half-round,e.g. ∩ 1.50 x 0.75 mm (REF 528-155-00), or round, e.g. ∅ 1.00 mm(REF 528-100-00), to form a lingual arch and then weld the wire to the bands. A palatal bar can be placed in the upper jaw between the first molars to provide stabilization, e.g. remanium® Goshgarian palatal bar or one made from round ∅ 1.0 mm remaloy® wire (REF 528-100-00). Weld the base of the retention hinges to the bands, e.g. using the spot welding unit Master 3000 (REF 079-000-00) or desktop Compact Laser (REF 090-578-50).

Alternatively, splints can be cast instead of using bands. The bases of the retention hinges should not be cast on but welded afterwards for reasons of accuracy and biocompatibility. Ensure that the screws are at the correct angle to the axis between the upper and lower jaw. It is important that the retention hinges operate without friction to ensure smooth functioning of the telescope mechanism. Friction can damage the appliance – especially the bands.

Fitting the telescopic tubes and rods:

First, shorten the telescopic tubes according to the construction bite and then fit the telescopic rods, which should be approx. 2 mm longer. If the telescopic rods are too short, they could come out of the telescopic tube if the patient’s mouth is opened too wide. If the telescopic rods are too long, they could cause mucosal irritation (risk of injury!). The correct length should be checked again intra-orally and adjusted if necessary.

Note: “Double bands” should be used to minimize the risk of the bands being damaged by shear forces. This involves welding a second band over the band already fitted (e. g. using the spot welding unit Master 3000 (REF 079-000-00) or the desktop Compact Laser (REF 090-578-50).

Generally the second band is 2 to 3 sizes larger than the first band fitted. Ready-welded “double bands” can be ordered from Dentaurum. The order form is available on request.

Fitting intra-orally:

The bands are fitted intra-orally using standard band cement. After the cement has set, the prepared retention hinges are fitted on the bases and secured with a lock washer, C-clip. When the hinges are fitted correctly on the base, it should be possible to insert the C-clip into the two slots of the ball unit so that the ends are clipped in at both sides. We recommend using Weingart universal pliers (REF 003-120-00) or How pliers (REF 003-111-00) to facilitate gripping the C-clips.

The first recall appointment should be one week following fitting of the appliance!

If necessary, the Herbst IV appliance can be activated using spacer rings; please refer to our product range.

4. Composition

Information about the composition is included in the material list, see www.dentaurum.com.

5. Safety instructions

Ensure that the C-clip is fitted parallel to the slot of the ball hinge! If the C-clip is fitted at an angle, it can be bent and become loose in the mouth. Always check that the C-clips are properly seated after fitting the Herbst IV appliance and at each recall appointment. If the C-clips are repeatedly fitted and removed, they lose their spring tension and should be replaced!

6. Information for single use products

The Herbst IV appliance is designed for single use only. Reconditioning of Herbst IV that has already been used (recycling) and its reuse on a patient are not permitted.

7. Delivery program

Herbst IV set with ball unit	REF 607-115-00	1 set
For welding to bands or cast splints. Large lateral movement. Fixed to base with C-clips. No screws required. High level of patient comfort. Content: <ul style="list-style-type: none">2 x telescopic tube with ball and socket hinge 2 x telescopic rod with ball unit 2 x coupling for ball unit 4 x base 10 x C-clip, lock washers 1 x Instructions for use	REF 607-116-00	10 pieces
Lock washer (C-clips) for Herbst IV	REF 607-118-00	10 pieces
Base for Herbst IV	REF 607-125-00	1 piece
Telescopic tube with ball and socket hinge for Herbst IV	REF 607-128-00	1 pièce
Telescopic rod with ball unit for Herbst IV		
Spacer rings for Herbst I/ II/IV		
Length: 1 mm	REF 607-103-00	10 pieces
Length: 2 mm	REF 607-104-00	10 pieces
Length: 3 mm	REF 607-105-00	10 pieces
Length: 4 mm	REF 607-107-00	10 pieces
Length: 5 mm	REF 607-112-00	10 pieces

8. Quality information

Dentaurum ensures faultless quality of its products. The recommendations in these instructions are based upon Dentaurum Implant’s own experiences. The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as we, Dentaurum, have no influence on the processing on site.

9. Explanation of symbols used on the label

 **Please refer to the label. Additional information can be found at www.dentaurum.com (Explanation of the symbols REF 989-313-00).**

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