Touch Up

Photo: © Christian Ferrario

Das innovative Touch Up Konzept – Weniger Aufwand, mehr Ästhetik	2 – 7
The innovative Touch Up concept – Minimal effort, enhanced aesthetics	8 – 13
L'innovant concept « Touch Up » – Moins de travail et plus d'esthétique	14 – 19
El innovador Concept Touch Up – Menos esfuerzo, más estética	20 – 25

Touch Up

ceraMotion[®]_{Ti}

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Touch Up

Minimal effort, enhanced aesthetics – The innovative Touch Up concept

The completely user-friendly Touch Up concept produces excellent aesthetics and can also be used for a further range of outstanding applications. The Touch Up special ceramic materials complement the basic programme of the main line ceraMotion® ceramic system. The reduced firing temperatures provide a whole new range of possibilities of use. The Touch Up Opaque, Base Dentines, Dentines and Incisal materials complement the ceraMotion® ceramic system perfectly and are exactly suited.

enhanced Aesthetics

- individual degrees of glaze possible on a single piece of work
- smooth and homogeneous surfaces
- optimal marginal areas, occlusal surfaces and contact points resembling that of natural dentition

time saving due to less effort

- easy to process correction and glaze firing in one
- safe build-up on objects which have been previously inserted
- easy, quick repair due to reduced firing temperature

ceraMotion®
Ti

From users for users.

Dental Technician Hubert Dieker

Master Dental Technician Waldemar Fritzler

Touch Up

Creative firing

- very natural marginal and incisal ridges
- customised degrees of glaze on a single restoration
- adjustment and glaze firing in one

Application of Touch Up material Dentin/Incisal, mixed with Modelling Liquid

Finished restoration

Note: Adjustments can be made with or without glaze material. When using glaze material, first cover the entire surface with glaze material, stain the restoration and then apply Touch Up material over the glaze material and staining.

	Start temp. (°C / °F)	Drying time (min)	Heat rate (°C / °F/min)	Vacuum start (°C / °F)	Vacuum end (°C / °F)	Final temp. (°C / °F)	Holding time *
Glaze and adjustment <u>with</u> glaze material	500 / 932	6	55 / 131	500 / 932	730 / 1346	730 / 1346	1 min
					730 / 1346	760 / 1400	20 s
Glaze and adjustment <u>without</u> glaze material	500 / 932	6	55 / 131	500 / 932	730 / 1346	730 / 1346	1 min
					730 / 1346	760 / 1400	20 s

* the required degree of glaze can be attained using a higher firing temperature and shorter holding time, or a lower firing temperature and longer holding time

Additional firing

- adding occlusal contact points at a later stage
- adding to the fitting surface

Application of Touch Up for missing occlusion

Adding to the fitting surface

Note: Adjustments can be made with or without glaze material. When using glaze material, first cover the entire surface with glaze material, stain the restoration and then apply Touch Up material over the glaze material and staining.

	Start temp. (°C / °F)	Drying time (min)	Heat rate (°C / °F/min)	Vacuum start (°C / °F)	Vacuum end (°C / °F)	Final temp. (°C / °F)	Holding time *
Glaze and adjustment <u>with</u> glaze material	500 / 932	6	55 / 131	500 / 932	730 / 1346	730 / 1346	1 min
					730 / 1346	760 / 1400	20 s
Glaze and adjustment <u>without</u> glaze material	500 / 932	6	55 / 131	500 / 932	730 / 1346	730 / 1346	1 min
					730 / 1346	760 / 1400	20 s

* the required degree of glaze can be attained using a higher firing temperature and shorter holding time, or a lower firing temperature and longer holding time

Touch Up

Repair

□ repair of a restoration worn by a patient

Note: Restorations that have been worn intraorally must be dried out in the preheat furnace. Clean the restoration, the surface must be roughened or sandblasted. Heat the restoration in the preheat furnace from room temperature to 400 °C / 752 °F at a rate of 5 °C / 41 °F/min. Hold time 4 hours minimum, allow to cool slowly.

Then apply the Touch Up materials Base Dentin, Dentin and Incisal, mixed with Modelling Liquid.

	Start temp. (°C / °F)	Drying time (min)	Heat rate (°C / °F/min)	Vacuum start (°C / °F)	Vacuum end (°C / °F)	Final temp. (°C / °F)	Holding time (min)
Dentin firing during repair	500 / 932	6	55 / 131	500 / 932	700 / 1292	700 / 1292	1
Glaze firing during repair <u>without</u> glaze material	500 / 932	4	75 / 167	-	-	700 / 1292	1
Glaze firing during repair <u>with</u> glaze material	500 / 932	6	55 / 131	500 / 932	690 / 1274	690 / 1274	1

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