Fig. 1 Fig. 2 Fig. 3

Fig. 4 Fig. 5


Fig. 6


Fig. 7


Fig. 8


Fig. 9 Fig. 10

Fig. 11 Fig. 12


Fig. 13 Fig. 14


Fig. 15 Fig. 16

Fig. 17 Fig. 18


Fig. 19


Fig. 20

Estimadas/os clientes,

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Dear Customer

Thank you for choosing a quality product from Dentaurum. It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative. As our products are regularly upgraded, we recommend that you always carefully read the current instructions for use supplied with the product and stored in the internet at www.dentaurum.com, even though you frequently use the same product.

1. Manufacturer

Dentaurum GmbH & Co. KG
Turnstr. 31 | 75228 Ispringen | Germany

2. Introduction

In addition to activator therapy, this system is being increasingly used in the treatment of distal occlusions. The fundamental principle of this type of system was developed by A. M. Schwarz.

Intermaxillary guide elements are integrated in the palatal plate and engage at a diagonal in the frontal area of the lingual plate, thereby creating a bimaxillary effect similar to that of the activator.

Besides shifting the occlusion, this system has the advantage of allowing the movement of individual teeth with the assistance of screws and spring components. This in turn permits variations in the formation of the upper and lower rows of teeth. Because of its fixed lingual and palatal position, the appliance can also be worn during the day, and there is little danger of it becoming dislodged during the night. In the course of time, the idea of A. M. Schwarz was developed into various systems. Industrially produced parts are extremely useful in fabricating this system.

3. Description

The flexible forward thrust double plate system acc. to Schaneng is a more advanced development and differs from other systems in the flexibility of the intermaxillary guide elements.

The guiding ridges connected with one another consist of a special 1.0 mm round spring wire. The double loop at the end of each creates the flexibility of the ridges. Due to its special seating device, this construction can be locked very accurately sideways onto the guide pins of the expansion screw (REF 600-302-30) (Figs. 1 – 3). The angle of the diagonal plane (60°) is set precisely in the lower jaw using a locating template (Fig. 4) as well as a narrow and a wide plastic adapter (Fig. 5).

4. Fabricating the flexible forward thrust double plate system

The first stage in making this flexible appliance begins by joining the two models with the functional occlusion in the fixator (Fig. 6). Retaining and moving elements as well as various screws (apart from the upper screw with guide ridges) can now be attached to the model in the same manner as with standard appliances (Figs. 7 + 8).

To set the angle of the diagonal level at 60° on the lower plate, the locating template and the plastic adapter are required.

The locating template is laid on the lingual model parallel to the occlusal plane. At the side, the 60° guide surfaces of the template are positioned adjacent to the screw retainer and tightly against the screw body. Notches in the front and rear areas of the template facilitate the alignment of the centre of the lower jaw (Fig. 9).

The plastic adapter is now pushed over the guide surfaces of the template and attached to the screw holder by means of adhesive wax (Fig. 10).

When the locating template has been removed, the position and distance of the plastic adapter in relation to the model and the screw can be seen and checked (Fig. 11).

The lower plate is made in the usual way (salt and pepper technique), polymerized in the autoclave and the plastic adapter is then removed (Fig. 12).

Before fitting the guide ridges of the upper plate onto the diagonal plane, the lower plate must first be fully finished and polished and the anterior region thinly covered in wax to avoid the possibility of the plates bonding together during polymerization (Fig. 13).

Note: The loops on the guide ridges of the upper screw are supplied waxed out. This ensures free movement/flexibility of the guide ridges in the acrylic. The layer of wax should be checked and adjusted if necessary before applying acrylic (Fig. 14).

The models are joined in the fixator before the palatal screw is positioned (Fig. 15). The guide ridges are now waxed into position parallel to the diagonal plane. The screw itself should be placed as close to the palate as possible (Fig. 16). Following this, the acrylic is applied in the palatal area with the fixator closed. The lingual plate should be coated lightly with vaseline beforehand to avoid any fusing of the plates during polymerization (Fig. 17).

The next stage is to finish the upper plate and to polish both plates. Ensure that the loops **are completely free of acrylic and that the guide ridges are not damaged** by the grinding tools (Fig. 18).

When the flexible forward thrust double plate is finished, it is again checked in the closed fixator to ensure that the guide ridges are in the correct position (Fig. 19). To ensure a satisfactory finished product (Fig. 20), it is essential to keep to the correct sequence of operations when fabricating the plates, as well as the correct adjustment of the angle of the diagonal level.

Important

The functional occlusion should always be determined on the patient in the following manner:

sagittal: **max. 5 mm** forward pressure
vertical: **3 mm** locking in front teeth area.

5. Summary

The development of the forward thrust double plate system into a "flexible forward thrust double plate system" has the following practical advantages:

- greater convenience for the patient through flexibility of the ridges while achieving the same effect
- no bending or breaking of the individual ridges due to monolithic design
- simple and accurate setting of diagonal level using locating template and plastic adapter
- guide ridges can be hooked onto a normal expansion screw
- no risk of injury through the guide ridges as these are connected with one another and rounded in shape.

6. Additional tips

Should the user and/or the patient gain knowledge of serious incidents arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the user and/or the patient is resident are informed.

7. Availability

Assortment (in wooden case) contains: 3 pcs. ridges with palatal screws 1 piece adapter, narrow 1 piece adapter, wide 1 piece locating template 1 piece instructions for use	600-000-00
Ridge for flexible forward thrust double plate with palatal screw 1 piece 10 pieces	600-001-00 600-001-30
Ridge for flexible forward thrust double plate without screw 1 piece 10 pieces	600-002-00 600-002-30
Adapter, narrow, for even expansion of upper and lower plates, 10 pieces	600-003-00
Adapter, wide, for normal expansion of upper plate and low expansion of lower plate, 10 pieces	600-004-00
Locating template, 1 piece	600-005-00

8. Quality information

Dentaurum ensures faultless quality of the products manufactured by us. These recommendations are based upon Dentaurum's own experiences. The user is responsible for the processing of the products. Responsibility for failures cannot be taken by Dentaurum as we have no influence on the processing on site.

9. Explanation of symbols used on the label

⚠ Refer to the label. Additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00).

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