Gebrauchsanweisung
Instructions for use
Mode d'emploi
Modo de empleo
Modalità d'uso

D
DENTAURUM

- Informationen zu Produkten finden Sie unter www.dentaurum.com
- For more information on our products, please visit www.dentaurum.com
- Vous trouverez toutes les informations sur nos produits sur www.dentaurum.com
- Descubra nuestros productos en www.dentaurum.com
- Informazioni sui prodotti sono disponibili nel sito www.dentaurum.com

Stand der Information | Date of information | Mise à jour | Fecha de la información | Data dell'informazione: 05/21

Turnstr. 31 | 75228 Ispringen | Germany | Tel. +49 7231/803-0 | Fax +49 7231/803-295
www.dentaurum.com | info@dentaurum.com

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© 2014 Dentaurum GmbH & Co. KG. All rights reserved.

EN

Dear Customer

Thank you for choosing a quality product from Dentaurum.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit.

Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current instructions for use supplied with the product and stored in the internet at www.dentaurum.com, even though you frequently use the same product.

1. Manufacturer

Dentaurum GmbH & Co. KG | Turnstr. 31 | 75228 Ispringen | Germany

2. General description

discovery® pearl is a ceramic bracket manufactured using the CIM (Ceramic Injection Molding) process. It consists of highly pure polycrystalline aluminum oxide. discovery® pearl ceramic brackets are supplied with an identification system (positioning guides + FDI marking on the base, see section 9) for all systems. This allows for easy identification and placement of the brackets.

3. Intended use

For the orthodontic treatment of malocclusions.

4. Working stages

4.1 Patient consultation

Inform the patient about the possible consequences and risks involved in treatment using brackets, in particular during debonding. This applies especially in the case of teeth with enamel that already exhibits signs of initial damage (cracks, crazing, splitting off, etc.).

4.2 Tooth conditioning

Adhere to the bracket adhesive manufacturer’s instructions for conditioning the teeth.

4.3 Selection of bracket adhesives

discovery® pearl brackets can be bonded using all commercially available bracket adhesives. We recommend CONTEC Ic (REF 163-111-00), ConTec Go! (REF 163-321-00) and SmartBond® (REF 163-506-00). Adhere to the manufacturer’s instructions for the bracket adhesive.

4.4 Removal and handling

Grip the discovery® pearl brackets using clamping tweezers (e.g. REF 025-277-00) mesio-distally and at the bar of the positioning guide to remove them from the packaging and handle according to the following steps. The base of the discovery® pearl brackets should not be touched with your fingers or contaminated in any way, as this could otherwise detrimentally affect the bond.

4.5 Application of the bracket adhesive

Teeth and bracket adhesive are prepared according to the manufacturer’s instructions. Apply a small amount of adhesive to the center of the bracket base, but enough to cover the entire bracket base. This guarantees an optimum bond.

4.6 Indirect bonding technique

discovery® pearl brackets are bonded on the plaster or the plastic model with the bracket adhesive. This adhesive layer should not be removed anymore! Not suitable are indirect bonding techniques using sugar solution or other water-soluble lab adhesives that later have to be completely removed. This process destroys the mechanical retention of discovery® pearl. A reliable fit of the bracket on the tooth is no longer guaranteed! Thermoformed splints or silicone transfers, e.g. Lingusil (REF 164-000-00) can be used as transfer matrix.

4.7 Placement and alignment of the bracket

After applying the adhesive, place, align and press the bracket in the usual way onto the tooth until the adhesive exudes from the sides. It is essential to avoid contact with the antagonist due to the risk of abrasion and/or premature bracket loss. If necessary, a palatal bite raise can be used e.g. Bite-Bumper™ (REF 705-000-00). Completely remove the exuded excess adhesive before it has cured. This is particularly important to ensure easy debonding at the end of the treatment.

4.8 Ligating the archwires

For the initial treatment phase we recommend the use of super elastic Tensic®, rematitan® LITE or translucent pearl archwires.

Do not use force to ligate the archwires. Do not torque the archwires when they are ligated.

Once the adhesive has cured (curing time according to the adhesive manufacturer’s instructions) and the placement aids have been removed, ligating can be completed in the usual way using wire or plastic ligatures (e.g. Dentalastics® pearl REF 774-563-00).

Note: Careful ligation avoids damage to the brackets!

4.9 Bracket debonding

First, remove the archwire. For debonding the discovery® pearl brackets we recommend the special discovery® pearl debonding instrument, REF 019-001-00, which grips the body of the bracket directly at the base. The bracket is gently removed from the bracket adhesive by slightly tipping it sideways in the direction of the slot and archwire.

Alternatively, the straight bracket removal pliers (REF 004-346-00) can be used for debonding the discovery® pearl brackets. Place the blades of the bracket removal pliers mesio-distally on the bracket base, **directly on the surface of the tooth**. The bracket is gently removed from the tooth by slightly twisting the pliers.

Note:

Despite their small, delicate shape there is no risk of fracture with discovery® pearl brackets provided they are used correctly. Use of other instruments can result in bracket fracture. If the bracket fractures during incorrect debonding, we recommend carefully removing the fractured segments and the adhesive from the enamel using bracket debonding pliers (e.g. REF 003-349-00) or a scaler (e.g. REF 027-349-00).

Use care when debonding the segments, as there is a risk of the enamel chipping. Use an aspirator and wear protective glasses.

5. Information about the composition

Information about the composition is included in the materials list, see catalog or www.dentaurum.com.

6. Safety information

The product should not be used if there is a known allergy to one or more of the ingredients.

7. Storage and shelf-life

Store discovery® pearl brackets in closed packaging.

discovery® pearl brackets do not have a limited shelf life.

8. Information for single-use products

discovery® pearl brackets are intended for single use only. Reconditioning of discovery® pearl brackets that have already been used (recycling) and reuse on patients is not permitted.

Micro cracks may result from loading during treatment and in particular bracket removal that could cause fractures in the ceramic during reuse!

9. Notes about the identification system (positioning guides)

9.1 The color of the positioning guides and the tooth symbol indicate the tooth / bracket combination.

9.2 Maxilla: right = red
left = green
Mandible: right = yellow
left = blue

9.3 The occlusal-gingival distance can be determined by means of the occlusal indicator. The distance between the occlusal end and the center of the slot is 6 mm.

9.4 The tooth symbol always points from occlusal to distal.

9.5 Removing positioning guides

Once the adhesive has cured (curing time according to the adhesive manufacturer’s instructions), grip the mid-section of the positioning guides using tweezers or Weingart pliers (e.g. REF 003-120-00) and remove.

The practical FDI notation on the base also aids clear identification.

10. Additional tips

- Technical information such as torque, angulation, size etc. can be found on the relevant product label or in the catalog information.
- There are no particular precautions to be taken when disposing of discovery® pearl brackets after treatment.
- Should the user and / or the patient gain knowledge of serious incidents arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the user and/or the patient is resident are informed.

11. Delivery program

For the complete range of discovery® pearl products please see our Orthodontics Catalog.

12. Quality information

Dentaurum ensures a faultless quality of the products manufactured by us. These recommendations are based upon our own experiences. The user himself is responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence on the processing on site.

13. Explanation of symbols used on the label

 Please refer to the label. Additional information can be found at www.dentaurum.com (Explanation of the symbols used on the label REF 989-313-00).

© 2014 Dentaurum GmbH & Co. KG. All rights reserved.

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