
D
DENTAURUM
1886

homas

homas PI
Palatal Implant

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Dentaurum GmbH & Co. KG · Turnstr. 31 · 75228 Ispringen · Germany

2. General device description

The torque ratchet is a precision instrument that can be disassembled, with a torque control of 10 Ncm – 50 Ncm and a blocking function for higher torques. The torque ratchet is a dental instrument for the manual preparation of the implant site, for tightening and loosening screws for implants and prosthetic components, and for inserting and removing dental implants and orthodontic implants.

3. Intended purpose

These torque ratchets are intended for temporary use when screwing screws in and out and for placement and removal of dental implants and orthodontic implants. And also for releasing connections within defined torque ranges for dental use in the fields of implant dentistry, osteosynthesis, surgery, prosthetic dentistry and orthodontics. The torque function can also be "blocked." In the blocking position, higher torques can be transmitted both for tightening and for loosening.

4. Intended users

The torque ratchet is supplied exclusively to qualified personnel such as surgeons, orthodontists, dentists and dental technicians. It may only be used by qualified personnel who are familiar with implant dentistry, including diagnosis, preoperative planning, surgical procedures and prosthetic restoration.

5. Medical indications/contraindications

5.1 Indications

- Use for manual thread cutting (for tomas® PI) in combination with the ISO shank/ratchet adapter with the respective defined torques
- Manual insertion or removal of implants (or tomas®-pin) into/from the jawbone with the respective defined torques
- Manual insertion or removal of fastening screws in combination with hex keys with the respective defined torques
- Use with system-internal instruments with a tool head specifically designed for the tool connection

5.2 Contraindications

- Use of third-party instruments or tool connections
- Exceeding the defined torques
- All indications not listed here are contraindicated

6. Patient target group

Orthodontic implants (tomas® PI or tomas®-pin) can be used for patients of any age, as long as they already have their permanent teeth. Most patients are adolescents. Please consider the indications/contraindications.

7. Information on use, processing, fitting and insertion

Prior to use, dental professionals should ensure that they have carefully read and understood all Instructions for use.

We recommend that, prior to using Dentaurum products, dental professionals attend a training course offered by Dentaurum to familiarize themselves with the Dentaurum products.

Carry out thorough radiological, physical and psychological examinations before use and inform patients of the process in detail. Dentaurum recommends full clinical, radiological, photographic and statistical documentation.

Important:

- The user should ensure the products cannot be aspirated during intra-oral use.
- The tomas® torque ratchet is fitted with an adapter for attaching the tomas® applicator (REF 302-004-20 or REF 302-004-70). The tomas® PI-torque ratchet is fitted with an adapter for attaching the tomas® PI-applicator (REF 303-320-00). Without the appropriate applicator, the torque ratchet cannot be used as intended with the tomas® pin or the tomas® PI 4.0 palatal implant.

7.1 Torque ratchet

- Immediately before each use, check the product for any signs of wear, functional impairment or loss, or signs of corrosion. Furthermore, the instrument must be assembled correctly.
- Immediately discard damaged products or those with the previously identified deficiencies and no longer use them in this condition!
- Products with damaged sterile packaging (after processing by the user) must not be used and must be reprocessed again in accordance with these instructions.

7.2 Possible default settings

7.2.1 Torque function:

The desired torque range can be adjusted continuously using the torque adjustment screw via the spring. The setting is visible on the marking of the scale sleeve.

7.2.2 Blocked function:

Turn the torque adjustment screw until it reaches the ∞ (infinity) mark. Do not turn too tightly.

CAUTION:

- Do NOT loosen the two screws (X) on the torque adjustment screw, as this will result in a loss of the factory default settings (see Figure 1).
- The appropriate applicator (see 7.) is inserted into the torque ratchet or the designated adapter, or pulled out for disassembly. The lock on the applicator can be released by applying firm pull or pressure.

Figure 1

7.3 Changing the tool (adapter)

Pull back the pin in the direction of the arrow (→) on both sides using your thumb and index finger and remove or insert the tool (adapter). (see Figure 1)

7.4 Correct handling of the torque release

- The pressure point for precise torque release is located exclusively on the handle of the torque adjustment screw (see arrow in Figure 2).
- Only activate with the touch of a finger.
- Do not touch the handle with your thumb and index finger for activation.
- When the adjusted torque has been reached, the scale sleeve will bend around the axis in the ratchet head. The release is audible and perceptible.

Figure 2

CAUTION:

After torque release, do not continue to apply more force. The ratchet or dental components could be damaged.

When letting go of the handle, the ratchet returns to its initial position.

8. Technical specification

Notes on composition

Product	Composition
Torque ratchet	1.4197 / PEEK

The SSCP is available at <https://ec.europa.eu/tools/eudamed> and www.dentaurum.com.

9. Safety instructions and precautions

- Each part belongs only to the instrument supplied. The exchange of components (even for identical instruments) is not permitted and requires renewed inspection by the manufacturer before the torque function can be used. To ensure maximum precision, the single components of the torque ratchet are coordinated and adjusted during manufacturing. During cleaning and subsequent assembly of the torque ratchet ensure that no components from different torque ratchets are used, as the tolerance values can then no longer be assured to be correct. The replacement parts cannot be purchased separately.
- Always store the torque ratchet with the spring relaxed at approx. 10 Ncm.
- Exceeding the torque specified by Dentauro can cause mechanical damage to components, to the implants, Mini pins and destruction of bone structures.
- The block function mode should be used with extreme caution. After use the value must be reset to standard torque to prevent subsequent faults.
- The precision of the instrument is $\pm 10.0\%$ on leaving the factory
- Do not use the product if there is a known allergic reaction to one or more of the material components.
- Due to the small size, the article could be swallowed and/or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. Therefore, secure all articles used intra-orally against swallowing and/or aspiration.
- Different types of alloy in the oral cavity can lead to galvanic reactions.
- If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.
- Observe the indications for the individual products. Consider all indications that are not listed explicitly to be contraindications.
- Disinfection, cleaning and sterilizing should not be done in the original packaging.
- Products must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.
- Technical advice (oral and written) is based on the scientific and technical knowledge available at the time the product is put on the market. It does not release the dental professional from the responsibility of personally checking the suitability of the products for the intended indications and applications. Advice is only given as a non-binding recommendation and cannot be assumed to provide any form of assurance, liability or guarantee.

Severe incidents:

Report all severe incidents arising from the use of the product to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.

Severe incidents are identified as malfunctions and a drop in performance or deterioration of the properties of a medical device which resulted or may have resulted in a direct or indirect temporary or permanent deterioration in a patient's state of health, in the death of a patient, a dental professional or other persons or as posing a serious danger to public health.

10. Storage and shelf life

Delivery state / storage of brand new products:

After sterilization, the products must be stored dust free and dry in the sterilization packaging. The delivery state is in a relaxed position at approx. 10 Ncm. This must also be maintained as a general storage status to ensure optimum functionality and related service life.

Product service life

The service life expires when the set torque has been attained 5000 times. Frequent reprocessing has little effect on these instruments – provided they are handled with care and are undamaged and fully functional. The end of a product's service life is usually dictated by wear and tear and damage due to use and depends on many factors – including the type, duration, and frequency of use, as well as how the instruments are handled, stored, and transported.

Damaged, blunt, or soiled instruments must not be used.

No liability is accepted if these instructions are disregarded. The same applies to damage caused by improper reprocessing or handling, such as excessive mechanical impact, dropping, overloading, etc.

11. Disposal

Adhere strictly to the rules that apply for your region and the details outlined in the safety data sheets, if available. Products to be disposed of should be treated as surgical waste and be decontaminated in compliance with the relevant regulations for medical waste disposal.

12. Information for multiple use products

The torque ratchet is supplied non-sterile and is intended for multiple use. The product must be disinfected, cleaned and sterilized before being used initially and before being used again. Torque ratchets must be replaced if they cannot be clearly identified or if their function is impaired due, for example, to poor readability of the markings and/or labels.

13. Information on disinfection/cleaning/sterilization

The torque ratchet is supplied non-sterile. It is thus compatible with the further handling described here for end users. The required necessary working methods have been validated at the supplier. To provide users with information about a functioning process for the individual cleaning steps, an efficiency check of reprocessing was conducted for each specified process in an accredited laboratory. The validation of the cleaning and sterilization procedures used is the responsibility of the user on site or of a central sterile supply department (CSSD).

CAUTION

Any change to the packaging or product also constitutes a change to the validated delivery state.

If several torque ratchets are in use, do not interchange the individual parts. Each individual part belongs to one instrument.

Do not use metal brushes or cleaning sponges.

- Reprocessing only by appropriately trained persons.
- The water used must at least be of drinking water quality (see information in the individual processing steps).
- This reprocessing instruction specifies the cleaning agents and disinfectants used in the efficiency verification. If an alternative cleaning agent and/or disinfectant is used, it must be listed by the RKI or the German Association for Applied Hygiene (VAH) and be compatible with the material. The pH value must lie between 4.5 and 10.
- The reprocessor is responsible for ensuring that the actual reprocessing performed with the equipment, materials, and personnel used in the processing facility achieves the desired results. As a rule, this requires validation and routine monitoring of the process and the equipment used.

When selecting water qualities for treating instruments, always prefer fully demineralized water so that corrosive deposits of salts and silicates do not accumulate and can be avoided or reduced to an absolute minimum.

13.1 Transport/place of use – reprocessing

The first steps of proper reprocessing begin immediately after use on the patient.

Coarse contamination, residues of e.g. filling material, disinfectants and other medicinal products must be removed before the instruments are set aside.

- Wherever and whenever possible, dry disposal (moistened, closed system) should be preferred. Otherwise, disposal is to be carried out in accordance with standard hospital procedures. The torque ratchet must be transported and disposed of in a closed container or a sealed protective cover.
- As a matter of principle, avoid allowing any residue from the application to dry!
- Avoid long waiting times for reprocessing, e.g. overnight or over the weekend, for both types of disposal (< 6 hours).

13.2 Cleaning and disinfection

Cleaning and disinfecting solutions with a pH value between 4.5 and 10 are recommended for cleaning – follow the manufacturer's instructions for these products (e.g. intended purpose, dosage, reaction time, etc.).

In general, when storing parts for the purpose of cleaning, make sure that they do not touch each other or lie on top of each other in order to avoid unwashed areas and to perform the cleaning process as efficiently as possible.

13.2.1 Basics

A mechanical procedure (washer disinfector) should be used if possible for cleaning and disinfection. Because of its inferior efficacy and reproducibility, a manual procedure – even with an ultrasonic cleaner – should only be used if a mechanical procedure is not an option.

The preparation and preconditioning described below are required in both cases.


13.2.2 Preparation for decontamination

Coarse impurities should be removed from the instruments immediately after use (within maximum 2 hours).

Before cleaning (regardless of the selected cleaning method), the torque ratchet must be dismantled into the individual parts. This can be done without tools. To do this, simply unscrew the torque adjustment screw ④ completely, remove the spring ③ and pull off the ratchet head with the threaded rod.

Take care not to lose the plastic washer as this would have a negative impact on the instrument's precision. The plastic washer needs only to be removed if there is visible contamination. The washer can be pulled off if necessary. After cleaning, press the washer back into place.

Removal of the adapter

Pull back the pin ② in the direction of the arrow on both sides using your thumb and index finger and remove the adapter.

13.2.3 Preconditioning

Procedure/Preconditioning

Regardless of the subsequent cleaning method, precleaning is always necessary. Hold the parts under cold mains water (drinking water quality, < 40 °C/104 °F) until all visible contamination has been removed. Use a soft brush to remove stubborn dirt. Rinse hollow sections and lumens intensively with a water pistol (or similar) using cold mains water (drinking water quality, < 40 °C/104 °F) (for > 30 seconds).

13.2.4 Machine process – thermal disinfection

Proof of basic suitability of the instruments for effective automatic cleaning and disinfection was provided by an independent, accredited test laboratory, recognized by the ZLG test laboratory (Central Authority of the Federal States) (§15 (5) MDR) using a G 7835 CD washer disinfector (thermal disinfection, Miele & Cie. KG, Gütersloh, Germany) and the cleaning agent neodisher® mediClean (Dr. Weigert GmbH & Co. KG, Hamburg, Germany). The procedure described was taken into account during the tests.

13.2.4.1 WD (washer disinfector) and media

When selecting a WD, please ensure

- that the efficacy of the washer disinfector has been certified (e.g. DGHM or FDA approved or CE marking according to DIN EN ISO 15883),
- which ideally includes using a certified program for thermal disinfection (A0 value > 3000 or – for older appliances – a minimum of 5 mins at 90 °C/194 °F), (with chemical disinfection there is the risk of disinfectant residue on the instruments),
- the program used is suitable for the instruments and has an adequate number of rinse cycles,
- that only demineralized water is used for rinsing,
- that the air used for drying is filtered (oil-free, low in particles and of low microbiological contamination) and that the WD is regularly serviced and inspected.

The material applications, concentrations, temperatures and reaction times given by the cleaning agent and disinfectant manufacturer along with instructions for rinsing should be strictly adhered to.

13.2.4.2 Machine cleaning/disinfection (→ RECOMMENDED)

The following program parameters were used to prove the suitability (program: Des-Var-TD / washer disinfector Miele G7835 CD):

- The parts must be placed in the mobile injection unit (E450/1) on a tray
- Preclean 1 minute (cold mains water in drinking water quality < 40 °C/104 °F) → Drain water → Preclean 3 minutes (cold mains water in drinking water quality < 40 °C/104 °F) → Drain water
- clean 10 minutes at 55 °C/131 °F ± 5 °C/41 °F with 0.2% alkaline cleaning agent (0.2% neodisher® MediClean) → Drain water
- Rinse 1 minute with deionized water < 40 °C/104 °F → Drain water → rinse 2 minutes with demineralized water < 40 °C/104 °F → Drain water
- Automatic disinfection >5 minutes at 92 °C/198 °F ± 2 °C/36 °F with deionized water.
- Automatic drying program in WD at 90 °C/194 °F ± 2 °C/36 °F for a minimum of 30 minutes (→ 60 °C/140 °F ± 5 °C/41 °F in rinsing area).
- Procedure during (repeated) reprocessing:
- Place the instruments in the WD. Ensure that the instruments do not come into contact with one another.
- Start the program.
- Remove the instruments immediately from the WD when the program has finished and ensure the instruments are dry enough before packing.
- Check and pack the instruments immediately after removal if possible.

13.2.4.3 Additional manual drying

Should it be necessary to subsequently dry the instruments manually in the case of residual moisture, use a lint-free cloth and/or blast the lumens with sterile, oil-free compressed air.

13.2.5 Manual process

Proof of basic suitability of the instruments for effective manual cleaning and disinfection was provided by an independent, accredited test laboratory, recognized by the ZLG (Central Authority of the Federal States) (§15 (5) MDR) using the cleaning and disinfection agents named below. The procedure described was taken into account during the tests.

13.2.5.1 Manual cleaning

1. Place the products in an alkaline detergent (e.g. 0.5% neodisher® MediClean) in an ultrasonic bath for an exposure time of 10 minutes. A maximum temperature of 40 °C/104 °F should not be exceeded. Follow the instructions of the manufacturer of the cleaning agent.
2. Clean the products with a soft brush. If present, rinse hollow sections and lumens intensively with a water pistol (or similar) (> 30 seconds).
3. Rinse the products under running mains water (drinking water quality) to remove the cleaning agent (> 15 seconds).

13.2.5.2 Manual disinfection

1. Dip the products into a disinfectant that is listed by the Robert Koch Institute (RKI) or the Association for Applied Hygiene (VAH). Follow the instructions of the manufacturer of the disinfectant. It must be assured that the disinfectant reaches all areas of the product (move the parts around in the disinfectant solution; if necessary, use a syringe - without cannula - to rinse obscure surfaces).
2. Rinse the products (full rinse on the inside, outside and in hollow areas) in demineralized water for >60 seconds.

13.2.5.3 Manual drying

1. Dry the products with a lint-free cloth for single use. To avoid moisture remaining in hollow areas as far as possible, blast the areas with sterile, oil-free compressed air.

13.3 Monitoring

The best method to identify and remove defective instruments is by examining them and carrying out functional tests before and after each use. Specific functional areas (e.g. connection for adapter, torque initiator) and moving parts require special attention. Allow the parts to cool down to room temperature. Parts with damaged surfaces, chippings, contamination, discoloration or corrosion must be removed. Discard any instruments that are deformed, worn in their functionality or otherwise damaged. Instruments that are still contaminated should be cleaned and sterilized again.

13.4 Maintenance (by user)

If several torque ratchets are in use, do not interchange the individual parts. Each individual part belongs to one instrument.

Assembly

To assemble the torque ratchet correctly, connect the components in the following order: first remove the pin ② as described above and insert the adapter ①.

Caution: to avoid confusion, the adapter a can only be inserted from one side.


◆ lightly lubricate marked areas with instrument maintenance oil.

Ensure that only instrument oils (paraffinic white oil without corrosion inhibitor or other additives) are used, which – depending on the maximum sterilization temperature used – are approved for steam sterilization and are certified as biocompatible. The oil should be used as sparingly as possible.

Reassemble the ratchet and perform a function test.

Slide the spring ③ back over the torque adjustment screw. Then pass the ratchet head ⑤ with the threaded rod through the scale sleeve ⑥ and screw to torque adjustment screw. After complete assembly and before each use, check the correct function. The instrument is ready for use when there is an audible regular ratchet noise and the mechanism of the torque limit functions.

After reassembly and prior to sterilization, the torque ratchet should be relaxed at max. 10 Ncm.

13.5 Packaging

The products should be sterilized in suitable sterilization packaging. The supplier used double sterilization packaging (standard for hospitals) for verification; it is therefore also possible to use suitable single sterilization packaging.

Flash sterilization or the sterilization of unpacked instruments is not permitted!

13.6 Sterilization

Proof of basic suitability of the instruments for effective sterilization was provided by an independent, accredited test laboratory, recognized by the ZLG (Central Authority of the German federal states) using an EHS3870 autoclave (Tuttnauer Europe B.V., Breda, NL), pre- and post-vacuum, and using sterilized packaging RB 51-3P and RB 52-3P (Steriking foil). The procedure described was taken into account during the tests. These specifications must be adhered to:

3 vacuum cycles | 132 °C/270 °F | ≥ 4 minutes holding time | Drying under vacuum for at least 20 minutes³

13.6.1 Sterilization process – fractionated vacuum process

Only the following sterilization procedures should be used for sterilization.

Other sterilization procedures are not permitted; the dental professional must provide their own proof of efficiency if another procedure is taken.

Fractionated vacuum process^{1,2} (with adequate product drying³)

- Steam sterilizer according to DIN EN 13060 / DIN EN 285 or ANSI AAMI ST 79 (for USA: FDA Clearance)
- Validated in accordance with DIN EN ISO 17665 (valid IQ/OQ (commissioning) and product-specific performance qualification (PQ))
- Maximum sterilization temperature 134 °C/273 °F plus tolerance in accordance with DIN EN ISO 17665
- Sterilization time

1 Min. 3 vacuum stages

2 Use of the less effective gravitation process is only permitted if a fractionated vacuum process is not available. The gravitation process requires much longer sterilization times which must be determined and validated for the specific instruments, equipment, procedure and parameters used at the responsibility of the dental professional in charge.

3 The actual drying time required for the product depends directly on the parameters used, for which the dental professional alone is responsible (loading configuration and density, state of sterilization) and which must therefore be determined by this person. Generally speaking, the drying time should not be less than 20 minutes.

13.7 Material resistance

- When selecting the cleaning agent and disinfectant please ensure that they do not contain the following components:
 - Organic, mineral, and oxidizing acids or strong alkalis
 - Organic solvents (e.g. alcohols, ethers, ketones, benzines)
 - Oxidation agents (e.g. hydrogen peroxides)
 - Halogens (chlorine, iodine, bromine)
 - Aromatic / halogenated hydrocarbons
 - Do not use acidic rinse aids or neutralizing agents!
- All instruments should only be exposed to temperatures not exceeding 138 °C/280 °F.

13.8 Repairs/Maintenance (by manufacturer)

- The default setting upon delivery is accurate to within ± 10% of the adjustable torque value.
- We recommend having the torque ratchet checked by the manufacturer once a year.
- No warranty shall apply to instruments that have been repaired by companies or persons who have not been expressly authorized to do so.
- The instrument must be visibly cleaned and sterilized before shipping in accordance with these instructions (e.g. using a wrap indicator on sealed sterilization foil). Otherwise, the product will be returned to the sender at his expense without further processing.

Further applicable guidelines for the reprocessing of medical devices can be found at www.rki.de or www.a-k-i.org

14. Product overview/Scope of delivery

For the product range, please see our Orthodontics Catalog. i.e., at www.dentaurum.com.

15. Quality information

Dentaurum guarantees the user that the products are of impeccable quality. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaurum has no influence on the processing on site.

16. Explanation of symbols

Please refer to the label. See section "Additional information".

17. Additional information

Additional information can be found at www.dentaurum.com/ifu

Instructions for use	REF
Explanation of symbols used on the label	989-313-00

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