tomas[®]-uprighting spring

- Informationen zu Produkten finden Sie unter www.dentaaurum.com
- For more information on our products, please visit www.dentaaurum.com
- Vous trouverez toutes les informations sur nos produits sur www.dentaaurum.com
- Descubra nuestros productos en www.dentaaurum.com
- Informazioni su prodotti sono disponibili nel sito www.dentaaurum.com

Stand der Information | Date of information | Mise à jour | Fecha de la información | Data dell'informazione: 02/16

Turnstr. 31 | 75228 Ispringen | Germany | Tel. +49 72 31/803-0 | Fax +49 72 31/803-295
www.dentaaurum.com | info@dentaaurum.com


Fig. 1:
Aufrichten ohne Intrusion und Extrusion
Uprighting without intrusion or extrusion
Redressement sans intrusion ni extrusion
Enderezamiento sin intrusión ni extrusión
Uprighting senza intrusione ed estrusione


Fig. 2:
Aufrichten ohne Intrusion und Extrusion
Uprighting without intrusion or extrusion
Redressement sans intrusion ni extrusion
Enderezamiento sin intrusión ni extrusión
Uprighting senza intrusione ed estrusione


Fig. 3a:
Aufrichten mit Intrusion, Variante 1
Uprighting with intrusion, version 1
Redressement avec intrusion, variante 1
Enderezamiento con intrusión, variante 1
Uprighting con intrusione, variante 1


Fig. 3b:
Aufrichten mit Intrusion, Variante 2
Uprighting with intrusion, version 2
Redressement avec intrusion, variante 2
Enderezamiento con intrusión, variante 2
Uprighting con intrusione, variante 2


Fig. 4a:
Aufrichten mit Extrusion, Variante 1
Uprighting with extrusion, version 1
Redressement avec extrusion, variante 1
Enderezamiento con extrusión, variante 1
Uprighting con estrusione, variante 1


Fig. 4b:
Aufrichten mit Extrusion Variante 2
Uprighting with extrusion, version 2
Redressement avec extrusion, variante 2
Enderezamiento con extrusión, variante 2
Uprighting con estrusione, variante 2


Fig. 5:
Aufrichten mit mesialem Wurzelortorque
Uprighting with mesial root torque
Redressement avec torque radicaire mérial
Enderezamiento con torque mesial de la raíz
Uprighting con torquemessiale della radice

tomas[®]-uprighting spring (Aufrichtefeder)

REF 302-009-00

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REF 302-009-00

Dear customer

Thank you for choosing a quality product from Dentaaurum. It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative. As our products are regularly upgraded, we recommend that you always carefully read the current instructions for use supplied with the product and stored in the internet at www.dentaaurum.com, even though you frequently use the same product.

- Manufacturer**
Dentaaurum GmbH & Co. KG | Turnstr. 31 | 75228 Ispringen | Germany
- General product description**
The tomas[®]-uprighting spring consists of a super-elastic nickel titanium spring (dimensions 0.42 mm x 0.62 mm; 16 x 24) and a steel wire (dimensions 0.46 mm x 0.64 mm; 18 x 25). Both wire elements are connected with a clamping tube. The nickel titanium spring is pushed into the clamping tube which is then crimped together to fix the wires in the desired position. If the tooth is tilted at approx. 40° the uprighting moment is between 10 and 20 Nmm. Depending on how the Alpha curve is formed, the tomas[®]-uprighting spring has a large plateau in the range of 8 to 15 Nmm and creates an intrusive force of approx. 0.5 to 1.0 N.
- Area of application**
The tomas[®]-uprighting spring is for uprighting teeth, particularly molars. The uprighting action can be combined with intrusion or extrusion, depending on how the spring is set. The tomas[®]-uprighting spring is skeletally anchored using the tomas[®]-pin SD. A dentoalveolar anchorage is also possible.
- Medical indications / contraindications**
4.1 Indication
▪ Uprighting of teeth, particularly molars.
4.2 Contraindications
▪ Teeth that cannot be uprighted due to the condition of crowns or the periodontium.
▪ Refer also to safety instructions.
- Application**
For uprighting lower molars, it is advisable to place the tomas[®]-pin SD between the first and second lower premolars. The pin head must be in close proximity to the attached gingiva. The tooth being uprighted can be fitted with a buccal tube for the adhesive technique or fitted conventionally with a band and welded buccal tube. A buccal tube may detach from the tooth, depending on the structure of the tooth, the fit of the buccal tube for the adhesive technique, the strength of the adhesive bond and the power of the transferring forces.
Step 1: Fitting the tomas[®]-uprighting spring
The effectiveness of the nickel titanium spring depends on its position in the clamping tube. If the nickel titanium spring is pushed as far as possible into the tube (refer to Fig. 1, distance C₀), the alpha and beta moments reach their maximum. The same applies to the intruding and extruding forces. Pulling the nickel titanium spring further out of the tube lowers the alpha and beta moments as well as the intruding and extruding forces (C < C₀). Once the desired position has been set, the tube is crimped with pliers to fix the nickel titanium spring in place. Slide the nickel titanium spring up to its first curve into the slot of the buccal tube. Hold the tomas[®]-uprighting spring parallel to the occlusal plane. Drop a perpendicular from the vertical slot in the tomas[®]-pin SD to the steel wire and mark the wire accordingly (Fig. 1). Depending on the objective of the treatment (refer to step 2), the steel wire is then curved at the mark and shortened as required.
Step 2: Setting the desired treatment objective
The tomas[®]-uprighting spring can be used for uprighting teeth and intruding, extruding or distalizing them at the same time. Teeth can also be uprighted by applying mesial root torque only. Distance A₀ is the horizontal distance between the vertical slot in the tomas[®]-pin SD and the mesial edge of the buccal tube. Distance B₀ is the vertical distance between the steel wire and the occlusal edge of the slot in the pin. Angle α is located occlusodistally from the intersection between the perpendicular of the pin slot and the steel wire. In the initial condition α = 90° (Fig. 1). The individual values can be varied to produce different results.
▪ **Uprighting**
When uprighting teeth, the above mentioned values should remain:

α = 90°	B = B ₀	A = A ₀	(Fig. 2)
---------	--------------------	--------------------	----------

▪ **Uprighting, with intrusion**
When uprighting and intruding teeth, the following parameters should be modified:

Version 1	α > 90° (max. 135°)	B = B ₀	A = A ₀	(Fig. 3a)
Version 2	α = 90°	B < B ₀	A = A ₀	(Fig. 3b)

▪ **Uprighting, with extrusion**
When uprighting and extruding teeth, the following parameters should be modified:

Version 1	α < 90° (min. 45°)	B = B ₀	A = A ₀	(Fig. 4a)
Version 2	α = 90°	B > B ₀	A = A ₀	(Fig. 4b)

▪ **Uprighting, with mesial root torque**
When uprighting teeth, the original values should remain:

α = 90°	B = B ₀	A = A ₀	(Fig. 5)
---------	--------------------	--------------------	----------

If an additional ligature is attached, the tooth can be uprighted with mesial root torque.
Step 3: Fitting the tomas[®]-uprighting spring
The nickel titanium spring must be fixed securely in the clamping tube. To facilitate care by the patient and depending on the setting, it is advisable to remove the nickel titanium wire overlapping the steel wire. Slide the nickel titanium spring into the buccal tube. Place the angled steel wire in the slot of the tomas[®]-pin SD and secure it in position with light-curing adhesive. Reactivation during treatment is usually not required.
Step 4: Removing the tomas[®]-uprighting spring
As soon as the desired treatment objective has been attained, the tomas[®]-uprighting spring is removed by carefully squeezing the adhesive on the head of the tomas[®]-pin SD with Weingart forceps (REF 003-120-00) until the adhesive releases. The steel wire is then removed from the slot and the nickel titanium spring pulled out of the slot in the buccal tube. If the tomas[®]-pin SD is not needed for further retention purposes, it can also be removed. Please refer to the instructions for use of the tomas[®]-pin SD (REF 989-534-00).

6. Composition

Information about the composition is included in the material list, see catalog or www.dentaaurum.com.

7. Safety instructions

The product should not be used if there is a known allergic reaction to one or more of the material components.

8. Information for single use products

The tomas[®]-uprighting spring is for single use only. Reconditioning of tomas[®]-uprighting springs that have been used previously (recycling) or reuse on patients is not permitted.

9. Quality information

Dentaaurum ensures flawless product quality. These recommendations are based upon Dentaaurum's own experiences. The user is personally responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence on the processing on site.

10. Explanation of symbols used on the label

Refer to the label. Additional information can be found at www.dentaaurum.com (Explanation of the symbols used on the label REF 989-313-00).

Subject to modifications