DE

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Abbreviated instructions for use

Dear customer,

Thank you for choosing a quality product from Dentaurm Implants.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit.

Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaurm.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurm Implants GmbH | Turnstr. 31 | 75228 Ispringen | Germany

2. General device description, intended purpose and characteristics, technical description

tiologic® ST implant types are designed for insertion in the endosteal region of the maxilla or mandible. Depending on the indication, appropriate transgingival abutments are secured on the implants and fitted with a prosthetic superstructure.

The tiologic® ST implant system contains specially coordinated instruments, abutments and accessories for placement of the implants and fabrication of the prosthetic restoration. Only original components of the tiologic® ST implant system should be combined in accordance with the Instructions for use / manuals.

3. Area of application

tiologic® ST implant types can be used in both the mandible and maxilla for surgical immediate implantation, delayed immediate implantation and delayed implantation using either the one-stage or two-stage technique.

4. Medical indications / contraindications

Indications

Indications for implant placement are small and large bounded saddles (one-tooth restorations, increasing the number of abutments), a shortened dentition and an edentulous jaw in the maxilla or mandible. Before deciding to use an implant, the risks and advantages / disadvantages of dental implants and alternative methods of treatment should be considered.

With all implantological cases, the implant diameter and length of the tiologic® ST implant types should be in proportion to the prosthetic restoration. If practical with the individual oral situation, implants with a minimum diameter of 4.2 mm should generally be used for restorations that subject the implant and superstructure to high mechanical loading

The tiologic® ST implant types S ø 3.3 mm are available for patients with narrow alveolar ridges. Due to the small diameter and low load capacity (compared to the tiologic® ST M ø 4.2 mm implants), these implants have a limited range of indications. In fully edentulous cases, four or more tiologic® ST implants type S ø 3.3 mm must be inserted with a splinted bar restoration. In partially edentulous cases, implant supported restorations must be combined with tiologic® ST ø 4.2 mm, ø 4.8 mm or ø 5.5 mm implants and a splinted fixed prosthetic restoration. In single restorations, tiologic® ST S ø 3.3 mm implants should only be used for the lower incisors or the upper lateral incisors and only with a minimum 11.0 mm implant length. Single restorations on tiologic® ST ø 3.7 mm, ø 4.2 mm, ø 4.8 mm or ø 5.5 mm implant types require a minimum 9.0 mm implant length. Care should be taken to avoid an excessive mechanical loading when using ball head abutments together with ø 3.3 mm implants.

Contraindications

General contraindications:

General contraindications for dental surgery procedures apply. These include: reduced immunodeficiency / steroid treatment / blood coagulation disorders / uncontrolled endocrine diseases / rheumatic disorders / bone system diseases / cirrhosis of the liver / drug, alcohol or tobacco abuse / depression, psychopathic disorders / poor patient compliance / chronic underlying inflammatory diseases / incomplete physical growth of patient.

Local / personal contraindications:

Osteomyelitis / radiotherapy in the head region / recurring mucosal diseases / temporomandibular joint dysfunctions / parafunctions / lack of vertical or horizontal bone availability, jaw defects, inadequate bone quality / poor oral hygiene.

5. Information on use, processing, initial implementation and handling

The tiologic® ST product range is supplied exclusively to doctors, dentists and dental technicians. It should only be used by doctors, dentists or dental technicians who are familiar with dental implantology, including diagnosis, preoperative planning, surgical techniques and prosthetic restorations.

Before use, dental professionals should ensure that they have carefully read and understood the full tiologic® ST Instructions for use / manuals. This quick reference guide does not replace the Instructions for use / manuals for the tiologic® ST implant system. As the instructions and manuals cannot provide all information for immediate use, we strongly recommend that, before using the system, implantologists attend a training course offered by Dentaurm to learn the correct techniques for the tiologic® ST implants.

- Please refer to the **Product Catalog** and the **Surgery Manuals** for information on precautions and the selection of components for the clinical procedure.
- Refer to the **Product Catalog** and the **Prosthetic Manuals** for information on precautions and the selection of components for the prosthetic procedure.

Before using this product, the implantologist must examine the patient thoroughly and explain the procedure in detail. Dentaurm recommends full clinical, radiological, photographic and statistical documentation.

The tiologic® ST implant system components can be documented, e.g. in the patient file or PatientPass, using the additional labels.

The implantologist should ensure the products cannot be aspirated during intra-oral use.

6. Composition

■ Titanium grade 4 (pure titanium) DIN EN ISO 5832-2

Chemical composition (% by mass)	O	0.4 % max.
	Fe	0.5 % max.
	C	0.08 % max.
	N	0.05 % max.
	H	0.012 % max.
	Ti	Residue
Physical and mechanical properties	0.2 % yield strength	520 MPa min.
	Tensile strength	680 MPa min.
	Elongation at rupture	10 % min.

■ Titanium grade 5 (titanium alloy) DIN EN ISO 5832-3

Chemical composition (% by mass)	Al	5.5 % – 6.75 %
	V	3.5 % – 4.5 %
	Fe	0.3 % max.
	C	0.08 % max.
	N	0.05 % max.
	H	0.015 % max.
Physical and mechanical properties	O	0.2 % max.
	Ti	Residue
	0.2 % yield strength	780 MPa
	Tensile strength	860 MPa
	Elongation at rupture	10 % min.

The SSCP is available at https://ec.europa.eu/tools/eudamed.

7. Adverse reactions

tiologic® ST implant types are not known to have any side effects or produce any interaction.

It cannot, however, be ruled out that in rare cases allergies to components used in the materials of the tiologic®ST implant system may occur or that there may be electrochemically-induced discomfort.

8. Safety information

■ The product should not be used if there is a known allergic reaction to one or more of the material components.

■ Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. All articles used intra-orally should therefore be secured against swallowing and/or aspiration.

■ All serious incidents arising from the use of the product should be reported to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.

■ Only types of alloy in the oral cavity can lead to galvanic reactions.

■ Only use tiologic® ST prosthetic components in combination with tiologic® ST implants.

9. Warnings

Though placement of dental implants has a high rate of success and implants have a long durability, successful treatment cannot be guaranteed. The implantologist should note and document any problematic cases and inform the manufacturer.

An inadequate number of implants, implants with insufficient length or diameter, unfavorable positioning of the implants or a statically-poor prosthetic restoration can cause fatigue fracture in implants, abutments, prosthetic restorations and screws under bio-mechanical loading. Placement of the implants and fabrication of the prosthetic restoration should take into account the individual oral situation of a patient to avoid overloading the components.

Using tiologic® ST implant system components in combinations contrary to the Instructions for use / manuals can also cause mechanical failure, damage to the tissue or unsatisfactory aesthetic results.

10. Storage and shelf life

tiologic® ST implants should always be stored in a dry, dark place at room temperature in the original packaging. The implant should not be used after the expiry date (see label).

11. Notes on sterile packaging

All tiologic® ST implants are supplied sterile with the appropriate closure screws in gamma-sterilized double packaging (foil and blister packaging).

The sterile packaging should not be removed until shortly before insertion.

Dentaurm guarantees the sterility of the tiologic® ST implants until the expiry date on the packaging, provided that the original packaging is undamaged. No guarantee can be given for the sterility after the expiration date and the tiologic® ST implants may no longer be used on patients after this date. If the sterile packaging is damaged, the product may not be used or re-sterilized. **Packaging that has been opened may not be re-sterilized.** A tiologic® ST implant that has been removed from the packaging but has not been inserted, may also no longer be used or re-sterilized.

All other components of the tiologic® ST implant system are supplied non-sterile.

12. Information for single use products

tiologic® ST implants are designed for single use only. Reconditioning of a tiologic® ST implant that has already been used (recycling) and its reuse on a patient are not permitted.

13. Information for multiple use products

Rotating instruments and many accessories are intended for multiple use. The product must be disinfected, cleaned and sterilized before being used again. This applies in particular to components that can be disassembled, e.g. the torque wrench. Rotary instruments – used with proper care and provided that they are not damaged or contaminated – can be reused in dense bone 15 to 20 times; any further reuse or the use of damaged and / or contaminated instruments must be avoided. The dental professional bears responsibility. Instruments and accessories must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels. No liability is accepted if these instructions are disregarded. Additional information can be found at www.dentaurm.com (Processing Instructions Instruments and Accessories REF 989-801-09).

14. Information on disinfection, cleaning and sterilizing

Please refer to the Surgery Manual tiologic® ST or www.dentaurm.com for information on disinfection.

15. Notes on compliance / Conformity declaration

Development, clinical testing, production and quality control of the tiologic® ST product range are completed in accordance with the Medical Device Directive 93/42/EEC. Sections 9 and 10 of our General Terms of Delivery and Payment apply with regard to warranty or liability – unless otherwise provided in this quick reference guide.

Warranty and liability are rendered void in particular if the products are not used by the dental professional or a third party in accordance with the Instructions for use; this also applies if the tiologic®ST product range is used in combination with products from other manufacturers which have not been specifically recommended for use by Dentaurm Implants. Dentaurm Implants has no control over processing and application of the product. These are the sole responsibility of the dental professional.

Technical advice (oral and written) is based on the scientific and technical knowledge available at the time the product is put on the market. It does not release the dental professional from the responsibility of personally checking the suitability of the products for the intended indications and applications. Advice is only given as a non-binding recommendation, which cannot be assumed to provide any form of assurance or guarantee.

All products are subject to continuous development based on current scientific knowledge and we reserve the right to make changes in the construction, design or material of the products.

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