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**5. Application**

All components of the TIOMESH program are supplied non-sterile and must be sterilized before use (e.g. steam sterilization: **121 °C / 250 °F, minimum holding time 15 minutes**).

**6. Availability, safety measures, guarantee**

The components of the TIOMESH program are supplied only to doctors or on their behalf. To avoid damage to anatomical structures, ensure the nutrition of soft tissue and safe anchoring of the titanium grids, they should be used only by qualified doctors who have substantial experience in bone and tissue surgery to ensure conservation of anatomical structures, nutrition of the soft tissue flap and secure fixation of the titanium mesh and microplate. Before using the TIOMESH system the user is advised to attend a training course on Dentaurum Implants since the instructions do not fully cover all possibilities for immediate application.

When using the TIOMESH product range make sure you only use the original components and follow the instructions at all times. Observe these instructions at all times. Any complications which may occur during treatment are attributable primarily to the procedures used and not to the product itself. Before treatment every patient has to undergo a thorough examination and must be fully informed.

**7. Legend**

-  Batch No. (CH.-B)
-  For single use only
-  See instructions for use
- REF** Order No.

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Dentaurum Implants GmbH | Turnstr. 31 | 75228 Ispringen | Germany

**2. Description**

The components of the TIOMESH bone augmentation system are for the regeneration of oral bone defects in the jaw area. The grid plates of pure titanium have a high degree of biological compatibility and mechanical stability, essential preconditions for bone regeneration.

**3. Indications**

The TIOMESH bone augmentation system is suitable for:

- ▼ reconstruction of defects in the alveolar ridge
- ▼ augmentation of circular and regional bone deficits
- ▼ reconstruction or augmentation in the pre-implant or post-traumatic phases
- ▼ treatment of bone defects in cases of periimplantitis
- ▼ simultaneous bone augmentation upon insertion of implants
- ▼ fixation of bone blocks

**4. Contraindications**

- ▼ General contraindications/acute illnesses: reduced immune defence, steroid treatment, disturbed blood coagulation, uncontrolled endocrine diseases, rheumatic illnesses, diseases of the osseous system, hepatic cirrhosis.
- ▼ Local contraindications/personal contraindications: Osteomyelitis, radiotherapy of the head region, recurrent diseases of the oral mucosa, inadequate oral hygiene, TMJ disorders, parafonctions, inadequate bone available vertically/horizontally, jaw defects, depressions, psychopathies.

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Dentaurum Implants GmbH | Turnstr. 31 | 75228 Ispringen | Germania