tomas® PI Gaumenimplantat

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tomas® PI palatal implant

Dear customer,

Thank you for choosing a quality product from Dentaurum.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaurum.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurum GmbH & Co. KG · Turnstr. 31 · 75228 Ispringen · Germany

2. General description

A temporary, reliable anchorage option is created for orthodontic treatment with the aid of tomas® PI (palatal implant). The tomas® PI has a self-tapping thread (Fig. 1) with a blasted, etched surface.

The tomas® PI-abutment multiple set (REF 303-400-00) serves as a closure cap for the implant during the healing phase. It can be found in the sterile packaging (Fig. 3 and 6). This abutment, alongside other tomas® PI-abutments, can be used later to connect the tomas® PI to the orthodontic appliance (Fig. 9).

The tomas® PI is made of the titanium alloy TiAl6V4.

3. Indications and contra-indications

The tomas® PI is used to provide temporary orthodontic anchorage and can be used for example for the following treatments:

- indirect anchorage to secure tooth positions.
- distalization, mesialization and intrusion of teeth.

3.1 General contraindications

The tomas® PI should not be used in the case of immunodeficiency, steroid therapy, blood coagulation disorders, uncontrolled endocrine disorders, rheumatic ailments, bone system diseases, cirrhosis of the liver or any other acute illnesses.

3.2 Specific contraindications

The tomas® PI should not be used in the case of osteomyelitis, cranial radiotherapy, recurring oral mucosa diseases or inadequate oral hygiene. Use is also contraindicated if there is no or insufficient bone available in the planned insertion site.

4. Precautions

The tomas® PI may only be inserted by orthodontists, dentists, oral surgeons or maxillofacial surgeons. Before use, operators should ensure that they have carefully read and observed the Instructions for use supplied.

Only original tomas® PI components may be used according to the Instructions for use.

All patients must be examined thoroughly and informed adequately before treatment with the tomas® PI.

5. Form of presentation

The sterile tomas® PI is protected by an acrylic sleeve in the blister (Fig. 3, ❶). The closure cap (tomas® PI-abutment multiple set; REF 303-400-00, Fig. 3, ❷) for the healing phase is in the small pocket of the packaging. Dentaurum guarantees sterility until the expiration date on the packaging, provided that the original packaging is

undamaged. Once the time limit has been exceeded, the sterility can no longer be guaranteed and the tomas® PI may no longer be used on a patient.

If the sterile packaging is damaged, do not re-sterilize. Packaging that has been opened may not be re-sterilized. A tomas® PI that has been removed from the packaging but has not been inserted, may also no longer be re-sterilized.

All instruments used for insertion (Fig. 2) are delivered non-sterile and must be cleaned, disinfected and sterilized before use on a patient. Please refer to the Instructions for the preparation of pliers and instruments (REF 989-590-20).

6. Insertion of tomas® PI

6.1 General information

Adhere to all essential hygienic measures for invasive surgery during insertion of the tomas® PI. Insertion is possible in the anterior palate, median or paramedian. To function correctly, the tomas® PI must sit firmly in the bone (primary stability).

6.2 Bone preparation

- Local anaesthesia.
- Removal of mucosa (Fig. 4, ❶).
- Measurement of mucosa thickness.
- A contra angle (500–800 Umin⁻¹) should be used intermittently for all rotating instruments. A drill extension (387-509-00) is available if the anatomical space is unfavorable. It is necessary to ensure adequate cooling using a sterile physiological saline solution. The markings on the drills (Fig. 4) or the depth gauge (303-240-54) help monitor the drill depth.

- Drillpoint marking (Fig. 4, ❷).
- Pilot drilling ø 2 mm (Fig. 4, ❸).
- Pilot drilling ø 3.4 mm (Fig. 4, ❹).

6.3 Insertion procedure

- Withdraw triangular sleeve (Fig. 5, ❶) from blister (Fig. 3).
- Remove yellow silicone plug (Fig. 5, ❷).
- **Note: palatal implant is loose in the sleeve. Do not tilt!**
- Remove the tomas® PI (Fig. 5, ❸).
- Screw the tomas® PI into the bone until the thread is fully inserted into the bone (Fig. 4, ❺). This can be done either with a hand-held instrument (max. 10 Umin⁻¹; Fig. 5, ❹) or manually (Fig. 5, ❺). Torque should not exceed 40 Ncm during insertion.
- Remove the closure cap from the grey silicone plug (Fig. 6, ❶) and screw this into the implant manually (Fig. 6, ❷). Torque should not exceed 15 Ncm.
- Use a counter ratchet to tighten or loosen the abutment screw (Fig. 7).

7. Fabrication of orthodontic appliance

- Remove the closure cap (Fig. 6, ❷).
- Screw the impression cap into the palatal implant (Fig. 8). The impression can be taken in analog (only using silicone or polyether) or digital form.
- Fabrication of model with the laboratory analog (Fig. 9).
- Wires, for example, can be welded onto the abutments (Fig. 9).
- Insert the appliance into the mouth and fix on the palatal implant using the abutment screw. Torque should not exceed 15 Ncm.

8. Removal of the tomas® PI

- Local anaesthesia.
- Loosen the abutment screw on the tomas® PI and remove the appliance.
- Alternative A (Fig. 10a): activate the blocking function of the tomas® PI-torque ratchet (set to ∞). Turn the implant to a maximum of 90° to the right using the

tomas® PI-driver 2.5 and ratchet (❶), then turn again in the opposite direction (❷) and remove.

- Alternative B (Fig. 10b): screw the tomas® PI-guide sleeve (303-500-10) into the palatal implant. Push the trephine drill (tomas® PI-trepan, 303-500-20) over the guide sleeve and mill the mucosa down to bone level. Rinse away soft tissue residues and lower the trephine drill until the top edge of the guide sleeve is in line with the second marking (10 mm). Use the contra angle (500–800 Umin⁻¹) intermittently. Ensure that adequate cooling is provided using a sterile physiological saline solution.

- Remove the palatal implant with the insertion instrument.

9. Composition

The material composition is included in the materials list, refer to the catalog or www.dentaurum.de.

10. Storage and shelf life

The tomas® PI should always be stored in a dry, dark place at room temperature in the original packaging. The palatal implant should not be used after the expiry date (see label).

11. Information for single use products

The tomas® PI is intended for single use only. Reconditioning of palatal implants that have been inserted previously (recycling) or reuse on patients is not permitted.

12. Further information

Should the user and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the user and/or the patient is resident is informed accordingly.

13. Quality information

Dentaurum ensures faultless quality of its products. These recommendations are based upon Dentaurum's own experiences. The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as Dentaurum has no influence on the processing on site.

14. Explanation of symbols

 Refer to the label. Additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00).

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