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Dentaurum GmbH & Co. KG · Turnstr. 31 · 75228 Ispringen · Germany

2. General description

A skeletal anchorage option is created for orthodontic treatment with the aid of the tomas® PI (palatal implant). The tomas® PI features a self-tapping thread with a blasted, etched surface.

The tomas® PI abutment serves as a closure cap for the implant during the healing stage and is included in the sterile packaging (Figs. 1). The tomas® PI abutments can be used later to connect the tomas® PI to the orthodontic appliance.

2.1 Materials used:

tomas® PI palatal implant	Titanium alloy according to ISO 5832-3
tomas® PI abutments	Stainless steel 1.4301
tomas® PI abutments with wire element	Stainless steel 1.4404 and 1.4301, and cobalt-nickel-chromium alloy according to ISO 5832-8

The material composition is included in the orthodontics materials list; please refer to the catalog or www.dentaurum.com. The SSCP is available at <https://ec.europa.eu/tools/eudamed> and www.dentaurum.com.

3. Intended purpose

The implantable products of the tomas® PI system have the following intended purpose:

- Orthodontic endosteal anchorage elements are inserted into the jawbone to serve as a temporary anchor.
- Abutments serve as a connection piece between an orthodontic endosteal anchorage element and orthodontic coupling element.
- The connection screws serve as a connection piece between the orthodontic endosteal anchorage element and the abutment.

3.1 Indications and contraindications

The tomas® PI serves as a temporary orthodontic anchor for the skeletal support of orthodontic appliances. It can also be used as part of pre-prosthetic treatment.

3.2 Intra-oral contraindications:

- Unfavorable anatomical jaw positions (e.g., reduced bone availability), inadequate bone quality, pathological findings in bone (e.g., cysts),
- Pathological and/or recurring findings in the mucosa (e.g., lichen),
- Inadequate compliance (e.g., oral hygiene, psychoses, drug, alcohol or tobacco abuse),
- Long-standing therapy-resistant functional disorders,
- Untreated periodontitis.

3.3 Temporary contraindications:

- Temporary medication intake (e.g., anticoagulants, immune suppressants, steroids), acute inflammatory diseases (e.g., maxillary sinusitis), pregnancy.

3.4 General medical contraindications

- Reduced immune defense, poorly managed diabetes, treatment with antiresorptive substances which may cause medication-related osteonecrosis of the jaw (MRONJ),
- Diseases relating to the bone structure,
- Radiotherapy in the head region, active tumors,
- Severe internal diseases (e.g., cirrhosis of the liver),
- Blood coagulation disorders,
- Rheumatic disorders,
- Titanium or nickel allergy.

Adverse reactions and residual risks

- Pain,
- Swelling,
- Excessive bleeding,
- Wound healing disorders,
- Sinus opening,
- Broken drill tip,
- Sensitization to one of the materials,
- Limitation of speech and chewing function,
- Soft tissue or bone infections,
- Nerve damage (paresthesia, dysesthesia, hypo- or hyperesthesia),
- Implant fracture or loss,
- Swallowing or aspiration.

4. Precautions

The tomas® PI should only be inserted by orthodontists, dentists, oral surgeons or maxillo-facial surgeons. Before use, dental professionals should ensure that they have carefully read and observed the Instructions for use supplied.

We recommend dental professionals complete a relevant training course prior to using the tomas® PI.

When using the tomas® PI system, only use the appropriate original components.

All patients must be examined thoroughly and informed adequately of the procedure, risks and possible complications before treatment with the tomas® PI!

5. Form of presentation

The sterile tomas® PI is protected by an acrylic sleeve in the blister. The closure cap (tomas® PI abutment) for the healing stage is in the small pocket of the packaging. Dentaurum guarantees sterility until the expiry date on the packaging, provided that the original packaging is undamaged. Once the time limit has been exceeded, sterility can no longer be guaranteed and the tomas® PI must no longer be used on a patient.

If the sterile packaging is damaged, do not re-sterilize. Packaging that has been opened may not be re-sterilized! A tomas® PI that has been removed from the packaging but has not been inserted, may also no longer be used or re-sterilized.

6. Insertion of the tomas® PI

6.1 General information

Adhere to all essential hygienic measures for invasive surgery during insertion of the tomas® PI. To function correctly, the tomas® PI must sit firmly in the bone (primary stability).

6.2 Procedure (Fig. 2)

- Rinse with disinfecting mouth rinse.
- Local anesthesia.
- Removal of the mucosa with the tomas® gingiva-punch ❶.
- Measurement of mucosal thickness.
- A contra angle (500 rpm – 800 rpm) should be used for all rotating instruments.
- It is necessary to ensure adequate cooling using a sterile physiological saline solution. The markings on the drills or the depth gauge (303-240-54) help monitor the drill depth.
- Drillpoint marking ❷.
- Pilot drilling Ø 2 mm ❸.
- Final drilling Ø 3.4 mm.


For insertion of the implant, proceed as follows (Fig. 3):

- Withdraw the triangular sleeve from blister ❶.
- Remove the plug ❷.
- Withdraw the tomas® PI with the insertion instrument ❸ ❹ ❺.

Screw the tomas® PI into the bone (Figs. 2 (last image) and 6) until the thread section (Figs. 3 and 5) is fully inserted into the bone. This can be done using a handpiece (max. 25 rpm) or manually. Torque of 15 to 25 Ncm is recommended during insertion.

Remove the abutment from the gray silicone plug and screw into the implant manually. Torque should not exceed 15 Ncm.

To tighten and loosen the abutment screw, the tomas® PI should be fixated with the tomas® PI-lock key. It is recommended to allow the tomas® PI to heal for at least 12 weeks after insertion.

7. Production of the orthodontic appliance

- Remove the closure cap.
- Screw the tomas® PI-impresion post into the tomas® PI.
- OR
- Attach the tomas® PI-transfer cap to the tomas® PI.
- Silicone or polyether are recommended for analog impression taking.
- Production of the model with the tomas® PI-laboratory analog.
- Wires can be welded onto the abutments.

Insert the appliance into the mouth and fixate on the tomas® PI using the abutment screw. Torque should not exceed 15 Ncm.


8. Removal of the tomas® PI

General information:

If the tomas® PI is no longer required for anchorage, explantation can be performed non-surgically (see Version A) or surgically (see Version B).


Procedure

- Local anesthesia.
- Loosen the abutment screw on the tomas® PI and remove the appliance.

Version A (Fig. 4):

Activate the blocking function of the tomas® PI-torque ratchet (set to ∞). Using the tomas® PI-driver 2.5 and tomas® PI-torque ratchet, turn the implant a quarter turn to the right and then unscrew in the opposite direction.

Version B (Fig. 5):

Screw the tomas® PI-guide sleeve into the tomas® PI implant. Slide the trephine drill (303-500-21) over the guide sleeve and lower until the top edge of the guide sleeve is in line with the second marking (10 mm). Use the contra angle (500 rpm – 800 rpm) intermittently. It is necessary to ensure adequate cooling using a sterile physiological saline solution. Remove the palatal implant with the insertion instrument.

9. Storage and shelf life

The tomas® PI should always be stored in a dry, dark place at room temperature in the original packaging. The tomas® PI should not be used after the expiry date (see label).

10. Product for single use

The tomas® PI is intended for single use only. Reconditioning of palatal implants that have been inserted previously or reuse on patients is not permitted.

11. Further information

- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. Therefore, secure all articles used intra-orally against swallowing and/or aspiration.
- Use of third-party components can negatively impact the function and safety. In such cases, Dentaureum GmbH & Co. KG accepts no liability for replacement products or guarantees.
- Instruments must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.
- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. Therefore, secure all articles used intra-orally against swallowing and/or aspiration.
- Use of third-party components can negatively impact the function and safety. In such cases, Dentaureum GmbH & Co. KG accepts no liability for replacement products or guarantees.
- Instruments must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.
- Should the dental professional and/or the patient become aware of serious problems arising from the use of the product, inform the manufacturer and the competent authority in the country in which the dental professional and/or the patient is resident.

12. Disposal

When disposing of the product, please adhere to the national regulations.

13. Quality information

Dentaureum ensures faultless quality of its products. The dental professional is solely responsible for the further processing of the products. Liability for failures cannot be taken, as Dentaureum has no influence on the further processing on site.

14. Explanation of symbols

Refer to the label. Additional information can be found on the internet at www.dentaureum.com (Explanation of symbols used on the label REF 989-313-00).

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