DE

rematitan® blank Ti5

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* Variations of ± 10% are possible depending on the batch.

Dear customer,

Thank you for choosing a quality product from Dentaurum. It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. Instructions for use cannot describe every eventuality and possible application. In case of questions or suggestions, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaurum.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurum GmbH & Co. KG | Turnstraße 31 | 75228 Ispringen | Germany

2. General product description

Non-precious material for crown and bridge technique.

3. Area of application

rematitan® milling blanks for subtractive manufacturing are dental materials made of pure titanium or based on titanium. rematitan® blank Ti5 is particularly suitable for the manufacture of constructions in the field of the crown and bridge technique, produced by milling, such as crowns, bridges, bar constructions, abutments in the form of fixed and removable restorations, and for metal ceramics.

4. Heat treatment (optional)

If the restoration produced with rematitan® blank Ti5 is to be given a ceramic veneering, we recommend the blank undergoes heat treatment before it is machined. This optimizes the material structure of the titanium alloy and increases its stability during the later process of veneering. Heat treatment is carried out (under ambient air) at a furnace temperature which corresponds to the bonder fire temperature (795 °C / 1463 °F when using ceraMotion® Ti). For this purpose, rematitan® blank Ti5 is placed in the preheated furnace, annealed for a period of 4 hours and then removed from the furnace to cool in ambient air.

5. Digital modeling

The dental products are modeled using a CAD software which has been approved for this purpose. Construction must be in line with dental regulations in order to fulfil mechanical requirements.

6. Milling

Machining can be carried out with most commercially available CAD/CAM milling machines or hand-guided copy milling machines which are designed for machining titanium. Suitable tools should be selected according to the machine manufacturer’s operating instructions. Before starting the milling process, ensure that the tools are clamped correctly and securely. Machining should be carried out with the milling parameters recommended by the machine manufacturer. The coolant supply depends on the milling strategy used and on the milling parameters used.

7. Post treatment

After completion of the milling operation, the framework must be removed from the blank. Support residues can be removed using finely toothed hard metal burs. Any particles or dust remaining on the surface must be removed with compressed air. It is important to observe the warnings and precautions to be taken.

8. Finishing

The finishing of the framework and the surfaces should be carried out by means of cross-toothed hard metal burs suitable for titanium. Work with low contact pressure to remove material at lower speeds. In this case too, it is important to observe the warnings and precautions to be taken.

9. Welding

Suitable welding processes such as laser welding produce mechanically strong and corrosion-resistant joints. It is important to observe the framework design, the surface structures, the welding sequence and the welding parameters recommended for each individual laser machine. Suitable welding wires are available as filler material for all titanium materials, e.g. rematitan® wire on coils 0.4 mm (REF 528-039-50) or 0.7 mm (REF 528-040-50) and rematitan® straight wire 1.0 mm (REF 528-041-00) or 1.2 mm (REF 528-042-00). It is important that the welding point is well covered with a shielding gas (argon).

10. Ceramic veneering

Titanium materials for the bonding technique can be veneered with all suitable ceramic materials, e.g. ceraMotion® Ti (Dentaurum), as long as the CTE is suitable for the material. In order to increase the process stability in the ceramic veneer, heat treatment of the unprocessed rematitan® blank Ti5 as described above is recommended. The highest recommended firing temperature is 880 °C / 1616 °F. The ceramic should cool down over a long period in the furnace unless the ceramic manufacturer has other specifications.

11. Contraindications and adverse reactions

Signs of intolerance to titanium and titanium alloys are extremely rare if the manufacturer’s Instructions for use are adhered to. If the patient has a proven allergy against any component within the material, it must not be used for safety reasons. There have been individual reportings of local irritations which were electrochemically induced. If various alloys have been used, it is possible that galvanic effects may occur. Exposure to dust and vapors may cause irritation of the eyes and/or respiratory tract. Please adhere to warnings and safety data sheets.

12. Warnings and precautions

Metal shavings and metal dust may be produced as the alloy is mechanically processed. The mechanical process will also create heat and possibly burs. Processed material may therefore be hot and/or have sharp edges. Due to the oxygen affinity of titanium, there is a

risk of burns and deflagration, particularly when the material is hot. Any type of ignition source should be kept at a distance. Metal fire powder extinguishers should always be positioned in close proximity to extinguish a fire if necessary. They should be in working order. The workpiece must be cooled during mechanical processing in the presence of oxygen. Processing at temperatures above the solidus temperature may generate vapors. The workplace should be fitted with a suitable suction or ventilation system. We generally recommend that gloves, protective clothing, safety glasses and respiratory protection equipment are worn. Should there be contact with the skin, please wash the affected area with soap and water for a few minutes. Should the product come into contact with the eye, please rinse the eye with a gentle stream of water or saline solution for at least 15 minutes.

Please observe the relevant safety data sheet (see www.dentaurum.com).

Interproximal or occlusal contact with other metallic elements and the use of different types of alloys in the same oral cavity should be avoided to prevent galvanic reactions or local irritations in individual cases. The material should not be ground or polished in the patient’s mouth. There is no scientific evidence on the safety or efficacy of treatment of pregnant women or nursing mothers or children.

13. Reuse

Fully-machined milling blanks must not be further processed for further use as a medical product, so that clear batch traceability is also possible.

14. Conditions for storing

Normal changes to the surrounding conditions (e.g. temperature, pressure, light) have no negative impact on the product characteristics of titanium and titanium alloys.

15. Disposal

Please adhere strictly to the rules that apply for your region, bearing in mind the details outlined in the safety data sheets.

16. Quality

Dentaurum ensures faultless quality of its products. These recommendations are based upon Dentaurum’s own experiences. The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as we, Dentaurum, have no influence on the processing on site.

17. Further information

Should the user and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the user and/or the patient is resident is informed accordingly.

The SSCP is available at https://ec.europa.eu/tools/eudamed and www.dentaurum.com.

18. General information

These Instructions for use contain the main steps for processing titanium and titanium alloys. You will find these Instructions for use and information on the current status at www.dentaurum.com/ifu.

Additional information can be found at www.dentaurum.com/ifu.

19. Explanation of symbols used on the label

 Please refer to the label. Additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00).

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
➔ For more information on our products, please visit www.dentaurum.com

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Dentaurum GmbH & Co. KG | Turnstraße 31 | 75228 Ispringen | Allemagne