**D**  
DENTALURUM  
1886

tiologic®  
TWINFIT

tiologic®ST

CITO mini®

---

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
## EN Preparation instructions · Implant system components

---

These instructions apply to all Dentaaurum, prosthetic components, instruments, rotary instruments and components that are approved for use on patients.



These instructions comply with the requirements of DIN EN ISO 17664 and RKI (Robert Koch Institute) guidelines.

## General requirements

All instruments and components must be cleaned, disinfected and if necessary, sterilized before each use. This also applies for initial use following delivery. Instruments and components are generally supplied non-sterile (clean and disinfect following removal of the protective transport packaging; sterilize after packing). Effective cleaning and disinfection is an indispensable requirement for effective sterilization. This includes rotary instruments (drills), hereafter referred to as instruments.

During use already ensure that contaminated instruments are either collected separately or replaced to their correct position in the Surgical Tray for tioLogic® TWINFIT or the Surgical Tray ADVANCED for tioLogic®. Clean / disinfect the contaminated instruments. Then replace them in the instrument tray and sterilize the fully loaded surgical tray if necessary. Single-use items may not be reprocessed. It is not possible to guarantee functionality and safe preparation.

## Manual preclean

Surface contamination must be removed from the instruments and components immediately after use, at the latest within one hour.

Disassemble the instruments and components as far as possible. (e.g. torque ratchet, silicone holder in the surgical tray for tioLogic® TWINFIT, in the Surgical Tray ADVANCED for tioLogic® and manual insertion keys). Preparation for manual precleaning, as well as disassembly and assembly of the torque ratchet, see Instructions for use torque ratchet (REF 989-719-00).

Rinse instruments and components under running water and immerse in a disinfectant solution.

Use only a soft brush (nylon brush) or a clean soft cloth intended for this purpose only to remove contamination manually.

Do not use metal brushes or steel wool.

If applicable: rinse all hollow sections of the instruments five times using a disposable syringe (minimum volume 5.0 ml).

Disinfectants should be aldehyde-free (otherwise fixation of blood residue), have a certified effectiveness (e.g. DGHM [German Society for Hygiene and Microbiology] or FDA approval and CE marking), be suitable for disinfecting instruments and compatible with the instruments (see chapter Material resistance).

Do not immerse osteotomes in a NaCl solution (danger of pitting corrosion or stress corrosion cracking).

Note that the disinfectant used for preconditioning is only for personal protection and cannot be regarded as a substitute for subsequent disinfection after cleaning.


## Machine-based cleaning and disinfection

- **When using a washer disinfectant, please make sure that:**

- the efficacy of the washer disinfectant has been certified (e.g. DGHM or FDA approved or CE marking according to DIN EN ISO 15883) and is validated specifically for the appliance and product,
- a certified program for thermal disinfection (minimum 5 mins at 90°C / 194 °F or an  $A_0 > 3000$ ) is used if possible (with chemical disinfection there is the risk of disinfectant residue on the instruments),
- the program used is suitable for the instruments and has an adequate number of rinse cycles,
- the water used for rinsing is sterile or has a low bacteria count (max. 10 bacteria / ml) and is low in endotoxins (max. 0.25 endotoxin units / ml) (e.g. purified water / highly purified water),
- the air used for drying is filtered,
- the disinfectant is regularly serviced and checked,
- the Instructions for use for the washer disinfectant are observed.
- **Important: To guarantee adequate cleaning conform with Robert Koch Institute regulations, the osteotomes must be cleaned separately and not in the osteotome tray. The osteotome tray serves only to hold the osteotome during sterilization and then to store the osteotome. For sterilization, the loaded osteotome tray is packed in disposable packaging dedicated to the sterilization process in accordance with ISO 11607.**

- **When choosing a cleaning agent system, make sure that:**

- it is suitable for cleaning metal and plastic instruments,
- an additional disinfectant with certified efficacy (e.g. DGHM or FDA approved and CE marking) is used – provided that thermal sterilization is not used – and that it is compatible with the cleaning agent used, and
- the chemicals used are compatible with the instruments (see chapter Material resistance).
- **Adhere to the concentrations given by the manufacturer of the cleaning agent and disinfectant.**

### Cleaning / disinfection procedure

- ❶ Disassemble the instruments and components as far as possible (e.g. torque ratchet, silicone holder in the Surgical Tray for tioLogic® TWINFIT, in the Surgical Tray ADVANCED for tioLogic® and manual insertion keys).
- ❷ Preparation of the torque ratchet for machine cleaning/disinfection, as well as disassembly and assembly of the torque ratchet, see Instructions for use torque ratchet (REF 989-719-00).
- ❸ Place the dismantled instruments in a lockable container in the washer disinfector. Instruments and components must not come into contact with one another. Place the insert of the Surgical Tray for tioLogic® TWINFIT and the Surgical Tray ADVANCED for tioLogic® in the upper rack of the washer disinfector.
- ❹ Start the program.
- ❺ Remove the instruments and components or the trays from the washer disinfector (WD) at the end of the program.
- ❻ Check, possibly reassemble and pack the instruments, components (torque ratchet and silicone holder for ratchet) and trays in a clean area as soon as possible after removal (see chapters Monitoring and maintenance, Packaging), if necessary after additional drying.

#### Note:

**If the washer disinfector does not have an automatic drying program, leave the door of the unit slightly open to dry the instruments.**

Proof of basic suitability for effective automatic cleaning and disinfecting was provided by an independent, accredited test laboratory using a G 7836 GD disinfector (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh, Germany) and the cleaning agent neodisher® Medizym (Dr. Weigert GmbH & Co. KG, Hamburg). GmbH & Co., Gütersloh) and the cleaning agent neodisher® Medizym (Dr. Weigert GmbH & Co. KG, Hamburg). The procedure described above was taken into account during the tests.

## Manual cleaning and disinfection.

### Observe the following points when choosing a cleaning agent and disinfectant:

- metal and / or plastic instruments and components must be suitable for cleaning and disinfecting,
- the cleaning agent must be suitable for ultrasonic cleaning (no foaming),
- use only disinfectants with certified effectiveness (e.g. DGHM or FDA approval and CE marking). Disinfectants should be compatible with the cleaning agent used.

Do not use combined cleaning / disinfecting agents.

The concentrations and reaction times given by the cleaning agent and disinfectant manufacturer should be strictly adhered to. Use only freshly mixed solutions, only water that is sterile or has a low bacteria count (max. 10 bacteria / ml) and is low in endotoxins (max. 0.25 endotoxin units / ml) (e.g. purified water / highly purified water), and only filtered air for drying.


**Important: To guarantee adequate cleaning conform with Robert Koch Institute regulations, the osteotomes must be cleaned separately and not in the osteotome tray. The osteotome tray serves only to hold the osteotome during sterilization and then to store the osteotome. For sterilization, the loaded osteotome tray is packed in disposable packaging dedicated to the sterilization process in accordance with ISO 11607.**

## Cleaning procedure

- 1 Disassemble the instruments and components as far as possible (e.g. torque ratchet, silicone holder in the Surgical Tray for tiologic® TWINFIT, in the Surgical Tray ADVANCED for tiologic® and manual insertion keys).
- 2 Preparation of the torque ratchet for manual cleaning/disinfection, as well as disassembly and assembly of the torque ratchet, see Instructions for use torque ratchet (REF 989-719-00).
- 3 Immerse the disassembled instruments and the components fully in the cleaning solution for the recommended reaction time (if required use an ultrasonic cleaner, or brush carefully with a soft brush).  
Instruments and components must not come into contact with one another. Rinse all hollow sections of the instruments using a disposable syringe (minimum volume 5.0 ml) before and after the reaction time. Special holders, e.g. bur blocks, should be used if necessary.
- 4 Remove the instruments and components from the disinfectant solution and rinse thoroughly at least three times with water.
- 5 If applicable: rinse all hollow sections of the instruments five times using a disposable syringe (minimum volume 5.0 ml).
- 6 **Check that the pliers and instruments have been properly cleaned**  
(see chapters Monitoring and maintenance, Packaging).

## Disinfection procedure

- 1 Place the disassembled instruments and components in the disinfectant solution according to the prescribed reaction time. Ensure that the instruments and components are covered and are not in contact with one another. Rinse all cavities in the instruments using a disposable syringe (minimum volume 5.0 ml) before and after the reaction time.
- 2 Remove the instruments and components from the disinfectant solution and rinse thoroughly at least three times with water.
- 3 If applicable: rinse all hollow sections of the instruments five times using a disposable syringe (minimum volume 5.0 ml).
- 4 Do not reassemble and pack the instruments and components until they are dry. Reassemble and pack dry instruments and components as soon as possible (see chapters Monitoring, Maintenance, Packaging).


Proof of basic suitability for effective manual cleaning and disinfecting was provided by an independent, accredited test laboratory using Bodedex® forte cleaning agent and Korsolex® plus disinfectant (Bode Chemie, Hamburg, Germany). The procedure described above was taken into account during the tests.

## Monitoring and maintenance

Check all instruments and components after cleaning or cleaning / disinfection for corrosion, damaged surfaces, chipped areas or contamination and separate damaged instruments and components. Instruments and components that are still contaminated must be cleaned and disinfected again.

Instrument oils should not be used if possible. If oil is to be used, ensure that only instrument oils (white oil) are used, which – depending on the maximum sterilization temperature used – are approved for steam sterilization and are certified as biocompatible.

## Packaging

**Arrange the cleaned and disinfected instruments and components as required in the appropriate sterilization tray. Pack the instruments, components or sterilization trays in disposable sterilization packaging (single or double packaging) and / or pack in sterilization containers. The following requirements must be met:**


- DIN EN ISO/ANSI AAMI ISO 11607-1/2 (formerly: DIN EN 868/ANSI AAMI ISO 11607)
- suitable for steam sterilization (temperature resistant to min. 134 °C / 273 °F and adequate steam permeability)
- adequate protection of the instruments and components and sterilization packaging against mechanical damage
- regularly maintained according to the manufacturer's instructions (sterilization containers)

## Sterilization

Only the following sterilization procedures should be used for sterilization. Other sterilization procedures are not suitable. Flash sterilization or gravitational method should not be used. Other procedures that should not be used are: hot-air sterilization, X-ray sterilization, formaldehyde or ethylene oxide sterilization or plasma sterilization.

### Steam sterilization


- fractional vacuum process (with adequate product drying)
- steam sterilizer in accordance with DIN EN 13060-2004 or DIN EN 285
- validated in accordance with DIN EN ISO / ANSI AAMI ISO 17665 (formerly: DIN EN 554 / ANSI AAMI ISO 11134) (valid commissioning and product-specific performance evaluation)
- maximum sterilization temperature 134 °C / 273 °F (plus tolerance in accordance with DIN EN ISO / ANSI AAMI ISO 17665 (formerly: DIN EN 554 / ANSI AAMI ISO 11134))
- sterilization time (exposure time at the sterilization temperature) minimum 5 mins at 134 °C / 273 °F
- max. pressure: 2.2 bar

Proof of basic suitability for effective steam sterilization was provided by an independent, accredited test laboratory using a EuroSelectomat steam sterilizer (MMM Münchener Medizin Mechanik GmbH, Planegg, Germany) and the fractional vacuum process as well as a Systec V-150 steam sterilizer (Systec GmbH Labor-Systemtechnik, Wettenberg, Germany) and the gravitation process. The procedure described above was taken into account during the tests.

**Note:** Preparation of the torque ratchet for sterilization, see Instructions for use torque ratchet (REF 989-719-00).

## Storage

Following sterilization, instruments and components should be stored dry and dust-free in the sterilization packaging.

## Documentation and approval

The preparation of medical products ends with a documented release for use.

## Material resistance

**When choosing the cleaning agent and disinfectant ensure that they do not contain the following components:**

- organic, mineral or oxidizing acids (maximum permitted pH 9.5, a neutral / enzymatic cleaner is recommended)
- strong alkali
- organic solvents (e.g. alcohols, ethers, ketones, benzines)
- oxidation agents (e.g. hydrogen peroxides)
- halogens (chlorine, iodine, bromine)
- aromatic/ halogenated hydrocarbons
- heavy metal salts

**Do not immerse osteotomes in a NaCl solution (danger of pitting corrosion or stress corrosion cracking).**

Do not clean instruments, components and sterilization trays with metal brushes or steel wool. Do not expose instruments, components and sterilization trays to temperatures above 134 °C!

## Reusability

All surgically invasive rotary instruments – when used with proper care and provided that they are not damaged or contaminated – can be reused in dense bone 15 to 20 times. The dental professional bears responsibility for any further reuse or use of damaged and/or contaminated surgically invasive rotary instruments. Instruments and components must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels. No liability is accepted if these instructions are disregarded.

Remove and, if necessary, dispose of damaged or corroded instruments and components.

**Important: To guarantee adequate cleaning conform with Robert Koch Institute regulations, the osteotomes must be cleaned separately and not in the osteotome tray. The osteotome tray serves only to hold the osteotome during sterilization and then to store the osteotome. For sterilization, the loaded osteotome tray is packed in disposable packaging dedicated to the sterilization process in accordance with ISO 11607.**

## Disposal

Products to be disposed of should be decontaminated observing the relevant regulations for medical waste disposal.

## References

Further tips on the correct conditioning of medical products in the practice and laboratory are provided by the German "Instrument Conditioning Research Group" (Arbeitskreis Instrumenten-Aufbereitung) in its yellow brochure.

[www.a-k-i.org](http://www.a-k-i.org)

[www.rki.de](http://www.rki.de)

Further information on the materials for implantology products can be found in the Materials list (989-801-06).

### Notes on the validations:

**Deviations from the listed procedures must be carefully evaluated for their effectiveness and possible adverse consequences. It is the responsibility of the dental professional performing the reconditioning to ensure that the equipment, materials and staff of the reconditioning facility are adequate to achieve the desired results. This requires validations and routine monitoring of the reconditioning process.**

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Dentaurum GmbH & Co. KG